Borehole reservoirs

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Instructions for use

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Mode d'emploi

Note for U.S. users:

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
Maßstab 1:1

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Indication

The borehole reservoir is used in cases of hydrocephalus to remove cerebrospinal fluid, administer medications and monitor pressure.

Description

One possible treatment in cases of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain to an appropriate region of the body. Our recommendation to use a reservoir in the shunt system is to drain to the abdominal cavity in accordance to the indication of use. This shunt system consists of a series of catheters and a differential pressure valve. In addition, the option of positioning the borehole reservoir in the borehole of the cranium also exists. It then becomes possible to measure intraventricular pressure, inject medication and check the status of the valve. Its solid titanium base is highly puncture-resistant. A puncture of the borehole reservoir or the SPRUNG RESERVOIR should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0,9 mm. 30 times of punctures are able without any restrictions.

With the SPRUNG RESERVOIR it is possible to flush CSF towards the valve due to a nonreturn valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The volume per pump procedure is approx. 0,08 ml. The opening pressure of the shunt system is not increased by the implantation of the SPRUNG RESERVOIR.

Caution: Frequent pumping of the SPRUNG RESERVOIR can lead to an excessive drainage and thus to unphysiological pressure situations. The patient should be informed about the risk.

Nothing may be passed through the borehole reservoir. A puncture should be performed as perpendicular to the reservoir surface as possible with a cannula whose maximum size is G 21. The borehole reservoir is available in various product configurations (individually or as a set) and as a borehole port. The borehole port doesn't have a lateral connector for a distal catheter. A pediatric borehole reservoir and pediatric borehole port adapted to the size of children's skulls are also available.

Implantation

The ventricle catheter is implanted with the aid of a introducing stylet. After being connected to the borehole reservoir the catheter is fixed with a ligature. The reservoirs are intended for use with catheters with an inner diameter of approx. 1,2 mm and an outer diameter of approx. 2,5 mm. When implanting a borehole reservoir without a distal catheter, care must be taken to ensure that the drainage catheter is also fixed in place with the aid of a ligature.

Borehole reservoir

It is recommended that a borehole with a diameter of 10 mm be performed. Owing to the fact that the ventricle catheter practically disappears in the borehole reservoir, the requisite ligature can be made in a very secure fashion (see Fig. 1).


Fig. 1: Borehole reservoir

Pediatric borehole reservoir

Unlike the standard *borehole reservoir*, the *pediatric borehole reservoir* does not have an additional reduction at the bottom of the reservoir. The outer diameter is much smaller.

SPRUNG RESERVOIR

As additional new feature of the *SPRUNG RESERVOIR* CSF can be flush towards the valve due to a nonreturn valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR* (see Fig. 2).


Fig. 2 *SPRUNG RESERVOIR*

MCLANAHAN RESERVOIR

The *MCLANAHAN RESERVOIR* functions similar to the *SPRUNG RESERVOIR*: By pressing down the membrane CSF-flow is produced in distal direction (s. Fig 3).


Fig. 3: *MCLANAHAN RESERVOIR*: Functions similar to the *SPRUNG RESERVOIR*

In addition to the above function, it is possible to inject specific drugs or fluids into either the ventricles or in a distal direction into the shunt system. To inject into the ventricles: A cannula is used to puncture the reservoir.

The cannula is then used to depress the floor of the reservoir. By depressing the floor of the reservoir completely, the distal outlet is closed and the proximal outlet is opened, thus allowing drainage into the ventricles (Fig 4).


Fig. 4: *MCLANAHAN RESERVOIR* after distal outlet is closed and proximal outlet is opened

When a cannula is used without depressing the reservoir floor, the proximal outlet is blocked and the distal outlet is opened. During the transition between the above mentioned conditions, both outlets are open.

Caution: The patency of the reservoir should be tested prior to implantation. In case of uncertainty, distilled water should be aspirated carefully from the distal connector and by immersing the proximal part in sterile liquid.

Reimplantation

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, as a validated decontamination process will compromise the functionality of the valve.

Compatibility with diagnostic procedures

Like the distal catheter (when applied), the *borehole reservoirs* consist of the non-magnetic substances titanium, sapphire and silicone. Consequently, the reservoir and catheters will not interfere with either nuclear magnetic resonance or computer tomographic procedures.

The *borehole reservoirs* are MRI conditional (ASTM-F2503-08). The provided catheters are MRI safe. Deflectors and connectors are MRI conditional.

Safety measures

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflammation and tension around the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation or vomiting often occur in cases of shunt dysfunction. In the event such symptoms occur or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

Sterilization

All products are carefully and thoroughly sterilized. Owing to the fact that the products are packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

Resterilization

The functional safety and reliability of resterilized products cannot be guaranteed.

Functional safety

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

Requirements of the MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

Note on the instructions for use

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

Medical products consultant

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
14469 Potsdam · Germany
Phone: +49(0) 7000 MIETHKE or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap-Platz
78532 Tuttlingen · Germany
Phone: +49 (0) 7461 95 0
Fax: +49 (0) 7461 95 26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042
AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

General information

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	see variants
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the shunt with its external dimensions:	
normal:	
pediatric:	

Variants

All reservoirs are available with distal connector or distal catheter.

Pediatric borehole reservoir

Borehole reservoir

SPRUNG RESERVOIR

MCLANAHAN RESERVOIR

Pediatric borehole port

Borehole port

ETOU RESERVOIR

connector can be adjusted to various positions

The reservoirs are also available as a set, including a ventricular catheter.

BOREHOLE RESERVOIR-SET

BOREHOLE RESERVOIR-SET for pediatric application

BOREHOLE PORT-SET for pediatric application

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland