## **TITANKONNEKTOREN/ TITANIUM CONNECTORS**

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Instructions for use

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Mode d'emploi

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Maßstab 1:1

**CAUTION**

**Federal law restricts this device to sale by or on order of a physician!**

**Indication**

One possible treatment in cases of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain. This shunt system is comprised of a series of catheters and a pressure relief valve. Connectors are used to link these components. Titanium connectors are biocompatible, very stable and generate hardly any artefacts during the performance of imaging procedures.

**Tube systems**

The *CONNECTORS* allow for the application of catheters with an inner diameter of 1.2 mm to 1.5 mm. The outer diameter of the catheter should be approximately twice that of the inner diameter. For spinal drainage special *TITANIUM CONNECTORS step down* are available. In every case, the catheter must be careful fastened to the *CONNECTOR* by means of a ligature. Any crimping of the catheter should be avoided.

**Reimplantation**

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

**Compatibility with diagnostic procedures**

Owing to the fact that the *CONNECTORS* are comprised of titanium, a nonmagnetic material, the product will not interfere with nuclear magnetic resonance and computer tomography procedures. The *CONNECTORS* are MRI conditional (ASTM F2503-08).

**Safety measures**

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflammation

and tension around the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation or vomiting often occur in cases of shunt dysfunction. In the event such symptoms occur or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

**Sterilization**

All products are carefully and thoroughly sterilized. Owing to the fact that the products are packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances. The functional safety and reliability of re-sterilized products cannot be guaranteed.

**Medical products consultant**

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke  
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG  
Ulanenweg 2  
14469 Potsdam · Germany  
Phone: +49(0) 7000 MIETHKE or  
Phone: +49(0) 331 620 83 0  
Fax: +49(0) 331 620 83 40  
e-mail: info@miethke.com



Service address in the US  
AESCULAP Inc.  
Attn. AESCULAP Technical Services  
615 Lambert Pointe Road  
Hazelwood, MO, 63042  
AESCULAP Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada  
AESCULAP Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
Phone: +1-800-282-9000  
www.aesculapusa.com

**Requirements of the MDD 93/42/EEC**

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

**General information**

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>CONNECTOR</i>
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the shunt with its external dimensions:	
normal:	
step down:	

**Variants**

**TITANIUM CONNECTOR straight**

Scale 1:1

**TITANIUM CONNECTOR straight, step down**

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
## Compatibility with diagnostic procedures

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland