## ***Bohrlochumlenker/ Deflector***

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**GB** Instructions for use

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**USA** Instructions for use

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## Indikation

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

*pädiatrischer Bohrlochumlenker*

**CAUTION**

**Federal law restricts this device to sale by or on order of a physician!**

**Indication**

The *deflector* is intended for use for CSF-drainage out of the ventricles into the peritoneum in cases of hydrocephalus.

**Description**

One possible treatment in cases of hydrocephalus is the implantation of a SHUNT-SYSTEM that achieves drainage of cerebrospinal fluid away from the ventricles of the brain into a suitable area of the body (in most cases the intestinal cavity).

This SHUNTSYSTEM consists of a valve, a proximal and distal catheter and a reservoir. If the borehole reservoir is placed into the borehole, there is no need for a deflector. If the reservoir is a prechamber, the reservoir is in line with the valve and the deflector is part of the SHUNTSYSTEM. The *deflector* is positioned in the borehole in the cranium and it is possible to select the length of the catheter to be inserted into the cranium prior to implantation. The *deflector* provides a turn (at a 90 degree angle) for the ventricular catheter in the borehole, which it seals off completely.

Due to the eyelets on both sides, it is possible to fix the *deflector* with ligatures in the tissue. Additional to the standard *deflector* a small *pediatric deflector* is also available.

Fig. 1 deflector with ventricular catheter

**Implantations and reimplantations**

The ventricular catheter is implanted with the aid of a mandrin. The length of catheter to be implanted in the cranium is selected in advance through the positioning of the *deflector* on the catheter. The *deflector* is designed and intended for use with a catheter with an inner diameter of ca. 1.2 mm and an outer diameter of ca. 2.5 mm. The lower crank of the standard variant has 9,5 mm in diameter. Therefore a borehole with a diameter of 10 mm or more should be prepared. Under no circumstances products that have already been implanted in one patient should be implanted in another.


Fig. 2 deflector

**Tolerability of diagnostic procedures**

Owing to the fact that the *deflector* is composed of titanium, which is a non-magnetic material, the *deflector* will not interfere with either nuclear magnetic resonance or computer tomography procedures. The *deflector* is MRI conditional (ASTM-F2503-05).

**Safety measures**

Following implantation, the patient's condition should be carefully and thoroughly monitored. Inflamed skin and tension in the area of the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation, or vomiting often occur in cases of shunt dysfunction. In the event such symptoms do occur, or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay. Particular care must be taken not to confuse the *deflector* with the *borehole reservoir*. The *borehole reservoir* is manufactured by Christoph Miethke GmbH & Co. KG. It allows both measurement of pressure and injection of medication to be carried out by means of a puncture.

**Sterilization**

All products are carefully and thoroughly steam sterilized. Owing to the fact that the product is packaged in two layers of sterile packaging, five years of sterility is guaranteed. The expiration date is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

**Resterilization**

The safe and reliable performance of resterilized products cannot be guaranteed.

**Requirements of the MDD 93/42/EEC**

The MDD for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted product (Charge/LOT/SN) are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

**Medical products consultant**

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke  
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG  
Ulanenweg 2  
14469 Potsdam · Germany  
Phone: +49(0) 7000 MIETHKE or  
Phone: +49(0) 331 620 83 0  
Fax: +49(0) 331 620 83 40  
e-mail: info@miethke.com

Please address any enquiries to:  
AESCULAP AG  
Am Aesculap-Platz  
78532 Tuttlingen · Germany  
Phone: +49 (0) 7461 95 0  
Fax: +49 (0) 7461 95 26 00  
e-mail: information@aesculap.de

Service address in the US  
AESCULAP Inc.  
Attn. AESCULAP Technical Services  
615 Lambert Pointe Road  
Hazelwood, MO, 63042

AESCULAP Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada  
AESCULAP Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
Phone: +1-800-282-9000  
www.aesculapusa.com

**General information**

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>Deflector</i>
Intended application	Treatment of hydrocephalus
Intended for one-time use only (disposable)	
Store in a clean and dry place	
Other configurations available:	
<i>deflector</i>	
<i>pediatric deflector</i>	

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland