## ***Prechamber*** **CONTROL RESERVOIR**

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Instructions for use

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Mode d'emploi

## Indikation

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Maßstab 1:1

**CAUTION**

Federal law restricts this device to sale by or on order of a physician!

**Indication**

The *prechamber* and the *CONTROL RESERVOIR* are used in cases of hydrocephalus to remove cerebrospinal fluid and monitor pressure.

**Description**

One possible treatment of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain to an appropriate region of the body (in most instances to the abdominal cavity). This shunt system consists of a series of catheters and a differential pressure valve. In addition, the option a *prechamber* also exists. It then becomes possible to measure intraventricular pressure, inject medication and check the status of the valve. With the *CONTROL RESERVOIR* it is possible to flush CSF towards the valve due to a nonreturn valve in the proximal inlet connector of the reservoir. By this mechanism a flow in the direction of the proximal catheter is avoided during the pumping procedure. The volume per pump procedure is approx. 0.25 ml. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

**Caution: Frequent pumping of the CONTROL RESERVOIR can lead to an excessive drainage and thus to unphysiological pressure situations. The patient should be informed about the risk.**

Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0,9 mm. 30 times of punctures are possible without any restrictions.


Fig. 1: CONTROL RESERVOIR

The reservoirs are available in various product configurations (individually or as a set) and as a *CSF-PORT*. The *CSF-PORT* has only one lateral connector for the ventricular catheter.


Fig. 2 CSF-PORT

A *pediatric prechamber* is available also, it is adapted to the size of children's skulls.

**Implantation**

When implanting a reservoir, care must be taken to ensure catheters are secured by means of a ligature, if this has not already been done by the catheter manufacturer. The reservoirs are intended for use with catheters with an inner diameter of approx. 1.2 mm and an outer diameter of approx. 2.5 mm. When implanting a reservoir set, the neurosurgeon should determine how the procedure is to be performed in accordance with experience. The ventricle catheter is implanted with the aid of a introducing stylet, after which it can be brought through the deflector into the 90°-angle. During this process (and prior to implantation of the ventricular catheter), the deflector can be precisely adjusted, and then used as a stop. The reservoir is used in conjunction, and is connected to a complete shunt system. A final check of the complete shunt system is recommended to ensure that it is not obstructed in any way.

**Reimplantation**

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

**Compatibility with diagnostic procedures**

Like other shunt components of the Christoph Miethke GmbH & Co. KG (when applied), the *prechamber* and the *CONTROL RESERVOIR* consist of the nonmagnetic substances titanium, sapphire and silicone. Consequently, the reservoirs and catheters will not interfere with either nuclear magnetic resonance (MRI) or computer tomographic (CT) procedures. The *prechamber* and the *CONTROL RESERVOIR* are MR conditional (ASTM-F2503-08). The provided catheters are MRI safe. All components are visible via X-ray.

**Safety measures**

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflammation and tension around the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation or vomiting often occur in cases of shunt dysfunction.

In the event such symptoms occur or if there is any leakage in the shunt system, either the individual shunt system components or the entire shunt system must be replaced without delay.

**Sterilization**

All products are carefully and thoroughly sterilised. Owing to the fact that the products are packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

**Resterilization**

The functional safety and reliability of resterilised products cannot be guaranteed.

**Functional safety**

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH<sub>2</sub>O during and after implantation.

**Requirements of the MDD 93/42/EEC**

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

**Note on the instructions for use**

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

**Medical products consultant**

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke  
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG  
Ulanenweg 2  
14469 Potsdam · Germany  
Phone: +49(0) 7000 MIETHKE or  
Phone: +49(0) 331 620 83 0  
Fax: +49(0) 331 620 83 40  
e-mail: info@miethke.com



Please address any enquiries to:  
AESCULAP AG  
Am Aesculap-Platz  
78532 Tuttlingen · Germany  
Phone: +49 (0) 7461 95 0  
Fax: +49 (0) 7461 95 26 00  
e-mail: information@aesculap.de

Service address in the US  
AESCULAP Inc.  
Attn. AESCULAP Technical Services  
615 Lambert Pointe Road  
Hazelwood, MO, 63042


AESCULAP Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada  
AESCULAP Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
Phone: +1-800-282-9000  
www.aesculapusa.com


**General information**

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	see variants
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the shunt with its external dimensions:	
normal:	
pediatric:	


## Variants


## Prechamber and CONTROL RESERVOIR-SET

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland