M.blue plus®

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www.miethke.com

ⓄUS This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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DIMENSIONEN


TABLE OF CONTENTS

INDICATION	18
TECHNICAL DESCRIPTION	18
FUNCTION OF THE VALVE	18
SELECTION OF APPROPRIATE PRESSURE LEVEL	19
PRESSURE RATING IDENTIFICATION IN X-RAY IMAGES	20
APPLICATION OF M.BLUE PLUS INSTRUMENTS	21
POSSIBLE SHUNT COMPONENTS	25
RESERVOIRS	25
TUBE SYSTEMS	25
IMPLANTATION	26
VALVE TEST	27
PRESSURE-FLOW CHARACTERISTICS	28
PRECAUTIONARY MEASURES AND CONTRAINDICATIONS	28
FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	29
ADVERSE REACTIONS AND INTERACTIONS	29
STERILISATION	29
REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)	29
MEDICAL DEVICE CONSULTANT	30
DIMENSIONS	30


Fig. 1: Cross section of M.blue plus

INDICATION

M.blue plus is used for cerebrospinal fluid (CSF) drainage in the treatment of hydrocephalus.

TECHNICAL DESCRIPTION

M.blue plus is a valve made from titanium. It consists of an adjustable differential pressure unit (hereinafter *proGAV 2.0*) and an adjustable gravitational unit (hereinafter *M.blue*) (Fig. 1). The adjustable differential pressure unit in the proximal part of the valve (1) consists of a stable titanium housing in which a ball-cone valve (2) is integrated in the front section. A torsion spring (3) determines the opening pressure of this unit. The pretension of the spring, and therefore the valve opening pressure, can be adjusted through the skin using a revolving rotor on a bearing (4).

The adjustable gravitational unit (5) contains a tantalum weight (7), which holds a sapphire ball in the ball seat via a lever (6). Depending on the body position of the patient, the influence of the tantalum weight on the sapphire ball and thus the valve opening pressure changes. Via

a rotor (8), the pretension of the torsion spring connected to the lever can be changed through the skin. This way, the influence of the tantalum weight on the sapphire ball can be influenced and thus the valve opening pressure can be adjusted.

FUNCTION OF THE VALVE

M.blue plus is a posture-dependent hydrocephalus valve. The opening pressure for *M.blue plus* is based on the opening pressures for the adjustable differential pressure unit and the adjustable gravitational unit.

Horizontal position

In the horizontal position, the gravitational unit is always open and does not present any resistance.

Consequently, the opening pressure of *M.blue plus* in the horizontal position is characterised by the adjustable differential pressure unit. The main operation method used by the adjustable differential pressure unit is shown in Fig. 2a and b.


Fig. 2: Operation of the differential pressure unit

In Fig. 2a the valve is closed, hence drainage is not possible.

If the patient's intraventricular pressure (IVP) exceeds the spring force which otherwise keeps the differential pressure unit closed, the sealing ball moves out of the cone, leaving a gap for CSF drainage (Fig. 2b).

Vertical position

As the patient's body moves from prone to vertical, the gravitational unit closes the discharge channel in the proximal part of the valve (Fig. 3a). Thus, the opening pressure of *M.blue plus* is increased in the upright position, because now the weight of the tantalum weight (opening pressure of the gravitational unit) must be overcome in addition to the opening pressure of the differential pressure unit. Drainage is only possible once the sum of IVP and hydrostatic suction is greater than the opening pressure of both units (Fig. 3b). For individual adaptation of the opening pressure to the patient, an opening pressure between 0 and 20 cm H₂O can be selected for the adjustable differential pressure unit, and additionally a valve opening pressure between 0 and 40 cm H₂O can be selected for the adjustable gravitational unit.


Fig. 3: Gravitational unit in vertical body position

During physical activity which is associated with shock (e.g. jogging) the opening pressure of *M.blue plus* may decrease temporarily according to laboratory results. Generally however, functionality is retained. At the end of physical activity, the opening pressure returns to its original level and remains stable.

SELECTION OF APPROPRIATE PRESSURE LEVEL

Horizontal position

In the lying position, the gravitational unit has no influence on the valve opening pressure. In this position, the valve opening pressure is thus determined exclusively by the differential pressure unit. In this case, the pressure level should be set in accordance with the clinical picture and indications. By default, a differential pressure unit with an opening pressure of 5 cm H₂O is recommended.

Vertical position

The *M.blue plus* opening pressure for the vertical body position is calculated from the sum of the opening pressure in the differential pressure unit and the gravitational unit. Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure for the gravitational unit (see pressure level recommendations at <https://www.miethke.com/produkte/downloads/>). This is a non-binding recommendation for the attending physician. According to

his diagnosis, the physician decides each case independently, without instructions and individually.

PRESSURE RATING IDENTIFICATION IN X-RAY IMAGES

Differential pressure unit

The selected pressure level of the differential pressure unit for *M.blue plus* should always be monitored using the *M.blue plus* compass or the *proGAV 2.0* compass, but it can also be checked using an x-ray image. The rotor setting is decisive in this case. The four magnets in the rotor can be seen on the x-ray image as white points and are located opposite each other in pairs. Two additional burr holes (right and left next to the magnet pairs) on one side of the rotor can be used as orientation. They can be seen as black points on the x-ray image. This side can be designated as the rotor rear side. The two front magnets are on the opposite side. The space between these two magnets can be considered as the triangle apex. The pressure stage can be read off using the orientation of this intermediate space. The triangle apex can take up any position except the space labelled as a non-adjustable area in Fig. 5. This means that the opening pressure of *proGAV 2.0* can be variably adjusted from 0 up to 20 cmH₂O. To ensure that the pressure stage is not read off as a mirror image, the valve is provided with a valve marking on one side which is visible as black in the x-ray image - on a plan view of the implanted valve as in Fig. 6 the recess on the right-hand side is visible.


Fig. 6: X-ray image of the adjustable differential pressure unit, position 14 cmH₂O

Gravitational unit

The selected pressure level for *M.blue plus* should always be monitored using the *M.blue plus* compass, but it can also be checked using an x-ray image. The rotor setting is decisive in this case. The four magnets in the rotor can be seen on the x-ray image as white points and are located opposite each other in pairs. Two additional burr holes (right and left next to the magnet pairs) on one side of the rotor can be used as orientation. They can be seen as black points on the x-ray image. This side can be designated as the rotor rear side. The two front magnets are on the opposite side. The space between these two magnets can be considered as the triangle apex. The pressure level can be read off using the orientation of this intermediate space (Fig. 7). The triangle apex can take up any position except the space labelled as a non-adjustable area in Fig. 7. This means that the opening pressure of *M.blue plus* can be variably adjusted from 0 up to 40 cm H₂O. In order not to read off the pressure level in reverse, the top view of the implanted valve in the housing ring has a recess with the tantalum weight to the right of the inlet nozzle (Fig. 7).


Fig. 5: Schematic representation of rotor in x-ray image

APPLICATION OF *M.blue plus* instruments

The *M.blue plus* Instruments may only be used by trained specialist personnel. The selected pressure level of *M.blue plus* can be determined, adjusted and monitored using the *M.blue plus* instruments.

The *M.blue plus* compass is used to locate and read the adjustable differential pressure unit (*proGAV 2.0*) of the *M.blue plus* as well as the adjustable gravitational unit (*M.blue*) of the *M.blue plus*.


Fig. 7: Schematic representation of rotor in x-ray image


Fig. 8: X-ray image (adjustable gravitational unit preset to 20 cmH₂O)


Fig. 9: *M.blue plus* compass

With the *M.blue plus* adjustment ring the opening pressure of the differential pressure unit of *M.blue plus* can be variably adjusted from 0 to 20 cmH₂O and the opening pressure of *M.blue plus* can be variably adjusted from 0 to 40 cmH₂O.

Fig. 10: *M.blue plus* adjustment ring

The opening pressure for the adjustable differential pressure unit and the adjustable gravitational unit can be changed before or after implantation. The opening pressure of the differential pressure unit is preset by the manufacturer to 5 cmH₂O. The opening pressure of the gravitational unit is preset by the manufacturer to 20 cmH₂O. The

following steps must be performed to adjust the opening pressures:

1. Localisation

When the *M.blue plus* compass is opened, a circular template becomes visible with which you can localise the valve in the patient's head using your index finger (Fig. 11).

Fig. 11: Localising the valve using the *M.blue* Compass

The direction markings on the template show the flow direction.

2. Test procedure

In order to determine the selected pressure level, the compass is then closed again. The float gauge should now be centred by moving the instrument in the designated circular marking (Fig. 10). Once the float gauge is centred, the currently set opening pressure of the differential pressure unit (*proGAV 2.0*) and/or of the gravitational unit *M.blue* can be read off via the line mark on the float gauge (Fig. 12).


Fig. 12: Determination of the pressure level with the *M.blue* Compass

There are two scales on the scale ring. For the opening pressure of the adjustable differential pressure unit (*proGAV 2.0*), the grey adjustment range of 0 - 20 cmH₂O on the outer scale applies. The opening pressure of the gravitational unit (*M.blue*) is set to the blue marked adjustment range of 0 - 40 cmH₂O of the inner scale.


Fig. 13: Outside: Scale of the differential pressure unit (*proGAV 2.0*) from 0-20 cmH₂O. Inside: Scale of the gravitational unit (*M.blue*) of 0-40 cmH₂O (opening pressure of the differential pressure unit (*proGAV 2.0*) in the picture example is 17 cmH₂O, opening pressure of the gravitational unit (*M.blue*) in the picture example is 16 cmH₂O)

CAUTION

The *M.blue plus* compass should be placed centrally over the valve, if possible; otherwise incorrect determination of the opening pressure may occur.

M.blue plus compass reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus* adjustment rings should not be in the immediate vicinity of the *M.blue plus* compass when the opening pressure is being determined. We recommend a minimum distance of 30 cm.

ADVICE

Possible air bubbles inside the Compass do not affect its functionality.

3. Adjustment process

In order to adjust the opening pressure, the compass is opened, but without changing the position of the scale ring. The adjustment ring is now inserted into the scale ring in such a manner that the line marking points to the desired value on the scale of the scale ring. (Fig. 14)


Fig. 14: Insertion of the adjustment ring

For the opening pressure of the differential pressure unit (*proGAV 2.0*), the grey adjustment range of 0 - 20 cmH₂O on the outer scale applies. The opening pressure of the *M.blue* gravitational unit uses the adjustment range of 0 - 40 cmH₂O of the blue inner scale.

Fig. 15: Setting the differential pressure unit (*proGAV 2.0*) in the picture example to 1 cmH₂O. Or rather Adjustment of the gravitational unit (*M.blue*) in the picture example to 32 cmH₂O.

By applying slight pressure with the index finger to the valve diaphragm located in the centre of the adjustment ring and under the skin, the rotor brake is released and the opening pressure of the differential or gravitational unit is changed to the desired value (Fig. 16).


Fig. 16: Adjustment using the *M.blue plus* adjustment ring

Both the differential pressure unit (*proGAV 2.0*) and the gravitational unit (*M.blue*) are equipped with a feedback mechanism. If targeted pressure is exerted on the valve, an acoustic signal (a clicking sound) is audible, or a resistance can be felt to give, as soon as the rotor brake has been released due to the valve housing design. In other words the valve shows both acoustically and haptically when the pressure is sufficient for uncoupling of the brake. Once this pressure has been released the rotor is adjustment-proof again. Although the click caused by releasing the rotor brake is easily audible before implantation, considerably reduced after implantation. This can be due to the CSF filling the valve, its position and the condition of the implants surroundings. Normally, however, it should be audible to the patients themselves or through the use of a stethoscope.

Adjustment with the *M.blue plus* Adjustment Assistant

Alternatively, the *M.blue plus* Adjustment Assistant can be used to adjust the opening pressure. To do this, insert the *M.blue plus* Adjustment Assistant into the adjustment ring aligned to the desired value and press it with your index finger (Fig. 17).

Fig. 17: *M.blue plus* Adjustment Assistant

WARNING

For adjustment of the differential pressure unit (proGAV 2.0), ensure that the opening pressure is changed by a maximum of 8 cmH₂O per adjustment process; otherwise errors can result.

Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in 2 stages: Initial adjustment from 3 to 11 and subsequently from 11 to 18 cmH₂O.

WARNING

When adjusting the gravitational unit (M.blue), please ensure that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process; otherwise errors may result.

Example: The opening pressure is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 6 to 22 and subsequently from 22 to 36 cmH₂O.

CAUTION

The M.blue plus adjustment ring emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.

4. Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the selected opening pressure. To do this, please proceed as in Points 1 and 2. If the measured pressure does not agree with the required pressure stage the adjustment process should be repeated. To do this, please start at Point 3 again. Adjustment can be made difficult by swelling of the skin for a few days after the operation! If checking the valve setting is not fully possible using the *M.blue plus compass*, we recommend carrying out a check using an imaging method.

proGAV Check-mate

The *proGAV Check-mate* is supplied sterile and can be re-sterilised. With the *proGAV Check-*

mate it is possible to carry out a pressure level change and the check before and during valve implantation directly on the adjustable differential pressure unit (*proGAV 2.0*) of *M.blue*. To determine the pressure level the *proGAV Check-mate* is placed centrally on the *proGAV 2.0*. The *proGAV Check-mate* automatically aligns itself over the valve. The pressure stage can be read off from the catheter in the proximal (leading to valve) direction. If the pressure level is supposed to be adjusted, the *proGAV Check-mate* is placed centrally on the *proGAV 2.0*. When doing so, the required pressure stage must point towards the proximal catheter (leading to valve). If the *proGAV Check-mate* is now pressed lightly onto the valve, the rotor brake in *proGAV 2.0* is released and the pressure level is set.

When making the adjustment, ensure that the opening pressure is changed by a maximum of 8 cmH₂O per adjustment process; otherwise errors may occur (see chapter "3. Adjustment process").

Fig. 18: proGAV Checkmate, Pressure levels 0-20 cmH₂O

M.blue Checkmate

The *M.blue checkmate* is supplied sterile and can be re-sterilised. With the *M.blue checkmate* it is possible to carry out a pressure level change and the check before and during valve implantation directly on the adjustable gravitational unit (*M.blue*) of *M.blue plus*. To determine the pressure level the *M.blue checkmate* is placed centrally on *M.blue*. The *M.blue checkmate* automatically aligns itself over the valve. The pressure stage can be read off from the catheter in the proximal (leading to valve) direction. If the pressure level is supposed to be adjusted, the *M.blue checkmate* is placed centrally on *M.blue*. When doing so, the required pressure stage must point towards the proximal catheter (leading to valve). If the *M.blue checkmate* is now pressed lightly onto the valve, the rotor brake in

M.blue is released and the pressure level is set. When adjusting, please ensure that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process otherwise errors can result (see Chapter "3. Adjustment process").

Fig. 19: M.blue checkmate, colour: blue, pressure levels 0-40 cmH₂O

CAUTION

Because of the magnets inside the M.blue plus instruments, M.blue plus instruments may not be used in the vicinity of active implants such as heart pacemakers. Furthermore, there is risk of damage to MRI instruments in their vicinity. Therefore, usage of M.blue plus instruments is forbidden in such locations!

It is absolutely essential to use only the *M.blue plus instruments* to determine, adjust and monitor the opening pressure of the gravitational unit in *M.blue plus*. The opening pressure of the differential pressure unit (*proGAV 2.0*) of *M.blue plus* can be determined, adjusted and monitored with the described *M.blue plus instruments* as well as with the *proGAV 2.0* instruments.

Fig. 20: M.blue plus instruments

POSSIBLE SHUNT COMPONENTS

M.blue plus can be ordered as a shunt system in a range of configurations. These configurations can be combined with the accessories presented below. In each case, versions for paediatric hydrocephalus and for hydrocephalus in adults are available.

Reservoirs

The use of a reservoir in combination with shunt systems provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to an integrated check valve in the *SPRUNG RESERVOIR* and the *CONTROL RESERVOIR*, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the drainage system as well as the valve catheter. During the pump action, access to the ventricular catheter is closed. The use of reservoirs does not increase the opening pressure of the shunt system. Puncturing the reservoir should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. A stable titanium floor prevents the bottom surface from being pierced. 30 punctures are possible without any restrictions.

WARNING

Frequent pumping can result in excessive drainage and thus lead to pressure conditions outside the normal physiological range. The patient should be informed about this risk.

Burr hole deflector

Because of the tight fit on the ventricular catheter, the burr hole deflector makes it possible to choose the length of catheter penetrating into the skull prior to implantation. The ventricular catheter is deflected at a right angle in the burr hole (see chapter "Implantation").

TUBE SYSTEMS

M.blue plus can be ordered as an individual valve unit or as a shunt system with integrated catheters (interior diameter 1.2 mm, exterior di-

ameter 2.5 mm). The supplied catheters do not fundamentally change the pressure-flow characteristics. If catheters by other manufacturers are used, a tight fit must be ensured. In any case, catheters have to be carefully fixed with a ligature to the valve's titanium connectors.

IMPLANTATION

Positioning of the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The required skin incision should be made in form of a lobule pedicled towards the draining catheter. If a burr hole deflector is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole.

M.blue plus is available in a range of different configurations: If a burr hole reservoir or a *SPRUNG RESERVOIR* is used, then the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the burr hole reservoir connected, with the connection secured with a ligature. When using a shunt system with a *CONTROL RESERVOIR*, a burr hole deflector is included. The deflector is used for adjusting the length of catheter to be implanted and for its positioning inside the ventricle. The ventricular catheter is deflected, and the *CONTROL RESERVOIR* is put into place. The position of the ventricular catheter should be inspected after the procedure by imaging (such as CRT or MRI).

Positioning of the valve

The adjustable differential pressure unit in *M.blue plus* is set to an opening pressure of 5 cmH₂O upon delivery. The adjustable gravitational unit in *M.blue plus* is set to an opening pressure of 20 cmH₂O upon delivery. The preset opening pressures can be changed to a different pressure before implantation. A location behind the ear is suitable as an implantation position, whereby the implantation height has no influence on the valve function. The adjustable units should be touching the

bone or the periosteum since pressure must be exerted on the valve during any later adjustment. A large arch-shaped or a small straight skin incision with a pocket for the valve should be made. The catheter is then pushed forward from the burr hole to the selected valve implantation location, shortened if necessary, and secured to the *M.blue plus* with a ligation. The valve should not be located directly under the skin incision. The valve unit has an arrow in the flow direction (arrow towards distal or downwards).

NOTE

The gravitational unit of *M.blue plus* is position-dependent. You must therefore ensure that the gravitational unit (*M.blue*) is implanted parallel to the body axis.

Therefore, if a shunt system in which the valve has been pre-fitted with a burr hole reservoir is being used, only the occipital access should be used.

WARNING

The adjustable valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e.g. underneath heavily scarred tissue).

CAUTION

The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

Positioning of the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. For example, it can be applied paraumbilical or at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. The recommendation is to pull the peritoneal catheter using a subcutaneous tunnelling tool from the valve to the intended position, if necessary with the aid an auxiliary incision. The peritoneal catheter that is usually securely attached to *M.blue plus* has an open distal end and no wall slits. Following the ex-

posure of the peritoneum or with the aid of a trocar, the peritoneal catheter (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

VALVE TEST

Preoperative valve test

M.blue plus should be vented before implantation and checked for permeability. The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological salt solution. The valve is continuous if saline solution can be extracted (Fig. 21).

WARNING

Contamination in the solution used for testing can impair the product's performance.


Fig. 21: Patency test

WARNING

Pressure admission through the single-use syringe should be avoided both at the proximal and the distal end (Fig. 22).


Fig. 22: Avoidance of pressurisation

Postoperative valve test

The *M.blue plus* has been constructed as a reliably functioning unit without pump or test function. The valve test can be performed by flushing, pressure measurement or pumping.

PRESSURE-FLOW CHARACTERISTICS

Horizontal position

Fixed differential pressure unit *M.blue 0*

The pressure-flow characteristic of the fixed differential pressure unit of *M.blue 0* is shown below.


Fig. 23: pressure-flow characteristic of *M.blue 0* in the horizontal body position

Adjustable differential pressure unit *proGAV 2.0*

The pressure flow characteristics for the adjustable differential pressure unit are shown below as an example for pressure levels 0, 10 and 20.


Fig. 24: Pressure-flow characteristics for selected pressure levels in the adjustable differential pressure unit

Vertical position

Adjustable gravitational unit *M.blue 0*

The pressure flow characteristics for various pressure level settings of *M.blue 0* in the vertical body position are shown below. Only values of the adjustable gravitational unit without consideration of the adjustable differential pressure unit are shown.


Fig. 25: Pressure-flow characteristic of *M.blue 0* in the vertical body position

The opening pressure refers to a reference flow of 20 ml/h.

PRECAUTIONARY MEASURES AND CONTRAINDICATIONS

Patients must be carefully monitored after implantation. Reddening of skin or tightness in the area of the drained tissue may be indications of infections at the shunt system. Symptoms such as headache, dizziness, confusion or vomiting often occur in conjunction with shunt dysfunction. These symptoms and a leakage within the shunt system require the immediate replacement of the affected shunt component or the entire shunt system. The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e.g. meningitis, ventriculitis, peritonitis, bacteraemia, septicaemia) in the region affected by the implantation.

FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. However, no guarantee can be given that these medical devices may not require replacement for medical or technical reasons. These medical devices are able to resist positive and negative pressures up to 200 cmH₂O during and after implantation. These medical devices have to be stored in a clean and dry environment at all times. Nuclear magnetic resonance (MRI) examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. The valve is MR compatible. The supplied catheters are MR safe. Reservoirs, deflectors and connectors are MR compatible.

WARNING

⚠ If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, it is not possible to rule out valve adjustment. In MRI imaging *M.blue plus* creates artefacts which are larger than the valve itself.

CAUTION

⚠ For people using cardiac pacemakers: It is possible that the function of the pacemaker is influenced by the implantation of a *M.blue plus*.

ADVERSE REACTIONS AND INTERACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): Infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in very rare cases noise development. Violent shocks for the outside (accident, fall) may put the integrity of the shunt system at risk.

STERILISATION

The products are sterilised with steam under strictly controlled conditions. The double wrapping in sterile bags ensures sterility for a five-year period. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used in any circumstances. No guarantee can be given for the functional safety and reliability of resterilised products.

REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the whereabouts of medical devices used in humans, especially for implants. The individual identification number of the implanted valve should therefore be recorded in the patient's medical records and patient passport to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website (<https://www.mietnke.com/produkte/downloads/>).

MEDICAL DEVICES CONSULTANT

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical product consultants as contacts for all product-related questions.

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz
Michaela Funk-Neubarth
Dipl.-Ing. Thoralf Knitter
Dr Andreas Bunge
Jan Mgel

Contact details can be found on the reverse of these instructions for use.

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany