# Reader Unit Set for *M.scio*<sup>®</sup>

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 [www.miethke.com](http://www.miethke.com)

Ⓔ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0  
info@miethke.com

**TABLE OF CONTENTS**

SCOPE OF THE SYSTEM	30
INDICATIONS	32
INTENDED PURPOSE	32
DEVICE DESCRIPTION	32
SAFE HANDLING	32
INTENDED USE	32
CONTRAINDICATIONS	32
OPERATING PRINCIPLE	33
SYMBOLS IN THE READER UNIT DISPLAY	34
WORKING WITH THE READER UNIT	35
SYSTEM START-UP	36
FUNCTION CHECKS	36
MEASUREMENT SETUP	37
OPERATING THE READER UNIT	38
GENERAL OPERATING INSTRUCTIONS	38
INFO MENU	38
TAKING MEASUREMENTS	39
DATA MANAGEMENT   DISPLAY OF MEASURED DATA	41
SETTINGS	44
CLEANING   DISINFECTION	48
CALIBRATING THE BAROMETRIC PRESSURE SENSOR   TECHNICAL SERVICE	49
SAFETY CHECKS (SC)	50
TROUBLESHOOTING	50
MALFUNCTIONS DISPLAYED WITH AN ERROR MESSAGE	51
OTHER MALFUNCTIONS   POSSIBLE FAULTS	54
TECHNICAL DATA	55
AMBIENT CONDITIONS	56
CE MARKING	56
STORAGE AND TRANSPORT	56
DISPOSAL	56
NOTES ON ELECTROMAGNETIC COMPATIBILITY	56
ELECTROMAGNETIC EMISSIONS	57
ELECTROMAGNETIC IMMUNITY	57
MEDICAL DEVICES CONSULTANTS	57

## SCOPE OF THE SYSTEM

Figure 1: Reader Unit Set

**Approved components of the Reader Unit Set:****1. Reader Unit**

- 1.1 ON/OFF switch
- 1.2 Display screen
- 1.3 Function key
- 1.4 Antenna socket
- 1.5 Connector for plug-in power unit
- 1.6 SD card slot with cover

**3. Plug-in power unit**

- 3.1 Plug
- 3.2 Control light
- 3.3 Attachment for EU/UK

**2. Antenna****Further scope of supply**

- a) Case for *Reader Unit Set*, incl. key
- b) Instructions for use

Symbol	Explanation
	Caution, general warning sign Caution, read the documentation supplied with the product
	Follow the instructions for use
	Non-ionising radiation: On devices fitted with RF emitters or that are intended to use electromagnetic RF energy for diagnostic or treatment purposes (e.g. RF surgical devices)
	ESD* sensitivity: On electric connectors of devices that require special protective measures against electrostatic discharge (* Electrostatic discharge)
	Labelling of electric and electronic devices in accordance with Directive 2002/96/EG (WEEE), see chapter on "Disposal"
	ON/OFF switch
	Function key: The 4 function keys are used to perform the options shown on the display screen
<b>IP44</b>	Reader 7502 0000
<b>IP44</b>	Antenna 7503 0000
<b>IP40</b>	Plug-in power unit 7504 0000
	Protection class II (protective insulation)
	Reader 7502 0000 Antenna 7503 0000
	Antenna symbol Antenna socket
	SD card SD card slot
	DCIN connector for plug-in power unit

## INDICATIONS

- Functional test of the shunt system implanted into the patient
- Technical functional diagnosis of the shunt system, e.g. detection and localisation of occlusions

## INTENDED PURPOSE

The *Reader Unit Set* is used for the pressure-based functional control of shunt systems.

## DEVICE DESCRIPTION

Designation	Article number
<i>Reader Unit Set</i>	75100000 (AAG ref. FV907X)

## SAFE HANDLING

### CAUTION

Risk of injury due to incorrect operation of the product! The product user must have participated in the product training before first starting-up the product. Please contact Christoph Miethke GmbH & Co. KG for information about product training.

The following general safety instructions must be followed when working with electrical medical devices:

- Prior to use, check product for proper condition and functioning.
- Clean factory-fresh product thoroughly after removing transport packaging (see chapter on "Manual Cleaning/ Disinfection").
- To prevent nosocomial infection and multiple resistances, the device should be disinfected after each use. Disinfection is preferably carried out by wiping with the disinfectant described in the section on "Cleaning/ Disinfection".
- To prevent damage as the result of improper assembly or operation and to avoid compromising warranty and liability, please note the following:
  - Use the product only as specified in these instructions for use.
  - Observe notes on electromagnetic compatibility (EMC).

- Observe safety information and maintenance instructions.
- Only combine products and use accessories by Christoph Miethke GmbH & Co. KG, see the "Scope of the system" chapter.

- Do not operate the product in the vicinity of flammable substances (e.g. anaesthetics).
- Device and accessories must only be operated and used by persons with the requisite training, knowledge or experience.
- Keep the instructions for use accessible for the user.
- The product must be set up so that the mains plug can easily be disconnected from the power supply.

## INTENDED USE

The *Reader Unit Set* (hereinafter referred to as *Reader Unit*) is used together with the *M.scio* for functional diagnostics. They are used for the reliable and non-invasive detection and localisation of occlusions within ventricular drainage systems (without intervention or radiological diagnostics) and the identification of mechanical damage to shunt valves.

## CONTRAINDICATIONS









- Intracranial pressure measurement
- Intracranial pressure measurement without a shunt system
- Intracranial pressure management based on measured values
- Functional diagnosis of extracorporeal fluid drainage system (cranio-cerebral trauma)

## OPERATING PRINCIPLE

The telemetric system serves to measure CSF pressure readings whose relative, time-dependent changes and characteristics may allow conclusions as to the function of an implanted ventricular drainage system. Diagnostic methods can be used to initiate a pressure change in the shunt system, which can then be interpreted. The *M.scio*'s measured

data can be read out and displayed with the *Reader Unit*. Measured data are automatically stored on the SD card and can thus be evaluated at a later stage. The relative behaviour of the CSF pressure in the shunt system can provide information about its function. This permits the non-invasive detection and localisation of an occlusion within the shunt system and also of a mechanical loss of function.

## SYMBOLS IN THE READER UNIT DISPLAY

	<Info>	For displaying additional information or starting with the information menu
	<Menu>	For starting with the selection menu
	<Fast measurement>	Fast measurement can be started without a further selection
	<Up arrow>	Upward navigation of cursor
	<Down arrow>	Downward navigation of cursor
	<Exit menu>	For quitting the current menu
	<OK>	Activates or confirms the selected function
	<Start>	Starts measurement
	<Stop>	Stops the current measurement
	<Delete>	For deleting a saved measurement in the [Data management] menu
	<Star>	For placing a marker during a [Continuous measurement]
	<Editing diagram>	For opening the diagram processing menu
	<Zoom menu>	For activating the zoom function
	<Cursor menu>	For activating the cursor function
	<Pressure axis values>	For adjusting the display of min. and max. values of the pressure axis
	<Exit cursor menu>	For returning to the diagram processing menu
	<Exit zoom menu>	For returning to the diagram processing menu
	<Right arrow>	Navigate forward
	<Left arrow>	Navigate backward
	<On>	On
	<Off>	Off
	<Confirmation>	Confirm
	<Zoom in>	For zooming into the display of a sequence in the [Data management] menu
	<Zoom out>	For zooming out of the display of a sequence in the [Data management] menu
	<Mute>	For deactivating the acoustic signal

## WORKING WITH THE READER UNIT

## CAUTION

Danger from use of MRI.

- Only operate the *Reader Unit* outside the MRI field of application.

## NOTE

Limited accuracy and no storage of measured data for measurements without SD card.

- For the *M.scio*, the SD card assembled by the manufacturer must be inserted into the *Reader Unit*. The implant (see patient data card) is associated with the respective SD card via the identification number (ID). Measured data can only be saved on the designated SD card.

## NOTE

Unexpected termination of measurement process if SD card is full.

- Check memory status prior to starting measurement and delete data no longer required.

## NOTE

The *Reader Unit* cannot be operated with a low battery status.

- Check battery status. Connect the plug-in power unit if battery status is low.

## NOTE

No or bad communication. The telemetric link between the *Reader Unit* and the implant may be disrupted if a further *Reader Unit* is used.

- Increase the distance between *Reader Units*.

## CAUTION

The telemetric link between the antenna and the implant may be influenced by a further implant within the vicinity of the antenna.

- Place the antenna in a central position above the implant to be read. Increase the distance between implants.

## NOTE

The telemetric link between the antenna and the implant may be disrupted by metal components within the vicinity of the implant.

- Increase the distance from metal components.

## NOTE

If the temperature of the *M.scio* falls below 20°C, use the procedure for "Measurement without SD card". The device will be operated with limited accuracy and measured data cannot be stored. The following error message can occur during measurement with SD card:

- During continuous or fast measurement: Pressure readings are invalid (Info: measurement will continue)
- During a single or fast measurement: Temperature readings are invalid (Info: measurement will not continue)

For a more detailed description of errors, see chapter on "Troubleshooting".

## SYSTEM START-UP

## NOTE

Allow the *Reader Unit* to acclimatize for approx. three hours at room temperature.

The *Reader Unit* (accessible part, classified as application part type BF) is equipped with a battery which, on a 100% charge, permits independent operation without external power supply for up to five hours. To that end, the battery needs to be charged when first using the device. The plug-in power unit is used to charge the Reader Unit battery. It takes approx. 6 hours for the battery to be completely charged. The *Reader Unit* can also be operated during the charging process (with connected plug-in power unit).

## NOTE

Charging is not possible at an ambient temperature of more than 35°C!

## Connecting the power supply

## WARNING

Incorrect accessories can damage the *Reader Unit* and pose a safety risk to user and patient!

- Only use an original plug-in power unit for the *Reader Unit*.

Mains voltage must correspond to the voltage range stated on the label of the *Reader Unit* plug-in power unit.

- Insert the connector of the plug-in power unit into the socket of the *Reader Unit*.
- Plug the plug-in power unit into a wall power socket (control lamp of the plug-in power unit must be on).

## Inserting the SD card

## NOTE

Improper use can damage the SD card!

- Do not touch the contact strips on the SD card

Insert the SD card designated for use with the respective *M.scio* (see patient data card) into the *Reader Unit* until it locks into place. SD cards can be ordered from the manufacturer if required.

Tap the SD card slightly to remove.

## NOTE

The device may be damaged by the ingress of liquids!


- After removing the SD card, the SD card slot connector must be inserted back into the device.

## FUNCTION CHECKS

- Before use, check that the entire system is in good condition and fully functional and check the system settings, such as the pressure unit (see the "Units" chapter).
- To check the battery status, perform a function check without the plug-in power unit.
- Ensure that neither the device nor the antenna used for the function test show any visible signs of damage. Check the function of the following components in the given order:
  1. Press the [On/Off] button.
  2. Automatic self-test after switching on including display and loudspeaker test
- Check the battery status; if necessary, plug in the *Reader Unit* plug-in power unit and charge the battery.


The following is displayed on the screen:

- Fade-in [Selftest ...]
- Fade-in [booting ...]

The following is displayed on the screen of the *Reader Unit* during operation with the plug-in power unit:

Display for operation with  
Reader Unit power supply ↓

#### NOTE

The control lamp of the plug-in power unit must be on!

- Check the system time, correct if necessary, see the [Settings] chapter.
- Activate the [On/Off] button.

#### NOTE

The operation of the *Reader Unit* can be safely shut down at any time by pressing the [On/Off] button.

## MEASUREMENT SETUP

The *Reader Unit* is intended for use by medical staff at professional healthcare facilities. For a pressure measurement, the reader must be placed on a base and the *Reader Unit*'s antenna is positioned about 10-30 mm away from the *M.scio*. The operation of the *Reader Unit* to start a measurement is described in the "Measurement" chapter.


Figure 2: The working distance for telemetric data transmission, i.e. the distance between the antenna and the *M.scio* is 10-30 mm

## OPERATING THE READER UNIT

### Key

- A Date
- B Time
- C SD card
- D Memory assignment (of the SD card)
- E Battery status
- F OK button
- G Down arrow
- H Exit menu

A B C D E

F G H

Activate the Reader Unit with [On/Off] switch. The start screen content is displayed.

## GENERAL OPERATING INSTRUCTIONS

The menu-controlled operation of the device uses the four function keys. The context-related function of these soft keys is indicated by the symbols shown on the display for each of these keys. The exact meaning of these symbols is explained in the list of symbols in the *Reader Unit* display in chapter "Symbols in the Reader Unit display".

The preselected submenu is highlighted by a frame. For improved navigation, the currently selected position in the submenu is highlighted by a scroll bar in addition to the frame. The <OK>, <Up arrow>, <Down arrow> and <Exit menu> buttons are used for standard operation.


When not in use, the *Reader Unit* switches to stand-by mode within 1-5 minutes, depending on the setting.

## INFO MENU

The following information can be called up from the info menus:

- [Hotline & Service]
- [System details]
- [Manufacturer info]
- [Only maintenance]

To access the [Info] menu, press the <Info> button. The following information is displayed on the screen:

The <Up arrow> and <Down arrow> buttons can be used to preselect the desired submenu, which is confirmed with the <OK> button. Return to the previous menu item with the button <Exit menu>.

### Submenu content

[Hotline & Service]

- [Hotline & Service]
- [Tel: +49 331 620 83-0]

### System details

- [Product name: *Reader Unit*]
- [Article number: 7510 0000]
- [Serial number: XXXXX]
- [Software version 2.0X]
- [Service date: dd/mm/yyyy]


### Manufacturer Info

- Christoph Miethke GmbH & Co. KG  
Ulanenweg 2 | 14469 Potsdam  
Germany


## TAKING MEASUREMENTS

Press the <Menu> button to open the [Menu]. The following information is displayed on the screen:

The <OK> button to open the [Measurement] submenu. The following information is displayed on the screen:


There are three types of measurement:

1. [Single measurement]: The mean of 8-10 measurements is displayed and saved to the SD card.
2. [Continuous measurement]: The sequence of single measurements within the selected time interval is shown as a time diagram and save to the SD card.
3. [Fast measurement]: The sequence of unaveraged single measurements at the maximum available measurement rate (max. 44 measurements per second) is shown as a time diagram and save to the SD card.

### Single measurement


Press the <OK> button to select the [Single measurement] menu.

The following information is displayed on the screen:


Press the <Start> button to start the [Single measurement].

The following screen content is displayed:


### Continuous measurement

In the [Measurement] menu, press the <Up arrow> or <Down arrow> button to select the [Continuous measurement] menu and confirm with <OK>. The following information is displayed on the screen:

Press the <Start> button to start the [Continuous measurement].

The following information is displayed on the screen:

Press the <Star> button to place a marker. Multiple markers can be placed during a measurement process. Press the <Stop> button to stop the measurement.

Markers permit the circumstantial analysis of the measured data.


**Symbols for reception quality**

Symbol	Explanation
●	Communication started
★	Distance from antenna to measuring cell: - good
↔	Distance from antenna to measuring cell too close: - increase distance
→←	Distance from antenna to measuring cell too large: - decrease distance

**Fast measurement**


There are two ways to start a [Fast measurement]:

On the start screen, press the <Fast measurement> button directly and then start the measurement with the <Start> button. Alternatively, you can press the <Menu> button to display the following submenu:


Press the <OK> button to open the [Measurement] menu. In the [Measurement] menu, press the <Up arrow> and <Down arrow> to select the [Fast measurement] menu.

The following is displayed on the screen:


Press the <OK> button to select the [Fast measurement].

Press the <Start> button to start the Fast measurement.


During fast measurement, the following is displayed on the screen:


Press the <Stop> button to stop the measurement. For an explanation of the symbols for reception quality: see section on “Continuous measurement”.

#### DATA MANAGEMENT | DISPLAY OF MEASURED DATA


In the start screen, press the <Menu> button to display the following submenu:

The <Up arrow> and <Down arrow> buttons can be used to select the [Data management] menu, which is confirmed with the <OK> button.

Measured data are listed chronologically (from the start of measurement) and shown as follows:

The <Up arrow> and <Down arrow> buttons can be used to select the desired file, which is confirmed with the <OK> button.

Measured data files are labelled as follows:


Symbol	Explanation
	Single measurement
	Continuous measurement
	Fast measurement

An hour glass is shown when loading the data for a continuous or fast measurement and the loading time displayed in seconds.

#### Single measurement

Press the <Delete> button to delete the measurement file.

Deletion must be confirmed by pressing the <Confirmation> button. Press the <Exit menu> button to return to the above view. Press this button again to return to the selection menu.

#### Continuous measurement


Press the <Delete> button to delete the measurement file. Deletion must be confirmed by pressing the <Confirmation> button. Press the <Exit menu> button to return to the above view. Press the <Editing diagram> button to open the diagram editing menu. The following information is displayed on the screen:


Press the <Zoom menu> button to open the zoom menu.

The following information is displayed on the screen:


Press the <Zoom in> or <Zoom out> buttons to display an expanded or compressed time diagram of the measured sequence. Press the <Exit zoom menu> button to return to the diagram editing menu.

Press the <Cursor menu> button in the diagram editing menu to open the following screen:


The <Exit cursor menu> button takes you back to the diagram editing menu. Press the <Right arrow> or <Left arrow> button to move the time axis to the left or right.

Press the <Pressure axis values> button in the diagram editing menu to open the following screen:


Use the <Up arrow> and <Down arrow> to manually set the maximum or minimum scale values on the cursor position. Make sure that the [Automatic] function is switched off. If the [Automatic] function is selected, the setting of maximum or minimum values is irrelevant. In this function, axis scaling is automatic according to the measured data. Press the <Right arrow> button to toggle the cursor between [Maximum], [Minimum] and [Automatic]. Press the <Exit menu> button to return to the measured value display.

**Fast measurement**


Press the <Delete> button to delete the measurement. Deletion must be confirmed by pressing the <Confirmation> button. Press the button to return to the previous view. Press the <Exit menu> button to return to the selection menu. Press the <Editing diagram> button to open the diagram editing menu.

The following information is displayed on the screen:

Press the <Zoom menu> button to open the <Zoom> menu.


Press the <Zoom in> or <Zoom out> buttons to display an expanded or compressed time diagram (zoom) of the measured sequence.


Press the <Exit zoom menu> button to return to the diagram editing menu.

Switching between displays can take a few seconds, depending on file size.

Press the <Cursor menu> button in the diagram editing menu to open the following screen:


The <Exit cursor menu> button takes you back to the diagram editing menu. Press the <Left arrow> or <Right arrow> button to move the time axis to the left or right. Press the <Pressure axis values> button in the diagram editing menu to open the following screen:


Use the <Up arrow> and <Down arrow> to manually set the maximum or minimum scale values on the cursor position. Make sure that the [Automatic] function is switched off. If the [Automatic] function is selected, the setting of maximum or minimum values is irrelevant. In this function, axis scaling is automatic according to the measured data. Press the <Right arrow> button to toggle the cursor between [Maximum], [Minimum] and [Automatic].

## SETTINGS


Press the <Menu> button in the start screen to open the following submenu:

Press the <Up arrow> and <Down arrow> buttons to select the [Settings] menu.


Press the <OK> button to open the [Settings] menu. The following is displayed on the screen:

### Date/Time

Press the <OK> button to open the [Date/Time] menu. The following information is displayed on the screen:


Press the <Right arrow> button to change the cursor position. Use the <Up arrow> and <Down arrow> buttons to change the values at the cursor position.


Modified values are saved immediately.

### Units

In the [Settings] menu, press the <Up arrow> and <Down arrow> buttons to select the [Units] menu.

Press the <OK> button to open the <Units> menu. The following is displayed on the screen:


The <Up arrow> and <Down arrow> buttons can be used to select the desired unit, which is confirmed with the <OK> button.

### Measurement interval


These settings only function in the continuous measurement mode.

In the [Settings] menu, press the <Up arrow> and <Down arrow> buttons to select the [Measurement interval] menu.

Press the <OK> button to open the [Measurement interval] menu.


The following information is displayed on the screen:


The <Down arrow> and <Up arrow> buttons can be used to select the desired measurement interval, which is confirmed with the <OK> button. Possible settings range from 1-300 s. The available memory space on the SD card is also shown.

### Hide graph

In the [Settings] menu, press the <Up arrow> and <Down arrow> buttons to select the [Hide graph] menu.


Press the <OK> button to open the [Hide graph] menu. The following is displayed on the screen:


[Display graph during continuous measurement] can be selected or deselected.

#### Language

In the [Settings] menu, press the <Down arrow> and <Up arrow> buttons to select the [Language] menu.


Press the <OK> button to open the [Language] menu. The following information is displayed on the screen:


The <Up arrow> and <Down arrow> buttons can be used to select the desired language, which is confirmed with the <OK> button.

#### Volume

In the [Settings] menu, press the <Up arrow> and <Down arrow> buttons to select the [Volume] menu.

Press the <OK> button to open the [Volume] menu. The following information is displayed on the screen:

The <Up arrow> and <Down arrow> buttons can be used to select the desired volume, which is confirmed with the <OK> button. Available settings are 1-5. When setting the volume, a corresponding acoustic signal is emitted.


**NOTE**


The setting determines the volume of acoustic alerts. Exception: Error alerts are generally emitted at volume level 5.

**Standby mode**

In the [Settings] menu, press the <Up arrow> and <Down arrow> buttons to select the [Standby mode] menu.


Press the <OK> button to open the [Standby mode] menu. The following information is displayed on the screen:


The <Up arrow> and <Down arrow> buttons can be used to select the desired time after which the *Reader Unit* is to switch to standby mode, which is confirmed with the <OK> button. Available settings are 1-5 minutes.

**Brightness**

In the [Settings] menu, press the <Up arrow> and <Down arrow> buttons to select the [Brightness] menu.


Press the <OK> button to open the [Brightness] menu. The following information is displayed on the screen:


The <Up arrow> and <Down arrow> buttons can be used to select the desired brightness of the background illumination, which is confirmed with the <OK> button. Possible settings are 000-100 (in steps of 5).

**Contrast**

In the [Settings] menu, press the <Up arrow> and <Down arrow> buttons to select the [Contrast] menu.

Press the <OK> button to open the [Contrast] menu. The following information is displayed on the screen:

The <Up arrow> and <Down arrow> buttons can be used to select the desired display contrast, which is confirmed with the <OK> button. Possible settings are 000-100 (in steps of 5).

## CLEANING | DISINFECTION

### WARNING

Risk of electric shock and fire!

- Disconnect from mains before cleaning.
- Do not use flammable or explosive detergents or disinfectants.
- Ensure that no fluids can penetrate into the product; insert the SD card slot connector into the *Reader Unit*.

### NOTE

Mechanical cleaning or disinfection as well as inappropriate detergents or disinfectants can damage or destroy the product!

- Only disinfect the product manually
- Never sterilize the product

### NOTE

Use detergents/disinfectants approved for surface cleaning only as specified by the manufacturer. Observe specifications regarding concentration, temperature and contact time.

## Manual cleaning/disinfection

Wipe disinfection for electrical devices without sterilization

Phase	I
Step	Wipe disinfection
T (°C/°F)	RT (room temperature)
t (min)	≥1
Conc. (%)	-
Water quality	-
Chemical	Meliseptol HBV wipes 50% Propan-1-ol

## Phase I

- Remove any visible residue with a disposable disinfectant wipe.
- Wipe all surfaces of the visually clean product with a fresh disposable disinfectant wipe.
- Observe the specified contact time (1 minute minimum).

## Inspection

- Check product for damages after each cleaning/disinfection.
- Immediately set aside product if damaged.

## Storage

- Store products protected from dust in a dry, dark and cool room with as little bacterial contamination as possible.

## CALIBRATING THE BAROMETRIC PRESSURE SENSOR | TECHNICAL SERVICE

The *Reader Unit* contains a barometric pressure sensor (BARD sensor). An annual calibration is required to ensure that predefined tolerance limits are maintained.

### NOTE

As part of the calibration, a comprehensive functional check of the device is also performed. Omitting to perform the annual calibration can result in a drift of the BARD sensor outside of the tolerance limits.

The device must therefore be sent to the technical service department on an annual basis. The requirements defined in the "Storage and transport" chapter as well as in the "Ambient conditions" chapter must be observed.

The due date for the next calibration can be found in the menu under [System details] > [Service date].

The battery can be replaced by the technical service department at the end of the associated operating life.

For calibration, technical service, repair and disposal, please contact Christoph Miethke GmbH & Co. KG.

### Technical Service:

Christoph Miethke GmbH & Co. KG  
Ulanenweg 2 14469 Potsdam, Germany  
Tel.: +49 331 62083-0  
Fax: +49 331 62083-40  
E-mail: info@miethke.com

### WARNING

Risk on injury and/or malfunction! The product must not be modified.

Modifications of medical devices can result in a loss of guarantee or warranty claims. Christoph Miethke GmbH & Co. KG only accepts responsibility for the safety, reliability and performance of the device if:

- the device is used in accordance with the instructions for use
- new settings, modifications and repairs have only been carried out by persons authorized by us
- the electrical installation of the respective room conforms to national standards (IEC definition).

## SAFETY CHECKS (SC)

The performance of safety checks (SC) is required by the Ordinance on the Operation of Medical Devices (MPBetreibV). No SC pursuant to the MPBetreibV is performed by the manufacturer as part of the annual calibration of the BARD sensor within the scope of the technical service. The operator is obliged to perform an SC after this type of maintenance measure and before start-up.

An annual inspection with the following scope is recommended:

1. Visual inspection (also refer to the scope of the visual inspection)
2. Check the functionality based on the operating manual
3. Check error messages on the display
4. Electrical safety – measurement of the leakage currents pursuant to IEC 62353:2014
5. Creation of a test report

### Scope of the visual inspection:


1. Is the maintenance book available?
2. Do the *Reader Unit* and the antenna as well as the antenna cable, power supply unit and connection points display mechanical defects?

3. Are all labels complete and legible?
4. Are all housing screws firmly tightened?
5. Is the antenna cable firmly connected to the reader?
6. Are parts within the housing (reader and antenna) loose? – carefully shake the device
7. Are all approved components in place?

## TROUBLESHOOTING

If an error occurs, it is shown on the screen.

Example of an error message:

Press the <Info> button for additional information.

## MALFUNCTIONS DISPLAYED WITH AN ERROR MESSAGE

Message displayed	Cause	Error detection / troubleshooting
Battery flat - Auto off	Battery capacity exhausted (0%)	All data are saved after 2 min. <i>Reader Unit</i> switches off automatically. Original plug-in power unit.
Battery voltage incorrect – use original power supply unit	Battery voltage of the <i>Reader Unit</i> is too low	<i>Reader Unit</i> switches off automatically after 20 seconds. Connect the original plug-in power unit.
Low battery voltage	Battery voltage too low	Backlight switches off after 3 seconds. Connect the original plug-in power unit. No interruption to ongoing measurements.
Flashing of crossed-out antenna symbol <Info> button: Antenna faulty	antenna faulty	Deactivate and reactivate device. If the error persists, contact Technical Service.
Flashing of crossed-out antenna symbol <Info> button: Antenna not plugged in	Antenna not connected at the start of measurement - or - Antenna was disconnected during measurement	Connect antenna: measurement starts again - or - Connect antenna: measurement continues.
Flashing of crossed-out antenna symbol <Info> button: No communication	Data recording interrupted during continuous measurement (interruption of telemetric coupling)	Measurement resumes automatically once communication is restored.
Flashing of crossed-out antenna symbol <Info> button: SD card has been removed. Measurement possible	SD card was removed during continuous measurement	Insert SD card. Start measurement again.
Dataset defective! <Info> button: File cannot be opened	Validation of file failed	File cannot be opened. Try again.
Continuous key activation Keypad error	Persistent pressing of a button > 60 seconds	Release the key.
Pressure readings out of range	Pressure data measured from implant not plausible - physiologically meaningless data	Measurement continues. If the error persists, contact Technical Service.
Problem with input voltage	The voltage of the plug-in power unit is too high	<i>Reader Unit</i> switches off automatically after 20 seconds. Use original plug-in power unit.
Wrong SD card inserted! Remove card! <Info> button: Measurement without data storage possible - or - Insert SD card with correct ID XXXXXXXXXX!	Measurement starts without SD card. During measurement, an SD card not designated for use with that implant is used during measurement	Insert the matching <i>M.scio</i> SD card into the <i>Reader Unit</i> . The implant (see patient data card) is associated with the respective SD card via the identification number (ID).

Message displayed	Cause	Error detection / troubleshooting
Wrong SD card inserted! Remove card! <Info> button: Measurement without data storage possible - or - Insert SD card with correct ID XXXXXXXXXX!	SD card with different ID to the implant has been used	Insert the matching <i>M.scio</i> SD card into the <i>Reader Unit</i> . The implant (see patient passport) is associated with the respective SD card via the identification number (ID).
Wrong implant - restart measurement! <Info> button: Switch to another implant during continuous measurement not possible!	After the start of a continuous measurement, data from another implant have been received	Increase the distance between the two implants.
Device temperature outside of range <Info> button: Device temperature range from 10°C to 40°C permissible!	Temperature in the <i>Reader Unit</i> exceeds calibrated interval	The <i>Reader Unit</i> can only be used at device temperatures from 10°C to 40°C. Current measurement is stopped.
Problem with internal voltage	Internal voltage of device too high/too low	<i>Reader Unit</i> switches off automatically after 20 seconds. Contact Technical Service.
SD card faulty! <Info> button: Measurement without data storage possible!	SD card not writeable or readable (dirt, corrosion, faulty contacts)	Use of implant's internal calibration data. Data are not saved.
SD card faulty!	SD card not readable (dirt, corrosion, faulty contacts)	Check SD card for damage or dirt.
SD card inserted! Restart measurement! <Info> button: Storage of measured data possible after restart of measurement!	Measurement starts without SD card. During measurement, the designated SD card for use with that implant is inserted	Start measurement again.
SD card missing! <Info> button: Insert SD card!	No SD card inserted in data management mode	Insert SD card.
SD card missing! <Info> button: Measurement without data storage possible - or - Insert SD card with correct ID!	SD card not inserted	Insert the matching <i>M.scio</i> SD card into the <i>Reader Unit</i> . The implant (see patient data card) is associated with the respective SD card via the identification number (ID).
SD card not readable! <Info> button: Measurement without data storage possible!	Insertion of wrongly formatted or unformatted SD card	Insert the matching <i>M.scio</i> SD card into the <i>Reader Unit</i> . The implant (see patient data card) is associated with the respective SD card via the identification number (ID).

Message displayed	Cause	Error detection / troubleshooting
SD card not readable! Info> button: Measurement without data storage possible!	SD card missing - or - SD card ID not readable - or - SD card contains no calibration data	Insert the matching <i>M.scio</i> SD card into the <i>Reader Unit</i> . Measurement with data storage possible, but with limited accuracy.
SD card memory full. Measurement without data storage possible!	SD card memory full (100%) during continuous measurement	Delete any data no longer required, measurement without data storage possible.
SD card memory full	SD card memory full (100%) during continuous measurement	Delete any data no longer required.
SD card memory almost full	SD card memory almost full (99%) during continuous or fast measurement	Delete any data no longer required.
System error <Info> button: Ambient pressure not readable	Barometric pressure (BARD) not readable	Device function locked. Deactivate and reactivate device. If the error persists, contact Technical Service.
System error - incompatibility	Device hardware and software versions not compatible	Device function locked. Deactivate and reactivate device. If the error persists, contact Technical Service.
System error - antenna incompatible	Hardware version of antenna and <i>Reader Unit</i> not compatible	Device function locked. Deactivate and reactivate device. If the error persists, replace antenna and/or contact Technical Service.
System error - ID data inadmissible	Implant identification data corrupt	Measurement is stopped. Deactivate and reactivate device. If the error persists, contact Technical Service.
System error - implant voltage out of range	Implant voltage outside of permissible range	Measurement is stopped. Deactivate and reactivate device. If the error persists, contact Technical Service.
System error - calibration data incorrect	Calibration data in implant corrupted or not readable (only applies if no SD card has been inserted)	Measurement is stopped. Deactivate and reactivate device. If the error persists, contact Technical Service. Measurement with inserted SD card possible.
System error - contact Technical Service	System test detected an error	Device function locked. Deactivate and reactivate device. If the error persists, contact Technical Service.
Keypad faulty	A pressed button is registered when switching on the <i>Reader Unit</i>	Release the key. Deactivate and reactivate device.

Message displayed	Cause	Error detection / troubleshooting
Temperature increase inadmissible	A spontaneous temperature increase in the implant of more than 2 K or above 39°C	Measurement is stopped. Leave to rest for 10 minutes. If the error persists, contact Technical Service.
Temperature readings out of range	Temperature data measured from implant not plausible - physiologically meaningless data	Measurement is stopped. Deactivate and reactivate device. If the error persists, contact Technical Service.
Time out! Restart measurement!	The time between starting measurement and successfully establishing communication was exceeded (60 seconds).	Restart measurement. Optimize the distance between antenna and the implant.
Ambient pressure out of range <Info> button: Permissible ambient pressure range from 800 to 1100 mbar	The permissible barometric pressure during measurement is too low or too high.	The <i>Reader Unit</i> can only be used at barometric pressures from 800 to 1100 hPa. Current measurement is stopped.

#### OTHER MALFUNCTIONS | POSSIBLE FAULTS

Malfunction	Cause	Error detection / troubleshooting
Device cannot be activated	Battery is completely discharged	Connect the plug-in power unit. It takes approx. 6 hours to charge the battery fully. The <i>Reader Unit</i> can also be operated during the charging process (with connected plug-in power unit). Note: Charging is not possible at an ambient temperature of more than 35°C.
Device switches off	Unfavourable operating conditions (e.g. low humidity or unsuitable floor coverings)	Switch device back on. If the error persists, contact Technical Service.

## TECHNICAL DATA

Description	Specifications and standards
Voltage range <i>Reader Unit</i> Plug-in power unit	6 V (DC) 100-240 V (50–60 Hz)
Power consumption <i>Reader Unit</i> Plug-in power unit	1.4 A (DC) 0.15-0.3A (50–60 Hz)
Operating frequency	133 kHz
Bandwidth	125 kHz – 135 kHz
Type of modulation	Amplitude modulation
Transmitting power of the antenna	max. 0.8W
Pressure measuring range	-66.67 hPa to +133.33 hPa
Precision	±5.32 hPa
Temperature range for pressure measurement: <i>Reader Unit</i> <i>M.scio</i>	10°C to 40°C 20°C to 39°C
Functional distance between the antenna and <i>M.scio</i>	10 to 30 mm
Battery life	5 years minimum
Battery charging cycle	at least 250 charging cycles
Battery self discharge	Remaining charge after 3 months (storage) > 70%
Flammability category of housing	UL 94 HB
Moisture protection or leak-tightness: <i>Reader Unit</i> Antenna plug-in power unit	IP44 IP44 IP40
Durability – impact test	pursuant to IEC 60601-1: 2012 15.3.2
Durability – drop test	pursuant to IEC 60601-1: 2012 15.3.4.1
Weight: <i>Reader Unit</i> Antenna Plug-in power unit	0.600 kg 0.215 kg 0.127 kg

Designation	Values and standards
Dimensions (W x H x D): <i>Reader Unit</i>	144 x 270 x 65 mm
Antenna (without cable)	100 x 250 x 25 mm
Plug-in power unit	77.5 x 31.5 x 41 mm
Protection class (plug-in power unit)	II
Conformity to standards	IEC 60601-1:2012
EMC	IEC 60601-1-2:2014

## AMBIENT CONDITIONS

Operating conditions	
Relative humidity	30 % to 75 %
Ambient temperature	10°C to 40°C, no condensation
Atmospheric pressure	800 to 1100 hPa

Transport conditions	
Relative humidity	15 % to 95 %
Ambient temperature	0°C to 50°C
Atmospheric pressure	500 hPa to 1100 hPa

Storage conditions	
Relative humidity	15 % to 95 %
Ambient temperature	10°C to 40°C
Atmospheric pressure	500 hPa to 1100 hPa

## CE MARKING

Approval for affixing the CE mark for active implantable medical devices in accordance with Directive 90/385/EEC) was initially granted in 2011.

## STORAGE AND TRANSPORT

- Store reprocessed products protected from dust in a dry, dark and cool room with as little bacterial contamination as possible.
- Ship product only in original packaging.
- For storage and transport conditions, see chapter on "Ambient conditions".

## NOTE

The product may be damaged if used too soon after storage/transport at temperatures below 10°C.

- Allow the *Reader Unit* to acclimatize for approx. three hours at room temperature.

## DISPOSAL

Adhere to national regulations when disposing of or recycling the product, its components or its packaging!

Products carrying this symbol have to be disposed of in the separate collection of electrical and electronic devices. Within the European Union, disposal is arranged free of charge by the manufacturer.

For recycling, please return the *Reader Unit* to the manufacturer.

If you have any questions regarding the disposal of the product, please contact Christoph Miethke GmbH & Co. KG; see also chapter on "Calibrating the barometric pressure sensor | Technical service".

## NOTES ON ELECTROMAGNETIC COMPATIBILITY

The device satisfies the requirements of IEC 60601-1-2:2014.

**ELECTROMAGNETIC EMISSIONS****CAUTION**

The device is only intended for use at professional healthcare facilities!

While the device complies with class B HF emissions in accordance with CISPR 11, it may only be used at professional healthcare facilities.

Guidelines and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in an electromagnetic environment as specified below. It is the user's responsibility to ensure that the device is operated in such an environment.		
Electromagnetic interference measurements	Compliance	Guidelines on the electromagnetic environment
RF emissions according to CISPR 11	Complies with Group 1	The device uses RF energy only for its internal function. Its RF emissions are thus very low and not likely to cause interference in nearby electronic devices.
RF emissions according to CISPR 11	Complies with Class B	The device is suitable for use in all establishments including domestic settings and areas immediately connected to a public power supply which also supplies buildings used for residential purposes.
Harmonic emissions according to IEC 61000-3-2	Complies with Class A (IEC 61000-3-2)	-
Voltage fluctuations/flicker emissions according to IEC 61000-3-3	Complies	-

**ELECTROMAGNETIC IMMUNITY****WARNING**

Creation of electromagnetic interferences!

Portable HF communication devices (radio devices) (including their accessories, such as antenna cables and external antennas) should not be used closer than a distance of 30 cm (or 12 inches) to the reader and the antenna. If electromagnetic interferences with spot frequencies of 385 MHz or 450 MHz arise, a distance of at least 80 cm must be ensured. Non-compliance may lead to a reduction of the device's performance characteristics. Electromagnetic interferences can lead to a device shutdown. In this case, the device must be restarted and the measurement repeated.

**MEDICAL DEVICES CONSULTANTS**

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical product consultants as contacts for all product-related questions.

You can reach our medical devices consultants at:

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Christoph Miethke GmbH & Co KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany