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

Our Strong Partner in Neurosurgery:

**B | BRAUN**  
SHARING EXPERTISE

www.braun.com

Adjustable DP (differential pressure) -unit		Gravitational unit	
Standard (children and NPH-patients)	5 cmH <sub>2</sub> O	Children up to 5 years	20 cmH <sub>2</sub> O
Defensive (e.g. patients with extremely wide ventricles and highly elevated ICP or aqueductal stenosis)	10 cmH <sub>2</sub> O	Children (over 5 years) & adults (up to 60 years)	25 cmH <sub>2</sub> O
Special (e.g. patients with pseudotumor cerebri)	15 cmH <sub>2</sub> O	Adults over 60 years	20 cmH <sub>2</sub> O

height in cm    Body-Mass-Index (BMI)    weight in kg

### Overweight

The peritoneal pressure inhibits drainage. Therefore the use of the gravitational unit should be considered for overweight patients as a function of body-mass-index (BMI):

25 - 29 BMI - 5 cmH<sub>2</sub>O below standard recommendation

30 - 34 BMI - 10 cmH<sub>2</sub>O below standard recommendation

35 - 39 BMI - 15 cmH<sub>2</sub>O below standard recommendation

> 40 BMI only apply the adjustable DP-unit (no gravitational unit necessary)

### Mobility

Standard pressure levels are suitable for active people. Bedridden patients should not be treated with a gravitational unit.

### Height

The hydrostatic suction effect normally depends on the height. Therefore we recommend the following corrections for the gravitational unit:

< 1,60 m height: - 5 cmH<sub>2</sub>O

> 1,80 m height: + 5 cmH<sub>2</sub>O

### Control

The adjustable DP-unit is adjusted to 5 cmH<sub>2</sub>O by the manufacturer. Please control the pressure setting before use.

### Body position

Only the adjustable DP-unit determines the pressure level in the horizontal position. The analogous pressure level for the vertical body position results from the addition of both - adjustable DP-unit and gravitational unit.

### Adjustment

Changing the adjustable DP-unit also influences the opening pressure level in the vertical body position.

All indicated pressure levels are in cmH<sub>2</sub>O. This is a non-binding recommendation for the attending physician. According to his diagnosis, the physician decides each case independently, without instructions and individually. The stated values consider the current scientific knowledge up to June 2018.

The MIETHKE App demonstrates the functionality of proGAV and the effect of patient posture on valve opening pressure. Now available in the App and Google Play Store.

Our Strong Partner in Neurosurgery: