M.blue[®]

THE BALANCED WAY OF LIFE
INSPIRED BY YOU

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 www.miethke.com

US This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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DIMENSIONEN

TABLE OF CONTENTS

INDICATION	16
TECHNICAL DESCRIPTION	16
FUNCTION OF THE VALVE	16
SELECTION OF APPROPRIATE PRESSURE LEVEL	17
PRESSURE RATING IDENTIFICATION IN X-RAY IMAGES	18
APPLICATION OF M.BLUE PLUS INSTRUMENTS	19
POSSIBLE SHUNT COMPONENTS	23
TUBE SYSTEMS	23
IMPLANTATION	23
VALVE TEST	24
PRESSURE-FLOW CHARACTERISTICS	25
PRECAUTIONARY MEASURES AND CONTRAINDICATIONS	26
FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	27
ADVERSE REACTIONS AND INTERACTIONS	27
STERILISATION	27
REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)	27
MEDICAL DEVICES CONSULTANTS	27
DIMENSIONS	28


Fig. 1: Cross section of *M.blue*

- 1. Adjustable gravitational unit
- 2. Sapphire ball
- 3. Tantalum weight
- 4. Rotor

- 5. Differential pressure unit
- 6. Sapphire ball
- 7. Micro-coil spring

INDICATION

M.blue is used for cerebrospinal fluid (CSF) drainage in the treatment of hydrocephalus.

TECHNICAL DESCRIPTION

M.blue is a valve made from titanium. It consists of an adjustable gravitational unit and a differential pressure unit (fig. 1). The adjustable gravitational unit (1) in the proximal part of the valve contains a tantalum weight (3), which holds a sapphire ball in the ball seat via a lever (2). Depending on the body position of the patient, the influence of the tantalum weight on the sapphire ball and thus the valve opening pressure changes. Via a rotor (4), the pretension of the torsion spring connected to the lever can be changed through the skin.

This way, the influence of the tantalum weight on the sapphire ball can be influenced and thus the valve opening pressure can be adjusted. In the distal part of the valve, a micro-coil spring (7) controls the opening pressure of the differential pressure unit (5). A sapphire ball (6) ensures precise opening and closing of the ball-cone unit.

FUNCTION OF THE VALVE

M.blue is a posture-dependent hydrocephalus valve. The opening pressure for *M.blue* consists of the opening pressures for the adjustable gravitational unit and the differential pressure unit combined.

Horizontal position

In the horizontal position, the gravitational unit is always open and does not offer any resistance.

Consequently, the opening pressure of *M.blue* in the horizontal position is characterised by the fixed differential pressure unit pre selected at implantation to be either 0, 5 or 10 cm H₂O. It is not adjustable. The main operation method used by the adjustable differential pressure unit is shown in Fig. 2a and b.


Fig. 2: Operation of the differential pressure unit
a) closed b) open

In fig. 2a the valve is closed, hence drainage is not possible.

If the patient's intraventricular pressure (IVP) exceeds the spring force of the micro-coil spring, the sealing ball moves out of the cone, leaving a gap for CSF drainage (fig. 2b).

Vertical position

As the patient's body moves from prone to vertical, the gravitational unit closes the discharge channel in the proximal part of the valve (Fig. 3a). Thus, the opening pressure of *M.blue* is increased in the upright position, because now the weight of the tantalum weight (opening pressure of the gravitational unit) must be overcome in addition to the opening pressure of the differential pressure unit. Drainage is only possible once the sum of IVP and hydrostatic suction is greater than the opening pressure of both units (fig. 3b). For individual adaptation of the opening pressure to the patient, a valve opening pressure between 0 and 40 cm H₂O can be selected for the adjustable gravitational unit.


Fig. 3: Gravitational unit in vertical body position
a) closed b) open

During physical activity which is associated with shock (e.g. jogging) the opening pressure of *M.blue* may decrease temporarily according to laboratory results. Generally however, functionality is retained. At the end of physical activity, the opening pressure returns to its original level and remains stable.

SELECTION OF APPROPRIATE PRESSURE LEVEL

Horizontal position

In the lying position, the gravitational unit has no influence on the valve opening pressure. In this position, the valve opening pressure is thus determined exclusively by the differential pressure unit. In this case, the pressure level should be set in accordance with the clinical picture and indications. By default, a differential pressure unit with an opening pressure of 5 cmH₂O is recommended.

Vertical position

The *M.blue* opening pressure for the vertical body position is calculated from the sum of the opening pressure in the differential pressure unit and the adjustable gravitational unit. Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure for the gravitational unit (see pressure level recommendations at <https://www.miethke.com/en/products/downloads>). This is a non-binding recommendation for the attending physician. According to his diagnosis, the physician decides each case independently, without instructions and individually.

PRESSURE RATING IDENTIFICATION IN X-RAY IMAGES

The selected pressure level for *M.blue* should always be monitored using the *M.blue plus Compass*, but it can also be checked using an x-ray image. The rotor setting is decisive in this case. The four magnets in the rotor can be seen on the x-ray image as white points and are located opposite each other in pairs. Two additional burrholes (right and left next to the magnet pairs) on one side of the rotor can be used as orientation. They can be seen as black points on the x-ray image. This side can be designated as the rotor rear side. The two front magnets are on the opposite side. The space between these two magnets can be considered as the triangle apex. The pressure level can be read off using the orientation of this intermediate space (fig. 4). The triangle apex can take up any position except the space labelled as a non-adjustable area in fig. 4. This means that the opening pressure of *M.blue* can be variably adjusted from 0 up to 40 cm H₂O. In order not to read off the pressure level in reverse, the top view of the implanted valve in the housing ring has a recess with the tantalum weight to the right of the inlet nozzle (fig. 4).


Fig. 4: Schematic representation of rotor in x-ray image

1. Coding burrhole gravitational unit
2. Coding burrhole differential pressure unit
3. Triangle tip
4. Nonadjustable range
5. proximal
6. distal


Fig. 5: X-ray image (adjustable gravitational unit preset to 20 cmH₂O; differential pressure unit: 5 cmH₂O)

The pressure level of the differential pressure unit can be recognised in the X-ray image by the encoding (fig. 5). The following pressure levels are possible for the differential pressure unit:


Fig. 6: Pressure ratings encodings of the differential pressure unit

1. Differential pressure unit
2. Coding burrhole

APPLICATION OF *M.blue plus* Instruments

The *M.blue plus* Instruments may only be used by trained specialist personnel. The selected pressure level of *M.blue* can be determined, adjusted and monitored using the *M.blue plus* Instruments. The *M.blue plus* Compass is used for localising and reading the adjustable gravitational unit of *M.blue*.


Fig. 7: *M.blue plus* Compass

a) open b) closed

1. Scale Ring 2. Float gauge compass

With the *M.blue plus* Adjustment Ring the opening pressure of *M.blue* can be variably adjusted from 0 to 40 cmH₂O.

Fig. 8: *M.blue plus* Adjustment Ring

The opening pressure for the adjustable gravitational unit can be changed before or after implantation. It is preset by the manufacturer to 20 cmH₂O. Please carry out the following steps to adjust the valve:

1. Localisation

When the *M.blue plus* Compass is opened, a circular template becomes visible with which you can localise the valve in the patient's head using your index finger (fig. 9).

Fig. 9: Localising the valve using the M.blue plus Compass

The direction markings on the template show the flow direction.

2. Test procedure

In order to determine the selected pressure level, the compass is then closed again. The float gauge should now be centred by moving the instrument in the designated circular marking (fig. 10). Once the float gauge is centred, the currently set opening pressure of the gravitational unit can be read off via the line mark on the float gauge (fig. 10).


Fig. 10: Determination of the pressure level with the M.blue plus Compass

There are two scales on the scale ring. The opening pressure of the M.blue gravitational unit is set to the blue marked adjustment range of 0 to 40 cmH₂O of the inner scale.


Fig. 11: 1. M.blue scale
Inside: M.blue scale of 0 to 40 cmH₂O (opening pressure of the gravitational unit in the picture example is 16 cmH₂O)

CAUTION

The M.blue plus Compass should be placed centrally over the valve, if possible; otherwise incorrect determination of the opening pressure may occur.

M.blue plus Compass reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the M.blue plus Adjustment Ring should not be in the immediate vicinity of the M.blue plus Compass when the opening pressure is being determined. We recommend a minimum distance of 30 cm.

NOTE

Possible air bubbles inside the Compass do not affect its functionality.

3. Adjustment process

In order to adjust the opening pressure, the compass is opened, but without changing the position of the scale ring. The adjustment ring is now inserted into the scale ring in such a manner that the line marking points to the desired value on the scale of the scale ring (fig. 12).


Fig. 12: Insertion of the Adjustment Ring
1. Adjustment Ring 2. Scale Ring

The opening pressure of the *M.blue* gravitational unit uses the adjustment range of 0 to 40 cmH₂O of the blue inner scale.


Fig. 13: Adjustment of the gravitational unit of *M.blue*
(in the picture example to 32 cmH₂O)

By applying slight pressure with the index finger to the valve diaphragm located in the centre of the adjustment ring and under the skin, the rotor brake is released and the opening pressure of *M.blue* is changed to the desired value (fig. 14).


Fig. 14: Adjustment using the
M.blue plus Adjustment Ring

M.blue is equipped with a feedback mechanism. If targeted pressure is exerted on the valve, an acoustic signal (a clicking sound) is audible, or a resistance can be felt to give, as soon as the rotor brake has been released due to the valve housing design. In other words the valve shows both acoustically and haptically when the pressure is sufficient for uncoupling of the brake. Once this pressure has been released the rotor is adjustment-proof again. Although the click caused by releasing the rotor brake is easily audible before implantation, this can be considerably reduced after implantation. This can be due to the CSF filling the valve, its position and the condition of the implants surroundings. Normally, however, it should be audible to the patients themselves or through the use of a stethoscope.

Adjustment with the *M.blue plus* Adjustment Assistant

Alternatively, the *M.blue plus* Adjustment Assistant can be used to adjust the opening pressure. To do this, insert the *M.blue plus* Adjustment Assistant into the adjustment ring aligned to the desired value and press it with your index finger (fig. 15).

Fig. 15: *M.blue plus* Adjustment Assistant

WARNING

When adjusting, please ensure that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process; otherwise errors may result.

Example: The opening pressure is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 6 to 22 cmH₂O and subsequently from 22 to 36 cmH₂O.

CAUTION

The *M.blue plus Adjustment Ring* emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.

4. Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the selected opening pressure. To do this, please proceed as in Points 1 and 2. If the measured pressure does not agree with the required pressure stage the adjustment process should be repeated. To do this, please start at Point 3 again. Adjustment can be made difficult by swelling of the skin for a few days after the operation! If checking the valve setting is not 100% possible using the *M.blue plus Compass*, we recommend carrying out a check using an imaging method.

M.blue Check-mate

The *M.blue Check-mate* is supplied sterile and can be re-sterilised. With the *M.blue Check-mate* it is possible to carry out a pressure level change and the check before and during valve implantation directly on *M.blue*. To determine the pressure level the *M.blue Check-mate* is placed centrally on *M.blue*. The *M.blue Check-mate* automatically aligns itself over the valve. The pressure stage can be read off from the catheter in the proximal (leading to valve) direction. If the pressure level is supposed to be adjusted, the *M.blue Check-mate* is placed centrally on *M.blue*. When doing so, the required pressure stage must point towards the proximal catheter (leading to valve). If the *M.blue* adjustment gyroscope is now pressed lightly onto the valve, the rotor brake in *M.blue* is released and the pressure level is set. When adjusting, please ensure that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process otherwise errors can result (see Chapter "3. Adjustment process").

Fig. 16: *M.blue Check-mate*, colour: blue, pressure levels 0 to 40 cmH₂O

CAUTION

Because of the magnets inside the *M.blue plus Instruments*, *M.blue plus Instruments* may not be used in the vicinity of active implants such as heart pacemakers. Furthermore, there is risk of damage to MRI instruments in their vicinity. Therefore, usage of *M.blue plus Instruments* is forbidden in such locations!

It is absolutely essential to use only the *M.blue plus Instruments* to determine, adjust and monitor the opening pressure of *M.blue*.


Fig. 17: *M.blue plus Instruments*

M.blue with proGAV 2.0 (M.blue plus)

The instruments can also be used for localising, reading and setting the adjustable differential pressure unit of *proGAV 2.0*. When combining the *M.blue* valve with the adjustable differential pressure unit of the *proGAV 2.0*, the *M.blue* valve is localised, checked and adjusted as described under points 1-4. The adjustable differential pressure unit (*proGAV 2.0*) can also be localised, checked and adjusted to a desired value between 0 and 20 cmH₂O with the *M.blue plus Instruments*, as described under points 1-4. For the opening pressure of the adjustable *proGAV 2.0* differential pressure unit, the grey adjustment range of 0 to 20 cmH₂O on the outer scale of the scale ring applies (fig. 18).


Fig. 18: 1. proGAV 2.0 scale

Outside: proGAV 2.0 scale of 0 to 20 cmH₂O (opening pressure of the adjustable differential pressure unit (proGAV 2.0) in the picture example is 17 cmH₂O)

POSSIBLE SHUNT COMPONENTS

M.blue can be ordered as a shunt system in a range of configurations. These configurations can be combined with the accessories presented below. In each case, versions for paediatric hydrocephalus and for hydrocephalus in adults are available.

Reservoirs

The use of a reservoir in combination with shunt systems provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to an integrated check valve in the *SPRUNG RESERVOIR* and the *CONTROL RESERVOIR*, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the drainage system as well as the valve catheter. During the pump action, access to the ventricular catheter is closed. The use of reservoirs does not increase the opening pressure of the shunt system. Puncturing the reservoir should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. A stable titanium floor prevents the bottom surface from being pierced. 30 punctures are possible without any restrictions.

WARNING

Frequent pumping can result in excessive drainage and thus lead to pressure conditions outside the normal physiological range. The patient should be informed about this risk.

Burrhole deflector

Because of the tight fit on the ventricular catheter, the burrhole deflector makes it possible to choose the length of catheter penetrating into the skull prior to implantation. The ventricular catheter is deflected at a right angle in the burrhole (see chapter "Implantation").

TUBE SYSTEMS

M.blue can be ordered as an individual valve unit or as a shunt system with integrated catheters (interior diameter 1.2 mm, exterior diameter 2.5 mm). The supplied catheters do not fundamentally change the pressure-flow characteristics. If catheters by other manufacturers are used, a tight fit must be ensured. In any case, catheters have to be carefully fixed with a ligature to the valve's titanium connectors.

IMPLANTATION

Positioning of the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The required skin incision should be made in form of a lobule pedicled towards the draining catheter. If a burrhole deflector is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burrhole.

M.blue is available in a range of different configurations: If a *Burrhole Reservoir* or a *SPRUNG RESERVOIR* is used, then the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the burrhole reservoir connected, with the connection secured with a ligature. When using a shunt system with a *CONTROL RESERVOIR*, a burrhole deflector is included. The deflector is used for adjusting the length of catheter to be implanted and for its positioning inside the ventricle. The ventricular catheter is deflected, and the *CONTROL RESERVOIR* is put into place. The position of the ventricular catheter should be inspected after the procedure by imaging (such as CRT or MRI).

Positioning of the valve

The adjustable gravitational unit in *M.blue* is set to an opening pressure of 20 cmH₂O upon delivery. This opening pressure can be changed to a different pressure before implantation. A location behind the ear is suitable as an implantation position, whereby the implantation height has no influence on the valve function. The adjustable valve should be touching the bone or the periosteum since pressure must be exerted on the valve during any later adjustment. A large arch-shaped or a small straight skin incision with a pocket for the valve should be made. The catheter is then pushed forward from the burrhole to the selected valve implantation location, shortened if necessary, and secured to the *M.blue* with a ligation. The valve should not be located directly under the skin incision. The valve unit has an arrow in the flow direction (arrow towards distal or downwards). The embossed blue surface of the valve with the arrow markings points to the outside.

NOTE

M.blue is position-dependent. For that reason, care must be taken to implant the valve parallel to the body axis.

Therefore, if a shunt system in which the valve has been pre-fitted with a *Burrhole Reservoir* is being used, only the occipital access should be used.

WARNING

The adjustable valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e.g. underneath heavily scarred tissue).

CAUTION

The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

Positioning of the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. For example, it can be applied paraumbilical or at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. The recommendation is to pull the peritoneal catheter using a subcutaneous tunnelling tool from the valve to the intended position, if necessary with the aid an auxiliary incision. The peritoneal catheter that is usually securely attached to *M.blue* has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the peritoneal catheter (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

VALVE TEST

Preoperative valve test

M.blue should be vented before implantation and checked for permeability. The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological salt solution. The valve is continuous if saline solution can be extracted (fig. 19).

WARNING

Contamination in the solution used for testing can impair the product's performance.


Fig. 19: Patency test

**WARNING**

Pressure admission through the single-use syringe should be avoided both at the proximal and the distal end (fig. 20).






Fig. 20: Avoidance of pressurisation

Postoperative valve test

The *M.blue* has been constructed as a reliably functioning unit without pump or test function. The valve test can be performed by flushing, pressure measurement or pumping.

PRESSURE-FLOW CHARACTERISTICS**Horizontal position**

The pressure flow characteristics for the *M.blue* differential pressure unit are shown below for pressure levels 0, 5, 10 and 15 in the horizontal valve position.

Fig. 21: Pressure in (cmH₂O); flow rate in (ml/h)

Pressure-flow characteristics of selected pressure levels of *M.blue* in the horizontal body position

Vertical position

In the vertical body position, the *M.blue* opening pressure is based on the setting of the differential pressure unit and the adjustable gravitational unit. The pressure flow characteristics for various pressure stage settings in the vertical body position are shown below.





Fig. 22: Pressure in (cmH₂O); flow rate in (ml/h)
Pressure-flow characteristics of selected pressure levels of *M.blue* in the vertical body position:

The opening pressure refers to a reference flow of 20 ml/h.

PRECAUTIONARY MEASURES AND CONTRAINDICATIONS

Patients must be carefully monitored after implantation. Reddening of skin or tightness in the area of the drained tissue may be indications of infections at the shunt system. Symptoms such as headache, dizziness, confusion or vomiting often occur in conjunction with shunt dysfunction. These symptoms and a leakage within the shunt system require the immediate replacement of the affected shunt component or the entire shunt system.

The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e.g. meningitis, ventriculitis, peritonitis, bacteriaemia, septicaemia) in the region affected by the implantation.

FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. However, no guarantee can be given that these medical devices may not require replacement for medical or technical reasons. These medical devices are able to resist positive and negative pressures up to 200 cmH₂O during and after implantation. These medical devices have to be stored in a clean and dry environment at all times. Nuclear magnetic resonance (MRI) examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. The valve is MR Conditional. The supplied catheters are MR Safe. Reservoirs, deflectors and connectors are MR Conditional.

WARNING

If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, it is not possible to rule out valve adjustment. In MRI imaging *M.blue* creates artefacts which are larger than the valve itself.

CAUTION

For people using cardiac pacemakers: It is possible that the function of the pacemaker is influenced by the implantation of a *M.blue*.

ADVERSE REACTIONS AND INTERACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): Infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in very rare cases noise development. Violent shocks for the outside (accident, fall) may put the integrity of the shunt system at risk.

STERILISATION

The products are sterilised with steam under strictly controlled conditions. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used in any circumstances. No guarantee can be given for the functional safety and reliability of resterilised products.

REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the whereabouts of medical devices used in humans, especially for implants. The individual identification number of the implanted valve should therefore be recorded in the patient's medical records and patient passport to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website (<https://www.miethke.com/en/products/downloads>).

MEDICAL DEVICES CONSULTANTS

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical products consultants as contacts for all product-related questions.

You can reach our medical devices consultants at:

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Christoph Miethke GmbH & Co KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany