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RECOMMENDED PRESSURE SETTINGS

Patient	Selection of pressure levels		Opening pressure in total	
	<p>adjustable differential pressure unit gravitational unit</p>			
Newborn	5	20	5	25
Children from 5 years	10	25	10	35
Adults	5	25	5	30
< 1.60 m	5	20	5	25
> 1.80 m	5	30	5	35
Adults from 65 years	5	20	5	25
< 1.60 m	5	15	5	20
> 1.80 m	5	25	5	30

The *proGAV 2.0* consists of an adjustable differential pressure unit and a gravitational unit.

The opening pressure of *proGAV 2.0* in the horizontal position is **only** determined by the differential pressure unit.

In the lying position the gravitational unit is not active.

The opening pressure in the vertical position is **the sum** of both: adjustable differential pressure unit and gravitational unit.

The Miethke App demonstrates the functionality of *proGAV 2.0* and the effect of patient posture on valve opening pressure.

Now available in the App and Google Play Store.

Our Strong Partner in Neurosurgery:

www.bbraun.com

All indicated pressure levels are in cmH₂O. This is a non-binding recommendation for the attending physician. According to his diagnosis, the physician decides each case independently, without instructions and individually. The stated values consider the current scientific knowledge up to June 2018.