Liegende Position

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
Horizontal position

In the horizontal position the *SHUNTASSISTANT* doesn't have any influence on the opening pressure. Only the chosen valve (fixed or programmable) determines the opening pressure of the whole shunt system.

$$\text{opening pressure of the whole shunt system} = \text{opening pressure of chosen valve} + 0 \text{ cmH}_2\text{O}$$

Vertical position

In the vertical position the opening pressure of the complete shunt system is the sum of the opening pressure of the chosen valve and the opening pressure of the *SHUNTASSISTANT*. In many cases this recommendation was useful:

$$\text{opening pressure of the whole shunt system} = \text{opening pressure of chosen valve} + \begin{cases} 20 \text{ cmH}_2\text{O} & \text{- for children up to 5 years} \\ 25 \text{ cmH}_2\text{O} & \text{- for children over 5 years} \\ & \text{adults up to 60 years} \\ 20 \text{ cmH}_2\text{O} & \text{- for adults over 60 years} \end{cases}$$

Overweight	The peritoneal pressure inhibits drainage. Therefore the use of the gravitational unit should be considered for overweight patients as a function of body-mass-index (BMI): 25 - 29 BMI - 5 cmH ₂ O below standard recommendation 30 - 34 BMI - 10 cmH ₂ O below standard recommendation 35 - 39 BMI - 15 cmH ₂ O below standard recommendation > 40 BMI no <i>SHUNTASSISTANT</i> necessary
Mobility	Standard pressure levels are suitable for active people. Bedridden patients should not be treated with a <i>SHUNTASSISTANT</i> .
Height	The hydrostatic suction effect normally depends on the height. Therefore we recommend the following corrections for the <i>SHUNTASSISTANT</i> : < 1,60 m height: - 5 cmH ₂ O > 1,80 m height: + 5 cmH ₂ O

The recommendations are based on common patient treatments, but can vary depending on the individual patient's condition.

CHRISTOPH MIETHKE GMBH & CO. KG

Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
 Phone +49 (0) 331 62083-0 | Fax: +49 (0) 331 62083-40 | www.miethke.com

Distributed by:

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
 Phone +49 (0) 74 61 95-0 | Fax: +49 (0) 74 61 95-26 00 | www.aesculap.com