M.scio[®]

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
Jan Mügel

Thoralf Knitter

"Everyone has a cross to bear - some more than others. The most important thing is not to give up and to live your life to its fullest."

Sarah, nurse and hydrocephalus patient

Foreword

Hydrocephalus is a disease that is frequently connected to the widest range of other existing problems. Some say that hydrocephalus should be considered a symptom of another underlying disorder.

Successful treatment of hydrocephalus only became possible in the 1950s. In a dramatic race for the life of his son Casey who had hydrocephalus, Philadelphia technician John D. Holder took only a few weeks to develop a silicone valve. The implantation of this valve in March 1956 marked a great step forward in the clinical treatment of this disease.

Christoph Miethke GmbH & Co. KG has built on the insights from over 50 years of valve treatment, and with the development of gravitation technology pioneered the development of valves that consistently consider the physical framework conditions in the drainage of cerebrospinal fluid, thus ensuring an intraventricular pressure in all body positions that matches that of a healthy person as closely as possible. These days, the primary aim of hydrocephalus therapy is no longer to prevent death or severe impairment as the result of the condition, but to enable patients to lead a life that is as normal as possible.

M.scio helps your physician perform a function check of this kind of valve or shunt system.

The purpose of this manual is to provide you and your family with a small insight into the treatment of hydrocephalus with the *M.scio*.

TABLE OF CONTENTS

Foreword	15
Basics	18
Hydrocephalus	18
Scope of application of the <i>M.scio</i>	20
Materials of M.scio	21
Patient data card and SD card	22
Treatment complications	22
Warnings	22
Conduct after the operation	23
About the company	25

Basics

The human brain is surrounded by a special fluid (scientific term: cerebrospinal fluid - CSF). Inside the head are several cerebral ventricles (scientific term: ventricles) where CSF is produced. These ventricles are interconnected with channels and constitute a complex drainage system.

CSF has the task of protecting the brain from mechanical damage, regulating intracranial pressure, maintaining moisture levels of the brain tissue and transporting metabolic products.

Fresh CSF is produced in the ventricles each day which is absorbed by the venous blood system. A healthy balance between the production and absorption of CSF is thus very important.

Hydrocephalus

Hydrocephalus means that the balance between production and absorption of CSF is impaired. If more CSF is being produced than can be drained, the cerebral ventricles expand and a hydrocephalus develops ("hydro": Greek for water; "cephale": Greek for head). The ventricles exert pressure on the surrounding brain tissue and can cause neurological damage, some of which irreversible.

Symptoms

Symptoms depend on the level of impairment and include headache, nausea, vomiting, coordination disorders, sleepiness up to the loss of consciousness. For children under the age of two, head circumference can increase significantly, because the bones in their skulls have not yet fused.


The causes of hydrocephalus are very varied and occur at any stage in life.

Treatment methods

In a hydrocephalus, intracranial pressure has to be reduced to within normal limits. In order to achieve controlled drainage of the cerebrospinal fluid, a valve and a tube system (together referred to as the shunt system) is implanted. This establishes a connection between the cerebral ventricles and another bodily cavity, generally the abdomen, in order to drain excess CSF from the ventricles.

There are three types of shunt implantation:


- 1 Ventriculo-atrial (from the ventricle to the right atrium of the heart)
- 2 Ventriculo-peritoneal (from the ventricle to the abdominal cavity)
- 3 Lumbo-peritoneal (from the spinal canal to the abdominal cavity)

Even though there have been consistent efforts to find alternative treatment options to the implantation of a valve, e.g. by drug-based treatment or recently by minimally invasive surgical procedures, to this day, in many cases, there is not really an alternative to the implantation of a drainage system, known as a "shunt system".

Scope of application of the *M.scio*


The *M.scio* enables the non-invasive control of the shunt function. To be able to draw conclusions about the proper functioning of the shunt system, the *M.scio* measures the pressure as well as pressure changes in the system. Measurement takes place with the help of the measuring cell, which is located within the *M.scio*. The treating physician can determine the relative pressure values using a *Reader Unit Set* and draw conclusions regarding the function. The *M.scio* is available in the following variants:

Besides pressure measurements and shunt function control, variants a) and c) also perform the functions of a conventional reservoir: Both provide the option of tapping the membrane to measure the pressure in the cerebral ventricles and injecting medication.

Function of the *M.scio*

The pressure in the *M.scio* can be measured by wireless communication with the antenna of the *Reader Unit Set* by the treating physician, which is completely painless.

Pressure measurement in the shunt system with M.scio and Reader Unit Set

The interpretation of the relative curve of the CSF pressure allows a shunt blockage as well as a mechanical loss of function of the shunt system to be identified, localised and evaluated non-invasively.

Materials of *M.scio*

The *M.scio* consists of high-quality materials that have been proven and standardised for use in implants; the main components are titanium and PEEK. The stable housing ensures a reduction of negative outside influences that could affect the function of the *M.scio*. This results in a high level of functional safety and thus a long operating life.

The catheters are made from silicone and are not made with natural rubber latex.

Patient data card and SD card

Each *M.scio* is accompanied by a patient data card. This is completed by the treating physician and thus contains important information for follow-up examinations. In addition to the patient data card, the patient is given an SD card on which all individual data about the implant are stored. If the SD card is lost, it can be reordered by the treating physician by specifying the serial number of the SD card or the ASIC ID.

Treatment complications

The treatment of hydrocephalus with a shunt system is not always free from complications. As with every surgical intervention, there is a risk of infection. Unfortunately, problems can occasionally develop that could be directly or indirectly linked to the implanted shunt system. Such complications include blockages of the shunt system or undesired excessive or insufficient drainage of cerebrospinal fluid (over- or underdrainage).

Immediately consult a physician if the patient suffers from skin rashes and tensions, severe headaches, dizzy spells, unnatural gait or the like.

Warnings

Violent shocks from the outside (accident, fall) may put the integrity of the shunt system and *M.scio* at risk.

The use of energy-emitting devices such as defibrillators and HF devices can cause the measuring cell to fail.

In cases where an electrical current from an external source comes into contact with the body, the measuring cell can sustain damage.

The use of radiation therapy and radionuclide imaging procedures can cause the measuring cell to fail.

During therapeutic ultrasound procedures, there is a risk of the concentration of the ultrasound field being unintentionally increased and thus injuring the patient.

Particular attention must be paid to the following for variants a) and c):

As too frequent pumping of this implant can result in excessive drainage (overdrainage) and thus unfavourable pressure conditions, this process should exclusively be performed by a doctor.

Conduct after the operation

Patients supplied with an *M.scio* are normally not restricted in their daily lives. You can safely travel by plane, do sports or go to work with the *M.scio*. Only for the period immediately after the surgery, you should discuss with your doctor to what extent you should avoid physical exertion and for how long that would be advisable in your specific case. However, high pressure and shock loads must be avoided.

The following applies for any MRI follow-up examinations:

The *M.scio* is MR Conditional up to 3 Tesla.

About the company

Christoph Miethke GmbH & Co. KG, based in Potsdam, has from 1992 concentrated on the development, production and distribution of innovative neurosurgical implants for the treatment of hydrocephalus.

Our daily aim is through the innovative development and production of high-quality reliable products to improve the life of hydrocephalus patients. One of our foremost task hereby in cooperation with hospitals, physicians and patients is primarily this: to understand. Again and again, we change perspectives and meet with our colleagues, partners and patients on an equal footing.

For detailed information about our company, please visit our website: www.miethke.com

More about the operating principle of our valves is available here: www.miethke.com/en/products

Or use our app:

Apple:

<https://itunes.apple.com/de/app/miethke/id450290015?mt=8>

Android:

<https://play.google.com/store/apps/details?id=com.miethke.graviton>


You are also welcome to contact our Medical Devices Consultants: +49 331 620 83-0 // info@miethke.com // www.miethke.com

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany