# *Burrhole Reservoir*<sup>®</sup> *SPRUNG RESERVOIR*<sup>®</sup>

ⓓ Gebrauchsanweisung | Ⓜ Instructions for use | ⓕ Mode d'emploi  
ⓔ Instrucciones de manejo | ⓘ Istruzioni per l'uso

**USA** This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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


Maßstab 1:1

## INDICATION

The burrhole reservoir is used in cases of hydrocephalus to remove cerebrospinal fluid, administer medications and monitor pressure.

## DESCRIPTION

One possible treatment in cases of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain to an appropriate region of the body. Our recommendation to use a reservoir in the shunt system is to drain to the abdominal cavity in accordance to the indication of use. This shuntsystem consists of a series of catheters and a differential pressure valve. In addition, the option of positioning the burrhole reservoir in the burrhole of the cranium also exists. It then becomes possible to measure intraventricular pressure, inject drugs and check the status of the valve. Its solid titanium base is highly puncture-resistant. A puncture of the burrhole reservoir or the *SPRUNG RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0,9 mm. 30 times of punctures are able without any restrictions. The volume per pump procedure is approx. 0.08 ml for the burrhole reservoir.

With the *SPRUNG RESERVOIR* it is possible to flush CSF towards the valve due to a non-return valve in the bottom of the reservoir. By this mechanism a flow in the direction of the Ventricular Catheter is avoided during the pumping procedure. The volume per pump procedure is approx. 0.2 ml for the *SPRUNG RESERVOIR*. The opening pressure of the shuntsystem is not increased by the implantation of the *SPRUNG RESERVOIR*.

**Caution: Frequent pumping of the *SPRUNG RESERVOIR* can lead to an excessive drainage and thus to unphysiological pressure situations. The patient should be informed about the risk.**

Nothing may be passed through the burrhole reservoir. A puncture should be performed as perpendicular to the reservoir surface as possible with a cannula whose maximum size is G 21. The burr-

hole reservoir is available in various product configurations (individually or as a set) and as a burrhole port. The burrhole port doesn't have a lateral connector for a distal catheter. A Pediatric burrhole reservoir and pediatric burrhole port adapted to the size of children's skulls are also available.

## IMPLANTATION

The Ventricular Catheter is implanted with the aid of a introducing stylet. After being connected to the burrhole reservoir the catheter is fixed with a ligature.

The reservoirs are intended for use with catheters with an inner diameter of approx. 1.2 mm and an outer diameter of approx. 2.5 mm. When implanting a burrhole reservoir without a distal catheter, care must be taken to ensure that the drainage catheter is also fixed in place with the aid of a ligature.

### *Burrhole Reservoir*

It is recommended that a burrhole with a diameter of 10 mm be performed. Owing to the fact that the ventricular catheter practically disappears in the burrhole reservoir, the requisite ligature can be made in a very secure fashion (see Fig. 1).


Fig. 1: Burrhole Reservoir

### *Pediatric Burrhole Reservoir*

Unlike the standard burrhole reservoir, the Pediatric burrhole reservoir does not have an additional reduction at the bottom of the reservoir. The outer diameter is much smaller.

### *SPRUNG RESERVOIR*

As additional new feature of the *SPRUNG RESERVOIR* CSF can be flushed towards the valve due to a non-return valve in the bottom of the reservoir. By this mechanism a flow in the

direction of the Ventricular Catheter is avoided during the pumping procedure. The opening pressure of the shuntsystem is not increased by the implantation of the *SPRUNG RESERVOIR* (see Fig. 2).


Fig. 2: *SPRUNG RESERVOIR*

#### *McLANAHAN RESERVOIR*

The *McLANAHAN RESERVOIR* functions similar to the *SPRUNG RESERVOIR*: By pressing down the membrane CSF-flow is produced in distal direction (see Fig 3).


Fig. 3: *McLANAHAN RESERVOIR*: function similar to the *SPRUNG RESERVOIR*

In addition to the above function, it is possible to inject specific drugs or fluids into either the ventricles or in a distal direction into the shunt system. To inject into the ventricles: A cannula is used to puncture the reservoir.

The cannula is then used to depress the floor of the reservoir. By depressing the floor of the reservoir completely, the distal outlet is closed and the proximal outlet is opened, thus allowing drainage into the ventricles (Fig 4).


Fig. 4: *McLANAHAN RESERVOIR* after distal outlet is closed and proximal outlet is opened

When a cannula is used without depressing the reservoir floor, the proximal outlet is blocked and the distal outlet is opened. During the transition between the above mentioned conditions, both outlets are open.

**Caution: The patency of the reservoir should be tested prior to implantation. In case of uncertainty, distilled water should be aspirated carefully from the distal connector and by immersing the proximal part in sterile liquid.**

#### REIMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, as a validated decontamination process will compromise the functionality of the valve.

#### COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

Like the distal catheter (when applied), the burrhole reservoirs consist of the non-magnetic substances titanium, sapphire and silicone. Consequently, the reservoir and catheters will not interfere with either nuclear magnetic resonance or computer tomographic procedures.

The burrhole reservoirs are MRI Conditional (ASTM F2503-08). The provided catheters are MRI Safe. Deflectors and connectors are MRI Conditional.

## **SAFETY MEASURES**

Following implantation of a shuntsystem, the patient's condition should be carefully and thoroughly monitored. Inflammation and tension around the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation or vomiting often occur in cases of shunt dysfunction. In the event such symptoms occur or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

## **STERILIZATION**

All products are carefully and thoroughly sterilized. Owing to the fact that the products are packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

## **RESTERILIZATION**

The functional safety and reliability of resterilized products cannot be guaranteed.

## **FUNCTIONAL SAFETY**

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH<sub>2</sub>O during and after implantation.

## **ADVERSE REACTION**

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

## **REQUIREMENTS OF THE MDD 93/42/EEC**

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

## **NOTE ON THE INSTRUCTIONS FOR USE**

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

**MEDICAL PRODUCTS CONSULTANT**

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke  
Dipl.-Ing. Roland Schulz  
Michaela Funk-Neubarth



Christoph Miethke GmbH & Co. KG  
Ulanenweg 2  
14469 Potsdam · Germany  
Phone: +49(0) 331 620 83 0  
Fax: +49(0) 331 620 83 40  
e-mail: info@miethke.com

Please address any enquiries to:  
AESCULAP AG  
Am Aesculap-Platz  
78532 Tuttlingen · Germany  
Phone: +49 (0) 7461 95 0  
Fax: +49 (0) 7461 95 26 00  
e-mail: information@aesculap.de

Service address in the US  
AESCULAP Inc.  
Attn. AESCULAP Technical Services  
615 Lambert Pointe Road  
Hazelwood, MO, 63042  
AESCULAP Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada  
AESCULAP Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
Phone: +1-800-282-9000  
www.aesculapusa.com

**GENERAL INFORMATION**

<b>Manufacturer</b>	<b>Christoph Miethke GmbH &amp; Co. KG</b>
Product name	see Variations
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the shunt with its external dimensions:	
standard:	
 <p>Standard shunt dimensions: diameter <math>\varnothing 20\text{ mm}</math>, height <math>5,65\text{ mm}</math>.</p>	
pediatric:	
 <p>Pediatric shunt dimensions: diameter <math>\varnothing 14\text{ mm}</math>, height <math>4,8\text{ mm}</math>.</p>	

## VARIATIONS

All reservoirs are available with distal connector or distal catheter.


*Pediatric Burrhole Reservoir*

*Pediatric Burrhole Port*

*Burrhole Reservoir*


*Burrhole Port*

*SPRUNG RESERVOIR*

*MCLANAHAN RESERVOIR*


*ETOU RESERVOIR*


connector can be adjusted to various positions

The reservoirs are also available as a set, including a Ventricular Catheter.

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany