# *Prechamber* **CONTROL RESERVOIR®**

ⓓ Gebrauchsanweisung | Ⓤ Instructions for use | ⓕ Mode d'emploi  
ⓔ Instrucciones de manejo | ⓘ Istruzioni per l'uso

Ⓤ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000.

A paper copy will be provided to you upon request at no additional cost.

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**CONTENT**

INDICATION	8
TECHNICAL DESCRIPTION	8
IMPLANTATION	8
REIMPLANTATIONS	9
CONTRAINdicATIONS	9
SAFETY MEASURES	9
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	9
FUNCTIONAL SAFETY	9
ADVERSE REACTION	9
STERILIZATION	9
RESTERILIZATION	9
MEDICAL PRODUCTS CONSULTANT	10
REQUIREMENTS OF THE MDD 93/42/EEC	10

## INDICATION

The *prechamber* and the *CONTROL RESERVOIR* are used in cases of hydrocephalus to remove cerebrospinal fluid and monitor pressure.

## TECHNICAL DESCRIPTION

One possible treatment of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain to an appropriate region of the body (in most instances to the abdominal cavity). This shunt system consists of a series of catheters and a differential pressure valve. In addition, the option a *prechamber* also exists. It then becomes possible to measure intraventricular pressure, inject medication and check the status of the valve. With the *CONTROL RESERVOIR* it is possible to flush CSF towards the valve due to a nonreturn valve in the proximal inlet connector of the reservoir. By this mechanism a flow in the direction of the proximal catheter is avoided during the pumping procedure. The volume per pump procedure is approx. 0.25 ml. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.


Fig. 1: *prechamber*

**Caution: Frequent pumping of the *CONTROL RESERVOIR* can lead to an excessive drainage and thus to unphysiological pressure situations. The patient should be informed about the risk.**

Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0.9 mm. 30 times of punctures are possible without any restrictions.

The reservoirs are available in various product configurations (individually or as a set) and as a *CSF Port*. The *CSF Port* has only one lateral connector for the ventricular catheter.


Fig. 2: *CONTROL RESERVOIR* with non return valve


Fig. 3: *CSF Port*

A *paediatric prechamber* is also available, it is adapted to the size of children's skulls.


Fig. 4: *paediatric prechamber*

## IMPLANTATION

When implanting a reservoir, care must be taken to ensure catheters are secured by means of a ligature, if this has not already been done by the catheter manufacturer. The reservoirs are intended for use with catheters with an inner diameter of approx. 1.2 mm and an outer diameter of approx. 2.5 mm. When implanting a reservoir set, the neurosurgeon should determine how the procedure is to be performed in accordance with experience. The ventricle catheter is implanted with the aid of an introducing stylet, after which it can be brought through the deflector into the 90°-angle. During this process (and prior to implantation of the ventricular catheter), the deflector can be precisely adjusted,

and then used as a stop. The reservoir is used in conjunction, and is connected to a complete shunt system. A final check of the complete shunt system is recommended to ensure that it is not obstructed in any way.

## REIMPLANTATIONS

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

## CONTRAINDICATIONS

The implantation of medical products is contraindicated if an infection (f.e. meningitis, verticillitis, peritonitis, bacteremia, septicemia) or the suspicion of an infection is diagnosed in a region of the patient's body where the shunt system will be implanted.

## SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible or of the entire shunt system.

## COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

Like other shunt components of the Christoph Miethke GmbH & Co. KG (when applied), the *prechamber* and the *CONTROL RESERVOIR* consist of the nonmagnetic substances titanium, sapphire and silicone. Consequently, the reservoirs and catheters will not interfere with either nuclear magnetic resonance (MR) or computer tomographic (CT) procedures. The *prechamber* and the *CONTROL RESERVOIR* are MR Conditional (ASTM-F2503-13). The

provided catheters are MR Safe. All components are visible via X-ray.

## FUNCTIONAL SAFETY

The medical products have been designed for long-term reliable and precise operation. Still, it cannot be excluded that the medical products need to be replaced for technical or medical reasons. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH<sub>2</sub>O during and after implantation. The medical products should be stored dry and clean.

## ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development. Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

## STERILIZATION

The products are sterilized with steam under closely monitored conditions. The double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

## RESTERILIZATION

The functional safety and reliability of resterilized products cannot be guaranteed, therefore resterilisation is not recommended.

## **MEDICAL PRODUCTS CONSULTANTS**

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical products consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke

Dipl.-Ing. Roland Schulz

Michaela Funk-Neubarth

Our contact details are written on the back page of this instructions for use.

## **REQUIREMENTS OF THE MDD 93/42/EEC**

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany