Burrhole Deflector

Ⓝ DE Gebrauchsanweisung | Ⓝ GB Instructions for use | Ⓝ FR Mode d'emploi
Ⓝ ES Instrucciones de manejo | Ⓝ IT Istruzioni per l'uso

Ⓝ US This Instructions for Use is NOT intended for United State users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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CAUTION

Federal law restricts this device to sale by or on order of a physician!

INDICATION

The Burrhole Deflector is part of a shuntsystem for the treatment of hydrocephalus. It serves as a right-angular guide for the ventricular catheter coming out of the cranial borehole of a hydrocephalus patient.

DESCRIPTION

One possible treatment in cases of hydrocephalus is the implantation of a shuntsystem that achieves drainage of cerebrospinal fluid away from the ventricles of the brain into a suitable area of the body (in most cases the intestinal cavity).

This shuntsystem consists of a valve, a proximal and distal catheter and a reservoir. If the Burrhole Reservoir is placed into the burrhole, there is no need for a Burrhole Deflector. If the reservoir is a prechamber, the reservoir is in line with the valve and the Burrhole Deflector is part of the shuntsystem. The Burrhole Deflector is positioned in the borehole in the cranium and it is possible to select the length of the catheter to be inserted into the cranium prior to implantation. The Burrhole Deflector provides a turn (at a 90 degree angle) for the ventricular catheter in the borehole, which it seals off completely. Additionally to the standard Burrhole Deflector a small Pediatric Burrhole Deflector is also available. Due to the eyelets on two sides of the standard Burrhole Deflector (and eyelets on three sides for the Pediatric Burrhole Deflector), it is now possible to fix the Burrhole Deflector with ligatures to the tissue.

Fig. 1: Burrhole Deflector with ventricular catheter

IMPLANTATION

The ventricular catheter is implanted with the aid of a mandrin. The length of catheter to be implanted in the cranium is selected in advance through the positioning of the Burrhole Deflector on the catheter. The Burrhole Deflector is designed and intended for use with a catheter with an inner diameter of ca. 1.2 mm and an outer diameter of ca. 2.5 mm. The lower crank of the standard variant has 9.5 mm in diameter (see fig. 2). There-fore a borehole with a diameter of 10 mm or more should be prepared.


Fig. 2: Burrhole Deflector

REIMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, as a validated decontamination process will compromise the functionality of the valve.

CONTRAINDICATIONS

The implantation of medical products is contraindicated if an infection (f.e. meningitis, ventriculitis, peritonitis, bacteremia, septicemia) or the suspicion of an infection is diagnosed in a region of the patient's body where the shunt system will be implanted.

INTERACTION WITH PRODUCTS FROM OTHER MANUFACTURERS

When using the Burrhole Deflector, the use of catheters of the Christoph Miethke GmbH & Co. KG is recommended.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

The Burrhole Deflectors and the catheters consist of the non-magnetic substances titanium and silicone. Consequently, the Burrhole Deflector and catheters will not interfere with either nuclear magnetic resonance or computer tomographic procedures. The Burrhole Deflectors are MR Conditional (ASTM F2503-13). The provided catheters are MR Safe.

SAFETY MEASURES

Following implantation, the patient's condition should be carefully and thoroughly monitored. Inflamed skin and tension in the area of the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation, or vomiting often occur in cases of shunt dysfunction. In the event such symptoms do occur, or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

Particular care must be taken not to confuse the Burrhole Deflector with the Burrhole Reservoir. The Burrhole Reservoir is manufactured by Christoph Miethke GmbH & Co. KG. It allows both measurement of pressure and injection of medication to be carried out by means of a puncture.

ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

STERILIZATION

All products are carefully and thoroughly steam sterilized. Owing to the fact that the product is packaged in two layers of sterile packaging, five years of sterility is guaranteed. The expiration date is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

RESTERILIZATION

The safe and reliable performance of resterilized products cannot be guaranteed.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted product (Charge/LOT/SN) are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.



MEDICAL PRODUCTS CONSULTANTS

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical products consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
 Dipl.-Ing. Roland Schulz
 Michaela Funk-Neubarth
 Dipl.-Ing. Thoralf Knitter
 Dr. Andreas Bunge
 Jan Mugel

Our contact details are written on the back page of this instructions for use.

GENERAL INFORMATION

| | |
|--|--|
| Manufacturer | Christoph Miethke GmbH & Co. KG |
| Product name | <i>Burrhole Deflector</i> |
| Intended use | Treatment of hydrocephalus |
| Intended for one-time use only (disposable) | |
| Store in a clean and dry place | |
| Schematic representation of the Burrhole Deflector with its external dimensions: | |
| <p><i>Burrhole Deflector</i></p>  <p><i>Pediatric Burrhole Deflector</i></p>  | |

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany