Valve Board

 Gebrauchsanweisung |  Instructions for use |  Mode d'emploi |  Instrucciones de manejo
 www.miethke.com

 This Instructions for Use is NOT intended for United States users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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TABLE OF CONTENTS

1 GENERAL INFORMATION	17
2 INDICATION	17
3 INTENDED PURPOSE	17
4 INTENDED USE	17
5 DESCRIPTION	17
6 POSSIBLE SHUNT COMPONENTS	18
7 IMPLANTATION	18
7.1 Lumboperitoneal implantation of the Valve Board	18
7.1.1 Placing of the Lumbar Catheter	18
7.1.2 Placing of the Valve Board	20
7.1.3 Placing of the Peritoneal Catheter	21
7.2 Thoracic implantation of the Valve Board for ventriculoperitoneal drainage	21
7.2.1 Mounting the valve on the Valve Board	22
7.2.2 Thoracic placement of valve and Valve Board	22
8 SINGLE-USE PRODUCT	22
9 PRECAUTIONARY MEASURES AND CONTRAINDICATIONS	23
10 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	23
11 IDENTIFICATION IN X-RAY IMAGES	23
12 ADVERSE REACTIONS AND INTERACTIONS	23
13 PATIENT EDUCATION	23
14 STERILISATION	24
15 DISPOSAL	24
16 GUARANTEE	24
17 REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)	24
18 MEDICAL DEVICES CONSULTANTS	24
19 GENERAL INFORMATION	24
20 VARIANTS AND DIMENSIONS	25

1 GENERAL INFORMATION

NOTE

Before use:

- ▶ **Read instructions for use carefully**
- ▶ **Check expiry date**
- ▶ **Check packaging for integrity**
- ▶ **Check product for intactness**

The product and the accompanying instructions for use are intended only for medical personnel who are qualified to use the product based on their training and experience.

The product may be used only in accordance with the defined indication and may only be adapted within the scope of the possibilities mentioned. Products with damaged packaging or products that are defective or expired must not be used.

2 INDICATION

Accessories for treating hydrocephalus.

3 INTENDED PURPOSE

The *Valve Board* is used for the positionally stable suture of Miethke Shunt Systems in the body.

4 INTENDED USE

The *Valve Board* is implanted together with Miethke Shunt Systems by neurosurgeons in the sterile operating field in patients who are treated with a CSF-draining shunt system on account of their clinical picture. Depending on the specified configuration, Miethke valves as well as reservoirs and catheters are inserted into the *Valve Board*. Adjustable Miethke valves inserted into the *Valve Board* can be adjusted postoperatively with the instruments approved for this purpose. The *Valve Board* and the entire shunt system are intended for single use only. The *Valve Board* and the entire shunt system are sterile packed and must not be resterilised.

5 DESCRIPTION

The *Valve Board* (Fig. 1) is made of silicone containing 9-20% BaSO₄ for improved X-ray visibility.


Fig. 1: The *Valve Board* with a proGAV 2.0 with CONTROL RESERVOIR, lumbar right

The *Valve Board* is available in different configurations for

- lumboperitoneal** drainage (Ch. 20, Tab. 1-3) or
- thoracic placement of valve and *Valve Board* for **ventriculoperitoneal** drainage (Ch. 20, Tab. 4).

The *Valve Board* for lumboperitoneal drainage is supplied with pre-assembled Miethke Shunt Systems. The *Valve Board* is available as a shunt system with the following valves:

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
Fig. 2: Left-sided placement of the *Valve Board*

For the lumboperitoneal application of Miethke valves (*M.blue*, *proGAV 2.0* and *GAV 2.0*), a variant is available for certain *Valve Board* configurations for a left-sided placement of the *Valve Board* and subsequent distal tun-

nelling into the peritoneum (Fig. 2) plus a *Valve Board* configuration for the right-sided placement of the *Valve Board* from the spinal column and subsequent distal tunnelling into the peritoneum (Fig. 3).


Fig. 3: Right-sided placement of the *Valve Board*

The *Valve Board* is pre-configured with Miethke Shunt Systems. For safe and effective use of the Miethke valves, please refer to the instructions for use of the relevant Miethke valve. Adjustable Miethke valves inserted into the *Valve Board* can be located, checked and adjusted with the approved *M.blue plus Instruments* (for the *M.blue* and *proGAV 2.0* valves) or the *proGAV 2.0 tools* or respectively *proGAV tools* (for the *proGAV 2.0* valve), with the exception of *Checkmates*. For safe and effective use, please refer to the instructions for use of the respective instruments.

6 POSSIBLE SHUNT COMPONENTS

The *Valve Board* can be implanted as a shunt system in various configurations with the *M.blue*, *proGAV 2.0* or *GAV 2.0* valves.

These configurations can be supplied with the accessories briefly presented below.

Reservoirs

The use of a reservoir in combination with shunt systems provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to an additional non-return valve in the *CONTROL RESERVOIR*, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the drainage system as well as the *Ventricular Catheter*. During the pump action, access

to the *Ventricular Catheter* is closed. The use of this reservoir does not increase the opening pressure of the shunt system. A puncture should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. 30 punctures are possible without any restrictions.

In *Valve Board* configurations with an adjustable valve and a *CONTROL RESERVOIR*, the distance between the centre of the adjustable valve to the centre of the reservoir is 22 mm (Ch. 20, Tab. 1-3).

WARNING

Frequent pumping can result in excessive drainage and thus lead to pressure conditions outside the normal physiological range. The patient should be informed about this risk.

7 IMPLANTATION

The *Valve Board* is implanted together with Miethke Shunt Systems by neurosurgeons in the sterile operating field in patients who are treated with a CSF-draining shunt system on account of their clinical picture, in either a lumbar or thoracic location.

7.1 LUMBOPERITONEAL IMPLANTATION OF THE *Valve Board*

Various surgical techniques are possible for lumboperitoneal implantation of the *Valve Board*. The surgical technique to be used is at the discretion of the surgeon according to his experience and the clinical condition of the patient.

The *Valve Board* configurations listed in Ch. 20, Table 1-3 are intended for lumboperitoneal use. The catheter attached to the proximal entrance of the *Valve Board* can be connected to a Lumbar Catheter (inner Ø 0.8 mm; outer Ø 1.6 mm) inserted into the spinal canal using the *step down Titanium Connector* provided.

7.1.1 PLACING OF THE *Lumbar Catheter*

For safe and effective use of the Miethke *Lumbar Catheter set* please refer to the instructions for use of the Miethke *Lumbar Catheter set*. The Miethke *Lumbar Catheter* is made of silicone with a 9-20% BaSO₄ content.

Various surgical techniques are possible for placing the *Lumbar Catheter*. It is up to the doc-

tor to choose the surgical technique according to his experience and the patient's clinical condition.

- ▶ Before introducing the Tuohy cannula, making a skin incision is left to the doctor's discretion.
- ▶ Introduce the Tuohy cannula with the inserted stylet and puncture the subarachnoid space in the spinal canal using a suitable technique. The puncture technique is left to the doctor. A possible approach would be to insert the bevelled cannula tip (cannula bevel) parallel to the longitudinal dural fibres when puncturing the dura.
- ▶ Check the correct insertion of the Tuohy cannula by carefully pulling back the stylet. If cerebrospinal fluid does not drip out, reinser the stylet into the Tuohy cannula and correct the placement of the Tuohy cannula.
- ▶ If cerebrospinal fluid drips out of the Tuohy cannula, remove the stylet from the Tuohy cannula and advance the *Lumbar Catheter*, with or without a guide wire depending on the variant, through the cannula to the desired lumbar area.

WARNING

When using the guide wire, make sure that it does not penetrate the drainage holes or another location on the *Lumbar Catheter*, as this can damage the catheter and/or tissue.

NOTE

When introducing a *Lumbar Catheter* with an open tip, it can be blocked with a sheathed clamp at the distal end in order to avoid the loss of cerebral fluid. A damaged catheter is trimmed after removing the clamp.

- ▶ After correct placement of the proximal *Lumbar Catheter*, leave the Tuohy cannula in place and expose a sufficiently large subcutaneous area on the Tuohy cannula to fix the catheter in the subcutaneous tissue at a later stage. The Tuohy cannula prevents damage to the catheter during the exposure of subcutaneous tissue. A smooth transition point is created to ensure only a slight bending of the catheter.
- ▶ Then Tuohy cannula is then carefully removed via the free end of the *Lumbar Catheter*. In that process, the *Lumbar Catheter* is held in place at the exit point to ensure that it remains at the implantation site. If used, the guide wire is removed in a subsequent step.

WARNING

To avoid damaging the *Lumbar Catheter*, the *Lumbar Catheter* must never be withdrawn individually but only in conjunction with the Tuohy cannula. When using the guide wire, first remove the Tuohy cannula and then the guide wire.

WARNING

Remove the guide wire slowly out of the *Lumbar Catheter* to avoid cutting the *Lumbar Catheter*.

- ▶ CSF is used to remove air from the catheter and the catheter then clamped.

CAUTION

Catheters should only ever be clamped using a sheathed clamp as they might be damaged otherwise.

- ▶ Fix the *Lumbar Catheter* in the subcutaneous tissue using the *Lumbar Catheter* suture wings.

WARNING

Fixing the *Lumbar Catheter* is particularly important to counteract possible migration.

- ▶ Then advance the distal end of the *Lumbar Catheter* to the selected valve implantation site.
- ▶ Shorten the *Lumbar Catheter* to the appropriate length and attach it to the valve or connector using ligature. The valve should not be located directly under the skin incision.
- ▶ The *Peritoneal Catheter* can be fixed in the subcutaneous tissue using the *Peritoneal Catheter* suture wings.

CAUTION

The catheter suture wings must not be confused; otherwise, the *Lumbar Catheter* will not be fixed and the *Peritoneal Catheter* can be disconnected.

WARNING

Catheters should not be fixed with the catheter suture wings in the immediate vicinity of the valve because of potential damage to the catheters.

7.1.2 PLACING OF THE Valve Board

The *Valve Board* must be implanted in a subcutaneous pocket so that the top of the valve or reservoir faces the skin surface and the arrow markings are visible (Fig. 4). The valves are provided with arrows in flow direction (arrows pointing distally).


Fig. 4: Valve board in subcutaneous pocket

The choice of placement is left to the surgeon. Possible implantation sites are:

- ▶ Dorsal at the side of the spinal column/paravertebrally at the level of the lumbar vertebrae on firm fascial tissue (fascia thoracolumbalis/lumbodorsal fascia)
- ▶ Lower costal arch/lower lateral thoracic region

The implantation height has no influence on the function of the Miethke valves.

The *Valve Board* with adjustable valves such as the *M.blue* or the *proGAV 2.0* should rest on a firm tissue base (bone or muscle fascia), as pressure must be applied to the valve during a subsequent adjustment. The adjustable Miethke valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e.g. underneath heavily scarred tissue).

WARNING

When adjusting the valves, the orientation of the adjustable unit must be considered, as otherwise faulty opening pressures might be read or set.

CAUTION

The operation of the gravitational valve mounted on the *Valve Board* is position-dependent. When placing the *Valve Board*, care must therefore be taken to ensure that the gravitational valves of the shunt system are implanted and fixed parallel to the body axis (Fig. 2 and 3).

Using one or more skin incisions, a suitable subcutaneous skin pocket is formed into which the *Valve Board* is inserted. The *Valve Board* should not be located directly under the skin incision.

After placement, to ensure permanent positional stability of the valves, the *Valve Board* must be fixed via the fixation loops to the underlying tissue with non-absorbable suture material in order to reduce the risk of migration of the *Valve Board*.

**CAUTION**

The ligatures should be tightened at the side of the **Valve Board** to prevent the catheters from being constricted (Fig. 5).


Fig. 5: Ligatures on the side of the **Valve Board**

**CAUTION**

During placement, it must be ensured that adjustable Miethke valves (*M.blue* or *proGAV 2.0*) are positioned a maximum of 10 mm below the skin surface to make it easier to locate the valve as well as to read and set the valve pressure level (Fig. 6). With a skin and tissue thickness greater than 10 mm, it may not be possible to locate the valve as well as read and set the valve pressure level.


Fig. 6: **Valve board** in subcutaneous pocket with a maximum distance of 10 mm to the skin surface

Using the step down Titanium Connector supplied with the Valve Board, the Lumbar Catheter can be connected to the shunt system (Fig. 7). It may be necessary to shorten the catheters to the appropriate length before connecting them. The Miethke Lumbar Catheter (inner Ø: 0.8 mm; outer Ø: 1.6 mm) is pushed onto the smaller nozzle of the step down Titanium Connector (Ø 1.4 mm). All connections must be secured by ligatures with non-absorbable suture material.


Fig. 7: Connecting the catheters to the step down Titanium Connector

**CAUTION**

To avoid tension and kinking of the catheters, it is important to provide sufficient room for trunk movements by the patient (such as bending or turning) when connecting a *Lumbar Catheter* to the **Valve Board** shunt system.

7.1.3 PLACING OF THE Peritoneal Catheter

The access site for the *Peritoneal Catheter* is left to the surgeon's discretion. For example, it can be used in a paraumbilical application or applied at the level of the epigastrium. Likewise, various surgical techniques can be employed for placing the *Peritoneal Catheter*. The recommendation is to pull the *Peritoneal Catheter* from the **Valve Board** to the intended position using a subcutaneous *Tunneller*, if necessary, with the aid of an auxiliary incision. The *Peritoneal Catheter* has an open distal end and no wall slits. Check the CSF outflow at the distal end of the *Peritoneal Catheter*. Following the exposure of the peritoneum or with the aid of a trocar, the *Peritoneal Catheter* (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

7.2 THORACIC IMPLANTATION OF THE Valve Board FOR VENTRICULOPERITONEAL DRAINAGE

Two **Valve Board** variants are available for thoracic implantation for ventriculoperitoneal drainage (Ch. 20, Tab. 4).

- Valve board* variant for the adjustable valves *M.blue* or *proGAV 2.0*
- Valve Board* variant for the non-adjustable valves *GAV 2.0* or *SA 2.0*

The Valve Board variants for thoracic placement with ventriculoperitoneal drainage are supplied without valve components. For safe and effective use of the Miethke valves, please refer to the instructions for use of the relevant Miethke valve.

7.2.1 MOUNTING THE VALVE ON THE Valve Board

The adjustable proGAV 2.0 or M.blue Miethke valves are inserted into the designated Valve Board (Tab. 4) so that the arrow markings on the top of the valve are visible.

After the valve with the connected distal catheter has been inserted into the Valve Board, the proximal catheter can be connected to the free valve sleeve (Fig. 8). The connection must be secured by a ligature with non-absorbable suture material. By gently pulling the proximal catheter, the valve and connector can be pulled into the Valve Board.


Fig. 8: Inserting an M.blue into the Valve Board

WARNING

The relevant Valve Board must only be used with the following valves: proGAV 2.0, M.blue, GAV 2.0 or SA 2.0.

7.2.2 THORACIC PLACEMENT OF VALVE AND Valve Board

The Valve Board must be implanted in a subcutaneous pocket so that the top of the valve faces the skin surface and the arrow markings are visible (Fig. 4). The valves are provided with arrows in flow direction (arrows pointing distally). The choice of placement is left to the surgeon. Possible implantation sites are:

- ▶ thoracic
- ▶ infraclavicular

The implantation height has no influence on the function of the Miethke valves.

The Valve Board with adjustable valves such as the M.blue or the proGAV 2.0 should rest on a firm tissue base (bone or muscle fascia), as pressure must be applied to the valve during a subsequent adjustment.

CAUTION

Miethke gravitational valves operate in a position-dependent manner. For that reason, care must be taken when placing the Valve Board to ensure that the gravitational valves of the shunt system are implanted and fixed parallel to the body axis.

CAUTION

Adjustable Miethke valves should not be implanted in an area that makes the detection or palpation of the valve difficult (e.g. underneath heavily scarred tissue).

Using one or more skin incisions, a suitable subcutaneous skin pocket is formed into which the Valve Board with mounted valve is inserted. The Valve Board should not be located directly under the skin incision.

To ensure permanent positional stability of the valves and reduce the risk of migration of the Valve Board, after placement the Valve Board must be fixed via the fixation loops to the underlying tissue with non-absorbable suture material.

CAUTION

During placement, it must be ensured that adjustable Miethke valves (M.blue or proGAV 2.0) are positioned a maximum of 10 mm below the skin surface to make it easier to locate the valve as well as to read and set the valve pressure level (Fig. 6). With a skin and tissue thickness greater than 10 mm, it may not be possible to locate the valve as well as read and set the valve pressure level.

8 SINGLE-USE PRODUCT

Products that have already been implanted must not be reimplanted in the same patient or in another patient, as valid cleaning without loss of function cannot be achieved.

9 PRECAUTIONARY MEASURES AND CONTRAINDICATIONS

CAUTION

- ▶ When using sharp instruments, care should be taken to avoid cuts and scratches in the silicone elastomer.
- ▶ The ligatures have to be correctly positioned and not tightened too tightly or too weakly to prevent blockage of the catheter lumen and tearing of the catheter. Damage may result in a loss of integrity of the shunt and therefore may require a revision.

Silicone is very electrostatic. Please ensure that neither the catheter nor the *Valve Board* come into contact with dry cloths, talcum powder or rough surfaces. Adherent particles can lead to tissue reactions.

Patients must be carefully monitored after implantation. Reddening of skin or tightness in the area of the drained tissue may be indications of infections at the shunt system. Symptoms such as headache, dizziness, confusion or vomiting often occur in conjunction with shunt dysfunction. These symptoms and a leakage within the shunt system require the immediate replacement of the affected shunt component or the entire shunt system.

The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e.g. meningitis, ventriculitis, peritonitis, bacteraemia, septicaemia) in the region affected by the implantation.

10 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are designed to ensure precise and reliable operation over long periods of time. However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons.

Nuclear magnetic resonance (MRI) examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function.

The *Valve Board* and the supplied catheters are MRI safe. Miethke valves, reservoirs, deflectors, connectors and the *Valve Board* with mounted valves are conditionally MRI safe.

The conditions for MRI compatibility of the products are found on our website:
<https://www.miethke.com/downloads/>

11 IDENTIFICATION IN X-RAY IMAGES

The *Valve Board* and the valves, reservoirs and catheters mounted on it can be seen in the X-ray image (Fig. 9).


Fig. 9: X-ray image of Valve Board with proGAV 2.0 with Prechamber, lumbar right
(pressure level of proGAV 2.0: 20 cmH₂O;
pressure level of SHUNTASSISTANT 2.0: 25 cmH₂O;)

The set pressure level of the adjustable Miethke valves *M.blue* and *proGAV 2.0* should always be checked with the instruments provided for this purpose (*M.blue plus Instruments* for *M.blue* or *proGAV 2.0 tools* or *proGAV tools* for *proGAV 2.0*), but can also be checked with the help of an X-ray image.

To read the set pressure level of the Miethke valve used with the help of the X-ray image, please refer to the instructions for use of the relevant Miethke valve.

12 ADVERSE REACTIONS AND INTERACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): Infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in very rare cases noise development. Violent external shocks (accident, fall, etc.) may put the integrity of the shunt system at risk. Disconnection of the shunt can lead to catheter migration.

13 PATIENT EDUCATION

The attending doctor is responsible for informing the patient and/or their proxy in advance.

This includes a comprehensive description of the course of the operation, the surgical technique, the medical devices to be implanted and the instruments to be used. The patient must be informed about complications associated with the operation and possible alternative treatment methods.

14 STERILISATION

The products are sterilised with steam under strictly controlled conditions. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used under any circumstances.

The *Valve Board* and the entire shunt system are intended for single use only. This product must not be reused or resterilised. No guarantee can be assumed for the functional safety of resterilised products.

15 DISPOSAL

Product components not used during implantation and, if necessary, product components removed during surgery must be disposed of professionally as clinical waste in accordance with the relevant regional regulations.

16 GUARANTEE

Christoph Miethke GmbH & Co. KG guarantees a faultless product for single use, which is free of material and manufacturing defects upon delivery. No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in the "Implantation" chapter, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose.

17 REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the fate of medical devices used in humans, especially for implants. The individual implant identification number should therefore be recorded in the patient's medical records and patient passport to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website: <https://www.miethke.com/downloads/>

18 MEDICAL DEVICES CONSULTANTS

In compliance with the European Directive concerning medical devices (Directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical devices consultants as contacts for all product-related questions.

You can contact our medical devices consultants via:



Tel. +49 331 62083-0

info@miethke.com

19 GENERAL INFORMATION



Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>Valve board</i>
Intended use / Specification of intended purpose	Positionally stable fixation of Miethke Shunt Systems in the body
Sterilisability	cannot be resterilised
Storage	Store in a clean and dry place
Material	Silicone with 9-20% BaSO ₄

20 VARIANTS AND DIMENSIONS**M.blue, lumbar**

M.blue, lumbar	 37 mm 25 mm	 7 mm
M.blue with CONTROL RESERVOIR, lumbar right	 37 mm 45 mm 22 mm	 7 mm
M.blue with CONTROL RESERVOIR, lumbar left	 37 mm 45 mm 22 mm	 7 mm



Tab. 1: Valve board with M.blue

proGAV 2.0, lumbar

proGAV 2.0, lumbar right	 A technical line drawing of the proGAV 2.0 valve board from a right-side perspective. It features a central circular valve component with two downward-pointing arrows. On either side of the valve are two curved, semi-circular cutouts. A vertical tube extends upwards from the top left, and another vertical tube extends downwards from the bottom right. Two small circular holes are located at the bottom left and bottom right. A dimension line indicates a width of 48 mm across the bottom. To the right of the board is a vertical dimension line indicating a height of 37 mm. 48 mm	 A technical line drawing of the handle for the proGAV 2.0 valve board. It is a long, rectangular component with a slight taper and a small protrusion at the top right. A dimension line to its right indicates a width of 7 mm. 7 mm
proGAV 2.0, lumbar left	 A technical line drawing of the proGAV 2.0 valve board from a left-side perspective. It features a central circular valve component with two downward-pointing arrows. On either side of the valve are two curved, semi-circular cutouts. A vertical tube extends upwards from the top right, and another vertical tube extends downwards from the bottom left. Two small circular holes are located at the top left and top right. A dimension line indicates a width of 48 mm across the top. To the right of the board is a vertical dimension line indicating a height of 37 mm. 48 mm	 A technical line drawing of the handle for the proGAV 2.0 valve board. It is a long, rectangular component with a slight taper and a small protrusion at the top left. A dimension line to its right indicates a width of 7 mm. 7 mm



Tab. 2a: Valve board with proGAV 2.0

proGAV 2.0, lumbar

proGAV 2.0 with CONTROL RESERVOIR, lumbar right		
proGAV 2.0 with CONTROL RESERVOIR, lumbar left		





Tab. 2b: Valve board with proGAV 2.0

GAV 2.0, lumbar

GAV 2.0, lumbar		
GAV 2.0 with CONTROL RESERVOIR, lumbar right		
GAV 2.0 with CONTROL RESERVOIR, lumbar left		

Tab. 3: Valve board with GAV 2.0

Valve board for thoracic implantation

Valve board for M.blue or proGAV 2.0	  37 mm 25 mm	 7 mm
Valve board for GAV 2.0 or SA 2.0	  37 mm 25 mm	 7 mm
Valve board for CONTROL RESERVOIR	  37 mm 25 mm	 7 mm

Tab. 4: Valve board variants for thoracic placement

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany