# MRT Sicherheitsinformation / MRI Safety Information

## Shuntkomponenten / Shunt Components

**US** This information is not intended for United States users.

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
\*ASTM F2503

Dipl.-Ing. Jörg Knebel  
Person Responsible for Regulatory Compliance (PRRC)

The following products (as „stand alone“ product and as part of a MIETHKE XABO® Shunt System, MIETHKE Shunt System) are **MR CONDITIONAL\***:

*DUALSWITCH VALVE  
SHUNTASSISTANT, SHUNTASSISTANT 2.0  
paediGAV  
GAV®, GAV® 2.0  
proGAV®, proGAV® 2.0  
miniNAV®  
proSA  
uniGAV  
M.blue®, M.blue plus®  
M.flow®  
Burrhole Reservoirs  
Prechambers  
Titanium Connectors  
Burrhole Deflectors  
M.scio®  
SENSOR RESERVOIR / SENSOR PRECHAMBER  
Ventricular Catheter with Burrhole Deflector  
Valve Board (pre-configured)*

The products mentioned above have been exposed to magnetic fields up to 3 Tesla in different MRI-scanners. The products were exposed to the magnetic field in different spatial positions. Non-clinical testing demonstrated that the devices are MR Conditional. A patient with one of these devices can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 1,400 Gauss/cm (14.0 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)

No product has changed the specification after MRI exposure. The products show artefacts. An adjustment of the *proGAV*, *proSA*, *proGAV 2.0*, *M.blue*, *M.blue plus* or *M.flow* by the MRI exposure could not be determined. A verification of the valve setting after MRI exposure is not required.

The following shunt components are **MR SAFE\***:

*Peritoneal Catheters, Ventricular Catheters, THOMALE KATHETER  
XABO® Peritoneal Catheters, XABO® Ventricular Catheters  
Lumbar Catheters  
Lumbar/ Peritoneal Catheter Suture Wings  
Valve Board*

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights. Detailed information can be requested from [info@miethke.com](mailto:info@miethke.com).