M.blue plus[®] Instruments

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 www.miethke.com

- ⒹS This Instructions for Use is NOT intended for United State users.
Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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TABLE OF CONTENTS

| | |
|---|----|
| INDICATION | 14 |
| APPLICATION OF <i>M.blue plus instruments</i> | 14 |
| SELECTION OF THE APPROPRIATE PRESSURE | 18 |
| CLEANING RECOMMENDATION FOR THE NON-STERILISABLE <i>M.blue plus instruments</i> | 19 |
| CLEANING AND DISINFEKTIION OF THE <i>M.blue adjustment gyroscope</i> | 19 |
| STERILISATION | 21 |
| RECOMMENDATION FOR STERILISATION OF THE <i>M.blue adjustment gyroscope</i> | 21 |
| MEDICAL PRODUCTS CONSULTANTS | 21 |
| GENERAL INFORMATION | 22 |

INDICATION

M.blue plus instruments can be used to determine, change and monitor the pressure level of the adjustable gravitational unit of the M.blue valve as well as the pressure level of the adjustable differential pressure unit of the proGAV 2.0 valve.

M.blue is a valve with location-dependent functioning for cerebrospinal fluid (CSF) drainage in the treatment of hydrocephalus. It consists of a differential pressure unit and an adjustable gravitational unit (Fig. 1).

Figure 1: M.blue is a combination comprising differential pressure unit and adjustable gravitational unit.

proGAV 2.0 is also a valve with location-dependent functioning for cerebrospinal fluid (CSF) drainage in the treatment of hydrocephalus. It consists of an adjustable differential pressure unit and a gravitational unit.

Figure 2: proGAV 2.0 is a combination comprising adjustable differential pressure unit and gravitational unit.

Figure 3: M.blue plus is a combination comprising differential pressure unit (proGAV 2.0) and adjustable gravitational unit (M.blue).

APPLICATION OF M.blue plus instruments

The M.blue plus Instruments may only be used by trained specialist personnel.

M.blue plus instruments can be used to determine, change and monitor the pressure level of the adjustable gravitational

unit of the M.blue valve as well as the pressure level of the adjustable differential pressure unit of the proGAV 2.0 valve.

M.blue plus compass serves for localising the valve and reading the set opening pressure of the gravitational unit of M.blue as well as the differential pressure unit of the proGAV 2.0.

Figure 4: M.blue plus compass

The M.blue plus adjustment ring allows the opening pressure of M.blue to be set from 0 to 40 cmH₂O. The M.blue plus adjustment ring also allows the opening pressure of the differential pressure unit of proGAV 2.0 to be set from 0 to 20 cmH₂O.

Figure 5: M.blue plus adjustment ring

The opening pressure for the adjustable gravitational unit of M.blue can be changed before or after implantation. The opening pressure of the gravitational unit is preset by the manufacturer to 20 cmH₂O.

The opening pressure for the adjustable differential pressure unit of *proGAV 2.0* can be changed before or after implantation. The opening pressure of the differential pressure unit is preset by the manufacturer to 5 cmH₂O.

The following steps must be performed to adjust the opening pressures:

1. Localisation

When the *M.blue plus compass* is opened, a circular cut-out becomes visible with which you can localise the valve in the patient's head using your index finger (Fig. 6).

Figure 6: Localising the valve using the *M.blue compass*

The direction markings show the flow direction.

2. Test procedure

In order to determine the selected pressure level, the compass is then closed again. The float gauge should now be centred by moving the instrument in the designated circular marking (Fig. 6). Once the float gauge is centred, the currently set opening pressure of the differential pressure unit (*proGAV 2.0*) and/or of the gravitational unit (*M.blue*) can be read off via the line mark on the float gauge (Fig. 7).

Figure 7: Determining the pressure level with the *M.blue compass*

There are two scales on the scale ring. The grey adjustment range from 0 to 20 cmH₂O on the outer scale applies for the opening pressure of the adjustable differential pressure unit of *proGAV 2.0*. The opening pressure of the gravitational unit of *M.blue* is set to the blue marked adjustment range from 0 to 40 cmH₂O of the inner scale.


Figure 8: Outer: Scale of the differential pressure unit of *proGAV 2.0* from 0 to 20 cmH₂O

Inner: Scale of the gravitational unit of *M.blue* from 0 to 40 cmH₂O (opening pressure of the differential pressure unit of *proGAV 2.0* in the example illustrated is 17 cmH₂O, opening pressure of the gravitational unit of *M.blue* in the example illustrated is 16 cmH₂O)

WARNING

The *M.blue plus* compass should be placed centrally over the valve, if possible; otherwise incorrect determination of the opening pressure may occur.

M.blue plus compass reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus* adjustment ring should not be in the immediate vicinity of the *M.blue plus* compass when determining the opening pressure. We recommend a minimum distance of 30 cm.

NOTE

Possible air bubbles inside the compass do not affect its functionality.

3. Adjustment process

In order to adjust the opening pressure, the compass is opened, but without changing the position of the scale ring. The adjustment ring is now inserted into the scale ring in such a manner that the line marking points to the desired value on the scale of the scale ring (Fig. 9).


Figure 9: Inserting the adjustment ring

The grey adjustment range from 0 to 20 cmH₂O on the outer scale applies for the opening pressure of the differential pressure unit of the *proGAV2.0*. The opening pressure of the *M.blue* gravitational unit uses the adjustment range from 0 to 40 cmH₂O of the blue inner scale.

Figure 10: Setting the differential pressure unit of *proGAV 2.0* in the example illustrated to 1 cmH₂O or adjusting the gravitational unit of *M.blue* in the example illustrated to 32 cmH₂O.

By applying slight pressure with the index finger to the valve diaphragm located in the centre of the adjustment ring and under the skin, the rotor brake is released and the opening pressure of the differential or gravitational unit is changed to the desired value (Fig. 11).

Figure 11: Adjustment using the *M.blue plus* adjustment ring

Both the differential pressure unit of *proGAV 2.0* and the gravitational unit of *M.blue* are equipped with a feedback mechanism. If targeted pressure is exerted on the valve, an acoustic signal (a clicking sound) is audible, or a resistance can be felt, as soon as the rotor brake has been released due to the valve housing design. In other words the valve shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released the rotor is adjustment-proof again. Although

the click caused by releasing the rotor brake is easily audible before implantation, this can be considerably reduced after implantation and filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patients themselves or through the use of a stethoscope.

Adjustment with the adjustment assistant

Alternatively, the *M.blue plus adjustment assistant* can be used to adjust the opening pressure. To do this, insert the *M.blue plus adjustment assistant* into the adjustment ring aligned to the desired value and press it with your index finger (Fig. 12).

Figure 12: *M.blue plus adjustment assistant*

WARNING

For adjustment of the differential pressure unit of *proGAV 2.0*, ensure that the opening pressure is changed by a maximum of 8 cmH₂O per adjustment process; otherwise errors can result.

Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 3 to 11 and subsequently from 11 to 18 cmH₂O.

WARNING

When adjusting the gravitational unit of *M.blue*, please ensure that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process; otherwise errors may result.

Example: The opening pressure is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 6 to 22 and subsequently from 22 to 36 cmH₂O.

CAUTION

The *M.blue plus adjustment ring* emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.

4. Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the set pressure level. To do this, please proceed as in Points 1 and 2. If the measured value does not agree with the required pressure stage the adjustment process should be repeated. To do this, please start at Point 3 again.

Adjustment can be made difficult by swelling of the skin for a few days after the operation. If checking the valve setting is not possible without ambiguity using the *M.blue plus compass*, we recommend carrying out a check using an imaging method.

M.blue adjustment gyroscope

The *M.blue adjustment gyroscope* is supplied sterile and can be re-sterilised. With the *M.blue adjustment gyroscope* it is possible to carry out a pressure level change and the check before and during valve implantation directly on *M.blue*. To determine the pressure level the *M.blue*

adjustment gyroscope is placed centrally on *M.blue*. The *M.blue adjustment gyroscope* automatically aligns itself over the valve. The pressure stage can be read off from the catheter in the proximal (leading to valve) direction. If the pressure level is supposed to be adjusted, the *M.blue adjustment gyroscope* is placed centrally on *M.blue*. When doing so, the required pressure stage must point towards the proximal catheter (leading to valve). If the *M.blue adjustment gyroscope* is now pressed lightly onto the valve, the rotor brake in *M.blue* is released and the pressure level is set. When adjusting, please ensure that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process otherwise errors can result (see Chapter "3. Adjustment process").

Figure 13: *M.blue adjustment gyroscope*, colour: blue, pressure levels 0 to 40 cmH₂O

CAUTION

Because of the magnets inside the *M.blue plus instruments*, *M.blue plus instruments* may not be used in the vicinity of active implants such as heart pacemakers. Furthermore, there is risk of damage to MRI instruments in their vicinity. Therefore, usage of *M.blue plus instruments* is forbidden in such locations!

It is absolutely essential to use only the *M.blue plus instruments* to determine, adjust and monitor the opening pressure of the gravitational unit in *M.blue*. The opening pressure of the differential pres-

sure unit in *proGAV 2.0* can be determined, adjusted and monitored with the described *M.blue plus instruments* as well as with the *proGAV 2.0 instruments*.

Figure 14: *M.blue plus instruments*

SELECTION OF THE APPROPRIATE PRESSURE

Horizontal position

In the lying position, the gravitational unit has no influence on the valve opening pressure. In this position, the valve opening pressure is thus determined exclusively by the differential pressure unit.

By default, a differential pressure unit with an opening pressure of 5 cmH₂O is recommended.

Vertical position

The opening pressure for the vertical body position is calculated from the sum of the opening pressure in the differential pressure unit and the gravitational unit for both *M.blue* as well as *proGAV 2.0*.

Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure for the gravitational unit (see pressure level recommendations at <https://www.miethke.com/produkte/downloads/>).

This is a non-binding recommendation. The physician must decide in each individual case.

CLEANING RECOMMENDATION FOR THE NON-STERILISABLE *M.blue plus instruments*

NOTE

M.blue plus instruments are made from thermolabile components that are sensitive to heat and moisture as well as chemically reactive. Do not immerse *M.blue plus instruments* in cleaning solutions or allow liquids to enter the housing, as this could impair the functioning as a result of moisture, corrosion and impurities.

Surface impurities on the *M.blue plus instruments* should always be removed immediately after use by wiping with alcohol-based detergents (minimum 75 % alcohol).

The exposure time should be at least 60 seconds and should be adapted to the degree of contamination. The instruments should then be wiped with a dry cloth.

The following cleaning methods are not suitable for preparation of the *M.blue plus instruments* (apart from *M.blue adjustment gyroscope*): Irradiation, ultrasound, sterilisation, mechanical preparation, immersion in cleaning fluids.

CLEANING AND DISINFEKTIION OF THE *M.blue adjustment gyroscope*

Avoid damaging the product due to unsuitable detergents/disinfectants and/or too high temperatures!

- Use detergents/disinfectants that are approved for surgical steels, according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and contact time.

- Do not exceed the maximum permissible cleaning temperature of 55°C.
- Perform ultrasonic cleaning:
 - as effective mechanical support for manual cleaning/disinfection
 - for precleaning products with dried residues before the mechanical cleaning/disinfection
 - as integrated mechanical support during mechanical cleaning/disinfection
 - for subsequently cleaning products with residues that have not been removed after mechanical cleaning/disinfection
- If the instruments cannot be fastened securely or safely for cleaning in machines or on storage aids, clean and disinfect the instruments mechanically.

Manual cleaning/disinfection

- Check visible surfaces for residues after manual cleaning/disinfection.
- If necessary, repeat the cleaning process.

Mechanical cleaning/disinfection

Place product on cleaning-compatible screen basket.

Mechanical cleaning/disinfection with manual precleaning

Manual precleaning with ultrasound

| Phase | I | II |
|---------------|---|----------------------|
| Step | Disinfected ultrasonic cleaning | Intermediate rinsing |
| T (°C/°F) | RT (cold) | RT (cold) |
| t (min) | 15 | 1 |
| Conc. (%) | 2 | - |
| Water quality | D-W | D-W |
| Chemical | B.Braun Stabimed®; free of aldehydes, phenols and QAV; pH = 9 | |

Mechanical alkaline cleaning and thermal disinfection

Phase I

- Clean product in the ultrasonic cleaning bath (frequency 35 kHz). When doing so, make sure that all accessible surfaces are wet.

Phase II

- Fully rinse/flush product (all accessible surfaces) under running water.

Device type: Single chamber cleaning/disinfection device without ultrasound

- Place product on cleaning-compatible screen basket (avoid missed areas)

D-W: Drinking water

RT: Room temperature

| Phase | I | II | III | IV | V | VI |
|---------------|-------------|--|---|-------------------------|----------------------|---------|
| Step | Pre-rinsing | Cleaning | Neutralisation | Intermediate inspection | Thermal disinfection | Drying |
| T (°C/°F) | <25/ 77 | 55/ 131 | 20/ 68 | 70/ 158 | 94/ 201 | 90/ 194 |
| t (min) | 3 | 10 | 2 | 1 | 10 | 40 |
| Water quality | D-W | DM-W | DM-W | DM-W | DM-W | - |
| Chemical | - | -Concentration alkaline: pH = 10.9 < 5 % anionic surfactants -1% solution pH = 10.5 | -Concentrate, acidic: pH = 2.6 base: Citric acid -1% solution pH = 3.0 | - | - | - |

D-W: Drinking water; DM-W: Demineralised water de-ionised

Monitoring, care, testing

- Allow product to cool down to room temperature.
- Check product after each cleaning and disinfection for: Cleanliness, functioning and damage, e.g. insulation, loose, bent, broken, cracked, worn and detached parts.
- Immediately set aside product if damaged. It is recommended that the instrument be packed in suitable double sterile packaging. This will ensure reuse in the best way possible.

STERILISATION

Apart from the *M.blue adjustment gyroscope* the *M.blue plus instruments* cannot be sterilised, only disinfected.

RECOMMENDATION FOR STERILISATION OF THE *M.blue adjustment gyroscope*

The *M.blue adjustment gyroscope* in its double sterile packaging must be sterilised in an autoclave (steam sterilisation, fractionated vacuum method) at 134°C and 5 minutes retention time.

MEDICAL PRODUCTS CONSULTANTS

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical products consultants as contacts for all product-related questions:

- Dipl.-Ing. Christoph Miethke
- Dipl.-Ing. Roland Schulz
- Michaela Funk-Neubarth
- Dipl.-Ing. Thoralf Knitter

- Dr. Andreas Bunge
- Jan Mugel

Contact details can be found on the reverse of these instructions for use.

GENERAL INFORMATION

| | |
|---------------------|---------------------------------|
| Manufacturer | Christoph Miethke GmbH & Co. KG |
| Product designation | <i>M.blue plus Instruments</i> |
| Intended use | Treatment of hydrocephalus |
| Storage | Store in a clean and dry place |

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany