
M.blue plus® Instruments

(DE) Gebrauchsanweisung | (EN) Instructions for use | (FR) Mode d'emploi | (ES) Instrucciones de manejo |

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 www.miethke.com

(US) This Instruction for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0
info@miethke.com

TABLE OF CONTENTS

0.00 PREFACE AND IMPORTANT INFORMATION	19
1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	19
1.01 Explanation of the warnings	19
1.02 Display conventions	19
1.03 Other accompanying documents and additional information	19
1.04 Feedback on the instructions for use	20
1.05 Copyright, disclaimer, warranty and other information	20
2.00 DESCRIPTION OF THE <i>M.blue plus Instruments</i>	20
2.01 Medical purpose	20
2.02 Clinical use	20
2.03 Indications	20
2.04 Contraindications	20
2.05 Intended patient groups	20
2.06 Intended users	20
2.07 Intended usage environment	21
2.08 Technical description	21
2.09 System components	22
3.00 PROPERTIES OF <i>M.blue plus Instruments</i>	22
3.01 Product description	22
3.02 Important safety information	23
3.03 Transport and storage	24
3.04 Use of the product	24
3.05 Disposal	29
3.06 Troubleshooting	30
3.07 Technical information	31
3.08 Symbols used for labelling	32
4.00 MEDICAL DEVICE CONSULTANTS	32

0.00 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing *M.blue plus Instruments*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product.

Your team at Christoph Miethke GmbH & Co. KG

Relevance of the instructions for use

WARNING

Improper handling and use of this product can cause risks and damage. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

The *M.blue plus Instruments* include the following components:

- ▶ *M.blue plus Compass*
- ▶ *M.blue plus Adjustment Ring*

Additional options:

- ▶ *M.blue plus Adjustment Assistant*
- ▶ *M.blue Checkmate*

1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

1.01 EXPLANATION OF THE WARNINGS

DANGER

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

WARNING

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

CAUTION

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

NOTE

Indicates a possibly harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for risk, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

1.02 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates product names

1.03 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

If you still require help after reading the instructions for use and the additional information, please contact us or your authorised distributor.

1.04 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

1.05 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

2.00 DESCRIPTION OF THE M.blue plus Instruments

2.01 MEDICAL PURPOSE

The *M.blue plus Instruments* are intended for adjusting and verifying the adjustable pressure levels of the *M.blue*, *M.blue plus* and *proGAV 2.0* valves.

2.02 CLINICAL USE

Treatment options

- ▶ Localising the adjustable valve unit using the instruments
- ▶ Identifying the pressure level of the gravitational unit of the *M.blue* valve and the differential pressure unit of the *proGAV 2.0* valve
- ▶ Checking the pressure level without diagnostic procedures such as X-rays
- ▶ Setting the pressure level without invasive procedures

2.03 INDICATIONS

M.blue plus Instruments are used for the following indications:

- ▶ Treatment of hydrocephalus

2.04 CONTRAINDICATIONS

The following contraindications apply to the *M.blue plus Instruments*:

- ▶ Intolerance to materials of the instruments

2.05 INTENDED PATIENT GROUPS

- ▶ Patients with an implanted, CSF-draining shunt system with an adjustable *M.blue*, *proGAV 2.0* or *M.blue plus* valve

2.06 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical specialists, e.g. neurosurgeons, nurses
- ▶ Knowledge about the function and proper use of the product
- ▶ Successful participation in product training

2.07 INTENDED USAGE ENVIRONMENT

Professional Healthcare Facilities

- Medical facilities, such as hospitals, doctor's surgeries and rehabilitation facilities

2.08 TECHNICAL DESCRIPTION

M.blue plus Instruments can be used to determine, change and monitor the pressure level of the adjustable gravitational unit of the *M.blue* valve as well as the pressure level of the adjustable differential pressure unit of the *proGAV 2.0* valve.

M.blue (Fig. 1) is a position-dependent valve to drain cerebrospinal fluid (CSF) in the treatment of hydrocephalus. It consists of a differential pressure unit and an adjustable gravitational unit.

Fig. 1: *M.blue*

proGAV 2.0 (Fig. 2) is also a position-dependent valve to drain cerebrospinal fluid (CSF) in the treatment of hydrocephalus. It consists of an adjustable differential pressure unit and a gravitational unit.

Fig. 2: *proGAV 2.0*

Fig. 3: *M.blue plus* combines an adjustable differential pressure unit (*proGAV 2.0*) and an adjustable gravitational unit (*M.blue*) and will be delivered in the presented order (*proGAV 2.0* proximal and *M.blue* distal).

M.blue plus Compass (Fig. 4) is used to locate the *M.blue* valve and to read the set opening pressure of its gravitational unit as well as the adjustable differential pressure unit of the *proGAV 2.0* valve.


Fig. 4: *M.blue plus Compass*

a) open

b) closed

1. Scale ring

2. Floating compass

The *M.blue plus Adjustment Ring* (Fig. 5) can be used to set the opening pressure of the gravitational unit of the *M.blue* from 0 to 40 cmH₂O. The *M.blue plus Adjustment Ring* can also be used to set the opening pressure of the differential pressure unit of the *proGAV 2.0* from 0 to 20 cmH₂O.


Fig. 5: *M.blue plus Adjustment Ring*

The opening pressure of the adjustable gravitational unit of the *M.blue* can be changed before or after implantation. It is preset by the manufacturer at 20 cmH₂O.

The opening pressure of the adjustable differential pressure unit of the *proGAV 2.0* can be changed before or after implantation. It is preset by the manufacturer at 5 cmH₂O.

2.09 SYSTEM COMPONENTS

M.blue plus Instruments can be used to set the opening pressures of the following valves:

Product name	Figure
M.blue	
M.blue plus	
proGAV 2.0	

3.00 PROPERTIES OF M.blue plus Instruments

3.01 PRODUCT DESCRIPTION


3.01.01 APPROVED COMPONENTS

The M.blue plus Instruments are comprised of the

- ▶ M.blue plus Compass (1) and the
- ▶ M.blue plus Adjustment Ring (2).

M.blue plus Instruments can also be supplemented with the

- ▶ M.blue plus Adjustment Assistant (3) and the
- ▶ M.blue Checkmate (4).

	Item no.	Figure
1	FX891T	
2	FX892T	
3	FX893T	
4	FX894T	

3.01.02 SCOPE OF DELIVERY

Box content (item no. FX890T)	Number
M.blue plus Compass	1
M.blue plus Adjustment Ring	1
M.blue plus Instruments instructions for use	1

Box content (item no. FX891T)	Number
M.blue plus Compass	1
M.blue plus Instruments instructions for use	1

Box content (item no. FX892T)	Number
<i>M.blue plus Adjustment Ring</i>	1
<i>M.blue plus Instruments</i> instructions for use	1

Box content (item no. FX893T)	Number
<i>M.blue plus Adjustment Assistant</i>	1
<i>M.blue plus Instruments</i> instructions for use	1

Box content (item no. FX894T)	Number
<i>M.blue Checkmate</i>	1
<i>M.blue plus Instruments</i> instructions for use	1

3.01.03 STERILITY

With the exception of the *M.blue Checkmate*, *M.blue plus Instruments* cannot be sterilised, only disinfected.

WARNING

If the packaging or the product is damaged, the products must not be used under any circumstances.

3.01.04 REPEATED USE AND RE-STERILISATION

The *M.blue Checkmate* must be sterilised in double sterile packaging in an autoclave (steam sterilisation, fractionated vacuum method) at 134 °C and 5 minutes retention time.

3.01.05 OPERATING CONDITIONS

M.blue plus Instruments are to be used at room temperature in medical facilities such as hospitals, doctor's surgeries and rehabilitation facilities.

3.01.06 PRODUCT LIFE

These medical devices are designed to operate accurately and reliably over long periods of time. The expected product life of *M.blue plus Instruments* is three years after their first use, provided that the product is exposed to normal usage conditions and is properly maintained (see chapter "3.04 Use of the product"). The instruments may be used beyond this time if they are in perfect working condition. However, no guarantee can be assumed that these

medical devices will not require replacement for medical or technical reasons.

3.01.07 PRODUCT CONFORMITY

The product meets current regulatory requirements.

Translations of these instructions for use into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

3.02 IMPORTANT SAFETY INFORMATION

3.02.01 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ If the sterile packaging of the *Checkmate* is damaged, the product must not be used in a sterile situation.
- ▶ After the expiry date the *Checkmate* must first be re-sterilised as described under chapter 3.01.04 before use in a sterile situation.

CAUTION

- ▶ Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be read carefully and understood before the product is used for the first time.
- ▶ Prior to use, it is essential to check that the *M.blue plus Instruments* are complete and intact.

3.02.02 COMPLICATIONS, SIDE EFFECTS AND RESIDUAL RISKS

Various side effects and complications can occur when using *M.blue plus Instruments*:

- ▶ Allergic reaction / intolerance to the materials of the instruments

Immediately consult a physician if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or the like.

3.02.03 REPORTING OBLIGATION

Report all serious incidents (damage, injuries, infections, etc.) that occur relating to the product to the manufacturer and the responsible state authority.

3.02.04 PATIENT INFORMATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product.

3.03 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place.

Storage conditions

Temperature range for storage	0 °C to +40 °C
-------------------------------	----------------

3.04 USE OF THE PRODUCT

3.04.01 INTRODUCTION

M.blue plus Instruments may only be used by trained specialists.

M.blue plus Instruments can be used to determine, change and monitor the pressure level of the adjustable gravitational unit of the *M.blue* valve as well as the pressure level of the adjustable differential pressure unit of the *proGAV 2.0* valve.

3.04.02 SAFETY AND WARNING NOTICES

CAUTION

- ▶ The *M.blue plus Adjustment Ring* emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.
- ▶ Because of the magnets inside *M.blue plus Instruments*, *M.blue plus Instruments* must not be used in the vicinity of active implants, such as cardiac pacemakers. Furthermore, there is risk of damage to MRI instruments in their vicinity. Therefore, *M.blue plus Instruments* must not be used in such locations!

CAUTION

- ▶ It is absolutely essential that only *M.blue plus Instruments* are used to determine, change and monitor the opening pressure of the gravitational unit of the *M.blue*. The opening pressure of the differential pressure unit of the *proGAV 2.0* can be determined, changed and monitored with the *M.blue plus Instruments* as described and also the *proGAV 2.0 Tools*.

3.04.03 REQUIRED MATERIALS

Adjustment of the adjustable gravitational unit of the *M.blue* or the adjustable differential pressure unit of the *proGAV 2.0* after implantation requires the *M.blue plus Compass*, the *M.blue plus Adjustment Ring* and possibly the *M.blue plus Adjustment Assistant*.

To adjust the adjustable gravitational unit of the *M.blue* before or during valve implantation, the resterilisable *M.blue Checkmate* can be used.

3.04.04 PREPARING FOR USE

Checking the sterile packaging

Before using the *M.blue Checkmate*, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Checking the integrity of the *M.blue plus Instruments*

Prior to use, the instruments must be checked for integrity. To this end, all instruments must be visually inspected.

3.04.05 ADJUSTING THE *M.blue*

In order to set the opening pressures, the following steps must be performed:

1. Localisation

WARNING

The *M.blue plus Compass* should be placed as centrally as possible over the valve, as otherwise the determined opening pressure may be incorrect.

NOTE

- ▶ Possible air pockets in the *M.blue plus Compass* do not affect its function.

NOTE

- ▶ The *M.blue plus Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus Adjustment Ring* should not be placed in the immediate vicinity of the *M.blue plus Compass* while determining the opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin can make adjustment difficult for a few days after the operation. If the valve setting cannot be checked conclusively using the *Compass*, we recommend carrying out a check using an imaging method.

When the *M.blue plus Compass* is opened, a circular cut-out becomes visible, which is used to localise the valve in the patient's head as centrally as possible using the index finger (Fig. 6).

Fig. 6: Localising the valve

The direction markings show the flow direction.

2. Test procedure

In order to determine the selected pressure level, the *Compass* is then closed again. The float gauge should now be centred by moving the instrument within the designated circular marking (Fig. 7). Once the float gauge is centred, the current setting of the opening pressure of the differential pressure unit (*proGAV 2.0*) and/or of the gravitational unit (*M.blue*) can be read from the line mark on the float gauge (Fig. 7).


Fig. 7: Determining the pressure level with the *M.blue plus Compass*

The scale ring is marked with two scales (Fig. 8). The grey adjustment range from 0 to 20 cmH₂O on the outer scale applies to the opening pressure of the adjustable differential pressure unit of the *proGAV 2.0*. The adjustment range from 0 to 40 cmH₂O of the inner scale, marked in blue, applies to the opening pressure of the gravitational unit of the *M.blue*.


Fig. 8: Scale ring of the *M.blue plus Compass*

1. Inner: Scale of the *M.blue* gravitational unit from 0 to 40 cmH₂O

2. Outer: Scale of the *proGAV 2.0* differential pressure unit from 0 to 20 cmH₂O

(The opening pressure of the *M.blue* gravitational unit in the illustrated example is 16 cmH₂O, the opening pressure of the *proGAV 2.0* differential pressure unit 17 cmH₂O)

3. Adjustment process

WARNING

- ▶ When adjusting the *proGAV 2.0* differential pressure unit, care must be taken that the opening pressure is changed by a maximum of 8 cmH₂O per adjustment process, otherwise errors can occur.

Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 3 to 11 cmH₂O and subsequently from 11 to 18 cmH₂O.

- ▶ When adjusting the *M.blue* gravitational unit, care must be taken that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process, otherwise errors can occur.

Example: The opening pressure is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 6 to 22 cmH₂O and subsequently from 22 to 36 cmH₂O.

In order to adjust the opening pressure, the *Compass* is opened, but without changing the position of the scale ring. The *Adjustment Ring* is now inserted into the scale ring so that the line marking points to the desired value on the scale of the scale ring (Fig. 9).


Fig. 9: Inserting the Adjustment Ring
1. Adjustment Ring
2. Scale ring

Fig. 10: Setting the *proGAV 2.0* differential pressure unit in the illustrated example to 1 cmH₂O or setting the *M.blue* gravitational unit to 32 cmH₂O.

By applying slight pressure with the index finger to the valve diaphragm located in the centre of the *Adjustment Ring* and under the skin, the rotor brake is released and the opening pressure of the differential pressure or gravitational unit is changed to the desired value (Fig. 11).


Fig. 11: Adjustment using *M.blue plus* Adjustment Ring

Both the differential pressure unit of the *proGAV 2.0* and the gravitational unit of the *M.blue* are equipped with a feedback mechanism. If targeted pressure is exerted on the valve, due to the valve housing design an acoustic signal (a clicking sound) is audible or a resistance can be felt, as soon as the rotor brake has been released. The valve thus shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released, the rotor is once again adjustment-proof to outside magnetic fields. Although the click caused by releasing the rotor brake is easily audible before implantation, it may be considerably reduced after implantation and the filling of the valve depending on its posi-

tion and the condition of the implant surroundings. Normally, however, it should be audible to the patient or by using a stethoscope.

Adjustment with the adjustment assistant

Alternatively, the *M.blue plus Adjustment Assistant* can be used to adjust the opening pressure. To do this, insert the *M.blue plus Adjustment Assistant* into the *Adjustment Ring* aligned to the desired value and press it with your index finger (Fig. 12).


Fig. 12: *M.blue plus Adjustment Assistant*

Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the set pressure level. To do this, proceed as in Points 1 and 2. Should the measured value not correspond with the required pressure level, the adjustment process should be repeated. To do this, start again at Point 3.

M.blue Checkmate CE 0297

The *M.blue Checkmate* is supplied sterile and can be re-sterilised. The *M.blue Checkmate* can be used to change and check the pressure level before and during valve implantation directly at the *M.blue* or through the packet containing the *M.blue*. To determine the pressure level, the *M.blue Checkmate* is placed centrally onto the *M.blue*. The *M.blue Checkmate* automatically aligns itself over the valve. The pressure level can be read from the direction of the proximal catheter I (leading towards the valve). If the pressure level is to be adjusted, the *M.blue Checkmate* is placed centrally onto the *M.blue*. When doing so, the required pressure level must point towards the proximal catheter (leading towards the valve). By slightly pressing the *M.blue Checkmate* onto the valve, the rotor brake in the *M.blue* is released and the pressure level set.

Fig. 13: *M.blue Checkmate*

Colour: blue

Pressure levels: 0 to 40 cmH₂O

3.04.06 CLEANING RECOMMENDATION FOR NON-STERILISABLE *M.blue plus Instruments*

NOTE

M.blue plus Instruments are made from thermolabile components that are sensitive to heat and moisture as well as chemically reactive. Do not immerse *M.blue plus Instruments* in cleaning solutions or allow liquids to enter the housing, as moisture, corrosion and impurities could impair proper functioning.

Surface impurities on the *M.blue plus Instruments* should always be removed immediately after use by wiping with alcohol-based cleaning agents (minimum 75% alcohol).

The exposure time should be at least 60 seconds and should be adapted to the degree of contamination. The instruments should then be wiped with a dry cloth.

The following cleaning methods are not suitable for processing the *M.blue plus Instruments* (apart from the *M.blue Checkmate*): Irradiation, ultrasound, sterilisation, mechanical processing, immersion in cleaning fluids.

3.04.07 CLEANING AND DISINFECTING THE *M.blue Checkmate*

Avoid damaging the product due to unsuitable cleaning agents/disinfectants and/or too high temperatures!

- ▶ Use cleaning agents/disinfectants approved for surgical steels in accordance with the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and contact time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 55°C.
- ▶ Perform ultrasonic cleaning:
 - ▶ as an effective mechanical support for manual cleaning/disinfection
 - ▶ for precleaning products with dried residues prior to mechanical cleaning/disinfection
 - ▶ as integral mechanical support during mechanical cleaning/disinfection
 - ▶ for the subsequent cleaning of products with residues still present after mechanical cleaning/disinfection
- ▶ Provided the instruments can be fastened securely or safely for cleaning in machines or on storage aids, clean and disinfect the instruments mechanically.

Manual cleaning/disinfection

- ▶ Check visible surfaces for residues after manual cleaning/disinfection.
- ▶ If necessary, repeat the cleaning process.

Mechanical cleaning/disinfection

- ▶ Place product on suitable screen basket.

Mechanical cleaning/disinfection with manual precleaning

- ▶ Manual precleaning with ultrasound

Phase	I	II
Step	Disinfected ultrasonic cleaning	Intermediate rinsing
T (°C/°F)	RT (cold)	RT (cold)
t (min)	15	1
Conc. (%)	2	-
Water quality	D-W	D-W
Chemical specifications	B. Braun Stabimed®; free of aldehyde, phenol and quaternary ammonium compounds; pH = 9	

D-W: Drinking water; RT: Room temperature

Mechanical alkaline cleaning and thermal disinfection

Phase I

- Clean product in the ultrasonic cleaning bath (frequency 35 kHz). Ensure that all accessible surfaces are wetted.

Phase II

- Fully rinse/flush product (all accessible surfaces) under running water.

Device type: Single chamber cleaning/disinfection device without ultrasound

- Place product on suitable screen basket (avoid missed areas)

Phase	I	II	III	IV	V	VI
Step	Pre-rinsing	Cleaning	Neutralisation	Inter-mediate inspection	Thermal disinfection	Drying
T (°C/°F)	< 25/77	55/131	20/68	70/158	94/201	90/194
t (min)	3	10	2	1	10	40
Water quality	D-W	P-W	P-W	P-W	P-W	-
Chemical specifications	-	Concentration, alkaline: pH = 10.9 < 5% anionic surfactants 1% solution pH = 10.5	Concentration, acidic: pH = 2.6; Base: Citric acid 1% solution pH = 3.0	-	-	-

D-W: Drinking water; P-W: Purified water demineralised

Checking, maintenance, testing

- Allow product to cool down to room temperature.
- Check product after each cleaning and disinfection for: Cleanliness, functioning and damage, e. g. insulation, loose, bent, broken, cracked, worn and detached parts.
- Discard damaged products immediately. After the cleaning process, it is recommended to place the instrument inside a suitable double sterile packaging. This will ensure reuse in the best way possible.

3.05 DISPOSAL

Adhere to national regulations when disposing of or recycling the product, its components or its packaging. *M.blue plus Instruments* are to be disposed of correctly as potentially infectious material in accordance with medical practice as well as the respective regional laws and regulations.

3.06 TROUBLESHOOTING

- ▶ The *M.blue plus Compass* should be placed as centrally as possible over the valve, as otherwise the determined opening pressure may be incorrect.
- ▶ It must be ensured that the arrow markings on the *COMPASS* housing (proximal, distal) correspond with the flow direction in the shunt system.
- ▶ To ensure a correct reading of the set pressure level, the float gauge in the *M.blue plus Compass* should be centred by slightly moving the instrument within the designated circular marking.
- ▶ The *M.blue plus Compass* is marked with two scale rings:
 - ▶ The grey adjustment range from 0 to 20 cmH₂O on the outer scale applies to the pressure levels of the *proGAV 2.0* adjustable differential pressure unit.
 - ▶ The blue marked adjustment range from 0 to 40 cmH₂O of the inner scale applies to the pressure levels of the *M.blue* gravitational unit.
- ▶ When adjusting the differential pressure unit of the *proGAV 2.0*, care must be taken that the opening pressure is changed by a maximum of 8 cmH₂O per adjustment process, otherwise errors can occur.
Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 3 to 11 cmH₂O and subsequently from 11 to 18 cmH₂O.
- ▶ When adjusting the *M.blue* gravitational unit, care must be taken that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process, otherwise errors can occur.
Example: The opening pressure is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 6 to 22 cmH₂O and subsequently from 22 to 36 cmH₂O.
- ▶ The *M.blue plus Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus Adjustment Ring* should not be placed in the immediate vicinity of the *M.blue plus Compass* while determining the opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin can make adjustment difficult for a few days after the operation. If the valve setting cannot be checked conclusively using the *M.blue plus Compass*, we recommend carrying out a check using an imaging method.
- ▶ When implanting the *M.blue* or *pro-GAV 2.0* valves, ensure that the distance between the surface of the valve and the skin is no more than 10 mm. If the skin and tissue is thicker than 10 mm, it may be difficult or impossible to locate the valve and read and set the valve pressure level.
- ▶ The adjustable Miethke valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e. g. underneath heavily scarred tissue).

3.07 TECHNICAL INFORMATION

3.07.01 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>M.blue plus Instruments</i>
Indication	Treatment of hydrocephalus
Sterilisability/disinfectability	<i>M.blue Checkmate</i> sterilisable <i>M.blue plus Compass</i> disinfectable <i>M.blue plus Adjustment Ring</i> disinfectable <i>M.blue plus Adjustment Assistant</i> disinfectable
Storage	Store in a clean and dry place
Dimensions	<i>M.blue plus Compass</i> : ø 58.7 mm, H: 16.4 mm <i>M.blue plus Adjustment Ring</i> : ø 43 mm; H: 11.5 mm <i>M.blue plus Adjustment Assistant</i> : ø 23.5 mm; H: 28.5 mm <i>M.blue Checkmate</i> : ø 20 mm; H: 27.5 mm

3.07.02 MATERIALS THAT COME INTO CONTACT WITH BODILY TISSUE/FLUIDS

<i>M.blue plus Compass</i>	Aluminium (AlMgSi0,5, AlMgSi1)
<i>M.blue plus Adjustment Ring</i>	Aluminium (AlMgSi0,5, AlMgSi1)
<i>M.blue plus Adjustment Assistant</i>	Aluminium (AlMgSi0,5, AlMgSi1), polyurethane
<i>M.blue Checkmate</i>	Titanium (TiAl6V4)

3.08 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxx indicates the number of the responsible Notified Body
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Production lot number, batch code
	Catalogue number
	UDI number (Unique Device Identifier)
	Sterilised using steam or dry heat
	Do not resterilise
	Do not reuse
	Do not use if package is damaged
	Store in a dry place
	Temperature limitation
	Consult instructions for use / electronic instructions for use
	Caution, read the documentation supplied with the product
	Pyrogen-free
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians.

Symbol	Explanation
	MR Unsafe
	Date
	Model number / European Medical Device Nomenclature Code

4.00 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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- *M.blue plus Compass*
- *M.blue plus Adjustment Ring*

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- *M.blue plus Adjustment Assistant*
- *M.blue Checkmate*

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany