proGAV[®] 2.0 Tools

IN TOUCH WITH YOU.

(DE) Gebrauchsanweisung | (EN) Instructions for use | (FR) Mode d'emploi | (ES) Instrucciones de manejo
(PT) Instruções de utilização | (IT) Istruzioni per l'uso | (NL) Gebruiksaanwijzing

 www.miethke.com

(US) This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0
info@miethke.com

TABLE OF CONTENTS

0.00 PREFACE AND IMPORTANT INFORMATION	18
1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	18
1.01 Explanation of the warnings	18
1.02 Display conventions	18
1.03 Other accompanying documents and additional information	19
1.04 Feedback on the instructions for use	19
1.05 Copyright, disclaimer, warranty and other information	19
2.00 DESCRIPTION OF THE <i>proGAV 2.0 Tools</i>	19
2.01 Medical purpose	19
2.02 Clinical use	19
2.03 Indications	19
2.04 Contraindications	19
2.05 Intended patient groups	19
2.06 Intended users	20
2.07 Intended usage environment	20
2.08 Technical description	20
2.09 System components	20
3.00 PROPERTIES OF THE <i>proGAV 2.0 Tools</i>	21
3.01 Product description	21
3.02 Important safety information	22
3.03 Transport and storage	22
3.04 Use of the product	23
3.05 Disposal	28
3.06 Troubleshooting	29
3.07 Technical information	30
3.08 Symbols used for labelling	31
4.00 MEDICAL DEVICE CONSULTANTS	31

0.00 PREFACE AND
IMPORTANT INFORMATION

Preface

Thank you for purchasing *proGAV 2.0 Tools*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product.

Your team at Christoph Miethke GmbH & Co. KG.

Relevance of the instructions for use

WARNING

Improper handling and use of this product can cause risks and damage. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

proGAV 2.0 Tools include the following components:

- ▶ *proGAV 2.0 Compass*
- ▶ *proGAV 2.0 Adjustment Tool*

Additional options:

- ▶ *proGAV Checkmate*
- ▶ *proGAV Adjustment Disc*
- ▶ *proGAV Compass*

1.00 INFORMATION ON HANDLING THESE
INSTRUCTIONS FOR USE

1.01 EXPLANATION OF THE WARNINGS

DANGER

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

WARNING

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

CAUTION

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

NOTE

Indicates a possibly harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for risk, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

1.02 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates product names

1.03 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

1.04 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

1.05 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

2.00 DESCRIPTION OF THE proGAV 2.0 Tools

2.01 MEDICAL PURPOSE

The *proGAV 2.0 Tools* are intended for adjusting and verifying the adjustable pressure levels of the *proGAV 2.0* valve.

2.02 CLINICAL USE

Treatment options

- ▶ Localising the adjustable valve unit using the instruments
- ▶ Identifying the pressure level of the differential pressure unit of the *proGAV 2.0* valve
- ▶ Checking the pressure level without diagnostic procedures such as X-rays
- ▶ Setting the pressure level without invasive procedures

2.03 INDICATIONS

proGAV 2.0 Tools are used for the following indications:

- ▶ Treatment of hydrocephalus

2.04 CONTRAINDICATIONS

The following contraindications apply to the *proGAV 2.0 Tools* :

- ▶ Intolerance to materials of the instruments

2.05 INTENDED PATIENT GROUPS

- ▶ Patients who, due to their clinical condition, are treated with a CSF draining shunt system with an adjustable *proGAV 2.0* valve

2.06 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical specialists, e. g. neurosurgeons
- ▶ Knowledge about the function and proper use of the product
- ▶ Successful participation in product training

2.07 INTENDED USAGE ENVIRONMENT

Professional Healthcare Facilities

- ▶ Healthcare facilities, such as hospitals, doctors' surgeries or rehabilitation facilities

2.08 TECHNICAL DESCRIPTION

proGAV 2.0 Tools can be used to determine, change and monitor the pressure level of the adjustable differential pressure unit of the proGAV 2.0 valve.

The proGAV 2.0 is a position-dependent valve to drain cerebrospinal fluid (CSF) in the treatment of hydrocephalus. It consists of an adjustable differential pressure unit and a gravitational unit (Fig. 1).

Fig. 1: proGAV 2.0

The proGAV 2.0 Compass (Fig. 2) is used for localising and reading out the adjustment unit in the proGAV 2.0.

Opening the instrument reveals a template that can be used to localise the valve in the patient's head using the index finger. After this, the proGAV 2.0 Compass template is adjusted to the flow direction of the CSF and placed over the valve. Folding the Compass down results in the automatic reading and display of the pressure level.

Fig. 2: proGAV 2.0 Compass

The proGAV 2.0 Adjustment Tool (Fig. 3) can also be used to set the opening pressure of the adjustment unit of the proGAV 2.0 from 0 to 20 cmH₂O.

Fig. 3: proGAV 2.0 Adjustment Tool

The opening pressure of the adjustable differential pressure unit of the proGAV 2.0 can be changed before or after implantation. It is pre-set by the manufacturer at 5 cmH₂O.

2.09 SYSTEM COMPONENTS

proGAV 2.0 Tools can be used to set the opening pressures of the following valves:

Product name	
proGAV 2.0	

3.00 PROPERTIES OF THE *proGAV 2.0 Tools*

3.01 PRODUCT DESCRIPTION

3.01.01 APPROVED COMPONENTS

The *proGAV 2.0 Tools* are comprised of the two components

- ▶ *proGAV 2.0 Adjustment Tool* (1) and
- ▶ *proGAV 2.0 Compass* (2).

proGAV 2.0 Tools can also be supplemented with the

- ▶ *proGAV Checkmate* (3),
- ▶ *proGAV Compass* (4) and the
- ▶ *proGAV Adjustment Disc* (5).

	Item no.	Figure
1	FX400T	
2	FX401T	
3	FV409T	
4	FV403T	
5	FV407T	

3.01.02 SCOPE OF DELIVERY

Box content (item no. FX400T)	Number
<i>proGAV 2.0 Adjustment Tool</i>	1
Instructions for use <i>proGAV 2.0 Tools</i>	1

Box content (item no. FX401T)	Number
<i>proGAV 2.0 Compass</i>	1
Instructions for use <i>proGAV 2.0 Tools</i>	1

Box content (item no. FX404T)	Number
<i>proGAV 2.0 Adjustment Tool</i>	1
<i>proGAV 2.0 Compass</i>	1
Instructions for use <i>proGAV 2.0 Tools</i>	1

Box content (item no. FV409T)	Number
<i>proGAV Checkmate</i>	1
Instructions for use <i>proGAV Tools</i>	1

Box content (item no. FX499T)	Number
<i>proGAV 2.0 Adjustment Tool</i>	1
<i>proGAV 2.0 Compass</i>	1
<i>proGAV Adjustment Disc</i>	1
<i>proGAV Compass</i>	1
Instructions for use <i>proGAV 2.0 Tools</i>	1

3.01.03 STERILITY

With the exception of the *proGAV Checkmate*, *proGAV 2.0 Tools* cannot be sterilised, only disinfected.

WARNING

If the packaging or the product is damaged, the products must not be used under any circumstances.

3.01.04 REPEATED USE AND RE-STERILISATION

The *proGAV Checkmate* must be sterilised in double sterile packaging in an autoclave (steam sterilisation, fractionated vacuum method) at 134 °C and 5 minutes retention time.

3.01.05 OPERATING CONDITIONS

proGAV 2.0 Tools are to be used at room temperature in healthcare facilities, such as hospitals, doctor's surgeries and rehabilitation facilities.

3.01.06 PRODUCT LIFE

These medical devices are designed to operate accurately and reliably over long periods of time. The expected product life of the *proGAV 2.0 Tools* is three years after first use, provided that the product is exposed to normal usage conditions and is properly maintained (see chapter "3.04 Use of the product"). The tools may be used beyond this time if they are in perfect working condition. However, no guarantee can be assumed that these medical devices will not require replacement for medical or technical reasons.

3.01.07 PRODUCT CONFORMITY

The product meets current regulatory requirements. Translations of these instructions for use into additional languages can be found on our website: (<https://www.miethke.com/downloads>).

3.02 IMPORTANT SAFETY INFORMATION

3.02.01 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ If the sterile packaging of the *Checkmate* is damaged, the product must not be used.

WARNING

- ▶ The *Checkmate* must not be used after the expiry date and must be re-sterilised as described under chapter 3.01.04.
- ▶ Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be read carefully and understood before the product is used for the first time.
- ▶ Prior to use, it is essential to check that the *proGAV 2.0 Tools* are complete and intact.

3.02.02 COMPLICATIONS, SIDE EFFECTS AND RESIDUAL RISKS

Various side effects and complications can occur when using *proGAV 2.0 Tools*:

- ▶ Allergic reaction / intolerance to the materials of the instruments

Immediately consult a physician if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or the like.

3.02.03 REPORTING OBLIGATION

Report all serious incidents (damage, injuries, infections, etc.) that occur relating to the product to the manufacturer and the responsible state authority.

3.02.04 PATIENT INFORMATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product.

3.03 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place.

Storage conditions

Temperature range for storage	0 °C to +40 °C
-------------------------------	----------------

3.04 USE OF THE PRODUCT

3.04.01 INTRODUCTION

proGAV 2.0 Tools may only be used by trained specialists.

proGAV 2.0 Tools can be used to determine, change and monitor the pressure level of the adjustable differential pressure unit of the *proGAV 2.0* valve.

3.04.02 SAFETY AND WARNING NOTICES

CAUTION

- The *proGAV 2.0* Adjustment Ring emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.
- Because of the magnets inside the *proGAV 2.0 Tools*, *proGAV 2.0 Tools* must not be used in the vicinity of active implants, such as cardiac pacemakers. Furthermore, there is risk of damage to MRI instruments in their vicinity. For that reason, the use of *proGAV 2.0 Tools* is not permitted there!
- It is essential only to use the *proGAV 2.0 Tools* or *M.blue plus Instruments* to adjust and check the pressure level of the differential pressure unit of the *proGAV 2.0*. If the *proGAV 2.0* is combined with an adjustable gravitational unit (*proSA* or *M.blue*), *proSA Tools* must be used to adjust and check the adjustable pressure level of the *proSA* valve. *M.blue plus Instruments* must be used to check and adjust the adjustable pressure level of the *M.blue*.

3.04.03 REQUIRED MATERIALS

Checking and adjusting the adjustable differential pressure unit of the *proGAV 2.0* after implantation requires the *proGAV 2.0 Compass* and the *proGAV 2.0 Adjustment Tool*.

To adjust the adjustable differential pressure unit of the *proGAV 2.0* before or during valve implantation, the resterilisable *proGAV Check-mate* can be used.

3.04.04 PREPARING FOR USE

Checking the sterile packaging

Immediately before using the *proGAV Check-mate*, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Checking the integrity of the *proGAV 2.0 Tools*

Prior to use, the instruments must be checked for integrity. To this end, all instruments must be visually inspected.

3.04.05 ADJUSTING THE *proGAV 2.0*

Adjusting the opening pressure of the *proGAV 2.0* requires the following steps:

1. Localisation

WARNING

The *proGAV 2.0 Compass* should be placed as centred as possible over the valve, as otherwise the determined opening pressure may be incorrect.

NOTE

- The *proGAV 2.0 Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *proGAV 2.0 Adjustment Tool* should not be placed in the immediate vicinity of the *proGAV 2.0 Compass* when determining the opening pressure. We recommend a distance of approx. 30 cm.
- Swelling of the skin can make adjustment difficult for a few days after the operation. If the valve setting cannot be checked conclusively using the *Compass*, we recommend carrying out a check using an imaging method.

Opening the *proGAV 2.0 Compass* reveals a template that can be used to localise the valve in the patient's head using the index finger (Fig. 4).


Fig. 4: Localising the valve

The *proGAV 2.0 Compass* template is then placed centred over the valve. The "proximal" and "distal" direction markings show the flow direction.

2. Test procedure

If you now fold the *Compass* down, the pressure stage is automatically displayed.


Fig. 5: Determination of the pressure level with the *proGAV 2.0 Compass*

proGAV Compass

In addition to the *proGAV 2.0 Compass*, there is a further option for checking the set opening pressure. The *proGAV Compass* (Fig. 6) is placed on the skin over the implanted valve and moved in small circles until the inner disc is firmly aligned. The opening pressure corresponds to the value displayed in the direction of the ventricular catheter.

NOTE

Small air bubbles in the compass chamber do not impair its function.

Fig. 6: *proGAV Compass*

3. Adjustment process

WARNING

- ▶ When adjusting the *proGAV 2.0* differential pressure unit, care must be taken that the opening pressure is changed by a maximum of 8 cmH₂O per adjustment process, otherwise errors can occur.

Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 3 to 11 cmH₂O and subsequently from 11 to 18 cmH₂O.

The *proGAV 2.0 Adjustment Tool* is centred over the valve. With the index finger, the valve is easy to identify by touch using the recess in the middle of the tool so that the tool is correctly placed. When doing so, the required pressure stage must be displayed on the scale in the direction of the valve inlet or the ventricle catheter. The rotor brake is released by lightly pressing the index finger on the adjustment unit and the pressure level of the *proGAV 2.0* is changed.

Fig. 7: Adjustment using the proGAV 2.0 Adjustment Tool


Fig. 8: Adjustment using the proGAV 2.0 Adjustment Tool

The differential pressure unit of the proGAV 2.0 is fitted with a feedback mechanism. If targeted pressure is exerted on the valve, due to the valve housing design an acoustic signal (a clicking sound) is audible or a resistance can be felt, as soon as the rotor brake has been released. The valve thus shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released, the rotor is once again adjustment-proof. Although the click caused by releasing the rotor brake is easily audible before implantation, it may be considerably reduced after implantation and the filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patient or by using a stethoscope.


Fig. 9: Rotor rotation during adjustment

- a) wrong direction
- b) correct direction

Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the set pressure level. To do this, proceed as in Points 1 and 2. Should the measured value not correspond with the required pressure level, the adjustment process should be repeated. To do this, start again at Point 3.

proGAV Adjustment Disc

The pressure level can also be changed using the proGAV Adjustment Disc. The proGAV Adjustment Disc is centred over the valve. When doing so, the desired pressure level must point in the direction of the ventricular catheter. Slightly pressing the proGAV Adjustment Disc onto the valve releases the rotor brake and changes the pressure level.

To use the Adjustment Disc for adjustment, it must be centred over the valve so that the desired pressure level (b) points in the direction of the inlet connector or the ventricular catheter (c) (Fig. 10).


Fig. 10: Adjustment with the proGAV Adjustment Disc

Press on the valve with the Adjustment Disc (Fig. 11) without changing its position. Then

remove the *Adjustment Disc* and check the pressure level as described in Points 1 and 2.


Fig. 11: Placing and pressing the *Adjustment Disc* onto the differential pressure unit

proGAV Checkmate CE 0297

The *proGAV Checkmate* is supplied sterile and can be re-sterilised. The *proGAV Checkmate* can be used to change and check the pressure level before and during valve implantation directly on the *proGAV 2.0*. To determine the pressure level, the *proGAV Checkmate* is centred over the *proGAV 2.0*. The *proGAV Checkmate* automatically aligns itself over the valve. The pressure level can be read from the direction of the proximal catheter (leading towards the valve). If the pressure level is to be adjusted, the *proGAV Checkmate* is centred over the *proGAV 2.0*. When doing so, the required pressure level must point towards the proximal catheter (leading towards the valve). By slightly pressing the *proGAV Checkmate* onto the valve, the rotor brake in the *proGAV 2.0* is released and the pressure level is set.

Fig. 12: *proGAV Checkmate*

3.04.06 CLEANING RECOMMENDATION FOR NON-STERILISABLE *proGAV 2.0 Tools*

NOTE

proGAV 2.0 Tools are made from thermolabile components that are sensitive to heat and moisture as well as chemically reactive. Do not immerse the *proGAV 2.0 Tools* in cleaning solutions or allow liquids to enter the housing, as moisture, corrosion and impurities could impair proper functioning.

Surface contamination on the *proGAV 2.0 Tools* should always be removed immediately after use by wiping with alcohol-based cleaning agents (minimum 75% alcohol).

The exposure time should be at least 60 seconds and should be adapted to the degree of contamination. The instruments should then be wiped with a dry cloth.

The following cleaning methods are not suitable for processing the *proGAV 2.0 Tools* (except the *proGAV Checkmate*): irradiation, ultrasound, sterilisation, mechanical processing, immersion in cleaning fluids.

3.04.07 CLEANING AND DISINFECTING THE *proGAV Checkmate*

Avoid damaging the product due to unsuitable cleaning agents/disinfectants and/or too high temperatures!

- ▶ Use cleaning agents/disinfectants approved for surgical steels in accordance with the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and contact time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 55°C.
- ▶ Perform ultrasonic cleaning:
 - ▶ as an effective mechanical support for manual cleaning/disinfection
 - ▶ for precleaning products with dried residues prior to mechanical cleaning/disinfection
 - ▶ as integral mechanical support during mechanical cleaning/disinfection
 - ▶ for the subsequent cleaning of products with residues still present after mechanical cleaning/disinfection

- Provided the instruments can be fastened securely or safely for cleaning in machines or on storage aids, clean and disinfect the instruments mechanically.

Manual cleaning/disinfection

- Check visible surfaces for residues after manual cleaning/disinfection.
- If necessary, repeat the cleaning process.

Mechanical cleaning/disinfection

- Place product on suitable screen basket.

Mechanical cleaning/disinfection with manual precleaning

- Manual precleaning with ultrasound

Phase	I	II
Step	Disinfected ultrasonic cleaning	Intermediate rinsing
T (°C/°F)	RT (cold)	RT (cold)
t (min)	15	1
Conc. (%)	2	-
Water quality	D-W	D-W
Chemical specifications	B. Braun Stabimed®; free of aldehydes, phenols and QAV; pH = 9	

D-W: Drinking water; RT: Room temperature

Mechanical alkaline cleaning and thermal disinfection

Phase I

- Clean product in the ultrasonic cleaning bath (frequency 35 kHz). Ensure that all accessible surfaces are wetted.

Phase II

- Fully rinse/flush product (all accessible surfaces) under running water.

Device type: Single chamber cleaning/disinfection device without ultrasound

- Place product on suitable screen basket (avoid missed areas)

Phase	I	II	III	IV	V	VI
Step	Pre-rinsing	Cleaning	Neutralisation	Inter-mediate inspection	Thermal disinfection	Drying
T (°C/°F)	< 25/77	55/131	20/68	70/158	94/201	90/194
t (min)	3	10	2	1	10	40
Water quality	D-W	P-W	P-W	P-W	P-W	-
Chemical specifications	-	Concentration, alkaline: pH = 10.9 < 5% anionic surfactants 1% solution pH = 10.5	Concentration, acidic: pH = 2.6; Base: Citric acid 1% solution pH = 3.0	-	-	-

D-W: Drinking water, P-W: Purified water demineralised

Checking, maintenance, testing

- Allow product to cool down to room temperature.
- Check product after each cleaning and disinfection for:
Cleanliness, functioning and damage, e. g. insulation, loose, bent, broken, cracked, worn and detached parts.
- Discard damaged products immediately.
After the cleaning process, it is recommended to place the instrument inside a suitable double sterile packaging. This will ensure reuse in the best way possible.

3.05 DISPOSAL

Adhere to national regulations when disposing of or recycling the product, its components or its packaging.

3.06 TROUBLESHOOTING

- ▶ The *proGAV 2.0 Compass* should be placed as centred as possible over the valve, as otherwise the determined opening pressure may be incorrect.
- ▶ It must be ensured that the arrow markings on the compass housing (proximal, distal) correspond with the flow direction in the shunt system.
- ▶ When adjusting the differential pressure unit of the *proGAV 2.0*, care must be taken that the opening pressure is changed by a maximum of 8 cmH₂O per adjustment process, otherwise errors can occur.
Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 3 to 11 cmH₂O and subsequently from 11 to 18 cmH₂O.
- ▶ The *proGAV 2.0 Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *proGAV 2.0 Adjustment Tool* should not be placed in the immediate vicinity of the *proGAV 2.0 Compass* when determining the opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin can make adjustment difficult for a few days after the operation. If the valve setting cannot be checked conclusively using the *proGAV 2.0 Compass*, we recommend carrying out a check using an imaging method.
- ▶ When implanting the *proGAV 2.0* valve, the distance between the surface of the valve and the skin should not exceed 10 mm. If the skin and tissue is thicker than 10 mm, it may be difficult or impossible to locate the valve and read and set the valve pressure level.
- ▶ The adjustable Miethke valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e.g. underneath heavily scarred tissue).

3.07 TECHNICAL INFORMATION

3.07.01 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>proGAV 2.0 Tools</i>
Indication	Treatment of hydrocephalus
Sterilisability/disinfectability	<i>proGAV 2.0 Compass</i> disinfectable <i>proGAV 2.0 Adjustment Tool</i> disinfectable <i>proGAV Checkmate</i> sterilisable <i>proGAV Adjustment Disc</i> disinfectable <i>proGAV Compass</i> disinfectable
Storage	Store in a clean and dry place
Dimensions	<i>proGAV 2.0 Compass</i> : ø 49.7 mm, H: 8.5 mm <i>proGAV 2.0 Adjustment Tool</i> : ø 65.7 mm, H: 20.3 mm <i>proGAV Checkmate</i> : ø 20.0 mm, H: 27.5 mm <i>proGAV Adjustment Disc</i> : ø 37.5 mm, H: 21.6 mm <i>proGAV Compass</i> : ø 57.8 mm, H: 14.7 mm

3.07.02 MATERIALS THAT COME INTO CONTACT WITH BODILY TISSUE/FLUIDS

<i>proGAV 2.0 Compass</i>	Aluminium (AlSi1MgMn), Makrolon 2458, sapphire
<i>proGAV 2.0 Adjustment Tool</i>	Makrolon 2458, 1.4021, aluminium
<i>proGAV Checkmate</i>	Titanium (TiAl6V4)
<i>proGAV Adjustment Disc</i>	PPSU, PUR
<i>proGAV Compass</i>	AlSi1MgMn, mineral glass MS 430

3.08 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxx indicates the number of the responsible Notified Body
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Production lot number, batch code
	Catalogue number
	UDI number (Unique Device Identifier)
	Sterilised using steam or dry heat
	Do not resterilise
	Do not reuse
	Do not use if package is damaged
	Store in a dry place
	Temperature limitation
	Consult instructions for use / electronic instructions for use
	Caution, read the documentation supplied with the product
	Pyrogen-free
	Free of natural rubber latex, latex-free

Symbol	Explanation
	Indicates that in the USA, the product may only be issued to physicians.
	MR Unsafe
	Date
	Model number / European Medical Device Nomenclature Code

4.00 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany