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Content

Indication	20
<i>proGAV VERIFICATION TOOL</i>	20
<i>proGAV MASTERDISC</i>	20
<i>proGAV VERIFICATION COMPASS</i>	20
<i>proGAV ADJUSTMENT TOOL</i>	21
<i>proGAV ADJUSTMENT DISC</i>	21
<i>proGAV CHECK-MATE</i>	22
Cleaning and disinfecting the	
<i>proGAV CHECK-MATE</i>	22
Cleaning recommendation for <i>proGAV</i>	
instruments which are not sterilisable	24
Recommendation for sterilisation	24
Recommendation for sterilisation of the	
<i>proGAV CHECK-MATE</i>	24
Technical description	25
Function of the <i>proGAV</i>	26
Physics background	27
Selecting the appropriate <i>proGAV</i>	28
Adjusting the <i>proGAV</i>	29
Reading the pressure setting from an X-ray image	30
Medical products consultant	31
General information	31
Recommendation of pressure levels	32
Variants	34

Indication

The *proGAV* is an implant for the treatment of hydrocephalus. With the *proGAV TOOL SET* the opening pressure of the valve can be verified, changed and controlled.

proGAV VERIFICATION TOOL

The *proGAV VERIFICATION TOOL* is used for reading the valve opening pressure setting. First, the *proGAV VERIFICATION TOOL* should be placed on the valve. The marking (3) on the *proGAV VERIFICATION TOOL* must be in line with the proximal (ventricular) catheter. The tool contains two magnets. As soon as the button on the instrument is pushed (2) the magnets in the tool align with the magnets in the valve. The opening pressure is shown on the scale (1).


Fig. 1: *proGAV VERIFICATION TOOL*

proGAV MASTERDISC

The *proGAV VERIFICATION TOOL* can be easily checked by using the *proGAV MASTERDISC* before measuring the opening pressure of the valve. On the disc the positions 0, 5, 10, 15 and 20 cmH₂O are indicated. If the *proGAV VERIFICATION TOOL* is placed on the disc the opening pressure shown by the instrument should be aligned to the value of the *proGAV MASTERDISC*. Example: The *proGAV VERIFICATION TOOL* is put on the *proGAV MASTERDISC* so that the marking on the instrument is in line with the value 10 cmH₂O on the *proGAV MASTERDISC*. The *proGAV VERIFICATION TOOL* should indicate the value of 10 cmH₂O.

Fig. 2: *proGAV MASTERDISC*

proGAV VERIFICATION COMPASS

Alongside the *proGAV VERIFICATION COMPASS* there is an additional device for measuring the adjusted opening pressure. The compass can be used to locate the valve when palpation is not possible. The measuring compass is set on the skin over the implanted valve and moved in a circling motion until the internal disc is fixed over the valve. The opening pressure corresponds to the value indicated towards the direction of the ventricular catheter.

Fig. 3: *proGAV VERIFICATION COMPASS*

proGAV ADJUSTMENT TOOL

The *proGAV ADJUSTMENT TOOL* is used for adjusting the valve opening pressure. First the intended pressure setting is selected at the knurled dial (1), the opening pressure is shown on a scale (2). Then the *proGAV ADJUSTMENT TOOL* is placed central on the valve. By pushing the button (3), the adjustment tip (4) appears, the brake is decoupled, the rotor turns and the adjusted pressure is achieved.

The marking (5) on the *proGAV ADJUSTMENT TOOL* has to point towards the proximal catheter (leading to the ventricle).

Warning note

If the implantation site is poorly selected or if the skin over the valve is too thick, adjustment of the *proGAV* can be difficult or sometimes impossible. The valve works then as gravitational valve with constant opening pressure. The patient should be informed about the risk.


Fig. 4: *proGAV ADJUSTMENT TOOL*

proGAV ADJUSTMENT DISC

The *proGAV ADJUSTMENT DISC* offers an other option to adjust the pressure setting. The *proGAV ADJUSTMENT DISC* is placed centrally over the valve. The desired pressure setting should be aligned with the proximal catheter (leading to the ventricle). By pressing down the *proGAV ADJUSTMENT DISC* on the valve, the brake is decoupled and the opening pressure of the *proGAV* is changed.


Fig. 5: *proGAV ADJUSTMENT DISC*

a) size S b) size M c) size L

Ensure that the pressure setting is changed no more than 8 cmH₂O. (see „Adjusting the *proGAV*“)

Use a longer tip (size L) for thicker skin. (Fig. 5)

proGAV CHECK-MATE

The *proGAV CHECK-MATE* is intended to be re-sterilised. It is possible to change and to verify an applied pressure setting on the valve directly. To verify the actual pressure setting the *proGAV CHECK-MATE* has to be put centrally over the valve. The *proGAV CHECK-MATE* will immediately start to move. If it remains stable, the pressure setting can be read in alignment to the inlet connector.

To adjust a new pressure setting, the *proGAV CHECK-MATE* has to be placed centrally over the valve. The new pressure setting has to point towards the proximal catheter (leading to the ventricle). By pressing down slightly the *proGAV CHECK-MATE*, the brake of the valve is decoupled, the rotor turns and the opening pressure of the *proGAV* is changed. Please be aware that the steps for changing the pressure setting should not be more than 8 cmH₂O per step (see page 30).

Fig. 6: *proGAV CHECK-MATE*

Cleaning and disinfecting the *proGAV CHECK-MATE*

Avoid damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for surgical steels according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55°C.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning for products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect instruments mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

Manual cleaning/disinfecting

- Check visible surfaces for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Mechanical cleaning/disinfecting

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).

Mechanical cleaning/disinfecting with manual pre-cleaning

Manual pre-cleaning with ultrasound

Stage	Step	T (°C/°F)	t (min)	Conc. (%)	Water quality	Chemical
I	Desinfecting ultrasound cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	

D-W: Drinking water, RT: Room temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots)

Stage	Step	T (°C/°F)	t (min)	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	-Concentrate, alkaline: pH = 10.9 <5% anionic tensides -1% solution: pH = 10.5
III	Neutralization	20/68	2	FD-W	- Concentrate, acid: pH = 2,6 Basis: Citric acid - 1% solution: pH = 3,0
IV	Intermediate rinse	70/158	1	FD-W	-
V	Thermal disinfecting	94/201	10	FD-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water, FD-W: Fully desalinated water (demineralized)

Control, care and inspection

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Set aside the product if it is damaged.

We recommend putting the instrument in a double sterile bag after cleaning. So the next use in the theatre (sterile area) will be prepared as its best.

Cleaning recommendation for proGAV instruments which are not sterilisable**Attention:**

proGAV instruments are made from thermal unstable components which are affectable by heat or humidity or chemical aggressive substances.

Do not steep the proGAV instrument in liquids and keep the inside of the instruments dry!

Remove surface pollution of the proGAV instruments after the use immediately with alcohol based cleaners (more than 75% alc.) by a wiping procedure.

The time of impact should be more than 60 sec. and should be depending on the level of pollution. For final cleaning use a dry wipe.

The following cleaning methods are **not allowed** for the cleaning of the proGAV instruments (except CHECK-MATE): Irradiation, Ultrasonic, Sterilization, Machine preparation, Inserting into liquids.

Recommendation for sterilisation

Except the proGAV CHECK-MATE, proGAV instruments can not be sterilised.

Recommendation for sterilisation of the proGAV CHECK-MATE

The instrument should be sterilised by steam at 134°C.

Technical description of the proGAV

The *proGAV* is a posture dependent hydrocephalus valve. It comprises a ball-cone valve with adjustable opening pressure and a fixed pressure gravitational unit. In this way, optimal CSF drainage is ensured for each individual patient in any body position.

The function of both units is described as follows. Fig. 8 shows a schematic cross section of the *proGAV*. The adjustable unit


Fig. 7: The proGAV is a combination of an adjustable unit and a gravitational unit.

is composed of a solid titanium body with a well-trying ball-cone valve (1) integrated in its proximal part. A bow spring (2) defines the opening pressure of the ball-cone valve. The pretensioning of the spring, and thus the valve opening pressure, can be adjusted by turning a rotor (3), with the valve implanted under the patient's skin. The gravitational unit contains a tantalum ball (4), which defines the opening pressure of this valve, and a sapphire ball (5), which ensures the precise closure of the valve. A connector under the silicone catheter (6) always allows a later disconnection of the catheter if necessary.

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- 1 sapphire ball
- 2 torsion bar
- 3 rotor
- 4 tantalum ball
- 5 sapphire ball
- 6 connector under the silicone catheter

Fig. 8: Schematic cross section of the proGAV

Function of the proGAV

The opening pressure of the *proGAV* is composed of the opening pressure of the adjustable unit and the opening pressure of the gravitational unit.


Fig. 9: Gravitational unit in horizontal position

Horizontal posture

When the patient is lying down, the gravitational unit is always open and therefore does not present any resistance to the fluid flow.

Hence, the opening pressure of the *proGAV* is defined by the adjustable unit. The operational principle of the adjustable unit is illustrated in figures 10a and b. In fig. 10a, the ball-cone valve is closed. The drainage is blocked. In fig. 10b, the adjustable unit is shown in the open condition. The patient's IVP is increased and the spring force, which otherwise keeps closed the ball-cone valve, is overcome. The closing ball moves out of the cone and a gap opens to allow drainage.


Fig. 10: Adjustable unit
a) closed b) open

Vertical position

As soon as the patient moves into an upright position, the gravitational unit closes, the opening pressure of the *proGAV* is significantly increased and CSF drainage is blocked (fig. 10a). Only when the sum of the IVP and the hydrostatic pressure exceeds the opening pressure of the *proGAV*, drainage will be possible again (fig. 11b). The opening pressure of the *proGAV* in the vertical position is the sum of the opening pressures of both the adjustable unit and the gravitational unit (weight of the tantalum ball).


Fig. 11: Gravitational unit in upright position a) closed, b) open

Physics background

The intraventricular pressure is positive in a healthy human in a horizontal position. To adjust this pressure through shunt drainage, one has to choose the appropriate pressure range, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig. 12).

The ventricular pressure in a healthy human in a vertical position becomes slightly negative. To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure.

- IVP Intraventricular pressure
- P_{Vli} Opening pressure in horizontal position (adjustable unit only)
- P_{Vst} Opening pressure in vertical position (adjustable + gravitational unit)
- P_B Pressure in the abdominal cavity
- P_{Hyd} Hydrostatic pressure

horizontal: $IVP = P_{Vli} + P_B$
standing: $IVP = P_{Hyd} - P_{Vst} - P_B$


Fig. 12: Pressure situation at the shunt for horizontal and vertical position

Selecting the appropriate *proGAV*

The *proGAV* is a position-dependent shunt, meaning the opening pressure changes depending on the position of the patient.

To choose the *proGAV* for an individual patient, one opening pressure is set for the horizontal position (patient lying down), and one for the vertical position (patient standing upright).


Horizontal position:

The opening pressure for the horizontal position is defined by the adjustable unit. The pressure level should be chosen according to the clinical situation and indication. The unit can be adjusted to a pressure setting between 0 cmH₂O and 20 cmH₂O.

Vertical position:

The opening pressure of the *proGAV* for the vertical position is calculated by the sum of the opening pressure of the adjustment and of the gravitational unit. The selection of the gravitational unit depends on the activity and the abdominal pressure (adiposity). The following opening pressure ranges for the gravitation unit are possible, the pressure range selected can be checked postoperatively on X-ray image:

Opening pressure for vertical posture	Coding of gravitational unit
10 cmH ₂ O	small, no ring
15 cmH ₂ O	large, no ring
20 cmH ₂ O	large, 1 ring
25 cmH ₂ O	large, 2 rings
30 cmH ₂ O	large, 3 rings
35 cmH ₂ O	large, 4 rings

Coding ring

a)

b)

Fig. 13: X-ray image of the gravitational unit a) large, 1 ring = 20 cmH₂O, b) small = 10 cmH₂O

Adjusting the proGAV

Each proGAV is calibrated under strict quality control procedures.

The presetting of the adjustable unit is 5 cmH₂O, but it must be checked before implantation. The setting is changed in the following steps:

1. The valve is located under the skin (fig. 14).


Fig. 14: Locating the valve

2. The proGAV VERIFICATION TOOL is positioned centrally on the valve. The marking on the instrument has to point towards the proximal catheter (leading to the ventricle). The button is pushed and the pressure setting is read (fig. 15).


Fig. 15: Measuring the pressure setting a) correctly, b) incorrectly

Caution: Placing the proGAV VERIFICATION TOOL in a non-central position on the valve can lead to erroneous readings!

3. The proGAV ADJUSTMENT TOOL is set to the required opening pressure by turning the knurled dial of the unit (fig. 16).

Abb. 16: Adjustment of the proGAV VERIFICATION TOOL

Caution: The new opening pressure setting of the valve must not differ from the measured opening pressure by more than 8 cmH₂O in any one setting (see page 30).

4. The proGAV ADJUSTMENT TOOL is positioned centrally on the valve. The marking on the instrument has to point towards the proximal catheter (leading to the ventricle).

5. As soon as the proGAV ADJUSTMENT TOOL has been positioned centrally on the valve, the button is pushed onto the valve and the adjustment tip appears. This triggers the mechanical decoupling and the valve setting is adjusted to the required opening pressure setting (fig. 17).

Caution: By pushing the button the instrument should stay very close to the valve.


Fig. 17: Adjusting the pressure setting

6. After the adjustment, the valve opening pressure is measured again, as described in step 2. If the pressure measured now differs from the intended pressure level, the adjustment procedure has to be repeated from step 3.

Caution: Due to postoperative swelling of the skin the adjustment of the valve setting may be difficult within the first few days.

During the adjustment the opening pressure of the adjustable unit should not be changed more than 8 cmH₂O per adjustment procedure. Example: Opening pressure is to be changed from 3 to 18 cmH₂O. With only one adjustment procedure the rotor would turn in a counter clockwise-

se direction (short way) and would stop at the position 0 cmH₂O (fig. 18a). The correct adjustment is in 2 steps: Adjustment from 3 to 11, and from 11 to 18 cmH₂O. The rotor turns now in a clockwise direction (fig. 18b).


Fig. 18: Rotor rotation during adjustment
a) false b) correct

Reading the pressure setting from an X-ray image

The pressure setting of the *proGAV* should be measured with the *proGAV VERIFICATION TOOL*, but if there is any discrepancy between the desired adjustment setting and the setting that is read by the *proGAV VERIFICATION TOOL*, then a radiographic confirmation (x-ray) can be performed to confirm the actual valve setting. The magnets are visible in the image as white dots. The direction of the rotor tip indicates the pressure setting. The rotor tip can take any position outside region (Fig. 19). Thus, the opening pressure of the *proGAV* can be adjusted in increments of 1 cmH₂O between 0 and 20 cmH₂O. In order to avoid misidentification of the adjusted opening pressure in the X-ray image, the valve is flattened on one side.


Fig. 19: Schematic X-ray image

Medical products consultant

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
D-14469 Potsdam
Phone: +49(0) 7000 6438453 or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap Platz
D-78532 Tuttlingen
Phone: +49 (0) 7461 95-9
Fax: +49 (0) 7461 95-26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

General information

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>TOOLS for proGAV</i>
Intended use	Treatment of Hydrocephalus
Store in a clean, dry place	


Recommendation of pressure levels

Adjustable unit

Standard (children and NPH-patients)	5 cmH₂O
Defensive (e.g. patients with extremely wide ventricles and highly elevated ICP or aqueductal stenosis)	10 cmH₂O
Special (e.g. patients with pseudotumor cerebri)	15 cmH₂O

Body-Mass-Index (BMI)

The recommendations are based on common patient treatments, but can vary depending on the individual patient's condition.

Gravitational unit

Children up to 5 years	20 cmH₂O
Children over 5 years and	25 cmH₂O
Adults up to 60 years	
Adults over 60 years	20 cmH₂O

Overweight	The peritoneal pressure inhibits drainage. Therefore the gravitational unit should consider the following adjustments concerning adipose patients as a function of body-mass-index (BMI): 25-29 BMI -5 cmH₂O below standard recommendation 30-34 BMI -10 cmH₂O below standard recommendation 35-39 BMI -15 cmH₂O below standard recommendation >40 BMI only apply the adjustable unit (no gravitational unit necessary)
Mobility	Standard pressure levels are suitable for active people. Patients who are bedridden should be treated with a maximum pressure of 20 cmH₂O in the gravitational unit.
Height	The hydrostatic suction effect normally depends on the height. Therefore we recommend the following corrections for the gravitational unit: < 1,60 m height: - 5 cmH₂O below standard recommendation > 1,80 m height: +5 cmH₂O above standard recommendation
Control	The adjustable unit is adjusted to 5 cmH₂O by the manufacturer. Please control the pressure setting before use.
Body position	Only the adjustable unit determines the pressure level in the horizontal body position. The analogous pressure level for the vertical body position results from the addition of both - adjustable and gravitational unit.
Adjustment	Changing the adjustable unit also influences the opening pressure level in the vertical body position.

Variants**proGAV VERIFICATION TOOL****proGAV ADJUSTMENT TOOL****proGAV VERIFICATION COMPASS****proGAV ADJUSTMENT DISC****proGAV MASTERDISC****proGAV CHECK-MATE**

CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
CE marking according to directive 93/42/EEC
Label CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE

Technische Änderungen vorbehalten
Technical alterations reserved
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche

Distributed by

B | BRAUN

AESCULAP AG
Am Aesculap-Platz
D-78532 Tuttlingen
Tel.: +49 (0) 7461 95-0
Fax: +49 (0) 7461 95-26 00
e-mail: information@aesculap.de
www.aesculap.de

Manufacturer acc. MDD 93/42/EEC

CHRISTOPH MIETHKE GMBH & CO. KG

Ulanenweg 2
D-14469 Potsdam
Tel.: +49 (0) 7000 MIETHKE or
Tel.: +49 (0) 331 62083-0
Fax: +49 (0) 331 62083-40
e-mail: info@miethke.com
www.miethke.com