proSA Tools

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Ⓜ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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CONTENT

INDICATION	20
<i>proSA Verification Tool</i>	20
<i>proSA Masterdisc</i>	20
<i>proSA Compass</i>	20
<i>proSA Adjustment Tool</i>	21
<i>proSA Adjustment Disc</i>	21
TECHNICAL DESCRIPTION	22
PHYSICS BACKGROUND	23
FUNCTION OF THE VALVE	24
SELECTING THE APPROPRIATE OPENING PRESSURE OF THE VALVE	26
ADJUSTING THE VALVE	28
READING THE PRESSURE SETTING FROM AN X-RAY IMAGE	30
<i>proSA Check-mate</i> 	31
CLEANING AND DISINFECTING THE <i>proSA Check-mate</i>	31
RECOMMENDATION FOR STERILISATION OF THE <i>proSA Check-mate</i>	32
CLEANING RECOMMENDATION FOR <i>proSA Tools</i> WHICH ARE NOT STERILISABLE	32
MEDICAL PRODUCTS CONSULTANT	32
GENERAL INFORMATION	32

INDICATION

The proSA is intended to divert cerebrospinal fluid (CSF) from the lateral ventricles of the brain into the peritoneum. It is the first hydrocephalus valve that is adjustable for the vertical body position, respectively for inclined body position.

proSA Verification Tool

The verification tool is used for reading the valve opening pressure setting (Fig. 1). Firstly, it is essential that the verification tool is placed centrally over the valve. The notch (3) on the verification tool must be in line with the proximal (ventricular) catheter. The tool contains two magnets. As soon as the button (2) on the instrument is pushed the magnets in the tool align with the magnets in the valve. The opening pressure is shown on the scale (1). The arrow on the underside of the verification tool should be aligned with the direction of flow.

Fig. 1: proSA Verification Tool

proSA Masterdisc

The verification tool can be easily checked by using the proSA Masterdisc before measuring the opening pressure of the valve. On the proSA Masterdisc the positions 0, 10, 20, 30 and 40 cmH₂O (Fig. 2) are indicated. If the verification tool is placed on the disc the opening pressure shown by the instrument should be aligned to the value of the proSA Masterdisc.

Example: The verification tool is put on the Masterdisc so that the marking on the instrument is in line with the value 10 cmH₂O on the Masterdisc. The verification tool should indicate the value of 10 cmH₂O.

Fig. 2: proSA Masterdisc

proSA Compass

Alongside the verification tool there is an additional device for checking the adjusted opening pressure. The compass can be used to locate the valve when palpation is not possible. The compass is set to the skin above the implanted valve. The opening pressure corresponds to the value indicated towards the direction of the ventricular catheter.

Fig. 3: proSA Compass

Caution: Airbubbles inside the compass do not affect its functionality.

proSA Adjustment Tool

The adjustment tool is used for adjusting the valve opening pressure. First the intended pressure setting is selected at the knurled dial (1), the opening pressure is shown on a scale (2). Then the adjustment tool is placed centrally on the valve. By pushing the button (3), the adjustment tip (4) appears, the brake is decoupled, the rotor turns and the adjusted pressure is set. The marking (5) on the adjustment tool has to point towards the proximal catheter (leading to the ventricle).

Fig. 4: proSA Adjustment Tool

proSA Adjustment Disc

The adjustment disc offers another option to adjust the pressure setting (Fig. 5). The adjustment disc is placed centrally over the valve. The desired pressure setting should be aligned with the proximal catheter (leading to the ventricle). By pressing down the adjustment disc on the valve, the brake is decoupled and the opening pressure of the valve is changed.

colour: red

Fig. 5: proSA Adjustment Disc one size only

Caution: Do not use the proSA Tools nearby MRI scanner, since there is a danger of damaging the MRI scanner due to the magnets inside the proSA Tools.

Warning note for carriers of pacemakers: The magnets inside the proSA Tools can alter the function of a pacemaker.

TECHNICAL DESCRIPTION

The proSA combines a fixed DP-unit or an adjustable DP-unit and an adjustable gravitational unit (Fig. 6). Each valve unit consists of a robust titanium housing and uses the reliable ball-cone principle.

A spiral spring (1) defines the fixed opening pressure of the fixed DP-unit. The sapphire ball (2) ensures the precise closure of the valve.


Fig. 6a: Schematic cross section of the proSA with DP-unit

A bow spring (3) defines the opening pressure of the adjustable DP-unit. The pre-tensioning of the spring and thus the opening pressure can be adjusted by turning the rotor (4), with the valve implanted under the patient's skin.


Fig. 6b: Schematic cross section of the adjustable DP-unit

The adjustable gravitational unit contains a tantalum weight (5), which is connected to the bow spring (3). The weight keeps the sapphire ball in its position. Depending on the body position, the influence of the tantalum weight changes and therefore the opening pressure. By turning the rotor (4) the pretension of the bow spring can be adjusted and the opening pressure can be changed with the valve implanted under the patient's skin.


Fig. 6c: Schematic cross section of the adjustable gravitational unit

PHYSICS BACKGROUND

The intraventricular pressure is positive in a healthy human in the horizontal position. To adjust this pressure through shunt drainage, the appropriate pressure range must be chosen, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (Fig. 7). The ventricular pressure in a healthy human when moving to the vertical position becomes slightly negative.

To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate for the change in pressure differential. This is equal to the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure.

Conventional shunts open immediately as soon as the patient stands up, which can lead to critical overdrainage.

- IVP Intraventricular pressure
- PV_{li} Opening pressure in horizontal position (DP-unit only)
- PV_{st} Opening pressure in vertical position (DP-unit + gravitational unit)
- PB Pressure in the abdominal cavity
- PH_{yd} Hydrostatic pressure

horizontal: $IVP = PV_{li} + PB$
 vertical: $IVP = PH_{yd} - PV_{st} - PB$


Fig. 7: Calculating the intraventricular pressure

FUNCTION OF THE VALVE

Horizontal body position

proSA with DP-unit

The operating principle of the DP-unit is illustrated in Fig. 8 a) and b). Fig. 8 a) shows the fixed DP-unit in the horizontal position. The ball-cone valve is closed and drainage is prevented. If the patient's IVP increases and continues to rise, the spring pressure of the ball-cone unit will be overcome. The sapphire ball will move away from the cone and a gap will open for fluid drainage.


Fig. 8: DP-unit
a) closed and b) open

proSA with adjustable DP-unit

In Fig. 9a the ball-cone unit of the adjustable DP-unit is closed, hence the drainage is blocked. Fig. 9b shows the adjustable DP-unit in an open condition. The patient's IVP is increased and the spring force, which otherwise keeps the ball-cone unit closed, is overcome. The sapphire ball moves out of the cone and a gap opens up to allow drainage.


Fig. 9: adjustable DP-unit
a) closed, b) open

With the *proSA* in the horizontal position (Fig. 10), the weight (1) does not affect the sapphire ball (2). Hence, in this configuration the *proSA* is open, and the force brought to bear by the weight does not counteract the fluid pressure. The ventricular pressure is now controlled solely by the differential-pressure valve installed.


Fig. 10: adjustable gravitational unit in horizontal position

Vertical body position

In the vertical position, the DP-unit and the adjustable gravitational unit of the *proSA* work together. The total opening pressure is the sum of the opening pressures of the DP-unit and the adjustable gravitational unit. If the IVP (intra-ventricular pressure) of the patient and the hydrostatic pressure exceed this opening pressure, the closing ball of the gravitational unit is pushed away from the cone seal and a gap opens for fluid drainage (Fig. 11).

The valve system prevents increases in intra-ventricular pressure above the physiological range.

To ensure optimal CSF drainage for each individual patient in any body position, the opening pressure of the adjustable gravitational unit can be adjusted from 0 to 40 cmH₂O.

The adjustment can be achieved non-invasively using the *proSA Tool Set*, preventing the need for further surgery.



Fig. 11: gravitational unit in vertical body position a) closed b) open

SELECTING THE APPROPRIATE OPENING PRESSURE OF THE VALVE

The *proSA* is a position-dependent valve, that means the opening pressure changes depending on the body position of the patient. Hence, one opening pressure for the horizontal position and one for the vertical position is set to fit the *proSA* individually to the patient.
(see recommendation for pressure levels: www.miethke.com)

Horizontal position

In the horizontal position the gravitational unit doesn't have any resistance. Only the differential pressure unit determines the opening pressure of the whole shunt system.
As standard configuration we recommend a differential pressure unit with opening pressure of 5 cmH₂O.

	opening pressure differential pressure unit	opening pressure gravitational unit
opening pressure of the whole shunt system	standard 5 cmH ₂ O	+ 0 cmH ₂ O
	defensive* 10 cmH ₂ O	
	special** 15 cmH ₂ O	

*e.g. patients with extremely wide ventricles and highly elevated ICP or aqueductal stenosis

**e.g. patients with pseudotumor cerebri

Vertical position

In the vertical position the opening pressure of the complete shunt system is the sum of the opening pressure of the differential unit and the opening pressure of the gravitational unit.
The opening pressure of the gravitational unit should be chosen depending on height, weight and age of the patient.

	opening pressure differential pressure unit	recommendation opening pressure gravitational unit
opening pressure of the whole shunt system	= opening pressure of the differential unit	children up to 5 years 20 cmH ₂ O
		adults up to 60 years 25 cmH ₂ O
		adults over 60 years 20 cmH ₂ O

The recommendations are based on common patient treatments, but can vary depending on the individual patient's condition.

ADJUSTING THE VALVE

Note: In order to locate, to verify and to adjust either the adjustable DP-unit or the adjustable gravitational unit, different verification and adjustment tools are used:

Tools used for the adjustable DP-unit:

proGAV Verification Tool
 proGAV Masterdisc
 proGAV Compass
 proGAV Adjustment Tool
 proGAV Adjustment Tool
 proGAV Check-mate

Tools used for the adjustable gravitational unit:

proSA Verification Tool
 proSA Masterdisc
 proSA Compass
 proSA Adjustment Tool
 proSA Adjustment Tool
 proSA Check-mate

Please verify specifically before using any tool for verifying or adjusting the opening pressure: for the adjustable DP-unit, use only proGAV Tools

and

for the adjustable gravitational unit, use only proSA Tools

The adjustment procedure for the adjustable gravitational unit (in the following named adjustable unit) is outlined in the steps below:

1. Localisation

The adjustable gravitational unit is located under the skin (Fig. 12).


Fig. 12: Locating the adjustable gravitational unit

2. Verifying:

The verification tool is positioned centrally over the valve. (Fig. 13). The notch on the instrument must point towards the proximal catheter (leading to the ventricle).


Fig. 13: Place centrally over the valve

The arrow on the bottom side of the verification tool indicates the direction of CSF-flow. The button is pushed and the pressure setting is read (Fig. 14).


Fig. 14: Measuring the pressure setting
 a) correctly, b) incorrectly

Caution: Failing to centre the verification tool over the valve can lead to erroneous readings!

3. Adjustment of the opening pressure

When adjusting the proSA preoperatively through the packaging, moderate force with the adjustment tool is sufficient. **DO NOT USE THE BUTTON.** Strong pressure can cause damage to the housing, which might affect valve function.

Caution: Due to postoperative swelling of the skin the adjustment of the valve may be difficult within the first few days!

Please ensure that the opening pressure is changed by no more than 16 cmH₂O per each adjustment for the adjustable gravitational unit.

Example: Opening pressure is to be changed from 6 to 36 cmH₂O. With only one adjustment step the rotor would turn in a counter clockwise direction (shortest path) and would stop at the position 0 cmH₂O (Fig. 15a).

A correct adjustment is performed in two steps: Adjustment from 6 to 22 cmH₂O, and then from 22 to 36 cmH₂O. The rotor turns now in a clockwise direction (Fig. 15b).


Fig. 15: Rotor rotation during adjustment
 a) false b) correct

3a. Adjustment with the adjustment tool

The adjustment tool is set to the required opening pressure by turning the knurled dial of the unit (Fig. 16).


Fig. 16: Setting the adjustment tool

The adjustment tool is positioned centrally on the valve. The notch on the instrument and the scale (a) must point towards the proximal catheter (b) (leading to the ventricle), see Fig. 17.


Fig. 17: Positioning the adjustment tool

As soon as the adjustment tool has been positioned centrally over the valve, the button is pushed and the adjustment tip appears to apply pressure to the valve. This triggers the mechanical decoupling of the rotor and the valve setting is adjusted to the required opening pressure setting (Fig. 18).


Fig. 18: Adjusting the pressure setting

Caution: Ensure that the instrument remains close to the valve during the adjustment procedure.

When treating patients who have a low tolerance to pain, local anaesthesia (e.g. applied through a plaster) should be considered, in cases where contraindication can be excluded.

3b. Adjustment with the adjustment disc

Center the adjustment disc over the gravitational unit of the and align the desired pressure setting (b) on top of the disc in direction of the ventricular catheter (c), see Fig. 19a.


Fig. 19a: Adjustment with the Adjustment Disc

For changing the opening pressure, press down the adjustment disc and release (Fig. 19b). Do not press and turn.


Fig. 19b: Press down slightly the Adjustment Disc and release

Finally, remove the adjustment disc and confirm the setting with the verification tool.

4. Verifying after adjustment

After the adjustment, the valve opening pressure has to be measured again, as described in step 2. If the pressure measured now differs from the intended pressure level, the adjustment procedure has to be repeated from step 3.

If the pressure configuration of the valve cannot be determined with complete certainty by the verification tool, the use of imaging techniques is recommended (excluding MRI: danger of artefacts).

MRI examinations must be performed at field strengths no greater than 3.0 tesla.

Caution: If the site of implantation is poorly selected or if the skin over the valve is too thick, an adjustment of the adjustable unit can be difficult or sometimes impossible. The adjustable gravitational unit then behaves like a gravitational unit with a fixed opening pressure for a given position.

READING THE PRESSURE SETTING FROM AN X-RAY IMAGE

Pressure settings for the adjustable gravitational unit

The pressure setting of the adjustable gravitational unit should be checked with the proSA verification tool. If there is any discrepancy between the desired adjustment setting and the setting that is read by the proSA verification tool, then a radiographic confirmation can be performed in addition to confirm the actual valve setting. The position of the rotor tip indicates the opening pressure. The rotor can take any position outside the region indicated shown below, see Fig. 20.


Fig. 20: Schematic X-ray image of adjustable gravitational unit

The opening pressure of the adjustable gravitational unit can be adjusted from 0 to 40 cmH₂O, in increments of 1 cmH₂O.

To further assist correct identification of the adjusted opening pressure, a cut-out in the casing ring is visible under X-ray imaging (Fig. 20b).


Fig. 20b: X-ray image adjustable gravitational unit: setting 31 cmH₂O

proSA Check-mate com

The proSA Check-mate is delivered sterile and is intended to be re-sterilised. It is possible to change and to verify an applied pressure setting on the adjustable DP-unit and on the adjustable gravitational unit directly in the OR.

To verify the actual pressure setting the Check-mate has to be put centrally over the adjustable gravitational unit. The Check-mate will immediately start to move. If it remains stable, the pressure setting can be read in alignment to the inlet connector.

To adjust a new pressure setting, the Check-mate has to be placed centrally over the adjustable DP-unit or the adjustable gravitational unit. The new pressure setting has to point towards the proximal catheter (leading to the ventricle). By pressing down slightly the Check-mate, the brake of the adjustable gravitational unit is decoupled, the rotor turns and the opening pressure of the valve is changed.

Fig. 21: proSA Check-mate, colour: gold pressure range 0-40 cmH₂O

CLEANING AND DISINFECTING THE proSA Check-mate

Avoid damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for surgical steels according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55°C.

Carry out ultrasound cleaning:

- as an effective mechanical supplement to manual cleaning/disinfecting.
- as a precleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
- as an integrated mechanical support measure for mechanical cleaning/disinfecting.
- for additional cleaning for products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect instruments mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

Manual cleaning/disinfecting

- Check visible surfaces for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Mechanical cleaning/disinfecting

- Place the product on a tray that is suitable for cleaning (avoid rinsing blindspots).

Note: proSA Tools can not be sterilised, except for the proSA Check-mate.

RECOMMENDATION FOR STERILISATION OF THE *proSA Check-mate*

Sterilise the *proSA Check-mate* by steam sterilization (fractionated vacuum process) at 134°C and 5 minutes cycle time.

Notice:

For use in the US, please check for further information concerning sterilisation of the *proSA Check-mate* (www.aesculapusa.com)

CLEANING RECOMMENDATION FOR *proSA Tools* WHICH ARE NOT STERILISABLE

Attention: *proSA Tools* are made from thermal unstable components which are affected by heat or humidity or chemical aggressive substances.

Do not steep the *proSA Tools* in liquids and keep the inside of the instruments dry!

Remove surface pollution of the *proSA Tools* after the use immediately with alcohol based cleaners (more than 75% alc.) by a wiping procedure.

The time of impact should be more than 60 sec. and should be depending on the level of pollution. For final cleaning use a dry wipe.

The following cleaning methods are not allowed for the cleaning of the *proSA Tools* (except *proSA Check-mate*):
Irradiation, Ultrasonic, Sterilization, Machine preparation, Inserting into liquids.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co.KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
D-14469 Potsdam
Phone: +49(0) 7000 6438453 or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap Platz
D-78532 Tuttlingen
Phone: +49 (0) 7461 95 0
Fax: +49 (0) 7461 95-26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042
AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>proSA Tools</i>
Intended use	Treatment of Hydrocephalus
Store in a clean, dry place	

CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
CE marking according to directive 93/42/EEC
Label CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE

Technische Änderungen vorbehalten
Technical alterations reserved
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche

Manufacturer acc. MDD 93/42/EEC:

■ CHRISTOPH MIETHKE GMBH & CO. KG

Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Phone +49 (0) 331 62 083-0 | Fax +49 (0) 331 62 083-40 | www.miethke.com

Distributed by:

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 74 61 95-26 00 | www.aesculap.com

Aesculap - a B. Braun company