Lesegerät Set für SENSOR RESERVOIR Reader Unit Set for SENSOR RESERVOIR

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CONTENT

READER UNIT SET FOR SENSOR RESERVOIR	34
SAFE HANDLING	6
PRODUCT DESCRIPTION	6
INTENDED USE	6
INDICATIONS	6
CONTRAINDICATIONS	7
OPERATING PRINCIPLE	7
SYMBOLS	8
USING THE READER UNIT SET FOR SENSOR RESERVOIR	9
SYSTEM SET-UP	10
CONNECTING THE POWER SUPPLY	10
CONNECTING THE ANTENNA	10
INSERTING THE SD-CARD	10
FUNCTION CHECKS	10
OPERATION	11
GENERAL OPERATING INSTRUCTIONS	12
INFO MENU	12
TAKING MEASUREMENTS	12
MANAGING DATA/DISPLAYING MEASUREMENTS	15
SETTINGS	19
CLEANING/DISINFECTION	23
MANUAL CLEANING/DISINFECTION	23
MAINTENANCE/CALIBRATION	24
CALIBRATION OF BARD SENSOR	24
TROUBLESHOOTING	24
MALFUNCTIONS DISPLAYED WITH AN ERROR MESSAGE	25
OTHER MALFUNCTIONS	29
TECHNICAL SERVICE	29
TECHNICAL DATA	30
AMBIENT CONDITIONS	31
CE MARKING	31
STORAGE AND TRANSPORT	31
DISPOSAL	31

Fig. 1 Reader Unit for SENSOR RESERVOIR

- 1 ON/OFF-switch
- 2 Display 3 Antenna socket 4 Antenna
- 5 SD-card slot

- 6 Function keys
- 7 Connector socket for power supply
 8 Connector (for Reader Unit)
 9 Control lamp

- 10 Power supply

Symbol	Explanation		
\wedge	Caution, see documentation supplied with the product		
F	Follow the instructions for use		
((v))	Non-ionized radiation		
X	Labeling of electric and electronic devices according to Directive 2002/96/EC (WEEE), see Chapter Disposal		
	ON/OFF-switch		
	Function key The options on the display can be selected by using the 4function keys		
IP44 IP44	Housing protection class to DIN EN 60529 Reader Unit 7502 0000 Power supply 7504 0000		
	Protection class II (protective insulation) (Power supply 7504 0000 only)		
Y	Antenna symbol Socket for antenna		
	SD-card SD-card slot		
DCIN	DCIN-jack for AC adapterl		

PURPOSE INTENDED

The *Reader Unit for SENSOR RESERVOIR* is used for pressure-based control function.

SAFE HANDLING

Caution: Risk of contamination if the Reader Unit or antenna is not used correctly in the surgical environment!

Caution with sterile and non-sterile areas.

Caution: Risk of injury due to incorrect operation of the product! Attend appropriate product training before using the product. For information regarding such training, please contact:

Christoph Miethke GmbH & Co. KG.

- Clean the new product thoroughly after removing from the transport packaging, see Chapter Manual Cleaning/Disinfection.
- To prevent nosocomial infection and multidrug resistance, the instrument should be disinfected after each use. Disinfection is carried out preferably by wiping. See chapter cleaning / disinfection section disinfectant
- Ensure that the product is in perfect condition and fully functional prior to use, see Chapter Inspection.
- To prevent damage as a result of incorrect installation or operation, and in order not to compromise warranty and manu-facturer liability:
- Use the product only as specified in these instructions for use.
- Observe all safety information and maintenance instructions.
- Use only Miethke products in combination, see Chapter Reader Unit System Components.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible to the user.

PRODUCT DESCRIPTION

System components: Reader Unit Set for SENSOR RESERVOIR die SENSOR PRECHAMBER

Description	Product no.
Reader Unit Set for SENSOR RESERVOIR SENSOR PRECHAMBER	75020000
Antenna Reader Unit Set for SENSOR RESERVOIR SENSOR PRECHAMBER	75030000
Power supply Reader Unit Set for SENSOR RESERVOIR SENSOR PRECHAMBER	75040000

INTENDED USE

The Reader Unit is used together with the telemetric SENSOR-RESERVOIR in shunt function diagnostic tests and permits reliable and non-invasive detection and localization of occlusions in ventricular drainage systems (with no intervention or X-ray diagnostics) as well as identification of mechanical damage to shunt valves.

INDICATIONS

- Functional diagnosis of implanted ventricular drainage systems (shunts)
- Detection and localization of occlusions in ventricular drainage sys-tems
- · Checking of shunt valve function

CONTRAINDICATIONS

- · Intracranial pressure measurement
- Intracranial pressure management based on measured values
- Intracranial pressure measurement without shunt
- Functional diagnosis of extracorporeal liquor drainage systems (trau-matic brain injury)

Fig. 1: The working distance for telemetric data transfer, i.e. the distance between antenna and SENSOR RESERVOIR is 10-30 mm, see Fig. 1.

OPERATING PRINCIPLE

The telemetric system is used to measure liquor pressure levels, the characteristics or relative time-dependent modification of which permit conclusions to be made regarding the function of an implanted ventricular drainage system. Diagnostic methods can be used to provoke and subsequently interpret a change in the level of pressure in the shunt. Using the Reader Unit, the data measured by the shunt sensor can be retrieved and displayed. The measured data are saved automatically on the SD-card and can thus be analyzed at a later point. Based on the relative course of theliquor pressure in the shunt, conclusions can be drawn as to its function. As a result, it is possible to identify and localize any occlusion in the shunt as well as a loss of mechanical function in the shunt valves using noninvasive methods.

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Info	For displaying additional information or starting with the info menu
Menü	For starting with the selection menu
●)))	Fast measurement - Fast measurement can be quickly started without having to select further options
	Upward navigation of the cursor
	Downward navigation of the cursor
↓	For quitting the current menu
ОК	Activates or confirms the selected function
Start	Starts measurement
Stop	Stops the present measurement
í	For deleting a saved measurement in the menu <data management=""></data>
*	For placing a marker during continuous measurement
	To call up the diagram processing menu
	To activate the zoom function
	To activate the cursor function
	For adjusting the presentation of the min. and max. values of the pressure axis
<u>⊥</u> ♀	For returning to the diagram processing menu
* *	For returning to the diagram processing menu
	Navigate forward
Ŧ	Navigate backward
on	On

off	Off
\checkmark	Confirm
~	For maximizing the display of a sequence in the menu <data management=""></data>
~	For minimizing the display of a sequence in the menu <data management=""></data>
×	For deactivating the acoustic signal

USING THE READER UNIT

Danger: Risk of contamination if the Reader Unit orantenna is not used correctly in the surgical environment!

- · Caution with sterile and non-sterile areas!
- Warning note: Limited accuracy and no storage of measured data if operated without the SD-card!
- · Insert the SD-card designated for use with the shunt sensor into the Reader Unit.
- The SD-card designated for use with the implant (see patient pass) is recognizable by its identification number (ID).
- Measured data can be saved only on the designated SD-card.

Warning note: Unexpected termination of measurement if the SD-card is full!

- Check the available memory prior to starting measurement.
- Delete any data that are no longer needed.

Warning note: Danger from use of MRI!

 Only operate the Reader Unit outside of the MRI field of use

Caution: The Reader Unit cannot be operated if the battery status is low!

- · Check the battery status.
- If the battery is low, connect to the power supply

Caution: Communication poor or non-existent!

The telemetric link between the Reader Unit and implant may be disrupted by the operation of a further scanning device.

- Increase the distance between the scanning devices.
- The telemetric link between the antenna and implant may be disrupted by a further implant in the vicinity of the antenna. Increase the distance between the implants.
- The telemetric link between the antenna and implant may be disrupted by metal components in the vicinity of the implant. Increase the distance from metal components!

Note: If the temperature of the telemetrical sensor is below 19°C, it is necessary to fall back on the "Measurement without SD-card". Operation will take place with less precision and it is not possible to save measurement data.

The following error messages occur during measurement with the SD-card:

- During a long-term measurement: Pressure values invalid, measurement will continue
- During an individual or collective measurement: Temperature values invalid, measurement will not continue

For a more detailed description of errors, see Troubleshooting

SYSTEM SET-UP

The *Reader Unit* is equipped with a battery which, when fully charged, permits wireless operation for up to 5 hours. The battery must be charged when first installing the device. The power supply is used for charging the Reader Unit battery. It takes about 6 hours to fully charge the battery. The Reader Unit can also be operated while it is being charged (with connected power supply).

Note: Charging is not possible at an ambient temperature of >35°C!!

Connecting the power supply

Danger: Damage to the Reader Unit from incorrect accessories and/or risk to the safety of user and patient!

 Use only the original power supply belonging to the Reader Unit SENSOR-RESER-VOIR (7504 0000).

The main voltage must be consistent with the voltage range stated on the label of the power supply for the Reader Unit shunt sensor.

- Insert the power supply connector of the device into the connector socket of the Reader Unit, see Fig. 2.
- Connect power supply by inserting into a wall socket.

Fig. 2

Connecting Antenna

Connect antenna to Reader Unit, see Fig. 3.

Inserting the SD-card

Danger: Damage to SD-card if not used correctly! Do not touch the contact strips of the SD-card.

Insert the SD-card designated for use with the relevant shunt sensor (see patient pass) into the Reader Unit until it locks in firmly, see Fig. 4. Note: Press the card lightly to remove.

Fig. 4

FUNCTION CHECKS

- Check that the entire unit is in good condition and fully functional prior to each use.
- To check the battery status, check the function without using the power supply.
- Ensure that the device and antenna used for the function check are notvisibly damaged.
- Check battery status, connect Reader Unit to mains power supply if necessary to charge the battery.
- Check the function of the following components in the given order:
 - 1. Activate the On/Off-switch (1).
 - 2. Automatic self-test after activation, including display and volume

The following is displayed on the screen when operating the reader until with main supply:

Display during power supply operation of Reader Unit

Fig. 5a

The following is displayed on the screen when operating the reader until with main supply:

Fig. 5b

Note: The control lamp of the Reader Unit power supply must be illuminated!

- Check system time, correct if necessary, see Chapter Settings.
- Activate the On/Off-switch (1).

OPERATION

Legend				
A Date	B Time			
C SD-card	D Memory map (on			
	SD-card)			
E Battery status	F Info key			
G Fast measurement	H Menu key			
A B	Ç D Ę			
* *	* * *			
38,85.2814 14:14h				
Deador	Unit for			
Reader Unit for				
SENSOR RESERVOIR				
Christoph Miethke GobH & Co. KG				
and the second se				
Info	AMO Monti			
Into	with the state			
Fig. 6 F	G H			

Activate Reader Unit with On/Off-Taste Start screen content is displayed, see fig. 7.

GENERAL OPERATING INSTRUCTIONS

enu-controlled operation of the device is facilitated by a total of four function keys.

The context-related function of these soft keys is indicated by the symbols shown on the display for each of the keys.

The exact meaning of the symbols is explained in the list "Display symbols". The pre-selected submenu is highlighted within a frame. For improved navigation, the position currently selected in the submenu is highlighted by a scroll bar in addition to the frame.

keys: OK The and are used for standard operation.

Note: When the Reader Unit is not in use it switches into standby-mode. Depending on the setting: 1-5 minutes.

INFO-MENU

The following information can be retrieved from the info menu:

- Hotline & Service
- System specifications
- Manufacturer specifications

To reach the Info menu, press Info

The following information is displayed on the screen:

pre-selected submenu

scroll bar

Fig. 7

Using the keys the desired submenu can be pre-selected and confirmed by pressing **OK**. Return to the previous menu by pressing

Content of 3 submenus

hotline & service hotline & service phone: +49 331 62083-0

System details

Product name: SENSOR RESERVOIR Reader Unit article number: 7502 0000 serial number:0XXXX software version:2.02 service date: tt/mm/jj

Manufacturer info

Manufacturer: Christoph Miethke GmbH & Co. KG Ulanenweg 2 14469 Potsdam Germany

TAKING MEASUREMENTS

To reach the submenu press Menü .

The following information is displayed on the screen:

Fig. 8

To reach the submenu push OK .

The following information is displayed on the screen:

There are three types of measurement:

- Single measurement: The mean of 8-10 measurements is displayed and saved to the SD-card.
- Continuous measurement: Here, the sequence of single measurements at the selected measurement intervals is displayed in a time diagram and saved to the SD-card.
- Fast measurement: Here, the sequence of unaveraged single measurements at the max. available measurement rate (approx. 40 measurements per second) is displayed in a time diagram and saved to the SD-card..

Single measurement

Use **OK** to select the single measurement menu.

The following information is displayed on the screen:

Start a single measurement by pressing **Start** Return to the Selection menu by pressing

The following information is displayed on the screen:

88.85.2914 14:17	hill	8:
Single	Measur	* 1
Date:	88.85.2014	
Tenperature:	14:17 h 25.22 °C	
10:	819C8F14	
Pres.:	3,43	mbar
Pressure Rong	e: DK	
Start		

Fig. 11

Continuous measurement

Select the continuous measurement menu by pressing the keys and and confirm by pressing **OK**.

The following information is displayed on the screen:

Fig. 12

Start continuous measurement by pressing **Start**. Return to the selection menu by pressing **Start**.

Fig. 10

The following information is displayed on the screen:

A marker can be placed by pressing -Multiple markers can be placed during one measurement. Stop measurement by pressing

Stop

Note: The markers permit circumstantial analysis of the measured data.

Symbols for reception quality

Symbol	Erklärung
	Communication started
*	Distance antenna to measuring unit good
↔	Distance antenna to measuring unit too small - increase distance
→←	Distance antenna to measuring unit too far - reduce distance

Fast measurement

There are two ways to start fast measurement:

Press the key: O)) on the start screen and then start measurement with the key: Start By pressing the key: Menü the following subnemu appears:

Fig. 14

To reach the measurement menu, press ок

Select the fast measurement menu by pressing						
the keys		and		and	confirm	by
pressing	ок					

The following information is displayed on the screen

Abb. 15

The following information is displayed on the screen:

Start fast measurement by pressing Start Return to the selection menu by pressin

During fast measurement, the following information is displayed on the screen:

Symbol for reception quality

Fig. 18

Stop measurement by pressing Stop .

Symbols for reception quality:

Please see Chapter Continuous measurement.

MANAGING DATA/DISPLAYING MEASUREMENTS

By pressing **Menü** on the start screen the following submenu appears:

Fig. 19

Select the data management menu by pressing the keys and and confirm by pressing **OK**.

Fig. 20

The data files are listed chronologically (time of starting measurement) and displayed as follows:

Fig. 21

The files are labeled as follows:

Symbol	Explanation	
Ť	Single measurement	
	Continuous measurement	
*** **	Fast measurement	

Note: When loading a continuous or fast measurement, an hour glass appears and the loading time is displayed in seconds.

Single measurement

Fig. 22

The data file can be deleted by pressing

Fig. 23

Deletion must be confirmed by pressing

Return to the selection menu by pressing

Continuous measurement

Fig. 23

The measurement can be deleted by pressing Deletion must be confirmed by pressing Return to the screen shown above by pressing Return to the selection menu by pressing Return.

By pressing **(C)** the diagram processing menu appears.

The following information is displayed on the screent:

Fig. 24

By pressing the zoom menu appears: 28.97.2018 99:22h 29.07.2018 99:22h 29.07.2018 99:18h cont. necsuren. MEM 29.97.2018-09:18h -0.34 nbar

Fig. 25

The and and keys enable the display of a measurement sequence to be expanded or compressed (zoom). Return to the diagram processing menu by pressing Return to the main display for continuous measurement by pressing .

By pressing the key in the diagram processing menu, the following screen appears:

Fig. 26

Return to the diagram processing menu by pressing **C**. The **C** or **C** bkeys enable the time axis to be shifted to the left and right, respectively. Return to the main display for continuous measurement by pressing **C**.

By pressing the **E** ikey in the diagram processing menu, the following screen appears:

Keys and can be used to adjust the max. and min. scale values at the cursor position. Using the **Key**, the cursor can be switched between maximum, minimum and automatic. If the "auto-matic" function is selected, the axes are scaled automatically inline with the measured data. If the "automatic" function is selected, the max. and min. value setting becomes irrelevant.

Return to the display of measured values by pressing

Fast measurement

Fig. 28

The measurement can be deleted by pressing Deletion must be confirmed by pressing Return to the screen shown above by pressing Return to the selection menu by pressing Return to the selection the diagram processing menu appears.

The following information is displayed on the screen:

Fig. 29

By pressing the zoom menu appears.

Fig. 30

The and and keys enable the display of a measurement sequence to be expanded or compressed (zoom). Return to the diagram processing menu by pressing to continuous measurement by pressing .

Note: Switching between displays can take a few seconds, depending on the file size.

By pressing the key in the diagram processing menu, the following screen appears:

Fig. 30

Return to the diagram processing menu by pressing . The constraints or constraints and right, respectively. Return to the main display for continuous measurement by pressing . By pressing the set way in the diagram processing menu, the following screen appears:

Fig. 31

Cursor

Keys and can be used to adjust the max. and min. scale values. Using the key, the cursor can be switched between maximum, minimum and automatic. If the "automatic" function is selected, the axes are scaled automatically in line with the measured data. If the "automatic" function is selected, the max. and min. value setting becomes irrelevant. Return to the display of measured values by pressing

SETTINGS

By pressing **Menü** on the start screen, the following submenu appears:

Fig. 32

Select the settings menu by pressing the keys and and a setting.

Fig. 33

To display the settings menu, press **OK** The following information is displayed on the screen:

Fig. 34

Date/Time

To display the date/time menu, press **OK** The following information is displayed on the screen:

The position of the cursor can be moved by pressing . Using the keys and , he values at the position of the cursorcan be amended. Return to the display of measured values by pressing .

Note: Modified values are saved immediately.

Units

ISelect the units menu in the settings menu using the keys and .

Fig. 36

To display the units menu, press **OK**. The following information is displayed on the screen:

13.07.2010 13:36h units	
□mmHg □cmH20 □hPa	□ Torr ⊠mbar
ок 🛧	

Fig. 37

Select the desired unit using the keys and , and confirm by pressing OK Return to the settings menu by pressing

Measurement interval

These settings function only in the continuous measurement mode.

Select the measurement interval menu in the settings menu using the keys and

Fig. 38

To reach the measurement interval menu, press **OK**. The following information is displayed on the screen:

Fig. 39

Select the desired measurement interval using the keys and , and confirm by pressing **OK**. The available settings are 1-300. The available memory map on the SDcard is also displayed. Return to the settings menu by pressing **C**.

Measurement display

Select the display data menu in the settings menu using the keys and the settings.

13.07.2010 13:35h	
date/tir	ne
units	
measur.	interval
display	data
	State State Sal
OK 1	

Fig. 40

To reach the display data menu press **OK** The following information is displayed on the screen:

Fig. 41

The measurement display can be selected by pressing \square or \square . Return to the settings menu by pressing \blacksquare .

Language

Select the language menu in the settings menu using the keys and .

Fig. 42

To reach the language menu press **OK**. The following information is displayed on the screen:

Volume

Select the volume menu in the settings menu using the keys and .

language volume	
volume	
standby-mode	I
brightness	

Fig. 44

To reach the volume menu press **OK**. The following information is displayed on the screen:

Select the desired volume using the keys and , and confirm by pressing ok. The available settings are 1-5. An acoustic signal will be emitted as the volume changed. Return to the settings menu by pressing

The setting determines the volume of the warning signals. Warning signals are, as a rule, emitted at volume level 5.

Standby-Mode

Select hte standby-mode menu in the settings menu using the keys and the settings.

Fig. 46

To reach the standby-mode menu press **OK**. The following information is displayed on the screen:

Fig. 47

Select the preferred time at which the *Reader Unit* should switch into standby-mode using the keys and and confirm by pressing **OK**. The available settings are 1-5 minutes. Return to the settings menu by pressing **CK**.

Brightness

Select the brightness menu using the keys and .

Fig. 48

To reach the brightness menu press **OK** betätigen. The following information is displayed on the screen:

Fig. 49

Select the desired brightness of the background illumination using the keys and and confirm by pressing **OK**. The available settings are 000–100 (in increments of 5). Return to the settings menu by pressing

Contrast

Select the contrast menu in the settings menu using the keys and the keys.

Fig. 50

To reach the contrast menu press **OK**. The following information is displayed on the screen:

Fig. 51

Select the desired contrast of the display using the keys and confirm by pressing **OK**. The available settings are 000-100 (in increments of 5). Return to the settings menu by pressing **CK**.

CLEANING/DISINFECTION

Danger: Risk of electric shock and fire!

- · Disconnect from the mains before cleaning.
- Do not use flammable or explosive detergents or disinfectants.
- Ensure that no fluids can infiltrate the product.

Caution: Damage to, or destruction of the product caused by mechanical cleaning/ disinfection!

- Only clean and disinfect the product manually.
- · Never sterilize the product.

Caution: Damage to the product due to inappropriate detergents/disinfectants!

Clean the surfaces only with approved detergents/disinfectants as specified by the manufacturer. Observe specifications regarding concentration, temperature and contact time.

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Phase	I
Step	wipe desinfection
T (°C/°F)	RT (room temperature)
t (min)	≥1
Conc. (%)	-
Water quality	-
Chemical	Meliseptol HBV wipes 50% Pro- pan-1-ol

Phase I

- Remove any visible residues with a disposable disinfectant wipe.
- Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- Observe the specified contact time (1 minute minimum).

Inspection

- After each complete cleaning/disinfection cycle, check the product for damage.
- · Set aside the product if it is damaged.

Storage

 Store reprocessed products protected from dust in a dry, dark and cool room which is exposed to as little bacterial contamination.

MAINTENANCE/CALIBRATION

The Reader Unit requires no maintenance.

CALIBRATION OF BARD SENSOR

The Reader Unit is equipped with a barometric pressure sensor (BARD sensor). To ensure that pre-defined tolerance limits are upheld, calibration is necessary on an annual basis.

For the purposes of calibration, the device should be sent each year to the following address:

 Service addresses:

 Christoph Miethke GmbH & Co. KG

 Ulanenweg 2

 14469 Potsdam / Germany

 Phone:
 +49 (0) 331 62083-0

 Fax:
 +49 (0) 331 62083-40

 E-Mail:
 info@miethke.com

Note: When undergoing calibration, the device will also be subjected to thorough functional and technical safety inspection.

TROUBLESHOOTING

If an error occurs, it will be displayed on the screen.

Example of an error message:

error message

Fig. 52

Additional information can be retrieved by pressing Infor Return to the previous menu by pressing Keen.

MALFUNCTIONS DISPLAYED WITH AN ERROR MESSAGE

Note: The table is sorted alphabetically by the initial letter of the error message appearing on the display

	Message displayed	Cause	Detection/remedy
A	Battery flat	Battery capacity exhausted (0%)	All data are saved after 2 min. The Reader Unit switches off automatically. Connect power supply.
	Battery voltage incompatible - use original power supply	Battery voltage of Reader Unit is too low	Reader Unit switches off automatically after 20 seconds.Connect original power supply.
	Low battery voltage	Battery voltage too low	The backlight switches off after 3 seconds. Connect power supply. Current measurements will not be inter- rupted.
В	Flashing, crossed-through an- tenna symbol Info : Antenna error	Antenna faulty	Deactivate and reactivate the device. If the error persists, please contact Technical Service.
	Flashing, crossed-through an- tenna symbol Info: Antenna not connected	Antenna not connected at start of measurement. -or- Antenna was disconnected during measurement.	Connect antenna. Measurement starts again. -or- Connect antenna. Measurement continues.
	Flashing, crossed-through an- tenna symbol Info: No communication	Data recording interrupted during continuous measure- ment (telemetric link dis- rupted)	Measurement resumes automatically once communi-cation is restored
	Flashing, crossed-through an- tenna symbol Info SD-card has been removed, measurement possible	SD-card was removed duri- ng continuous measurement.	Insert the SD-card. Start measurement again
D	Data error! Info: File opening not possible	Validation of the file failed	File cannot be opened
	Continues key activation key pad error	Permanent actuation of a key >60 se- conds	Release the key
	Pressure values out of range	Pressure data measured from implant not plausible - physiologically meanin- gless data	Measurement continues. If the error persists, please contact Technical Service

	Message displayed	Cause	Detection/remedy
E	Input voltage incompatible	Voltage of power supply is too high.	Reader Unit switches off automatically after 20 seconds. Use original power supply
F	Wrong SD-card inserted! SD-card! Info: Measurement w/o storage of data possible - or - Insert SD-card with correct ID XXXXXXXXXX	Measurement starts wit- hout SD-card. SD-card not designated for use with the implant is employed during measurement	Insert the SD-card designated for use with the shunt sensor into the Reader Unit. Use the identification number (ID) to assign the cor- rect SD-card to the implant (see patient pass).
	Wrong SD-card inserted! Info: Measurement w/o storage of data possible -or - Insert SD-card with correct ID XXXXXXXXX!	SD-card with ID which dif- fers to that of the implant	Insert the SD-card designated for use with the shunt sensor into the Reader Unit. Use the identification number (ID) to assign the cor- rect SD-card to the implant (see patient pass).
	Wrong implant -restart measurement! Info: Switch to another implant du- ring continuous measurement not possible!	Once continuous measure- ment has started, data are received from a different implant	Increase the distance between the two implants
G	Unit-temperature incorrect Info: Unit temperature from 10°C to 40°C valid!	Temperature in Reader Unit exceeds calibrated interval	Reader Unit can only be used at a device tem- perature of 10°C to 40°C. Current measurement is interrupted
I	Internal voltage incompatible	Internal voltage of device too high/too low	Reader Unit switches off automatically after 20 seconds. Contact Technical Service
S	SD-card faulty! Info: Measurement w/o storage of data possible!	SD-card not writable or readable (dirt, corrosion, faulty contacts)	Use of implant's internal calibration data. Data are not saved
	SD-card faulty!	SD-card not readable (dirt, corrosion, faulty contacts)	Check SD-card for damage or dirt
	SD-card inserted! Restart measurement! Info: Storage of reading after restart of measurement possible!	Measurement starts without SD-card. SD-card belonging to the implant is inserted during measurement	Start measurement again

Message displayed	Cause	Detection/remedy
SD-card inserted! Restart measurement! Info: Storage of reading after restart of measurement possible!	Measurement starts without SD-card. SD-card belonging to the implant is inserted during measurement	Start measurement again
SD-card missing! Info: Insert SD-card!	No SD-card inserted in the data manage-ment mode	Insert SD-card
SD-card missing! Info: Measurement w/o storage of data possible - or - Insert SD-card with correctID XXXXXXXXX!	SD-card not inserted	Insert the SD-card designated for use with the shunt sensor into the Reader Unit. Use the identification number (ID) to assign the correct SD-card to the implant (see patient pass).
SD-card not readable! : Measurement w/o storage of data possible!	Insertion of wrongly formatted or unfor- matted SD-card	Insert the SD-card designated for use with the shunt sensor into the Reader Unit. Use the identification number (ID) to assign the correct SD-card to the implant (see patient pass)
SD-card not readable! Info Measurement with storage of data possible!	No SD-card - or - SD-card ID not readable - or - SD-card contains no cali- bration data	Insert the SD-card designated for use with the shunt sensor into the Reader Unit. Measurement with data storage possible, but with limited accuracy. Measurement with data storage possible, but with limited accuracy
SD-card memory full Option to measure w/o saving!	Memory map on SD-card exhausted (100%) during continuous measurement	Delete any data that are no longer needed, measurement possible without storage
SD-card memory full	Memory map on SD-card exhausted (100%) during continuous measurement	Delete any data that are no longer needed
SD-card memory almost full	Memory map on SD-card almost exhausted (99%) during continuous measure- ment or fast measurement	Delete any data that are no longer needed
System error Info Ambient pressure not readable	Barometric (BARD) pressure not readable	Device function locked. Deactivate and reac- tivate the device. If the error persists, please contact Technical Service

	Message displayed	Cause	Detection/remedy
	System error -incompatibility	Device hardware and soft- ware not compatible	Device function locked.Deactivate and reac- tivate the device. If the error persists, please contact Technical Service
	System error -antenna incompatible	Hardware status of antenna and Reader Unit incompatible	Device function locked.Deactivate and reac- tivate the device. If the error persists, replace the antenna and/ or contact Technical Service
	System error -ID-data inadmissible	Identification data of implant corrupted	Measurement is stopped.Deactivate and re- activate the device. If the error persists, please contact Technical Service
	System error -implant voltage out of range	Implant voltage outside permissible range	Measurement is stopped.Deactivate and re- activate the device. If the error persists, please contact Technical Service
	System error -calibration data incorrect	Calibration data in implant corrupted or not readable (applies only if no SD-card has been inserted)	Measurement is stopped. Deactivate and re- activate the device. If the error persists, please contact Technical Service. Measurement possible with inserted SD-card
	System error -contact Service	The system test detected an error	Device function locked.Deactivate and reac- tivate the device. If the error persists, please contact Technical Service
T	Keypad faulty	A pressed key is registered when activating the Reader Unit	Release the key. Deactivate and reactivate the device.
	Temperature increase out of range	Spontaneous temperature increase in implant above 2 K or beyond 39.0°C	Measurement is stopped.Leave to rest for 10min. If the error persists, please contact Technical Service
	Temperature value out of range	Temperature data measured from implant not plausible - physiologically meanin- gless data	Measurement is stopped.Deactivate and re- activate the device. If the error persists, please contact Technical Service
	Timeout! Restart measurement!	The time between starting measurement and suc- cessfully establishing com- munication was exceeded (60sec)	Start measurement again. Optimize distance between antenna and implant
	Ambient pressure incorrect Info Ambient pressure from 800 to 1.100 mbar valid	The permissible barometric pressure during measure- ment is under - or- over-cut	Reader Unit can be used only at ambient pressure levels of 800–1.100 mbar. Current measurement is interrupted

OTHER MALFUNCTIONS

Malfunction	Cause	Remedy
Device cannot be activated	Battery is completely empty	Connect power supply.It takes about 6 hours to fully charge the battery. The Reader Unit can also be operated while it is being charged (connected to mains supply). Note: Charging is not possible at an ambient temperature of >35°C.
System turns off	Unfavorable operating con- ditions (eg, low humidity or unsuitable flooring)	Turn system on. If the error occurs again, please contact the Technical Service

TECHNICAL SERVICE

Warning: Risk of injury and/or malfunction! Do no modify the product.

For service and repairs, please contact Christoph Miethke GmbH & Co. KG Ulanenweg 2 14469 Potsdam/ Germany Phone: +49 (0) 331 62083-0 Fax: +49 (0) 331 62083-40 E-Mail: info@miethke.com

Modifications carried out on medical technical equipment may result in a loss of guarantee/ warranty rights and forfeiture of applicable licenses.

Christoph Miethke GmbH & Co. KG is liable for the safety, reliability and performance of the device only if: the device is used in accordance with the instructions for use resetting, modification or repair ensures only by persons authorized by Christoph Miethke GmbH & Co. KG the electrical installations in the appropriate facilities comply with national standards (IEC Standard).

Christoph Miethke GmbH & Co. KG will gladly provide technical documentation upon request. For service, maintenance and repairs, please contact Christoph Miethke GmbH & Co. KG. Service address: Christoph Miethke GmbH & Co. KG Ulanenweg 2 14469 Potsdam/ Germany Phone: +49 (0) 331 62083-0 Fax: +49 (0) 331 62083-40 E-Mail: info@miethke.com

TECHNICAL DATA

Product no.	Descriptiion
Voltage range Reader Unit Power supply	6 V (DC 100-240 V (50–60Hz)
Power consumption Reader Unit Power supply	1,4A (DC) 0,25A (50–60Hz)
Operating frequency	133kHz
Pressure measurement range	-66,67 hPa to 3.33 hPa ± 5.32 hPa
Temperature range for pressure measurement	20°C to 45°C
Working distance between antenna and shunt sensor	10 to 30mm
Battery lifetime	min. 5 years
Battery charging cycles	min. 250 charge cycles
Battery self-discharge	Residual charge after 3 months (storage) > 70%
Flammability class of housing	UL 94 HB
Damp protection or impermeability Reader Unit Antenna Power supply	IP44 IP44 IP40
Durability - shock testing	to IEC 60601-1 (3rd edition) 15.3.2
Durability - drop testing	to IEC 60601-1 (3rd edition) 15.3.4.1
Weight: Reader Unit Antenna Power supply	0.600 kg 0.215 kg 0.115 kg
Dimensions (W x H x D) Reader Unit Antenna (without cable) Power supply	144 x 270 x 65mm 100 x 250 x 0.98 55 x 65 x 2.76
Safety rating (Power supply)	11
Conforming to standard	IEC/DIN EN60601-1/ EN 45502-1
EMC	IEC/DIN EN 60601-1-2

AMBIENT CONDITIONS

Operating conditions	
Relative humidity	30% to 75%
Ambient temperature	10°C to 40°C, no con- densation
Atmospheric pressure	800 hPa to 1.100 hPa
Transport conditions	
Relative humidity	15% to 95%
Ambient temperature	0°C to 50°C
Atmospheric pressure	500 hPa to 1.100 hPa
Storage conditions	
Relative humidity	15% to 95%
Ambient temperature	10°C to 40°C
Atmospheric pressure	500 hPa to 1 100 hPa

CE-MARKING

Approval to affix the CE marking as defined by Directive 90/385/EEC for active implantable medical devices was granted in 2011.

STORAGE AND TRANSPORT

- Store reprocessed products protected from dust in a dry, dark and cool room which is exposed to as little bacterial contamination as possible.
- Ship the product only in its original packaging.
- For storage and transport conditions, see Chapter Technical Data.
- Caution: Damage to the product if used too soon after storage/shipment at temperatures below +10°C.
- · Leave the Reader Unit to acclimatize to room temperature for about three hours.

DISPOSAL

 Adhere to national regulations when disposing of or recycling the product, its components and packaging

Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is arranged free of charge by the manufacturer.

Please contact Christoph Miethke GmbH & Co. KG if you have any questions concerning the disposal of the product, see Chapter Technical Service.

For recycling of the product please send the *Reader Unit Set for SENSOR RESERVOI*R back to the manufacturer.

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