*M.scio*®

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**TABLE OF CONTENTS**

|   |    |
|---|----|
| 00.00 GENERAL NOTE  | 17 |
| 01.00 PRODUCT IDENTIFICATION AND GENERAL INFORMATION  | 17 |
| 02.00 INTENDED USE OF THE DEVICE  | 18 |
| 02.01 Intended purpose  | 18 |
| 02.02 Indication and intended patient groups  | 18 |
| 02.03 Contraindications   | 18 |
| 03.00 PRODUCT DESCRIPTION   | 18 |
| 03.01 Product description incl. materials and substances<br>that come into contact with the body or tissue                      | 18 |
| 03.02 Description of the intended mode of operation   | 19 |
| 03.03 References to previous generations or<br>variants and a description of the differences                                    | 19 |
| 03.04 Description of accessories, if present  | 20 |
| 03.05 Information about contained medicines and medicinal products  | 20 |
| 03.06 Description of the products with which the device is supposed<br>to be combined in order to achieve the intended function | 21 |
| 04.00 RISKS, WARNINGS AND PRECAUTIONS   | 22 |
| 04.01 Description of risk control measures  | 22 |
| 04.02 Residual risks  | 22 |
| 04.03 Undesirable side effects  | 22 |
| 04.04 Warnings and precautions  | 23 |
| 04.05 Summary of all field safety measures<br>(Field Safety Note / FSN, Field Safety Corrective Action / FSCA), if applicable   | 25 |
| 05.00 SUMMARY OF THE POST-MARKETING CLINICAL EVALUATION AND<br>CLINICAL FOLLOW-UP MONITORING                                    | 25 |
| 05.01 Summary of clinical data related to the equivalent product  | 25 |
| 05.02 Clinical evidence of CE marking   | 27 |
| 05.03 Summary of clinical data from other sources   | 27 |
| 05.04 Summary of clinical performance and safety  | 27 |
| 05.05 Planned continued Post Market Clinical Follow-up (PMCF)   | 28 |
| 06.00 POTENTIAL DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES  | 28 |
| 07.00 USER TRAINING   | 28 |
| 08.00 REFERENCE TO HARMONISED STANDARDS AND COMMON SPECIFICATIONS   | 29 |
| 09.00 CHANGE HISTORY  | 29 |
| 10.00 APPENDIX  | 29 |

## 00.00 GENERAL NOTE

This Summary of Safety and Clinical Performance Report is intended to provide public access to the updated summary of key aspects of the safety and clinical performance of the medical device. The summary report is not intended to replace the instructions for use as a document describing the safe use of the product, nor is it intended to provide diagnostic or therapeutic suggestions/advice to the user or patient. This SSCP is intended to provide information for medical specialists, e.g. neurosurgeons.

## 01.00 PRODUCT IDENTIFICATION AND GENERAL INFORMATION

|  |  |
|--|--|
| <b>Manufacturer</b>  | Christoph Miethke GmbH & Co. KG (CMG)<br>Ulanenweg 2, 14469 Potsdam, Germany             |
| <b>SRN</b>   | DE-MF-000010822  |
| <b>Name of the medical device</b>                                | <i>M.scio</i> (incl. associated SD-card)   |
| <b>Basic UDI-DI</b>  | 40419060000000000000002RY  |
| <b>Medical device nomenclature</b>                               | MDA 0101/A   |
| <b>Risk class</b>  | Rule 8.7; Class III  |
| <b>Year of first placing on the market of the medical device</b> | 07/04/2020 (Predecessor <i>SENSOR RESERVOIR</i> and <i>SENSOR PRECHAMBER</i> 19/12/2011) |
| <b>Name and identification number of the notified body</b>       | DQS, Reg.-No ZLG-BS-241.10.11,<br>identification number: 0297                            |

## 02.00 INTENDED USE OF THE DEVICE

### 02.01 INTENDED PURPOSE

The *M.scio* System is used for the diagnostic measuring of the intracranial pressure in the cerebrospinal fluid.

Thanks to the silicone membrane, the *M.scio* “dome” variants fulfil the characteristics of a conventional reservoir due to their puncturability and pumpability. In other words, they provide the facility for therapeutic pressure release through extraction of cerebrospinal fluid, as well as diagnostic extraction of cerebrospinal fluid, administration of fluids and verification of pressure values.

### 02.02 INDICATION AND INTENDED PATIENT GROUPS

The indications for the *M.scio* System are hydrocephalus and subarachnoid haemorrhage. Extended indications are shunt dependency, shunt dysfunctions and treatment optimisation. At implantation of the *M.scio*, patients must weigh more than 10 kg. Apart from that, there are no patient group limitations as to the use of the *M.scio* System.

### 02.03 CONTRAINDICATIONS

The use of *M.scio* is contraindicated in cases of coagulation disorders (risk of post-operative bleeding), blood in the cerebrospinal fluid, infections or a suspected infection in the regions of the body affected by the implantation (e.g. skin infection, meningitis, ventriculitis, bacteraemia, septicaemia, also peritonitis if the *M.scio* is used in the shunt).

Relative contraindications include high pressure and shock loads caused by patient actions (including diving, boxing, football), aggressive/autoaggressive patient behaviour that may limit patient compliance during aftercare and complicate the readout using the *Reader Unit Set*. Such behaviour may damage the *M.scio* and increase the risk of complications from wound conditions.

## 03.00 PRODUCT DESCRIPTION

### 03.01 PRODUCT DESCRIPTION INCL. MATERIALS AND SUBSTANCES THAT COME INTO CONTACT WITH THE BODY OR TISSUE

#### ***M.scio***

The housing design of the *M.scio*, dome-angled, is based on the reservoir shapes of (*Burrhole Reservoir*) already approved by Christoph Miethke GmbH & Co. KG. The *M.scio*, dome-inline is based on the already approved shapes for the *Prechambers*.

After implantation of the ‘dome’ variants, the following materials are in contact with body tissue and fluids:

- ▶ Measuring cell housing: titanium alloy
- ▶ Housing of the *M.scio*: plastic (PEEK), silicone
- ▶ Connectors: titanium alloy

The housing of the flat-angled and flat-inline variants of the *M.scio* was designed to match the dome-angled and dome-inline variants of the *M.scio*. However, the silicone membrane was removed with the aim of reducing implant height, while all other components are retained. The housing design thus includes the following components:

- ▶ Measuring cell housing: titanium alloy
- ▶ Housing of the *M.scio*: plastic (PEEK)
- ▶ Connectors: titanium alloy

The *M.scio* dome-angled, flat-angled, dome-inline and flat-inline variants are available with and without a preconnected distal silicone catheter.

The measuring cell is enclosed by a titanium alloy housing, hermetically encapsulating the electronics inside.

### Accessories Reader Unit Set and SD-card

The measured values of the *M.scio* can be read and displayed telemetrically and non-invasively using the *Reader Unit Set*.

The scope of delivery of the *M.scio* includes an SD-card. The measured data are automatically stored by the *Reader Unit Set* on the SD-card associated with the *M.scio*.

## 03.02 DESCRIPTION OF THE INTENDED MODE OF OPERATION

The medical device presented here is a system with both a diagnostic and a therapeutic function. The *M.scio* is implanted to measure the pressure as well as dynamic pressure changes within the cerebrospinal fluid. Intracranial pressure can be determined via connection to a *Ventricular Catheter*. In addition, the *M.scio* can be integrated into a shunt in order to determine the intracranial pressure within the shunt and to conduct a shunt function diagnosis, for example. The *M.scio* constitutes an addition to the shunt without impacting the drainage function.

In both cases, the *Reader Unit Set* is used for the telemetric reading and display of the pressure values.

## 03.03 REFERENCES TO PREVIOUS GENERATIONS OR VARIANTS AND A DESCRIPTION OF THE DIFFERENCES

Overview of previous generations of the product made by the manufacturer:

- ▶ Combined approval of OEM measuring cell and approval of PEEK *burrhole reservoir*
- ▶ *SENSOR RESERVOIR + SENSOR PRECHAMBER* (with and without distal catheter) with the purpose of pressure-based functional monitoring of shunt systems (August 2013; CE575455 and CE575454)
- ▶ Name change to *M.scio* and portfolio expansion to include flat variants

From a technological perspective, the *M.scio* system described here is based on the already AIMDD-approved system consisting of *M.scio*, SD-card for *Reader Unit* and *Reader Unit Set*. No technological changes have been made.

Overview of the product variants made by the manufacturer:

|               |                                  |               |
|---------------|----------------------------------|---------------|
| <i>M.scio</i> | dome-angled                      | 4041906505011 |
| <i>M.scio</i> | dome-angled with distal catheter | 4041906505028 |
| <i>M.scio</i> | dome-inline                      | 4041906505059 |
| <i>M.scio</i> | dome-inline with distal catheter | 4041906505066 |
| <i>M.scio</i> | flat-angled                      | 4041906505035 |
| <i>M.scio</i> | flat-angled with distal catheter | 4041906505042 |
| <i>M.scio</i> | flat-inline                      | 4041906505073 |
| <i>M.scio</i> | flat-inline with distal catheter | 4041906505080 |

### 03.04 DESCRIPTION OF ACCESSORIES, IF PRESENT

#### **Reader Unit Set (Basic UDI-DI: 404190600000000000001RW)**

The measured values of the *M.scio* can be read and displayed telemetrically and non-invasively using the *Reader Unit Set*. The *M.scio* does not have a battery. Its energy is supplied telemetrically, thus without a cord from outside of the body through the *Reader Unit Set*. For a pressure measurement, the antenna of the *Reader Unit Set* must be positioned about 10 to 30 mm away from the *M.scio*. The *Reader Unit Set* consists of the Reader Unit, antenna and plug-in power unit and be battery or mains operated.

#### **SD-card (Basic UDI-DI: 40419060000000000000002RY)**

The scope of delivery of the *M.scio* includes an SD-card, on which all individual information of the *M.scio* (ID and calibration data) is stored during production. To carry out measurements, the SD-card is inserted into the card slot of the *Reader Unit Set*. The measured data are automatically stored by the *Reader Unit Set* on the SD-card associated with the *M.scio*, allowing their evaluation at a later time.

At the start of a measurement, the *M.scioReader Unit Set* compares the IDs stored on the SD-card to ensure that the measured values are exclusively stored on the SD-card associated with the *M.scio*.

### 03.05 INFORMATION ABOUT CONTAINED MEDICINES AND MEDICINAL PRODUCTS

The system does not contain any medicinal products, radiopharmaceuticals, nanoparticles, biologically degradable substances, biological tissues of human or animal origin, or substances manufactured from biological tissue.

### 03.06 DESCRIPTION OF THE PRODUCTS WITH WHICH THE DEVICE IS SUPPOSED TO BE COMBINED IN ORDER TO ACHIEVE THE INTENDED FUNCTION

We recommend using products by Christoph Miethke GmbH & Co. KG in combination with the *M.scio*. The *M.scio* System can be safely combined with our range of implantable shunt components. The following products are particularly useful for intracranial pressure measurement in combination with the *M.scio*:

| Product name   | Order no. |
|--|-----------|
| <i>Ventricular catheter</i> (with mandrin, length 250 mm)  | FV077P    |
| <i>Ventricular catheter</i> (with mandrin, length 180 mm) with deflector (small, diameter 13 mm) | FV076P    |
| <i>Ventricular Catheter</i> (with mandrin, length 250 mm) with deflector (small, diameter 16 mm) | FV078P    |
| <i>Prechamber</i> (small, diameter 14 mm)  | FV035T    |
| <i>Prechamber</i> (large, diameter 20 mm)  | FV033T    |
| <i>Paediatric CONTROL RESERVOIR</i> (small, diameter 14 mm)                                      | FV066T    |
| <i>CONTROL RESERVOIR</i> (large, diameter 20 mm)   | FV047T    |
| <i>Burrhole reservoir</i> (small, diameter 14 mm)  | FV039T    |
| <i>Burrhole reservoir</i> (large, diameter 20 mm)  | FV028T    |
| <i>Titanium shutting plug</i>  | FV024T    |

To measure the pressure within a shunt system, the *M.scio* can be safely combined with further components, such as valves and *peritoneal catheters*.

If a (*Paediatric*) *CONTROL RESERVOIR* is used, care must be taken not to place this product between the ventricle and the *M.scio* when connecting to shunt components. Otherwise, the dynamics of the pressure signal may be distorted. For that reason, the combination of the *M.scio* with a *SPRUNG RESERVOIR* is ruled out. The *M.scio* must be placed between the ventricle and the valve in order to determine the intracranial pressure.

## 04.00 RISKS, WARNINGS AND PRECAUTIONS

### 04.01 DESCRIPTION OF RISK CONTROL MEASURES

In general, the risks identified through a risk analysis have been reduced to an acceptable level through risk control measures such as inherently safe designs, control measures (e.g. 100% final inspections or 100% in-process controls) and process specifications in production and/or downstream phases, product configurations based on norms and harmonised standards, user training as well as instructions in the product documentation.

### 04.02 RESIDUAL RISKS

Residual risks are defined as remaining physical injuries or damage to the health of patients, users or goods, as well as to the environment, which may occur with a certain probability of occurrence and a certain extent of damage despite risk minimisation measures.

The following residual risks are possible when using the *M.scio* system for a 5 year product life:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma (SDH)
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation
- ▶ Local shunt irritation/allergic reaction

### 04.03 UNDESIRABLE SIDE EFFECTS

The following undesirable side effects/complications may occur in connection with the *M.scio* system for a 5 year product life:

- ▶ Headache, dizzy spells, mental confusion, vomiting in cases of possible leakage from the *M.scio*/shunt and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid
- ▶ Wound healing disorders because of the height of the dome-angled *M.scio*

The mentioned residual risks and undesirable side effects are monitored via reactive and proactive post-marketing surveillance activities and evaluated qualitatively and, whenever possible, quantitatively.

The hazard assessment, risk evaluation, effectiveness verification of the risk control measures over the entire life cycle and consideration of the recognised state of the art confirm the suitability of the *M.scio* System for its intended purpose.

The *M.scio* is designed to ensure that for use under the intended conditions and for the intended purposes, all risks are acceptable when weighed against the benefits to the patient.

## 04.04 WARNINGS AND PRECAUTIONS

### M.scio incl. SD-card

#### WARNING

- ▶ Improper handling and use of this product can cause hazards and damages. For that reason, please read and closely follow the instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.
- ▶ If the sterile packaging is damaged, the products must not be used under any circumstances.
- ▶ The product must not be implanted after its expiry date.
- ▶ Products that have already been used must not be reused on the same or another patient.

#### CAUTION

- ▶ Because of the risk of injury from incorrect use, the product user must have participated in product training before using the product for the first time. Please contact Christoph Miethke GmbH & Co. KG for information about product training.
- ▶ Silicone is very electrostatic. Please ensure that the catheter does not come into contact with dry cloths, talcum powder or rough surfaces. Adherent particles can lead to tissue reactions.
- ▶ The functionality of the *M.scio* must be tested prior to implantation.
- ▶ Patency should be checked by aspiration of a sterile fluid from the distal end and air should be removed from the implant prior to implantation.
- ▶ If there are any doubts as to the precision or general appropriateness of the measured intracranial pressure reading, e.g. because of the absence of detectable pulse-synchronous rhythmic oscillations, the plug of the *Ventricular Catheter* may be occluded or included air cause dampening of the pressure dynamics. The electronic function to measure the ICP may also be impaired. When in doubt, other procedures should be employed for verification of the obtained pressure readings in these and other cases, such as transcutaneous reference pressure measurement by puncture, and the patient's overall clinical picture should be assessed.
- ▶ Therapeutic ultrasound procedures are associated with a risk of unintentional concentration of the ultrasound field and thus the risk of injury to the patient.
- ▶ Violent external shocks (accident, fall, etc.) may put the integrity of the *M.scio* and the shunt system at risk.
- ▶ Frequent pumping of the "dome" variants can result in excessive drainage, thus leading to pressure conditions outside of the normal physiological range. The patient should be informed about this risk.

#### NOTE

The following situation can cause the measuring cell of the *M.scio* to fail:

- ▶ Use of energy-emitting devices such as defibrillators
- ▶ Passing of an electrical current from an external source through the body (monopolar forceps)
- ▶ Use of radiation therapy and radionuclide patient imaging procedures

## Reader Unit Set

### WARNING

- ▶ Improper handling and use of this product can cause hazards and damages. For that reason, please read and closely follow the instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.
- ▶ The use of non-approved components poses a risk to users and patients and may damage the *Reader Unit Set*. The product must not be modified. Only the original plug-in power unit may be used.
- ▶ To avoid electric shocks and damage to the device due to penetrating liquids, the connector must be inserted back into the device after removing the SD-card.
- ▶ Electric shock and fire hazard: Disconnect from mains before cleaning. Ensure that no fluids can penetrate into the product; insert the connector into the SD-card slot of the Reader Unit. Do not use flammable or explosive detergents or disinfectants.
- ▶ Creation of electromagnetic interferences! Portable HF communication devices (radio) (including their accessories, such as antenna cables and external antennas) should not be used at a distance closer than 30 cm (or 12 inches) to the Reader Unit and antenna. A distance of at least 80 cm must be ensured in the presence of electromagnetic interference with spot frequencies of 385 MHz or 450 MHz. Non-compliance may lead to a reduction of the device's performance characteristics. Electromagnetic interference can lead to a device shut-down. In this case, the device must be restarted and the measurement repeated.

### CAUTION

- ▶ Because of the risk of injury from incorrect use, the product user must have participated in product training before using the product for the first time. Please contact Christoph Miethke GmbH & Co. KG for information about product training.
- ▶ Observe notes on electromagnetic compatibility (EMC).
- ▶ Observe the maintenance instructions.
- ▶ Prior to use, check the product for proper condition and functioning.
- ▶ Do not operate the product in the vicinity of flammable substances (e.g. anaesthetics).
- ▶ The product must be set up so that the mains plug can easily be disconnected from the power supply.
- ▶ Only operate the *Reader Unit Set* outside of the operating range of an MRI.
- ▶ Thoroughly clean the factory-new product after removing the transport packaging according to the manufacturer's instructions.
- ▶ To prevent nosocomial infection and multiple resistances, the device should be disinfected before and after each use.
- ▶ Before each use, check that the *Reader Unit Set* is in good condition and fully functional and check the system settings, such as the pressure unit.
- ▶ Failure to perform the annual calibration can result in the drifting of the barometric pressure sensor outside of the tolerance limits.

**NOTE**

- ▶ **Allow the *Reader Unit Set* to acclimatise at room temperature for approx. 3 hours.**
- ▶ **To prevent possible damage during transport, the *Reader Unit Set* must be shipped in its original packaging.**
- ▶ **To prevent damage due to improper use, do not touch the contact points of the SD-card.**
- ▶ **Mechanical cleaning or disinfection as well as inappropriate detergents or disinfectants can damage or destroy the product: Only clean/disinfect the product manually; never sterilise the product; only use cleaning/disinfecting agents approved for surface cleaning according to the manufacturer's instructions; follow the instructions for concentration, temperature and contact time.**

#### **04.05 SUMMARY OF ALL FIELD SAFETY MEASURES (FIELD SAFETY NOTE / FSN, FIELD SAFETY CORRECTIVE ACTION / FSCA), IF APPLICABLE**

Until recently, a Field Safety Notice for the *M.scio* incl. accessories has been published on the website of the Federal Institute for Drugs and Medical Devices: An internal audit revealed an error in the instructions for use (IFU) on the part of the supplier. The audit found a few pages missing in some IFUs. However, the majority of the audited IFUs were free from errors. Nevertheless, it could not be ruled out that instructions for use with missing pages were enclosed with already delivered products. Moreover, the audit showed that the symbols/pictograms for the notices and warnings were not present on all pages of the printed copies. However, the notices and warnings were included in the IFUs as intended. In this case, the printing company was provided with an inadequate print template by Christoph Miethke GmbH & Co. KG.

There was no risk to patients or users at any time, as the products themselves were not defective.

#### **05.00 SUMMARY OF THE POST-MARKETING CLINICAL EVALUATION AND CLINICAL FOLLOW-UP MONITORING**

##### **05.01 SUMMARY OF CLINICAL DATA RELATED TO THE EQUIVALENT PRODUCT**

All data used for this clinical evaluation come from the predecessor product of the *M.scio* System, the *SENSOR RESERVOIR*, which from 2020 has also been approved under the trade name *M.scio*. They originate mainly from the scientific literature and from scientific conference contributions. Service reports, our in-house complaints database and the relevant public databases were also evaluated. Additional data come from our in-house feedback database, which is maintained by the product management department and which registers every piece of user feedback from conversations, phone calls, emails, etc.

## Literature data overview

- ▶ Clinical benefits
  - ▶ Optimisation of diagnosis and therapy thanks to telemetric measurement of intracranial pressure values
    - ▶ Use of a long-term implant
    - ▶ Simple and easy reading of pressure values
    - ▶ Identification of pathological pressure situations
    - ▶ Low risk thanks to non-invasive measuring method
    - ▶ Conditionally MRI-compatible for magnetic fields up to 3 tesla
    - ▶ Increased sense of security in concerned patients and their families thanks to easily accessible readings
    - ▶ Different variants for individual patient-specific treatment requirements
    - ▶ Optional expandability of the *M.scio* System to a shunt system
  - ▶ Optimised patient management in shunt patients
    - ▶ Improved patient outcomes
      - ▶ Optimised valve settings on the basis of the determined pressure readings
      - ▶ Reduction of over- and underdrainage
    - ▶ Reduced strain on patients
      - ▶ Avoidance of unnecessary clinical diagnostic procedures and their associated risks (such as exposure to radiation for imaging procedures and the use of invasive diagnostic techniques)
      - ▶ Avoidance of unnecessary revisions thanks to functional monitoring of the shunt and ruling out of occlusions and shunt malfunctions
    - ▶ Cost savings
      - ▶ Avoidance of unnecessary clinical procedures (such as medical imaging, pressure measurement and revisions)
  - ▶ Optimised diagnosis and treatment options by using the *M.scio* “dome” variants
    - ▶ Increased options thanks to puncture
      - ▶ CSF removal for manual pressure release and laboratory analyses
      - ▶ Option of external reference pressure measurement
      - ▶ Administration of fluids
    - ▶ Reduced strain on patients
      - ▶ Pump tests for functional monitoring of shunts
    - ▶ Cost savings
      - ▶ Avoidance of unnecessary clinical procedures (such as medical imaging and revisions)

► Biological compatibility

Biological complications associated with the *SENSOR RESERVOIR*, such as allergic, inflammatory and rejection reactions, encapsulation, frustrated phagocytosis and similar, are not mentioned or addressed anywhere in either market feedback or literature.

► Risk/safety

The only product-specific risk reflected in the available clinical data arises from the height of the conventional *SENSOR RESERVOIR*, as in rare cases (2/158 = 1%) this can lead to wound healing disorders due to mechanical tension of the skin.

Regarding the use of the *SENSOR RESERVOIR / M.scio* as an ICP sensor, no particular risk emerged from the literature, since its clinically required accuracy is ensured and no measures in the field of shunt therapy associated with a particular diagnostic or therapeutic decision are ever carried out based on a single criterion.

Even a complete failure of the clinical function of "ICP measurement within the CSF" does not lead to any new identifiable risks, as it does not affect the function of the shunt.

Accordingly, the *M.scio* does not have to be explanted in such a case.

► Usability

There is no doubt in the literature about the overall good usability of the *SENSOR RESERVOIR* and the *Reader Unit Set* as far as the clinical benefit of the short-term ICP measurement and the optimisation of the shunt setting described above are concerned.

► PMS/PMCF activities

From the review of the reactive and proactive market feedback, it can be summarised that for the *SENSOR RESERVOIR*, *SENSOR PRECHAMBER* and the new *M.scio* products (both "dome" and "flat" variants), there was only one new product risk during the assessment period and no corrective action had to be taken. The product risk of "wound healing complications" due to height has already been addressed and new flatter *M.scio* variants have been made available on the market.

## 05.02 CLINICAL EVIDENCE OF CE MARKING

No clinical trials were conducted prior to market approval.

## 05.03 SUMMARY OF CLINICAL DATA FROM OTHER SOURCES

Clinical data of a similar medical device from the scientific literature and the MAUDE database query were used for comparison. The clinical data do not show any unknown risks.

## 05.04 SUMMARY OF CLINICAL PERFORMANCE AND SAFETY

The performance characteristics specified for the *M.scio* System have been shown to meet the performance requirements based on the clinical background of the state of the art. The available clinical data demonstrate that the *SENSOR RESERVOIR/M.scio* has become established as a second telemetric system alongside the similar older telemetric system. For the diagnostic part of the intended purpose ("ICP measurement within the CSF"), neither known nor previously unknown risks are evident. The clinical data show no evidence of impairment of the relevant performance characteristics such as insufficient accuracy, clinically relevant zero drift, failure of the "resolvability of ICP pulse waves", telemetry, long-term robustness of the measuring cell or its CSF conductivity. Hypothetical but potential cases are addressed in the general warnings. A potential total failure of the measurement function, to date only documented for a single complaint but not at all in the literature, does not lead to any new risk.

For the therapeutic part of the intended use ("puncturability and pumpability" for the "dome" variants), there are also neither known nor previously unknown risks. The warning addresses the currently hypothetical but potential risk of excessive drainage due to the (too) frequent use of the pumping function.

The clinical data include no systematic significant evidence to indicate impairment of the "long-term implant" performance characteristic or any risks and hazards associated with such impairment. There are no documented known or as yet unknown risks and hazards resulting from interactions of the *M.scio* system with other diagnostic and therapeutic procedures, such as MRI or ultrasound. Hypothetical but potential risks that could result from such interactions are addressed in the warnings.

Active implants are usually sterilised with ethylene oxide (ETO) to protect the electronics. No cases of allergies are known so far for the *SENSOR RESERVOIR*. Based on the requirements of the harmonised standard DIN EN ISO 10993-17:2009-08, a limit is set for the body weight of a patient to limit the risk of such an ETO allergy. We have made generous allowances for this by specifying a minimum patient weight of 10 kg.

#### **05.05 PLANNED CONTINUED POST MARKET CLINICAL FOLLOW-UP (PMCF)**

The following PMCF activities are planned for post-market clinical follow-up:

- ▶ Literature research in selected literature databases
- ▶ Database research on potential incidents in Germany, the EU and the USA
- ▶ Monitoring of similar or identical medical devices for possible unknown risks

#### **06.00 POTENTIAL DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES**

In the field of bedside (non-mobile) ICP measurement (e.g. in trauma patients), conventional partially implantable catheters have been used up to now. Measuring probes are placed inside the skull and connected by cable to an external device to display measured values and determine pressure values. However, these ICP catheters are unsuitable for long-term applications because of the risk of infection.

In the field of mobile (non-bedside) ICP measurement, a pressure probe exists with the pressure data read out telemetrically. The implantation period of three months must not be exceeded. If hydrocephalus is diagnosed based on the ICP measurement, the pressure probe must be explanted in another surgical procedure, or a further burr hole must be made in the skull of the patient for the insertion of a shunt.

The implantable products on the market only offer a purely diagnostic function. A direct extension of the implant with a therapeutic option, such as the possibility of administering fluid or removing CSF via a shunt, is not possible.

#### **07.00 USER TRAINING**

A concept for suitable product training will be developed in preparation for the market launch of the product. The training is intended for users and multipliers (i.e. employees of the distribution and trade partners). Implementation of the training is documented.

## 08.00 REFERENCE TO HARMONISED STANDARDS AND COMMON SPECIFICATIONS

The complete list of applied standards for the system component *M.scio* and *Reader Unit Set* can be found in the Appendix under 10.00.

## 09.00 CHANGE HISTORY

| SSCP revision number | Issue date | Description of change   | Revision validated by the notified body   |
|----------------------|------------|---|---|
| 01                   | 05.02.2021 | Initial preparation as part of the conformity assessment according to Annex X for initial approval under (EU) Regulation 2017/745 | <ul style="list-style-type: none"> <li>▶ Yes</li> <li>▶ Validated language: German</li> </ul> |
| 02                   | 04.03.2022 | New layout and content-related adaptation   | Not yet validated by the notified body  |

## 10.00 APPENDIX

### Appendix 1 - Applied harmonised standards

| Number | Complete identifier         |
|--------|-----------------------------|
| N22    | EN ISO 11135:2014 + A1:2019 |
| N33    | EN ISO 11737-2:2020         |

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany