M.scio®

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TABLE OF CONTENTS

00.00 GENERAL NOTE	11
01.00 PRODUCT IDENTIFICATION AND GENERAL INFORMATION	11
02.00 INTENDED USE OF THE DEVICE	12
02.01 Medical basics	12
02.02 Intended purpose	12
02.03 Indication and intended patient group	12
02.04 Contraindications	12
03.00 PRODUCT DESCRIPTION	13
03.01 Product description incl. materials and substances that come into contact with the body or tissue	13
03.02 Information about contained medicines and drugs	13
03.03 Description of the intended operating principle	13
03.04 Description of accessories, if present	13
03.05 Description of the products with which the device is supposed to be combined in order to achieve the intended function	13
04.00 RISKS, WARNINGS AND PRECAUTIONS	14
04.01 Description of risk control measures	14
04.02 Adverse side effects/complications	14
04.03 Residual risks	14
04.04 Warnings and precautions	15
04.05 Summary of all field safety measures (FSN/FSCA), if applicable	15
05.00 SUMMARY OF THE POST-MARKETING CLINICAL EVALUATION AND CLINICAL FOLLOW-UP MONITORING	16
05.01 Clinical background of the product	16
05.02 Clinical evidence of CE marking	16
05.03 Description of the PMCF activities to ensure continued safety and performance of the product	16
06.00 POTENTIAL DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES	16
07.00 USER TRAINING	16
08.00 CHANGE HISTORY	17

00.00 GENERAL NOTE

This summary report is intended to provide public access to the updated summary of key aspects of the safety and clinical performance of the medical device. The SSCP is not intended to provide general advice for the treatment of hydrocephalus. Please contact your doctor if you have any questions about your medical condition or use of the medical device.

This summary report is no substitute for an implant card or the instructions for use for information on the safe use of the medical device. This report is intended for the information of patients and laypersons.

01.00 PRODUCT IDENTIFICATION AND GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG (CMG) Ulanenweg 2, 14469 Potsdam, Germany
SRN	DE-MF-000010822
Name of the medical device	<i>M.scio</i> (incl. associated SD-card)
Basic UDI-DI	404190600000000000000002RY
Year of first placing on the market of the medical device	07.04.2020 (Predecessor <i>SENSOR RESERVOIR</i> and <i>SENSOR PRECHAMBER</i> 19/12/2011)

02.00 INTENDED USE OF THE DEVICE

02.01 MEDICAL BASICS

The human brain is surrounded by a special fluid (scientific term: cerebrospinal fluid – CSF). The cerebrospinal fluid protects the brain from mechanical damage, regulates the pressure inside the skull (the intracranial pressure), keeps the brain tissue moist and transports metabolic products. Every day, new cerebrospinal fluid is produced in the cerebral ventricles and reabsorbed by the venous blood system. A healthy balance between the production and absorption of CSF is thus very important.

The body has mechanisms for keeping the intracranial pressure comparatively constant. However, various diseases and injuries, such as haemorrhages and hydrocephalus, can lead to increased intracranial pressure and life-threatening situations. The symptoms depend on the degree of impairment and include nausea, headache, vomiting, coordination disorders, drowsiness, and even unconsciousness.

In certain situations, it is beneficial to measure the intracranial pressure in order to gather important diagnostic information and ultimately to tailor the required treatment effectively to the patient's needs. Measuring intracranial pressure provides an important basis for selecting the appropriate treatment option.

In the case of hydrocephalus, normalisation of the intracranial pressure can be achieved through controlled drainage of the cerebrospinal fluid – usually into the abdominal cavity. To accomplish this task, a valve and a tube system (in conjunction known as a shunt system) are implanted.

02.02 INTENDED PURPOSE

The implantable *M.scio* makes it possible for your doctor to measure the intracranial pressure within the cerebrospinal fluid for diagnostic purposes, to perform a functional check of the shunt system and to evaluate the result of therapeutic measures. The wireless readout of the pressure values of the *M.scio* is carried out externally through the *Reader Unit Set*.

The “dome” variant of the *M.scio* system additionally offers the function of puncturability and pumpability as in a conventional reservoir thanks to a silicone membrane. In other words, it enables therapeutic pressure reduction through the extraction of cerebrospinal fluid as well as the diagnostic extraction of cerebrospinal fluid, administration of fluids and also the checking of pressure values.

02.03 INDICATION AND INTENDED PATIENT GROUP

The *M.scio* is used in patients with hydrocephalus and subarachnoid haemorrhages whose weight is greater than 10 kg. Other reasons for using the *M.scio* are shunt dependency, shunt dysfunction and treatment optimisation.

02.04 CONTRAINDICATIONS

The use of *M.scio* is not indicated in cases of coagulation disorders (risk of post-operative bleeding), blood in the cerebrospinal fluid, infections or a suspected infection in the regions of the body affected by the implantation (e.g. skin infection, meningitis, ventriculitis, bacteraemia, septicæmia, also peritonitis if the *M.scio* is used in the shunt).

The *M.scio* must not be implanted in cases of expected high pressure and shock loads (including diving, boxing, football). For aggressive/auto-aggressive patient behaviour with the potential of limiting the patient's willingness and cooperation in follow-up care and complicate the readout process with the *Reader Unit Set*, the use of the *M.scio* is only indicated after stringent consideration in individual cases. Such behaviour may damage the *M.scio* and increase the risk of complications from wound conditions.

03.00 PRODUCT DESCRIPTION

03.01 PRODUCT DESCRIPTION INCL. MATERIALS AND SUBSTANCES THAT COME INTO CONTACT WITH THE BODY OR TISSUE

Materials of the *M.scio*

The *M.scio* consists of high-quality materials that have been tested and standardised for use as implant materials. The following materials of the *M.scio* may come into contact with the body and tissues: Titanium and PEEK; in addition silicone for the "dome" versions. The catheters are also made from silicone. The silicone is latex-free.

03.02 INFORMATION ABOUT CONTAINED MEDICINES AND DRUGS

The system does not contain any medicinal products, radiopharmaceuticals, nanoparticles, biologically degradable substances, biological tissues of human or animal origin, or substances manufactured from biological tissue.

03.03 DESCRIPTION OF THE INTENDED OPERATING PRINCIPLE

Thanks to the wireless communication of the *M.scio* with the antenna of the *Reader Unit Set*, the intracranial pressure can be measured completely painlessly by the attending physician. The *M.scio* contains a measuring cell for this purpose.

Product life of the *M.scio*

The *M.scio* has been designed to work accurately and reliably over long periods of time. The maximum product life of the *M.scio* is five years. In the event of a measuring cell failure, the *M.scio* "dome" operates like a conventional burrhole reservoir or prechamber without restrictions, while the "flat" variants act like a deflector or connector. The integrated measuring cell poses no risk when the *M.scio* remains in the body.

03.04 DESCRIPTION OF ACCESSORIES, IF PRESENT

SD-card

For each *M.scio*, the patient receives an SD-card, which is used to store all individual data about the implant and save the measured values. If the SD-card is lost, it can be reordered by the attending physician, stating the serial number of the *M.scio* or the associated identification number (ID).

03.05 DESCRIPTION OF THE PRODUCTS WITH WHICH THE DEVICE IS SUPPOSED TO BE COMBINED IN ORDER TO ACHIEVE THE INTENDED FUNCTION

The *M.scio* can be safely combined with our implantable shunt components during implantation by the physician. For example, this includes *Ventricular Catheters*, *Peritoneal Catheters*, *Prechambers*, reservoirs, valves and *Titanium Shutting Plugs*. We recommend using the Christoph Miethke GmbH & Co. KG products in combination with the *M.scio*.

04.00 RISKS, WARNINGS AND PRECAUTIONS

Consult your doctor if you think you may experience side effects related to the product or its use, or if you are concerned about any potential risks. This document is no substitute for advice from your treating physician.

04.01 DESCRIPTION OF RISK CONTROL MEASURES

The risks identified by means of a risk analysis have been reduced to an acceptable level through various risk control measures. For example, they include safe design, 100% final inspections, process specifications in production, configuration of products based on norms and harmonised standards, user training and instructions in the product documentation.

04.02 ADVERSE SIDE EFFECTS/COMPLICATIONS

Measuring intracranial pressure with the *M.scio* and treating hydrocephalus with a shunt system is not always without complications. As with every surgical intervention, there is a risk of infection. Unfortunately, problems can occasionally develop that could be directly or indirectly linked to the implanted shunt system. Such complications include:

- ▶ Headache, dizzy spells, mental confusion, vomiting in cases of potential leakage from the *M.s-cio*/shunt and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Blockages due to protein and/or blood in the cerebrospinal fluid
- ▶ Wound healing disorders because of the height of the dome-angled *M.scio*

A physician must be consulted immediately in case of skin rashes and tightness, severe headaches, dizzy spells or similar.

04.03 RESIDUAL RISKS

The following residual risks are possible with the use of the *M.scio*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis)/allergic shock
- ▶ Acute hygroma / subdural haematoma (SDH)
- ▶ Cerebrospinal fluid accumulations (cerebrospinal fluid under the skin)
- ▶ Tissue damage/puncture
- ▶ Chronic hygroma
- ▶ Skin irritation
- ▶ Local shunt irritation/allergic reaction

The mentioned residual risks and undesirable side effects are monitored by us via numerous activities and evaluated qualitatively and, whenever possible, quantitatively.

The suitability of the *M.scio* for the described purpose is confirmed by the hazard assessment, risk evaluation, effectiveness verification of the risk control measures over the entire life cycle and consideration of the recognised state of the art. The *M.scio* was designed to ensure that, when it is used under the intended conditions and for the intended purposes, all risks are acceptable when weighed against the benefits to the patient.

04.04 WARNINGS AND PRECAUTIONS

M.scio incl. SD-card

WARNING

- ▶ Improper handling and use of this product can cause hazards and damages.

CAUTION

- ▶ Therapeutic ultrasound procedures are associated with a risk of unintentional concentration of the ultrasound field and thus the risk of injury to the patient.
- ▶ Violent external shocks (accident, fall, etc.) may put the integrity of the *M.scio* and the shunt system at risk.
- ▶ Frequent pumping of the *M.scio* "dome" variants can result in excessive drainage, thus leading to pressure conditions outside of the normal physiological range. The patient should be informed about this risk.

NOTE

The following situations can cause the measuring cell of the *M.scio* to fail:

- ▶ Use of energy-emitting devices such as defibrillators
- ▶ Passing of an electrical current from an external source through the body (monopolar forceps)
- ▶ Use of radiation therapy and radionuclide patient imaging procedures

04.05 SUMMARY OF ALL FIELD SAFETY MEASURES (FSN/FSCA), IF APPLICABLE

Until recently, a Field Safety Notice for the *M.scio* incl. accessories has been published on the website of the Federal Institute for Drugs and Medical Devices: An internal audit revealed an error in the instructions for use (IFU) on the part of the supplier. The audit found a few pages missing in some IFUs. However, the majority of the audited IFUs were free from errors. Nevertheless, it could not be ruled out that instructions for use with missing pages were enclosed with already delivered products. Moreover, the audit showed that the symbols/pictograms for the notices and warnings were not present on all pages of the printed copies. However, the notices and warnings were included in the IFUs as intended. In this case, the printing company was provided with an inadequate print template by Christoph Miethke GmbH & Co. KG.

There was no risk to patients or users at any time, as the products themselves were not defective.

05.00 SUMMARY OF THE POST-MARKETING CLINICAL EVALUATION AND CLINICAL FOLLOW-UP MONITORING

05.01 CLINICAL BACKGROUND OF THE PRODUCT

The *M.scio* system is primarily used for diagnostic intracranial pressure measurement.

05.02 CLINICAL EVIDENCE OF CE MARKING

Clinical evidence of CE marking was obtained during the collection of clinical and non-clinical data of the equivalent product *SENSOR RESERVOIR* as the predecessor product.

05.03 DESCRIPTION OF THE PMCF ACTIVITIES TO ENSURE CONTINUED SAFETY AND PERFORMANCE OF THE PRODUCT

The following PMCF activities are planned for post-market clinical follow-up:

- ▶ Literature research in selected literature databases
- ▶ Database research on potential incidents in Germany, the EU and the USA
- ▶ Monitoring of similar or identical medical devices for possible unknown risks

06.00 POTENTIAL DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES

If you are considering alternative treatment methods, please contact your treating physician who can take your individual situation into account.

In the field of bedside (non-mobile) intracranial pressure measurement, conventional partially implantable systems have been used up to now. Measuring probes are placed inside the skull and connected by cable to an external device to display measured values and determine pressure values. However, these measuring probes are unsuitable for long-term applications because of the risk of infection.

In the field of mobile (non-bedside) intracranial pressure measurement, there is a sensor with the pressure data read out wirelessly. The implantation period of three months must not be exceeded. If hydrocephalus is diagnosed based on intracranial pressure measurement, the sensor must be explanted in another surgical procedure, or another burrhole must be made in the skull of the patient for insertion of a shunt.

The implantable products on the market only offer a purely diagnostic function. A direct extension of the implant with a therapeutic option, such as the possibility of administering fluid or removing CSF via a shunt, is not an option with the alternative treatment methods.

07.00 USER TRAINING

A concept for suitable product training will be developed in preparation for the market launch of the product. Training is provided to users and staff of distribution and retail partners. Implementation of the training is documented.

08.00 CHANGE HISTORY

SSCP revision number	Issue date	Description of change	Revision validated by the notified body
01	05.02.2021	Initial preparation as part of the conformity assessment according to Annex X for initial approval under (EU) Regulation 2017/745	► Yes ► Validated language: German
02	07.09.2021	Basic UDI-DI SRN Risk class Rule 8.7; Class III	► Yes ► Validated language: German
03	04.03.2022	New layout and content-related adaptation	Not yet validated by the notified body

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany