Reader Unit Set

FOR THE USE WITH *M.scio*® IMPLANTS

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 www.miethke.com

(US) This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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info@miethke.com

TABLE OF CONTENTS

0.00 PREFACE AND IMPORTANT INFORMATION	36
1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	36
1.01 Explanation of the warnings	36
1.02 Display conventions	36
1.03 Other accompanying documentation and supplementary Information material	37
1.04 Feedback on the instructions for use	37
1.05 Copyright, disclaimer, warranty and other information	37
2.00 DESCRIPTION OF THE <i>M.scio</i> SYSTEM	37
2.01 Intended use	37
2.02 Clinical benefits	37
2.03 Indications	38
2.04 Contraindications	38
2.05 Intended patient groups	38
2.06 Intended users	38
2.07 Intended use environment	39
2.08 Functional principle	39
2.09 System components	39
2.10 Summary of safety and clinical performance	40
2.11 Supplementary product information	40
3.00 SYSTEM COMPONENT <i>Reader Unit Set</i>	41
3.01 Product description	41
3.02 Important safety information	42
3.03 Transport and storage	43
3.04 Use of the product	44
3.05 Manual cleaning and disinfection of the <i>Reader Unit Set</i>	56
3.06 Technical support	57
3.07 Disposal	58
3.08 Troubleshooting	58
3.09 Technical data and performance data	62
3.10 Electromagnetic compatibility	63
3.11 Symbols used on the product and labelling	64
4.00 MEDICAL DEVICE CONSULTANTS	64

0.00 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing the *Reader Unit Set*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product.
Your team at Christoph Miethke GmbH & Co. KG.

Relevance of the instructions for use

**WARNING**

Improper handling and use of this product can cause risks and damage. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

The *Reader Unit Set* forms a part of the *M.scio* system, which consists of the following components:

- *M.scio* with associated SD card
- *Reader Unit Set*

The *M.scio* system can be safely combined with our range of shunt components. These instructions for use apply to the *Reader Unit Set*

- FV907X
- from software version 2.04

including the approved components of the Reader Unit, antenna and plug-in power unit.

Operation of the *M.scio* and shunt components is described in the respective instructions for use.

Basic UDI-DI

The Basic UDI-DI of the *Reader Unit Set* and associated approved components is: 404190600000000000000001RW.

1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

1.01 EXPLANATION OF THE WARNINGS

**DANGER**

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

**WARNING**

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

**CAUTION**

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

**NOTE**

Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for risk, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

1.02 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates product names
[...]	Square brackets indicate selectable menu items or information that is shown on the <i>Reader Unit Set</i> display.
<...>	Pointed brackets indicate context-dependent symbols on the <i>Reader Unit Set</i> display

1.03 OTHER ACCOMPANYING DOCUMENTATION AND SUPPLEMENTARY INFORMATION MATERIAL

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

1.04 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests or remarks about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use.

1.05 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

2.00 DESCRIPTION OF THE *M.scio* SYSTEM

2.01 INTENDED USE

The *M.scio* system is used for diagnostic intracranial pressure measurement within the cerebrospinal fluid.

Due to the silicone membrane, the dome variants of the system component *M.scio* have the pumpability and puncturability of a conventional reservoir, i. e. they offer the possibility of therapeutic pressure reduction by removing cerebrospinal fluid, diagnostic removal of cerebrospinal fluid, administration of fluids as well as verification of the pressure values.

2.02 CLINICAL BENEFITS

Optimisation of diagnosis and therapy thanks to telemetric measurement of intracranial pressure values

- ▶ Use of a long-term implant
- ▶ Simple and easy reading of pressure values
- ▶ Identification of pathological pressure situations
- ▶ Low risk thanks to non-invasive measuring method
- ▶ MR Conditional for magnetic fields up to 3 tesla
- ▶ Increased sense of security in concerned patients and their relatives thanks to easily accessible readings
- ▶ Different variants for individual patient-specific treatment requirements
- ▶ Optional expandability of the *M.scio* system to a shunt system

Optimised patient management in shunt patients

Improved patient outcomes

- ▶ Optimised valve settings on the basis of the determined pressure readings
- ▶ Reduction of over- and underdrainage

Reduced strain on patients

- ▶ Avoidance of unnecessary clinical diagnostic procedures and their associated risks (such as exposure to radiation for imaging procedures and the use of invasive diagnostic techniques)
- ▶ Avoidance of unnecessary revisions thanks to monitoring of the shunt function and ruling out of occlusions and shunt malfunctions

Cost savings

- ▶ Avoidance of unnecessary clinical procedures (such as medical imaging, pressure measurement and revisions)

Optimised diagnosis and treatment options by using the *M.scio* “dome” variants

Increased options thanks to puncture

- ▶ CSF removal for manual pressure release and laboratory analyses
- ▶ Possibility of external reference pressure measurement
- ▶ Administration of fluids

Reduced strain on patients

- ▶ Pump tests for functional monitoring of shunts

Cost savings

- ▶ Avoidance of unnecessary clinical procedures (such as medical imaging and revisions)

2.03 INDICATIONS

The following indications apply to the *M.scio* system:

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Extended indications

- ▶ Shunt dependency
- ▶ Shunt dysfunction
- ▶ Treatment optimisation

2.04 CONTRAINDICATIONS

The following contraindications apply to the *M.scio* system:

Contraindications

- ▶ Coagulation disorders (risk of postoperative bleeding)
- ▶ Blood in the cerebrospinal fluid
- ▶ Infections or a suspected infection in the regions of the body affected by the implantation (e. g. skin infection, meningitis, ventriculitis, bacteraemia, septicaemia, also peritonitis if the *M.scio* is used in the shunt)

Relative contraindications

- ▶ High pressure and shock load caused by patient actions (including diving, boxing, football)
- ▶ Aggressive/autoaggressive patient behaviour may limit patient compliance during aftercare and complicate the readout using the *Reader Unit Set*. The *M.scio* may sustain damage from such behaviour, increasing the risk of complications from the wound conditions.

2.05 INTENDED PATIENT GROUPS

When implanting the *M.scio*, patients must weigh more than 10 kg. Apart from that, there are no patient group limitations as to the use of the *M.scio* System.

2.06 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical specialists, e. g. neurosurgeons
- ▶ Knowledge about the function and proper use of the product
- ▶ Successful participation in product training

2.07 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions
- ▶ Reading and evaluation of intracranial pressure values
- ▶ Use of the *M.scio* “dome” variant’s pump and puncture function

2.08 FUNCTIONAL PRINCIPLE

The *M.scio* is implanted to measure the pressure as well as dynamic pressure changes within the cerebrospinal fluid. The intracranial pressure can be determined via connection to a *Ventricular Catheter*.

In addition, the *M.scio* can be integrated into a shunt in order to determine the intracranial pressure in the shunt and to conduct a shunt function diagnosis. The *M.scio* is an addition to the shunt without impacting the drainage function.

The pressure is measured using a measuring cell located inside the *M.scio*. The measured values can be read and displayed telemetrically and non-invasively using the *Reader Unit Set*. The *M.scio* does not have a battery for this purpose. Its energy is supplied telemetrically and thus cordlessly from outside of the body through the *Reader Unit Set*. For a pressure measurement, the antenna of the *Reader Unit Set* must be positioned about 10 to 30 mm away from the *M.scio* (Fig. 1).

Fig. 1: The ideal working distance for telemetric data transmission, i.e. the distance between the Reader Unit Set antenna and the *M.scio*, is 10-30 mm.

The *Reader Unit Set* automatically stores the measured data on the SD card associated with the *M.scio*, allowing them to be evaluated at a later time.

2.09 SYSTEM COMPONENTS

Reader Unit Set

The measured data of the *M.scio* may only be read out using the *Reader Unit Set* (FV905X / FV907X). The operation of the *M.scio* is described in the associated instructions for use.

M.scio and SD card

The scope of delivery of the *M.scio* includes an SD card, on which all individual information of the *M.scio* (ID and calibration data) is stored during production. To carry out measurements, the SD card is inserted into the card slot of the *Reader Unit Set*. When a measurement is started, the *Reader Unit Set* compares the IDs stored on the SD card and the *M.scio* to ensure that measured values are exclusively stored on the SD card associated with the *M.scio*.

If the SD card is lost, it can be reordered by specifying the *M.scio* serial number or the associated ID. The ID can be read with [Single measurement] from the *M.scio* and shown on the *Reader Unit Set* display. The use of a standard SD card is not possible.

Combination with shunt components

The implantable shunt components from Christoph Miethke GmbH & Co. KG can be safely combined with the *M.scio*. We recommend using implantable products from our company in combination with the *M.scio*.

The following products are particularly useful for intracranial pressure measurement in combination with the *M.scio*:

Product name	Item no.
<i>Ventricular Catheter</i> (with mandrin, length 250 mm)	FV077P
<i>Ventricular Catheter</i> (with mandrin, length 180 mm) with deflector (small, diameter 13 mm)	FV076P
<i>Ventricular Catheter</i> (with mandrin, length 250 mm) with deflector (large, diameter 16 mm)	FV078P

Product name	Item no.
<i>Prechamber</i> (small, diameter 14 mm)	FV035T
<i>Prechamber</i> (large, diameter 20 mm)	FV033T
<i>Pediatric CONTROL RESERVOIR</i> (small, diameter 14 mm)	FV066T
<i>CONTROL RESERVOIR</i> (large, diameter 20 mm)	FV047T
<i>Burrhole Reservoir</i> (small, diameter 14 mm)	FV039T
<i>Burrhole Reservoir</i> (large, diameter 20 mm)	FV028T
<i>Titanium Shutting Plug</i>	FV024T

For pressure measuring within a shunt system, the *M.scio* can be safely combined with further components, such as valves and *Peritoneal Catheters*.

If a (*Pediatric*) *CONTROL RESERVOIR* is used, care must be taken not to place this product between the ventricle and the *M.scio* when connecting to shunt components. Otherwise, the dynamics of the pressure signal may be distorted. For that reason, the *M.scio* cannot be combined with a *SPRUNG RESERVOIR*. The *M.scio* must be placed between the ventricle and the valve in order to determine the intracranial pressure.

The handling of the shunt components is described in the corresponding instructions for use.

2.10 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance can be downloaded from the following address:

<https://www.miethke.com/downloads/>

2.11 SUPPLEMENTARY PRODUCT INFORMATION

In accordance with EN 45502

- The *M.scio* identification number (ID) can be read with a [Single measurement] from the *Reader Unit Set* display, making it possible to identify the implant. The assignment of the ID to the serial number (SN) of the *M.scio* can be found on the SD card label that forms part of each *M.scio* scope of delivery.
- Approval for applying the CE mark for active implantable medical devices (in accordance with Directive 90/385/EEC) was first granted in 2011.

In accordance with ISO 7197

- Nuclear magnetic resonance imaging (MRI) up to a field strength of 3 tesla or computed tomography (CT) examinations can be performed without risk or impairment to the function of the *M.scio*. The *M.scio* is MRI Conditional. Artefacts may occur in MRI scans. The supplied catheters are MRI Safe. The documents on MRI safety are available from this website: <https://miethke.com/downloads/>
- The *M.scio* as well as the entire shunt system are safely able to resist positive and negative pressures up to 100 mmHg during and after the procedure.

3.00 SYSTEM COMPONENT *Reader Unit Set*

3.01 PRODUCT DESCRIPTION

3.01.01 APPROVED COMPONENTS

The *Reader Unit Set* consists of the Reader Unit, antenna and plug-in power unit components. No additional accessories are required.

1. Reader Unit

- 1.1 ON/OFF switch
- 1.2 Display screen
- 1.3 Function key
- 1.4 Antenna socket
- 1.5 Connector for plug-in power unit
- 1.6 SD card slot with sealing plug

2. Antenna

3. Plug-in power unit

- 3.1 Plug
- 3.2 Control light
- 3.3 Attachment for EU/UK

The Reader Unit and antenna are type BF application parts.

3.01.02 SCOPE OF DELIVERY

Box content	Number
<i>Reader Unit Set</i> (including EU/UK attachment for plug-in power unit)	1
Instructions for use for the <i>Reader Unit Set</i>	1
Case (including key)	1
Original packaging incl. mechanical damping	1

3.01.03 CALIBRATION

The *Reader Unit Set* contains a barometric pressure sensor for measuring the air pressure. Annual calibration of this sensor is necessary (see Section 3.06 Technical support). The user does not need to calibrate the product.

3.01.04 OPERATING CONDITIONS

Operating conditions of the <i>Reader Unit Set</i>	
Relative humidity	30 % to 75 %, non-condensing
Ambient temperature	10°C to 40°C
Atmospheric pressure	800 to 1100 hPa

3.01.05 PRODUCT LIFE

These medical devices are designed to ensure their precise and reliable operation over long periods of time. The expected service life of the *Reader Unit Set* is 5 years after first use, provided that the product is exposed to normal usage conditions and is properly maintained (see Section 3.06 Technical support). However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons.

3.01.06 PRODUCT CONFORMITY

The product meets current regulatory requirements as amended:

- ▶ (EU) 2017/745 (MDR)
- ▶ IEC 60601-1
- ▶ IEC 60601-1-2
- ▶ EN 45502-1
- ▶ ANSI/AAMI NS28

3.02 IMPORTANT SAFETY INFORMATION

3.02.01 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions to avoid injuries and life-threatening situations and to avoid jeopardising the warranty and liability.

WARNING

- ▶ The use of non-approved components poses a risk to users and patients and may damage the *Reader Unit Set* (see Section 3.01.01 Approved components). The product must not be modified.
- ▶ To avoid electric shocks and damage to the device due to penetrating liquids, the sealing plug must be inserted back into the device after removing the SD card.

CAUTION

- ▶ Because of the risk of injury from incorrect operation, the product user must have participated in the product training before using the product for the first time. Please contact Christoph Miethke GmbH & Co. KG for information about product training.
- ▶ Observe notes on electromagnetic compatibility (EMC)
- ▶ Observe maintenance instructions
- ▶ Prior to use, check the proper condition and functioning of the product

**CAUTION**

- ▶ **Do not operate the product in the vicinity of flammable substances (e.g. anaesthetics)**
- ▶ **The product must be set up so that the mains plug can easily be disconnected from the power supply**
- ▶ **Only operate the *Reader Unit Set* outside the operating range of an MRI**
- ▶ **Thoroughly clean the factory-new product after removing the transport packaging according to the manufacturer's instructions**
- ▶ **To prevent nosocomial infection and multiple resistances, the device should be disinfected before and after each use.**

3.02.02 COMPLICATIONS AND RESIDUAL RISKS

The following complications can occur in conjunction with the *M.scio* system:

- ▶ Headaches, dizzy spells, mental confusion, vomiting in cases of possible leakage from the *M.scio*/shunt and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid
- ▶ Wound healing disorders because of the height of the *M.scio*, dome-angled

A physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the *M.scio*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis)/allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation
- ▶ Local shunt irritation/allergic reaction

3.02.03 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible authority of the EU member state in which you work.

In Germany, the relevant authority is the Federal Institute for Drugs and Medical Devices (BfArM). For current contact information, please visit the BfArM website: <https://www.bfarm.de>.

3.02.04 PATIENT INFORMATION

The attending doctor is responsible for informing the patient and/or his/her proxy in advance. This includes a comprehensive description of the course of the operation, the surgical technique and the medical devices to be used. With regard to implanted medical devices, the patient must be informed about

- ▶ Warnings, precautionary measures to be taken; restrictions on use in relation to the medical device; information that ensures the safe use of the medical device; contraindications
- ▶ General information on quality and quantity of the materials and substances that the patient may be exposed to
- ▶ Anticipated service life of the medical device together with all necessary follow-up measures.

3.03 TRANSPORT AND STORAGE

3.03.01 TRANSPORT

Transport conditions for the *Reader Unit Set*

Ambient temperature	0 °C ... 50 °C
Atmospheric pressure	596 hPa ... 1100 hPa
Relative humidity	15 % ... 95 %

**NOTE**

To prevent possible damage during transport, the *Reader Unit Set* must be shipped in its original packaging.

3.03.02 STORAGE

The medical devices must always be stored in a clean and dry place.

Storage conditions for the *Reader Unit Set*

Ambient temperature	10 °C ... 40 °C
Atmospheric pressure	800 hPa ... 1100 hPa
Relative humidity	15 % ... 95 %

3.04 USE OF THE PRODUCT

3.04.01 INTRODUCTION

The *M.scio* system can be used in two scenarios in order to determine the intracranial pressure:

- ▶ *M.scio* implanted without shunt system
- ▶ *M.scio* integrated into a shunt system

In both scenarios, the *Reader Unit Set* is used for the telemetric reading and display of the pressure values.

3.04.02 SAFETY NOTICES AND WARNINGS

**WARNING**

The use of non-approved components poses a risk to users and patients and may damage the *Reader Unit Set*. Only the original plug-in power unit may be used.

3.04.03 INITIAL START-UP

**NOTE**

Allow the *Reader Unit Set* to acclimatize for approx. 3 hours at room temperature.

The *Reader Unit Set* is equipped with a battery which permits independent operation without external power supply for up to 5 hours on a 100 % charge. The battery needs to be charged before first using the device. The plug-in power unit can be used for this purpose. It takes approx. 3 hours for the battery to be fully charged.

The *Reader Unit Set* can only be operated if the battery has a sufficient charge level. If the *Reader Unit Set* has an insufficient charge, the device switches off. The plug-in power unit must be connected for charging. The *Reader Unit Set* can also be operated during the charging process (with connected plug-in power unit).

Charging is not possible at an ambient temperature of more than 35 °C.

Connecting the power supply

The mains voltage must correspond to the voltage range stated on the label of the plug-in power unit for the *Reader Unit Set*.

1. Insert the connector of the plug-in power unit into the socket of the *Reader Unit*.
2. Plug the plug-in power unit into a wall power socket.

Insert SD card

NOTE

To prevent damage due to improper use, do not touch the contact points of the SD card.

1. Remove sealing plug from the SD card slot
2. Insert the SD card assigned to the respective *M.scio* into the SD card slot of the Reader Unit until it locks into place.

Tap the SD card slightly to remove.

WARNING

To avoid electric shocks and damage to the device due to penetrating liquids, the sealing plug must be inserted back into the device after removing the SD card.

3.04.04 FUNCTION TEST

Preparation

CAUTION


Before use, check that the *Reader Unit Set* is in good condition and fully functional and check the system settings, such as the pressure unit (see the "Units" section in section 3.04.09 Settings).

- To be able to determine the battery charge level, the function test can be carried out without a plug-in power unit. We recommend regularly checking the charge level of the Reader Unit battery.
- If the charge level is not sufficient, the battery can be charged by plugging in the plug-in power unit of the *Reader Unit Set*. The control lamp of the plug-in power unit must light up as soon as it is plugged into the socket of the on-site installation.

- Make sure that the *Reader Unit Set* has no visible damage, such as to the housing, keyboard, display and power unit.
- The function of the following components must be checked in the given order:


Switch on

1. Press On/Off button
2. Automatic self-test of the device after activation including display and speaker volume test
3. The following is displayed on the screen:
 - [Selftest ...]
 - [booting ...]
4. The following content is then shown on the display:

The following is displayed on the screen of the *Reader Unit Set* during operation with the plug-in power unit:

Display for operation
with plug-in power unit

If necessary, the system time can be corrected (see "Settings" section).

Switch off

- Press On/Off button


Operation of the *Reader Unit Set* can be safely shut down at any time by pressing the on/off button.

3.04.05 GENERAL OPERATING INFORMATION

The menu-controlled operation of the device uses the four function keys. The context-related function of these soft keys is indicated by the symbols shown on the display for each of these keys.

Display

- A Date
- B Time
- C SD card
- D Assignment of memory space on SD card
- E Battery status
- F OK button
- G Down arrow
- H Exit menu

The preselected submenu is highlighted by a border. For improved navigation, the currently selected position in the submenu is highlighted by a scroll bar in addition to the border. Standard operation involves using the <OK>, <Up arrow>, <Down arrow> buttons and the <Exit menu> button.

When not in use, the *Reader Unit Set* switches to Standby Mode. Depending on the setting selected, this switches the monitor off after 1 to 5 minutes. The device is reactivated by pressing one of the four function keys.

The following context-related symbols are shown on the display of the *Reader Unit Set*:

	<Info>	For displaying additional information or starting with the information menu
	<Menu>	For starting with the selection menu
	<Fast measurement>	Fast measurement can be started without further selection
	<Up arrow>	Upward navigation of cursor
	<Down arrow>	Downward navigation of cursor
	<Exit menu>	For exiting the current menu
	<OK>	Activates or confirms the selected function
	<Start>	Starts measurement
	<Stop>	Stops the current measurement
	<Delete>	For deleting a saved measurement in the [Data management] menu
	<Star>	For placing a marker during a [Continuous measurement]
	<Editing diagram>	For accessing the diagram processing menu
	<Zoom menu>	For accessing the zoom function
	<Cursor menu>	For accessing the cursor function
	<Pressure axis values>	For adjusting the display of min. and max. values of the pressure axis
	<Exit cursor menu>	For returning to the diagram processing menu
	<Exit zoom menu>	For returning to the diagram processing menu
	<Right arrow>	Navigate forward
	<Left arrow>	Navigate backward
	<On>	On
	<Off>	Off


	<Confirmation>	Confirm
	<Zoom in>	For zooming into the display of a sequence in the [Data management] menu
	<Zoom out>	For zooming out of the display of a sequence in the [Data management] menu
	<Mute>	For deactivating the acoustic signal

3.04.06 INFO MENU

The following information can be accessed from the info menu:

- ▶ [Hotline & Service]
- ▶ [System details]
- ▶ [Manufacturer info]
- ▶ [Only maintenance]

To access the [Info] menu, press the <Info> button. The following information is displayed on the screen:

Preselect the desired submenu with the <Up arrow> and <Down arrow> keys, and confirm with <OK>. Return to the previous menu item with the <Exit menu> button.

The submenus contain the following information:

[Hotline & Service]

- ▶ [Hotline & Service]
- ▶ [Tel: +49 331 620 83-0]

[System details]


- ▶ [Product Name: Reader Unit]
- ▶ [Article Number: 7510 0000]
- ▶ [Serial Number: XXXXX]
- ▶ [Software Version: 2.XX]
- ▶ [Service Date: dd/mm/yyyy]

[Manufacturer info]


- ▶ [Christoph Miethke GmbH & Co. KG]
[Ulanenweg 2 | 14469 Potsdam]
[Germany]

3.04.07 PERFORMING MEASUREMENTS

Open the [Menu] by pressing <Menu>. The following information is displayed on the screen:

Open the [Measurement] submenu by pressing <OK>. The following information is displayed on the screen:

There are three measuring modes:

1. [Single measurement]: The pressure value measured at a point is displayed as a single value and saved on the SD card. Manual averaging of the measured values from 8 to 10 repeated single measurements is recommended.
2. [Continuous measurement]: Sequential single measurements are carried out, the recorded measured values are displayed as a curve and stored on the SD card.
3. [Fast measurement]: This records sequential single measurements at a high sampling rate (approx. 44 measurements per second), which are then displayed as a curve and stored on the SD card.


The telemetric link between the antenna of the *Reader Unit Set* and the *M.scio* may be disrupted by metal parts or operation of another *Reader Unit Set* in the vicinity of the implant. In this case, increase the distance to the metal parts or to the additional *Reader Unit Set*. After that, a measurement can be started.

A raised patient body temperature may cause a malfunction. In readout mode, the temperature in the *M.scio* may rise. An integrated temperature safety device automatically stops the measurement at 39 °C in the implant.

Measurement may be interrupted if the free memory capacity of the SD card is exhausted. The memory status must be checked before starting the measurement. Measurement data that is no longer required can be deleted. Measured data can only be saved on the SD card associated with the *M.scio*.

Single measurement

Press <OK> to select the [Single measurement] menu. The following information is displayed on the screen:


Press the <Start> button to start the [Single measurement]. The following information is displayed on the screen:

In addition to the measurement data, the identification number (ID) of the *M.scio* is also displayed. If the SD card is lost, it can be reordered by specifying the ID or serial number of the *M.scio*.


Continuous measurement

In the [Measurement] menu, press <Up arrow> or <Down arrow> to select the [Continuous measurement] menu and confirm with <OK>. The following information is displayed on the screen:

Press the <Start> button to start the [Continuous measurement].

The following information is displayed on the screen:

Press <Star> to place a marker. Multiple markers can be added during a measurement. Press the <Stop> button to stop the measurement. Markers permit the circumstantial analysis of the measured data.


Symbols for reception quality

Symbol	Explanation
	Communication started
	Distance from antenna to measuring cell: - OK
	Distance from antenna to measuring cell too close: - increase distance
	Distance from antenna to measuring cell too large: - decrease distance


Fast measurement

There are two ways to start a [Fast measurement]:


On the start screen, press the <Fast measurement> button directly and then start the measurement with the <Start> button. Alternatively, you can press the <Menu> button to display the following submenu:


Go to the [Measurement] menu by pressing the <OK> button. In the [Measurement] menu, use the <Up arrow> and <Down arrow> buttons to select the [Fast measurement] menu. The following is displayed on the screen:

Select the [Fast measurement] with the <OK> button. The following is displayed on the screen:


Press the <Start> button to start the fast measurement. During fast measurement, the following is displayed on the screen:


Press the <Stop> button to stop the measurement. For an explanation of the symbols for reception quality: see the “Continuous measurement” section.

3.04.08 ADMINISTRATION OF THE MEASURED DATA


On the start screen, press <Menu> to open the following submenu:

Press <Up arrow> and <Down arrow> to select the [Data management] menu and confirm with <OK>.

Measured data are listed chronologically (from the start of measurement) and shown as follows:

Press <Up arrow> and <Down arrow> to select the desired file and confirm with <OK>. Measured data files are labelled as follows:


Symbol	Explanation
	Single measurement
	Continuous measurement
	Fast measurement

An hour glass is shown when loading the data for a continuous or fast measurement and the loading time is displayed in seconds.

Single measurement


Press <Delete> to delete the measurement data file.

Deletion has to be confirmed by pressing <Confirmation>. Press <Exit menu> to return to the previous view. Press this button again to return to the selection menu.

Continuous measurement

Press <Delete> to delete the measurement data file. Deletion has to be confirmed by pressing <Confirmation>. Press <Exit menu> to return to the previous view. Press <Editing diagram> to open the <Diagram processing> menu. The following information is displayed on the screen:

Press <Zoom menu> to open the <Zoom> menu. The following information is displayed on the screen:


Press <Zoom in> or <Zoom out> to display an expanded or compressed time diagram of the measured sequence. Press <Exit zoom menu> to return to the <Diagram processing> menu.

Press <Cursor menu> in the <Diagram processing> menu to open the following screen:


Press <Exit cursor menu> to return to the <Diagram processing> menu. Press <Right arrow> or <Left arrow> to move the time axis to the left or right.

Press <Pressure axis values> in the <Diagram processing> to open the following screen:


Use the <Up arrow> and <Down arrow> keys to manually set the maximum or minimum scale values at the cursor position. Make sure that the [Automatic] function is switched off. If the [Automatic] function is selected, the setting of maximum or minimum values is irrelevant. In this function, axis scaling is automatic according to the measured data. Press <Right arrow> to toggle the cursor between [Maximum], [Minimum] and [Automatic]. Press <Exit menu> to return to the measured value display.

Fast measurement


Press <Delete> to delete the measurement. Deletion has to be confirmed by pressing <Confirmation>. Press the button to return to the previous view. Press <Exit menu> to return to the selection menu. Press <Editing diagram> to open the <Diagram processing> menu.

The following information is displayed on the screen:


Press <Zoom menu> to open the <Zoom> menu.


Press <Zoom in> or <Zoom out> to display an expanded or compressed time diagram (zoom) of the measured sequence. Press <Exit zoom menu> to return to the <Diagram processing> menu. Switching between displays can take a few seconds, depending on file size.

Press <Cursor menu> in the <Diagram processing> menu to open the following screen:


Press <Exit cursor menu> to return to the <Diagram processing> menu. Press <Left arrow> or <Right arrow> to move the time axis to the left or right. Press <Pressure axis values> in the <Diagram processing> menu to open the following screen:


Use the <Up arrow> and <Down arrow> keys to manually set the maximum or minimum scale values at the cursor position. Make sure that the [Automatic] function is switched off. If the [Automatic] function is selected, the setting of maximum or minimum values is irrelevant. In this function, axis scaling is automatic according to the measured data. Press <Right arrow> to toggle the cursor between [Maximum], [Minimum] and [Automatic].

3.04.09 SETTINGS


On the start screen, press <Menu> to open the following submenu:

Use the <Up arrow> and <Down arrow> buttons to select the [Settings] menu.


Go to the [Settings] menu by pressing the <OK> button. The following is displayed on the screen:

Date/Time


Go to the [Date/Time] menu by pressing the <OK> button. The following information is displayed on the screen:


Press <Right arrow> to change the cursor position. Use the <Up arrow> and <Down arrow> keys to change the values at the cursor position. Modified values are saved immediately.

Units

In the [Settings] menu, use the <Up arrow> and <Down arrow> buttons to select the [Units] menu.

Go to the [Units] menu by pressing the <OK> button. The following is displayed on the screen:


Press <Up arrow> and <Down arrow> to select the desired unit and confirm with <OK>.


Measurement interval

These settings only work in Continuous Measurement mode.

In the [Settings] menu, use the <Down arrow> and <Up arrow> buttons to select the [Measurement interval] menu.


Go to the [Measurement interval] menu by pressing the <OK> button. The following information is displayed on the screen:


Press <Down arrow> and <Up arrow> to select the desired measurement interval and confirm with <OK>. Possible settings range from 1–300 s. The available memory space on the SD card is also shown.

Hide graph

In the [Settings] menu, use the <Down arrow> and <Up arrow> buttons to select the [Hide graph] menu.


Go to the [Hide graph] menu by pressing the <OK> button. The following is displayed on the screen:


[Display graph during continuous measurement] can be selected or deselected.

Language

In the [Settings] menu, use the <Down arrow> and <Up arrow> buttons to select the [Language] menu.


Go to the [Language] menu by pressing the <OK> button. The following information is displayed on the screen:


Press <Up arrow> and <Down arrow> to select the desired language and confirm with <OK>.

Volume

In the [Settings] menu, use the <Up arrow> and <Down arrow> buttons to select the [Volume] menu.

Go to the [Volume] menu by pressing the <OK> button. The following information is displayed on the screen:


Press <Up arrow> and <Down arrow> to select the desired volume and confirm with <OK>. Available settings are 1–5. When setting the volume, a corresponding acoustic signal is emitted.


The setting determines the volume of acoustic alerts. Exception: Error alerts are generally emitted at volume level 5.

Standby mode

In the [Settings] menu, use the <Up arrow> and <Down arrow> buttons to select the [Standby mode] menu.


Go to the [Standby mode] menu by pressing the <OK> button. The following information is displayed on the screen:


Press the <Up arrow> and <Down arrow> buttons to select the desired time after which the Reader Unit switches to Standby Mode and confirm with <OK>. Available settings are 1-5 minutes.

Brightness

In the [Settings] menu, use the <Up arrow> and <Down arrow> buttons to select the [Brightness] menu.


Go to the [Brightness] menu by pressing the <OK> button. The following information is displayed on the screen:


Press <Up arrow> and <Down arrow> to select the desired background brightness and confirm with <OK>. Available settings are 000–100 (in steps of 5).

Contrast

In the [Settings] menu, use the <Up arrow> and <Down arrow> buttons to select the [Contrast] menu.

Go to the [Contrast] menu by pressing the <OK> button. The following information is displayed on the screen:

Press <Up arrow> and <Down arrow> to select the desired contrast and confirm with <OK>. Available settings are 000–100 (in steps of 5).

3.05 MANUAL CLEANING AND DISINFECTION OF THE *Reader Unit Set*

WARNING

Risk of electric shock and fire

- ▶ Disconnect from mains before cleaning.
- ▶ Ensure that no fluids can penetrate into the product; insert the sealing plug into the SD card slot of the Reader Unit.
- ▶ Do not use flammable or explosive detergents or disinfectants.

NOTE

Mechanical cleaning or disinfection as well as inappropriate detergents or disinfectants can damage or destroy the product

- ▶ Only disinfect the product manually
- ▶ Never sterilize the product
- ▶ Use detergents/disinfectants approved for surface cleaning only as specified by the manufacturer.
- ▶ Observe specifications regarding concentration, temperature and contact time.

Procedure

Wipe disinfection for electrical devices without sterilization before and after each use:

Phase I

- ▶ Remove any visible residue with a disposable disinfectant wipe.
- ▶ Wipe all surfaces of the visually clean product with a fresh disposable disinfectant wipe.
- ▶ Observe the prescribed contact time.

Parameter	Description
Step	Wipe disinfection
T (°C/°F)	RT (room temperature)
t (min)	≥ 1
Conc. (%)	-
Water quality	-
Chemical specifications	Meliseptol HBV wipes 50 % Propan-1-ol

Inspection

- ▶ Always check the product for damage after cleaning/disinfection.
- ▶ Immediately discard damaged products.

Storage

- ▶ Pack *Reader Unit Set* in case

3.06 TECHNICAL SUPPORT

3.06.01 CALIBRATION, MAINTENANCE AND REPAIR

The *Reader Unit Set* contains a barometric pressure sensor for measuring the air pressure. Annual calibration is required to ensure that predefined tolerance limits are maintained. If the calibration of the barometric pressure sensor has expired, a corresponding message appears on the display of the *Reader Unit Set*.

CAUTION

Failure to perform the annual calibration can result in the drifting of the barometric pressure sensor outside of the tolerance limits.

The device must therefore be sent to our Technical Support department annually. The requirements defined in the "Storage and transport" section as well as in the "Operating conditions" section must be observed.

As part of the calibration, a comprehensive functional check of the device is also performed.

The due date for the next calibration can be found in the menu under [Menu info] > [System details] > [Service Date].

The battery can be replaced by Technical Support at the end of the associated operating life. For calibration, maintenance and repair, please contact Technical Support:

Technical Support:

Christoph Miethke GmbH & Co. KG

Technical Support

Ulanenweg 2

14469 Potsdam

Tel.: +49 331 62083-0

Fax: +49 331 62083-40

E-mail: technicalsupport@miethke.com

Modifications of medical devices can result in a loss of guarantee or warranty claims. Christoph Miethke GmbH & Co. KG only accepts responsibility for the safety, reliability and performance of the device if:

- ▶ The device is used in accordance with the instructions for use.
- ▶ New settings, modifications and repairs have only been carried out by persons authorized by us.
- ▶ The electrical installation of the respective room conforms to national standards (IEC definition).

3.06.02 SAFETY CHECK

The Ordinance on the Operation of Medical Devices (MPBetreibV) requires safety checks to be carried out. No corresponding check pursuant to the MPBetreibV is carried out by the manufacturer as part of the annual calibration of the barometric pressure sensor within the scope of the technical service. The operator is obliged to conduct a safety check after this type of maintenance measure and before start-up.

An annual inspection with the following scope is recommended:

1. Visual inspection (also refer to the scope of the visual inspection)
2. Check the functionality based on the operating manual
3. Check error messages on the display
4. Electrical safety – Measurement of the leakage currents pursuant to IEC 62353 as amended
5. Creation of a test report

Scope of the visual inspection:

1. Is the maintenance book available?
2. Do the reader and the antenna as well as the antenna cable, power supply unit and connection points display mechanical defects?
3. Are all labels complete and legible?
4. Are all housing screws firmly tightened?
5. Is the antenna cable firmly connected to the Reader Unit?
6. Are parts within the housing (reader and antenna) loose? The device must be shaken carefully for the test.
7. Are all approved components in place?

3.07 DISPOSAL


Adhere to national regulations when disposing of or recycling the product, its components or its packaging!

The disposal of electrical appliances is regulated by the German Electrical and Electronic Equipment Act (ElektroG), which is derived from the European WEEE Directive. According to the requirements of the applicable ElektroG, a product marked with the symbol shown above must be taken to the separate collection of electrical and electronic equipment by handing it over to an appropriately certified disposal company. Alternatively, the disposal of the *Reader Unit Set* within the European Union can be taken over by the manufacturer free of charge. If you have any questions regarding the disposal of the product, please contact Christoph Miethke GmbH & Co. KG; also refer to section 3.06 Technical support.

3.08 TROUBLESHOOTING

If an error occurs, it is shown on the display of the *Reader Unit Set*.

Example of an error message:

Press <Info> for additional information.

Faults with error text on the display of the *Reader Unit Set*

Indication on the display	Cause	Fault detection/rectification
Battery flat - Auto off	Battery capacity exhausted (0 %)	All data are saved after 2 min. Reader Unit switches off automatically. Connect original plug-in power unit.
Battery voltage incorrect - use original power supply unit	Battery voltage of the <i>Reader Unit Set</i> is too low	Reader Unit switches off automatically after 20 seconds. Connect original plug-in power unit.
Low battery voltage	Battery voltage too low	Backlight switches off after 3 seconds. Connect original plug-in power unit. No disruption to ongoing measurements.
Flashing of crossed-out antenna symbol <Info> button: Antenna faulty	Antenna defective	Switch the device off and back on. If the error persists, contact Technical Support.
Flashing of crossed-out antenna symbol <Info> button: Antenna not plugged in	Antenna not connected at the start of measurement - or - Antenna was disconnected during measurement	Connect antenna: Measurement restarts - or - Connect antenna: measurement continues.
Flashing of crossed-out antenna symbol <Info> button: No communication	Data recording interrupted during continuous measurement (interruption of telemetric link)	Measurement resumes automatically once communication is restored.
Flashing of crossed-out antenna symbol <Info> button: SD card has been removed. Measurement possible	SD card was removed during continuous measurement	Insert SD card. Restart measurement.
Dataset defective! <Info> button: File cannot be opened	File validation failed	File cannot be opened. Try again.
Continuous key activation Keypad error	Persistent activation of a key > 60 seconds	Release the key.
Pressure readings out of range	Pressure data measured from implant not plausible - physiologically meaningless data	Measurement continues. If the error persists, contact Technical Support.
Problem with input voltage	The voltage of the plug-in power unit is too high	Reader Unit switches off automatically after 20 seconds. Use original plug-in power unit.
Wrong SD card inserted! Remove card! <Info> button: Measurement without data storage possible - or - Insert SD card with correct ID XXXXXXXXX!	Measurement starts without SD card. An SD card not designated for use with the implant is used during measurement	Insert the SD card appropriate to the <i>M.scio</i> into the <i>Reader Unit Set</i> . The implant is assigned to the respective SD card via the identification number (ID).
Wrong SD card inserted! Remove card! <Info> button: Measurement without data storage possible - or - Insert SD card with correct ID XXXXXXXXX!	SD card with different ID to the implant has been used	Insert the SD card appropriate to the <i>M.scio</i> into the <i>Reader Unit Set</i> . The implant is assigned to the respective SD card via the identification number (ID).

Indication on the display	Cause	Fault detection/rectification
Wrong implant - restart measurement! <Info> button: Switch to another implant during continuous measurement not possible!	Data from another implant received after starting a continuous measurement	Increase the distance between the two implants.
Device temperature outside of range <Info> button: Device temperature range from 10 °C to 40 °C permissible!	Temperature in the <i>Reader Unit Set</i> exceeds calibrated interval	The <i>Reader Unit Set</i> can only be used at device temperatures from 10 °C to 40 °C. Current measurement stopped.
Problem with internal voltage	Internal voltage of device too high/ too low	Reader Unit switches off automatically after 20 seconds. Contact Technical Support.
SD card faulty! <Info> button: Measurement without data storage possible!	SD card not writeable or readable (dirt, corrosion, faulty contacts)	Use of implant's internal calibration data. Data are not saved.
SD card faulty!	SD card not readable (dirt, corrosion, faulty contacts)	Check SD card for damage or dirt.
SD card inserted! Restart measurement! <Info> button: Storage of measured data possible after restart of measurement!	Measurement starts without SD card. The designated SD card for use with that implant is inserted during measurement	Restart measurement.
SD card missing! <Info> button: Insert SD card!	No SD card inserted in Data Management mode	Insert SD card.
SD card missing! <Info> button: Measurement without data storage possible - or - Insert SD card with correct ID!	SD card not inserted	Insert the SD card appropriate to the <i>M.scio</i> into the <i>Reader Unit Set</i> . The implant is assigned to the respective SD card via the identification number (ID).
SD card not readable! <Info> button: Measurement without data storage possible!	Insertion of incorrectly formatted or unformatted SD card	Insert the SD card appropriate to the <i>M.scio</i> into the <i>Reader Unit Set</i> . The implant is assigned to the respective SD card via the identification number (ID).
SD card not readable! <Info> button: Measurement without data storage possible!	SD card missing - or - SD card ID not readable - or - SD card contains no calibration data	Insert the SD card appropriate to the <i>M.scio</i> into the <i>Reader Unit Set</i> . Measurement with data storage possible.
SD card memory full. Measurement without data storage possible!	SD card memory exhausted (100 %) during continuous measurement	Delete any data no longer required, measurement without data storage possible.
SD card memory full	SD card memory exhausted (100 %) during continuous measurement	Delete any data no longer required.
SD card memory almost full	SD card memory almost exhausted (99 %) during continuous or fast measurement	Delete any data no longer required.
System error <Info> button: Ambient pressure not readable	Barometric pressure not readable	Device function locked. Switch the device off and back on. If the error persists, contact Technical Support.

Indication on the display	Cause	Fault detection/rectification
System error - incompatibility	Hardware and software versions not compatible	Device function locked. Switch the device off and back on. If the error persists, contact Technical Support.
System error - antenna incompatible	Hardware version of antenna and Reader Unit not compatible	Device function locked. Switch the device off and back on. If the error persists, replace antenna and/or contact Technical Support.
System error - ID data inadmissible	Implant identification data corrupted	Measurement stopped. Switch the device off and back on. If the error persists, contact Technical Support.
System error - implant voltage out of range	Implant voltage outside permissible range	Measurement stopped. Switch the device off and back on. If the error persists, contact Technical Support.
System error - calibration data incorrect	Calibration data in implant corrupted or not readable (only applies if no SD card has been inserted)	Measurement stopped. Switch the device off and back on. If the error persists, contact Technical Support. Measurement with inserted SD card possible.
System error - contact Technical Service	System test detected an error	Device function locked. Switch the device off and back on. If the error persists, contact Technical Support.
Keypad faulty	A pressed key was registered when switching on the <i>Reader Unit Set</i>	Release the key. Switch the device off and back on.
Temperature increase inadmissible	A spontaneous temperature increase in the implant above 39 °C	Measurement stopped. Leave to rest for 10 minutes. If the error persists, contact Technical Support.
Temperature readings out of range	Temperature data measured from implant not plausible - physiologically meaningless data	Measurement stopped. Switch the device off and back on. If the error persists, contact Technical Support.
Time out! Restart measurement!	The time between starting measurement and successfully establishing communication was exceeded (60 seconds).	Restart measurement. Optimize the distance between the antenna and the implant.
Ambient pressure out of range <Info> button: Permissible ambient pressure range from 800 to 1100 mbar	The permissible barometric pressure during measurement is too low or too high.	The <i>Reader Unit Set</i> can only be used at barometric pressures from 800 to 1100 hPa. Current measurement stopped.

Other faults

Indication on the display	Cause	Fault detection/rectification
Device cannot be switched on	Battery is completely flat	Connect the plug-in power unit. It takes approx. 6 hours to fully charge the battery. The <i>Reader Unit Set</i> can be operated during the charging process (with connected plug-in power unit). Note: Charging is not possible at an ambient temperature of more than 35°C.
Device switches off	Unfavourable operating conditions (e.g. low humidity or unsuitable floor coverings)	Switch device back on. If the error persists, contact Technical Support.

3.09 TECHNICAL DATA AND PERFORMANCE DATA

Designation	Values and standards
Voltage range: Reader unit Plug-in power unit	6 V (DC) 100–240 V (50–60 Hz)
Power consumption: Reader unit Plug-in power unit	1.4 A (DC) 0.15–0.3 A (50–60 Hz)
Telemetry: Operating frequency Bandwidth Type of modulation Transmission power of the antenna Magnetic field strength at 10 m Functional distance from antenna to <i>M.scio</i>	133 kHz 125 kHz – 135 kHz Amplitude modulation max. 0.8 W <30 dB μ A/m 10 to 30 mm
Measurement accuracy of the <i>M.scio</i> system Measurement accuracy within the first 10 days of the pressure range	Pressure range (relative to atmospheric pressure): -50 mmHg ... +100 mmHg -50 mmHg ... -20 mmHg: +/- 10 % -20 mmHg ... +20 mmHg: +/- 2 mmHg +20 mmHg ... +100 mmHg: +/- 10 %
Reading accuracy	Depending on the scaling of the y axis
Battery: Type Service life Charging cycle Self-discharge	Lithium polymer min. 5 years at least 250 charging cycles Remaining charge after 3 months (storage) > 70%
Flammability category of housing	UL 94 HB
Moisture protection or protection against ingress of liquids: Reader unit Antenna Plug-in power unit	 IP44 IP44 IP40
Strength: Impact test Drop test	pursuant to IEC 60601-1: 2012 15.3.2 pursuant to IEC 60601-1: 2012 15.3.4.1
Weight: Reader unit Antenna Plug-in power unit	0.600 kg 0.215 kg 0.127 kg
Dimensions (W x H x D): Reader unit Antenna (without cable) Plug-in power unit	144 x 270 x 65 mm 100 x 250 x 25 mm 77.5 x 31.5 x 41 mm
Protection class (plug-in power unit)	II

3.10 ELECTROMAGNETIC COMPATIBILITY

The *Reader Unit Set* fulfils the requirements of IEC 60601-1-2 in the current version.

3.10.01 ELECTROMAGNETIC EMISSION

Although the device complies with class B HF emissions in accordance with CISPR 11, it may only be used for professional healthcare facilities.

Guidelines and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in an electromagnetic environment as specified below. It is the user's responsibility to ensure that the device is operated in such an environment.		
Electromagnetic interference measurements	Compliance	Guidelines on the electromagnetic environment
HF emissions according to CISPR 11	Complies with Group 1	The device uses HF energy only for its internal function. Its HF emissions are thus very low and not likely to cause interference in nearby electronic devices.
HF emissions according to CISPR 11	Complies with Class B	The device is suitable for use in all establishments including domestic settings and areas immediately connected to a public power supply which also supplies buildings used for residential purposes.
Harmonic emissions according to IEC 61000-3-2	Complies with Class A (IEC 61000-3-2)	-
Voltage fluctuations/flicker emissions according to IEC 61000-3-3	Complies	-

3.10.02 ELECTROMAGNETIC IMMUNITY

WARNING

Risk of electromagnetic interference!

Portable HF communication devices (radio devices) (including their accessories, such as antenna cables and external antennas) should not be used at a distance closer than 30 cm (or 12 inches) to the Reader Unit and antenna. A distance of at least 80 cm must be ensured in the presence of electromagnetic interference with spot frequencies of 385 MHz or 450 MHz. Non-compliance may lead to a reduction of the device's performance characteristics. Electromagnetic interference can lead to a device shut-down. In this case, the device must be restarted and the measurement repeated.

3.11 SYMBOLS USED ON THE PRODUCT AND LABELLING

Symbol	Explanation
	EU conformity marking, xxxx indicates the identifier of the responsible notified body
	Medical device
	Manufacturer
	Date of manufacture
	Production lot number, batch
	Article number
	Serial number
	UDI number (Unique Device Identifier)
	UDI-DI number
	Do not use if package is damaged
	Store in a dry place
	Temperature limitation
	Atmospheric pressure limitation
	Humidity limitation
	Consult instructions for use/electronic instructions for use
	Electronic instructions for use/electronic accompanying documents
	Caution, read the documentation supplied with the product
	Free of natural rubber latex, latex-free

Symbol	Explanation
	Indicates that, in the USA, the product may only be issued to physicians
	Not MRI safe
	Non-ionising electromagnetic radiation
	Electrostatic sensitive devices
	Electronic devices: Dispose of properly, do not dispose of electronic products with household waste
	ON/OFF switch
	Function key: The 4 function keys are used to perform the options shown on the display screen
	IP code, degree of protection of the housing against foreign bodies and water
	IP4X - protection against granular foreign bodies, IPX0 - no protection against water, IPX4 - splash protection
	Protection class II
	Application part of type BF
	Antenna symbol
	Antenna socket
	SD card
	SD card slot
	DCIN connector for plug-in power unit

4.00 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany