M.scio®

READING INNER VALUES. FOR THE BIG PICTURE.

(DE) Gebrauchsanweisung | (EN) Instructions for use | (FR) Mode d'emploi | (ES) Instrucciones de manejo
(PT) Instruções de utilização | (IT) Istruzioni per l'uso | (NL) Gebruiksaanwijzing

 www.miethke.com

(US) This Instructions for Use is NOT intended for United States users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0
info@miethke.com

TABLE OF CONTENTS

0.00 PREFACE AND IMPORTANT INFORMATION	20
1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	21
1.01 Explanation of the warnings	21
1.02 Display conventions	21
1.03 Other accompanying documents and additional information	21
1.04 Feedback on the instructions for use	21
1.05 Copyright, disclaimer, warranty and other information	21
2.00 DESCRIPTION OF THE <i>M.scio</i> SYSTEM	22
2.01 Intended use	22
2.02 Clinical benefits	22
2.03 Indications	22
2.04 Contraindications	23
2.05 Intended patient groups	23
2.06 Intended users	23
2.07 Intended use environment	23
2.08 Functional principle	23
2.09 System components	24
2.10 Brief report on safety and clinical performance	24
2.11 Supplementary product information	25
3.00 SYSTEM COMPONENTS OF THE <i>M.scio</i> AND SD CARD	25
3.01 Product description	25
3.02 Important safety information	26
3.03 Transport and storage	28
3.04 Use of the product	28
3.05 Pressure measurement in exceptional situations	31
3.06 Explantation and disposal	31
3.07 Troubleshooting	31
3.08 Technical information	32
3.09 Symbols used for labelling	34
4.00 MEDICAL DEVICE CONSULTANTS	34

0.00 PREFACE AND
IMPORTANT INFORMATION

Preface

Thank you for purchasing the *M.scio*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product.
Your team at Christoph Miethke GmbH & Co. KG.

Relevance of the instructions for use

WARNING

Improper handling and use of this product can cause risks and damage. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

The *M.scio* forms a part of the *M.scio* System, which consists of the following components:

- *M.scio* with associated SD card
- *Reader Unit Set*

The *M.scio* System can be safely combined with our range of shunt components.

These instructions for use apply for the following product variants of the *M.scio* together with the respective associated SD card:

Product	Variant
<i>M.scio</i>	dome-angled
<i>M.scio</i>	dome-angled, with distal catheter
<i>M.scio</i>	flat-angled
<i>M.scio</i>	flat-angled, with distal catheter
<i>M.scio</i>	dome-inline
<i>M.scio</i>	dome-inline, with distal catheter
<i>M.scio</i>	flat-inline
<i>M.scio</i>	flat-inline, with distal catheter

Operation of the *Reader Unit Set* and the shunt components is described in the respective instructions for use.

Basic UDI-DI

The Basic UDI-DI of the *M.scio* product variants including the associated SD cards is:
404190600000000000000002RY

1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

1.01 EXPLANATION OF THE WARNINGS

DANGER

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

WARNING

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

CAUTION

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

NOTICE

Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for risk, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

1.02 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates product names
[...]	Square brackets indicate selectable menu items or information that are shown on the <i>Reader Unit Set</i> display.
<...>	Pointed brackets indicate context-dependent symbols on the <i>Reader Unit Set</i> display.

1.03 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

The latest version of the instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

The delivery includes a Patient Implant Card (IC) with IC label containing product information. The respective healthcare facility must record the information about the patient, the implant as well as the physician in charge on the IC,

and the IC label must be attached. The Patient Implant Card is thus intended to make all of the important information available to the patient in a compact form. The mentioned website also contains information about the symbols used on the Patient Implant Card together with a description on how the Patient Implant Card is to be completed by the healthcare institution. Patient labels containing information about the product are enclosed with the delivery. The patient labels provide the treating physician with all the product information in a compact form for the patient record. For individually purchased shunt components, the patient label must be added to the existing Implant Card. Helpful information for patients, particularly concerning the symbols on the Implant Card and the labels, is available on this website:

<https://www.miethke.com/ic/>

If you still need additional information despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

1.04 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests or remarks about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use.

1.05 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

2.00 DESCRIPTION OF THE M.scio SYSTEM

2.01 INTENDED USE

The M.scio system is used for diagnostic intracranial pressure measurement within the cerebrospinal fluid.

Due to the silicone membrane, the dome variants of the system component M.scio have the pumpability and puncturability of a conventional reservoir, i. e. they offer the possibility of therapeutic pressure reduction by removing cerebrospinal fluid, diagnostic removal of cerebrospinal fluid, administration of fluids and verification of the pressure values.

2.02 CLINICAL BENEFITS

Optimisation of diagnosis and therapy thanks to telemetric measurement of intracranial pressure values

- ▶ Use of a long-term implant
- ▶ Simple and easy reading of pressure values
- ▶ Identification of pathological pressure situations
- ▶ Low risk thanks to non-invasive measuring method
- ▶ MR Conditional for magnetic fields up to 3 tesla
- ▶ Increased sense of security in concerned patients and their relatives thanks to easily accessible readings
- ▶ Different variants for individual patient-specific treatment requirements
- ▶ Optional expandability of the M.scio system to a shunt system

Optimised patient management in shunt patients

Improved patient outcomes

- ▶ Optimised valve settings on the basis of the determined pressure readings
- ▶ Reduction of over- and underdrainage

Reduced strain on patients

- ▶ Avoidance of unnecessary clinical diagnostic procedures and their associated risks (such as exposure to radiation for imaging procedures and the use of invasive diagnostic techniques)
- ▶ Avoidance of unnecessary revisions thanks to functional monitoring of the shunt and ruling out of occlusions and shunt malfunctions

Cost savings

- ▶ Avoidance of unnecessary clinical procedures (such as medical imaging, pressure measurement and revisions)

Optimised diagnosis and treatment options thanks to using the M.scio “dome” variants

Increased options thanks to puncture

- ▶ CSF removal for manual pressure release and laboratory analyses
- ▶ Possibility of external reference pressure measurement
- ▶ Administration of fluids

Reduced strain on patients

- ▶ Pump tests for functional monitoring of shunts

Cost savings

- ▶ Avoidance of unnecessary clinical procedures (such as medical imaging and revisions)

2.03 INDICATIONS

The following indications apply to the M.scio System:

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Extended indications

- ▶ Shunt dependency
- ▶ Shunt dysfunction
- ▶ Treatment optimisation

2.04 CONTRAINDICATIONS

The following contraindications apply to the *M.scio* system:

Contraindications

- ▶ Coagulation disorders (risk of post-operative bleeding)
- ▶ Blood in the cerebrospinal fluid
- ▶ Infections or a suspected infection in the regions of the body affected by the implantation (e.g. skin infection, meningitis, ventriculitis, bacteraemia, septicaemia, also peritonitis if the *M.scio* is integrated into a Shunt System)

Relative contraindications

- ▶ High pressure and shock load caused by patient actions (including diving, boxing, football)
- ▶ Aggressive/autoaggressive patient behaviour may limit patient compliance during aftercare and complicate the readout using the *Reader Unit Set*. The *M.scio* may sustain damage from such behaviour, increasing the risk of complications from the wound conditions.

2.05 INTENDED PATIENT GROUPS

At implantation of the *M.scio*, patients must weigh more than 10 kg. Apart from that, there are no patient group limitations as to the use of the *M.scio* System.

2.06 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical specialists, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product. These can be acquired, e. g., in training offers organised by the Christoph Miethke GmbH & Co. KG (see chap. 4.00 Medical device consultants).

2.07 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions
- ▶ Reading and evaluation of intracranial pressure values
- ▶ Use of the *M.scio* "dome" variant's pump and puncture function

2.08 FUNCTIONAL PRINCIPLE

The *M.scio* is implanted to measure the pressure as well as dynamic pressure changes within the cerebrospinal fluid. The intracranial pressure can be determined through connection to a *Ventricular Catheter*.

In addition, the *M.scio* can be integrated into a shunt in order to determine the intracranial pressure in the shunt and to conduct a shunt function diagnosis. The *M.scio* is an addition to the shunt without impacting the drainage function.

The pressure is measured using a measuring cell located inside the *M.scio*. The measured values can be read and displayed telemetrically and non-invasively using the *Reader Unit Set*. The *M.scio* does not have a battery for this purpose. Its energy is supplied telemetrically and thus cordlessly from outside of the body through the *Reader Unit Set*. For a pressure measurement, the antenna of the *Reader Unit Set* must be positioned about 10 to 30 mm away from the *M.scio* (Fig. 1).

Fig. 1: Ideally, the working distance for telemetric data transmission, i.e. the distance between the *Reader Unit Set* antenna and the *M.scio*, is 10 to 30 mm.

The measured data are automatically stored by the *Reader Unit Set* on the SD card associated with the *M.scio*, allowing them to be evaluated at a later time.

2.09 SYSTEM COMPONENTS

Reader Unit Set

The measured data of the *M.scio* may only be read out using the *Reader Unit Set* (FV905X / FV907X). The operation of the *Reader Unit Set* is described in the associated instructions for use.

M.scio and SD card

The scope of delivery of the *M.scio* includes an SD card, on which all individual information of the *M.scio* (ID and calibration data) are stored during production. To carry out measurements, the SD card is inserted into the card slot of the *Reader Unit Set*. When a measurement is started, the *Reader Unit Set* compares the IDs stored on the SD card and the *M.scio* to ensure that measured values are exclusively stored on the SD card associated with the *M.scio*.
If the SD card is lost, it can be reordered by specifying the *M.scio* serial number or the associated ID. The ID can be read with [Single measurement] from the *M.scio* and shown on the *Reader Unit Set* display. The use of a standard SD card is not possible.

Combination with shunt components

The *M.scio* System can be safely combined with our range of implantable shunt components. We recommend using products from Christoph Miethke GmbH & Co. KG in combination with the *M.scio*.
The following products are particularly useful for intracranial pressure measurement in combination with the *M.scio*:

Product name	Item no.
<i>Ventricular Catheter</i> (with stylet, length 250 mm)	FV077P
<i>Ventricular Catheter</i> (with stylet, length 180 mm) with deflector (small, diameter 13 mm)	FV076P
<i>Ventricular Catheter</i> (with stylet, length 250 mm) with deflector (large, diameter 16 mm)	FV078P

Product name	Item no.
<i>Prechamber</i> (small, diameter 14 mm)	FV035T
<i>Prechamber</i> (large, diameter 20 mm)	FV033T
<i>Pediatric CONTROL RESERVOIR</i> (small, diameter 14 mm)	FV066T
<i>CONTROL RESERVOIR</i> (large, diameter 20 mm)	FV047T
<i>Burrhole Reservoir</i> (small, diameter 14 mm)	FV039T
<i>Burrhole Reservoir</i> (large, diameter 20 mm)	FV028T
<i>Titanium Shutting Plug</i>	FV024T

For pressure measuring within a shunt system, the *M.scio* can be safely combined with further components, such as valves and *Peritoneal Catheters*.
If a (*Pediatric*) *CONTROL RESERVOIR* is used, care must be taken not to place this product between the ventricle and the *M.scio* when connected to shunt components. Otherwise, the dynamics of the pressure signal may be distorted. For that reason, the combination of the *M.scio* with a *SPRUNG RESERVOIR* is ruled out. The *M.scio* must be placed between the ventricle and the valve in order to determine the intracranial pressure.

2.10 BRIEF REPORT ON SAFETY AND CLINICAL PERFORMANCE

The brief report on safety and clinical performance can be downloaded from the following address:
<https://www.miethke.com/downloads/>

2.11 SUPPLEMENTARY PRODUCT INFORMATION

In accordance with EN 45502-1

- The *M.scio* identification number (ID) can be read with [Single measurement] from the *Reader Unit Set* display, making it possible to identify the implant. The assignment of the ID to the serial number of the *M.scio* can be found on the SD card label that forms part of each *M.scio* scope of delivery.
- Approval for marking with the CE mark for active implantable medical devices (in accordance with Directive 90/385/EEC) was first granted in 2011.

In accordance with ISO 7197

- Nuclear magnetic resonance (MRI) examinations up to a field strength of 3 tesla or computed tomography (CT) examinations can be performed without risk or impairment to the function of the *M.scio*. The *M.scio* is MRI Conditional. Artefacts may occur in MRI scans. The supplied catheters are MR Safe. The documents on MRI safety are available from this website:
<https://miethke.com/downloads/>
- The *M.scio* as well as the entire shunt system are safely able to resist positive and negative pressures up to 75 mmHg (100 cmH₂O) during and after the procedure.

3.00 SYSTEM COMPONENTS OF THE *M.scio* AND SD CARD

3.01 PRODUCT DESCRIPTION

3.01.01 *M.scio* VARIANTS

The *M.scio* is available in four variants with and without a preconnected 60 cm distal catheter.


Fig. 2: *M.scio* with integrated measuring cell (X) in the following variants:

- a) dome-angled
- b) flat-angled
- c) dome-inline
- d) flat-inline

3.01.02 SCOPE OF DELIVERY

Box content	Number
Sterile packaging with <i>M.scio</i>	1
Pierce protection within the sterile bags (for <i>M.scio</i> angled)	1
SD card in protective cover	1
Instructions for use of the <i>M.scio</i>	1
Patient Implant Card	1
Package insert with information regarding the Patient Implant Card	1

3.01.03 STERILITY

The products are sterilised with ethylene oxide under strictly controlled conditions. The double wrapping in sterile bags ensures sterility for a three-year period. The expiry date is printed on the package.

If the sterile packaging is damaged, the products must not be used under any circumstances.

3.01.04 CALIBRATION

M.scio are calibrated during production. The associated calibration data are stored within the measuring cell as well as on the SD-card included in the scope of delivery. When performing a measurement with the *Reader Unit Set*, the data are automatically read and processed. A calibration of the product by the user is not required.

3.01.05 OPERATING CONDITIONS

Operating conditions of the <i>M.scio</i>	
Body temperature	35 °C ... 39 °C, function checks also possible at room temperature with reduced precision
Atmospheric pressure	800 to 1100 hPa

3.01.06 PRODUCT LIFE

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. The anticipated service life of the *M.scio* after implantation as well as of the SD card after first use is 5 years.

However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons. If the measuring cell fails, the *M.scio* “dome” variants continue to function as a conventional *Burrhole Reservoir* or a conventional *Prechamber* without any limitations, while the “flat” variants function as a connector or deflector. The integrated measuring cell poses no risk when the *M.scio* remains inside the patient.

3.01.07 DISPOSABLE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the *M.scio*. No guarantee can be assumed for the functional safety of resterilised products.

3.01.08 PRODUCT CONFORMITY

- The product meets current regulatory requirements as amended:
- ▶ (EU) 2017/745 (MDR)
 - ▶ EN 45502-1
 - ▶ ANSI/AAMI NS28
 - ▶ EN ISO 7197

3.02 IMPORTANT SAFETY INFORMATION

3.02.01 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ **If the sterile packaging is damaged, the products must not be used under any circumstances.**
- ▶ **The product must not be implanted after its expiry date.**
- ▶ **Products that have already been used must not be reused on the same or another patient.**

**CAUTION**

- ▶ **Silicone is very electrostatic. Please ensure that the catheter does not come into contact with dry cloths, talcum powder or rough surfaces. Clinging particles can lead to tissue reactions.**
- ▶ **The functionality of the *M.scio* must be tested prior to implantation.**
- ▶ **If there are any doubts as to the precision or general appropriateness of the measured intracranial pressure reading, e.g. because of the absence of detectable pulse-synchronous rhythmic oscillations, the *Ventricular Catheter* may be occluded or the air inside could cause dampening of the pressure dynamics. There may also be an impairment to the electronic function of measuring the ICP. When in doubt, for verification of the obtained pressure readings in these and other cases, other procedures should be employed, such as transcutaneous reference pressure measurement by puncture, and the patient's overall clinical picture should be assessed.**
- ▶ **During therapeutic ultrasound procedures, there is a risk of unintentionally increasing the concentration of the ultrasound field and thus the risk of injury to the patient.**
- ▶ **Violent external shocks (accident, fall, etc.) may put the integrity of the *M.scio* and the shunt system at risk.**
- ▶ **Frequent pumping of the "dome" variants can result in excessive drainage, thus leading to pressure conditions outside of the normal physiological range. The patient should be informed about this risk.**

**NOTICE**

- ▶ **The following situations can cause the measuring cell of the *M.scio* to fail:**
 - ▶ **Use of energy-emitting devices such as defibrillators**
 - ▶ **Passing of an electrical current from an external source through the body (monopolar forceps)**
 - ▶ **Use of radiation therapy and radionuclide patient imaging procedures**

3.02.02 COMPLICATIONS AND RESIDUAL RISKS

The following complications can occur in conjunction with the *M.scio* system:

- ▶ Headache, dizzy spells, mental confusion, vomiting in cases of possible leakage from the *M.scio* / shunt and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid
- ▶ Wound healing disorders because of the height of the dome-angled *M.scio*

A physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the *M.scio*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation
- ▶ Local shunt irritation/allergic reaction

3.02.03 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible authority of the EU member state in which you work.

In Germany, the relevant authority is the Federal Institute for Drugs and Medical Devices (BfArM). For current contact information, please visit the BfArM website: <https://www.bfarm.de>

3.02.04 PATIENT INFORMATION

The attending doctor is responsible for informing the patient and/or his/her proxy in advance. This includes a comprehensive description of the course of the operation, the surgical technique and the medical devices to be used. With regard to implanted medical devices, the patient must be informed about

- ▶ Warnings, precautionary measures to be taken; restrictions on use in relation to the product that ensure safe use of the medical device; contraindications (particularly Section 3.02 Important safety information)
- ▶ General information on quality and quantity of the materials and substances that the patient may be exposed to (Section 3.08.02 Materials that come into contact with bodily tissue/fluids)
- ▶ Anticipated service life of the medical device together with all necessary follow-up measures (Section 3.01.06 Product life).

3.03 TRANSPORT AND STORAGE

3.03.01 TRANSPORT

Transport conditions

Ambient temperature	-30 °C ... +50 °C
Atmospheric pressure	596 hPa ... 1100 hPa

3.03.02 STORAGE

The medical devices must always be stored in a clean and dry place.

Storage conditions

Temperature range for storage	0 °C ... +50 °C
Pressure range for storage	800 hPa ... 1100 hPa

3.04 USE OF THE PRODUCT

3.04.01 INTRODUCTION

The M.scio system can be used in two scenarios in order to determine the intracranial pressure:

- ▶ M.scio implanted without shunt system
- ▶ M.scio integrated into a shunt system

In both scenarios, the Reader Unit Set is used for the telemetric reading and display of the pressure values.

3.04.02 SAFETY AND WARNING NOTICES

CAUTION

- ▶ Testing the functionality of the implant prior to implantation is an essential requirement.
- ▶ Patency should be checked by aspiration of a sterile fluid from the distal end and air should be removed from the implant prior to implantation.
- ▶ If there are any doubts as to the precision or general appropriateness of the measured intracranial pressure reading, e.g. because of the absence of detectable pulse-synchronous rhythmic oscillations, the plug of the Ventricular Catheter may be occluded or the air inside could cause dampening of the pressure dynamics. There may also be an impairment to the electronic function of measuring the ICP. When in doubt, for verification of the obtained pressure readings in these and other cases, other procedures should be employed, such as transcutaneous reference pressure measurement by puncture, and the patient's overall clinical picture should be assessed.

3.04.03 REQUIRED MATERIALS

The *M.scio* is designed for safe use with the shunt components described in Section “2.09 System components”. We advise using a catheter with an inner diameter of 1.2 mm and an outer diameter of approx. 2.5 mm for connection. In any case, catheters have to be carefully fixed with a ligature to the connectors of the *M.scio*. Any kinks in the catheter must be avoided.

3.04.04 PREPARING FOR IMPLANTATION

Checking the sterile packaging

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Removing the pierce protection

The *M.scio* „angled“ variants are stored inside the sterile packaging in a transparent pierce protection. It must be removed prior to implantation of the *M.scio*.

Checking for proper functionality of the *M.scio*

The functionality and telemetric readability of the *M.scio* must be verified prior to implantation. To do so, the SD-card associated with the *M.scio* must be inserted into the SD-card slot and the antenna of the *Reader Unit Set* must be placed close to the sterile-packed *M.scio*. The display of readings after the start of measurement (e.g., [Fast measurement] – see the *Reader Unit Set* instructions for use) confirms the readability of the *M.scio*. Readings taken at room temperature do not permit conclusions to be drawn about accuracy at body temperature.

Testing the patency and removal of air from the implant

Prior to implantation, the *M.scio* should be checked for patency, and the air should be removed from the implant. The most careful way of filling the *M.scio* is by aspiration through an attached sterile single-use syringe. To that end, the *M.scio* is connected and immersed in a sterile physiological salt solution. In the “dome” variants, the removal of air bubbles can

be supported by pressing the silicone membrane in the direction of the measuring cell.

3.04.05 IMPLANTATION OF THE *M.scio* WITHOUT SHUNT SYSTEM

Implantation of the *Ventricular Catheter*

Several surgical techniques are available for placement of the *Ventricular Catheter*. The required skin incision should be made in form of a lobule pedicled towards the draining catheter. The skin incision must not be located right above the planned burr hole. For use of an “angled” *M.scio* variant, the burr hole should ideally have a diameter of 10 mm (max. 14 mm). To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole.

M.scio is available in a range of different variants: If an “angled” variant is used, the *Ventricular Catheter* is implanted first. Once the stylet has been removed, the patency of the *Ventricular Catheter* can be tested by checking if cerebrospinal fluid is dripping out.

The implantation of the “inline” variants additionally requires a burr hole deflector to orient the *Ventricular Catheter* at the required 90° angle. The deflector is used for adjusting the length of the catheter to be implanted and for its positioning inside the ventricle. As an alternative, a *Burrhole Reservoir* can be placed at the burr hole that is then connected to the *Ventricular Catheter*.

Implantation of the *M.scio*

The *M.scio* must always be placed outside of the skullcap.

If the *M.scio* is implanted for the “subarachnoid hemorrhage” indication, because of the contraindication “blood in the cerebrospinal fluid”, it must be ensured that the blood has already been removed from the cerebrospinal fluid.

Placing the *M.scio* requires a sufficiently sized pocket between skin and bone. The orientation has to be chosen in such a way as to enable the connection of a shunt system in a subsequent surgery.

Prior to implantation of a “dome” variant of the *M.scio*, the previously placed catheter is shortened. The *Ventricular Catheter* and the *M.scio* are then connected.

In the implantation of an “inline” variant, the *Ventricular Catheter* brought into 90° orien-

tation with the deflector is connected to the *M.scio*.

Alternatively, if a *Burrhole Reservoir* is used, a short catheter segment is used to for the connection between the burr hole reservoir and the *M.scio*.

The *M.scio* should not be located directly underneath the skin incision.

Distal closure of the *M.scio*

A short catheter segment can be used to connect the *Titanium Shutting Plug* to the distal part of the *M.scio*. Prior to the distal closure, it is recommended to check the patency of the connected components and to remove any air from the system.

When using sharp instruments, care should be taken to avoid cuts and scratches of the silicone catheter.

The individual connection of the implanted components should be secured with a ligature. The connection of the components must be stable, but excessive tightening of the ligature must be avoided.

The position of the *Ventricular Catheter* should be inspected after the procedure by imaging (such as CT or MRI).

3.04.06 INTEGRATION OF THE *M.scio* INTO A SHUNT SYSTEM

The only difference for integration of the *M.scio* into a shunt system from the implantation procedure described in the previous section is that additional shunt components (such as valves) are connected to the distal end instead of a plug. Here, too, care must be taken to secure the individual connections by a ligature.

3.04.07 PRESSURE MEASUREMENT AFTER IMPLANTATION

After implantation, the functionality of the implanted products should be checked once again using the *Reader Unit Set* and the associated SD card (see the *Reader Unit Set* instructions for use). To this end, a plausibility test can be conducted. Using the [Fast measurement] mode, pulse-synchronous rhythmic oscillations (pulse waves) and possibly also ventilatory fluctuations of the intracranial pressure signal can be observed when everything functions properly.


Fig. 3: Typical measurement curves on the display of the *Reader Unit Set* with pulses and ventilatory fluctuations in [Fast measurement] mode after implantation.

After implantation, the “dome” *M.scio* variants provide the additional option to observe a pressure change in the measurement curves by applying gentle pressure to the silicone membrane of this implant with the patient in a stable position.

The absence of any pulse-synchronous rhythmic variations in the intracranial pressure signal may indicate an occlusion of the *Ventricular Catheter*.

The telemetric link between the antenna of the *Reader Unit Set* and the *M.scio* may be disrupted by metal parts or operation of another *Reader Unit Set* within the vicinity of the implant. In this case, the distance to the metal parts or to the additional *Reader Unit Set* must be increased.

If the patient has a raised temperature, functional impairment of the *Reader Unit Set* may occur (see the *Reader Unit Set* instructions for use). In readout mode, the temperature within the *M.scio* may rise. An integrated temperature safety device stops the measurement automatically at 39 °C in the implant. Furthermore, measurements are impossible if the body temperature exceeds 39 °C.

If the measurement is interrupted, the process can be restarted while observing the functional distance between the *M.scio* and the antenna of the *Reader Unit Set*.

In order to ensure the comparability of the measured values from the *M.scio* with those of other pressure probes, the measured values should, for hydrostatic reasons, refer to a standard reference point (e.g. the interventricular foramen). Otherwise, particularly if the patient is in an upright position, the determined pressure values may differ by several cmH₂O or rather mmHg.

3.04.08 USE OF THE PUMP AND PUNCTURE FUNCTION

Thanks to the silicone membrane, the “dome” variants of the *M.scio* also have the pumping and puncturing capability of a standard reservoir. Any puncturing of the silicone membrane should be performed perpendicular to the reservoir surface with a cannula of max. Ø 0.9 mm. 30 punctures are possible without any restrictions. The robust titanium housing prevents any accidental damage to the measuring cell. The volume per pump cycle is approx. 0.062 ml.

3.05 PRESSURE MEASUREMENT IN EXCEPTIONAL SITUATIONS

In exceptional situations, a pressure measurement of the *M.scio* can be taken even without the associated SD card, using the *Reader Unit Set*. However, it is then impossible to store the measured data.

3.06 EXPLANATION AND DISPOSAL

3.06.01 EXPLANATION

The explanation of the *M.scio* should be performed according to the state of the art and in compliance with medical practice.

3.06.02 DISPOSAL

M.scio and shunt components

Products and product parts not used in the implantation or surgically removed must be disposed of correctly as potentially infectious material in accordance with medical practice as well as respective regional laws and regulations.

Explanted medical devices must not be reused.

SD card

Adhere to national regulations when disposing of or recycling the product, its components or its packaging!

Because of the stored data, the SD card must be disposed of in a safe and data-protection compliant manner. Care must be taken to ensure that the restoration of the information on the card is either rendered impossible or only achievable with great effort.

3.07 TROUBLESHOOTING

Potential faults occurring at the reading of the measured values from the *M.scio* together with error messages on the display of the *Reader Unit Set* and the procedures for troubleshooting are described in detail in the *Reader Unit Set* instructions for use.

For faults such as a loss of communication between the *Reader Unit Set* and the *M.scio* and for verification of the determined pressure values, recourse can be taken to a transcutaneous reference pressure measurement.

3.08 TECHNICAL INFORMATION

3.08.01 TECHNICAL DATA

Designation	Values
Pressure measurement sampling rate	max. 44 Hz in the [Fast measurement] mode
Maximum frequency response including ICP display at peak pressures of 10, 20 and 50 mmHg	+10 mmHg: 20 Hz +20 mmHg: 19 Hz +50 mmHg: 18 Hz
Maximum frequency response at peak pressures of 10, 20 and 50 mmHg	Pressure increase: +10 mmHg: 408 mmHg/s +20 mmHg: 759 mmHg/s +50 mmHg: 1777 mmHg/s Pressure decrease: -10 mmHg: -462 mmHg/s -20 mmHg: -786 mmHg/s -50 mmHg: -1922 mmHg/s
Time constant for full alignment	Pressure increase: +10 mmHg: 18 ms +20 mmHg: 18 ms +50 mmHg: 18 ms Pressure decrease: -10 mmHg: 15 ms -20 mmHg: 18 ms -50 mmHg: 19 ms
Measurement accuracy of the <i>M.scio</i> system (The information refers to a mean ambient pressure of 1013 hPa.)	Pressure range (relative to atmospheric pressure): -50 mmHg ... +100 mmHg Measurement accuracy determined over 10 days in the pressure range: -50 mmHg ... -20 mmHg: +/- 10% -20 mmHg ... +20 mmHg: +/- 2 mmHg +20 mmHg ... +100 mmHg: +/- 10%
Expected drift of the zero point value of the <i>M.scio</i> system at 37 °C (stated values each apply from the date of production of the <i>M.scio</i>)	Long-term (4 years): < 2 mmHg
Stability of pressure measurements from the <i>M.scio</i> system within the temperature range from 20 °C to 39 °C	The temperature range is covered by the calibration; the greatest accuracy is achieved in the temperature range from 35 °C to 39 °C and after completing a single temperature-related running-in period of at least 120 hours after implantation.
Telemetry	see the <i>Reader Unit Set</i> instructions for use
Dimension of the intracranial part of the <i>M.scio</i> system	Depending on the size of the <i>Ventricular Catheter</i> (VC) used: VC external diameter: 2.5 mm VC length: patient-specific
Size of the required burr hole for implantation of the <i>M.scio</i> "dome" variant	ideally 10 mm, max. 14 mm

Designation	Values
Total fluid volume within the <i>M.scio</i> and the <i>Ventricular Catheter</i> (VC)	“dome” variants: $450 \text{ mm}^3 + (1.13 \times \text{VC length}) \text{ mm}^3$ “flat” variants: $250 \text{ mm}^3 + (1.13 \times \text{VC length}) \text{ mm}^3$ VC length: patient-specific
Weight of <i>M.scio</i> variants:	
dome-angled	3.5 g
dome-angled, with distal catheter	6.0 g
flat-angled	3.2 g
flat-angled, with distal catheter	5.6 g
dome-inline	3.5 g
dome-inline, with distal catheter	6.0 g
flat-inline	3.2 g
flat-inline, with distal catheter	5.6 g

3.08.02 MATERIALS THAT COME INTO CONTACT WITH BODILY TISSUE/FLUIDS

- ▶ Measuring cell housing: titanium alloy
- ▶ *M.scio* housing: plastic (PEEK), for “dome” variants also silicone
- ▶ Connectors: titanium alloy
- ▶ Catheter (optional): silicone

Quantitative material properties

Material	<i>M.scio</i> dome	<i>M.scio</i> flat
PEEK:		
Volume [mm ³]	1145	1049
Weight [g]	1.50	1.37
Contact surface [mm ²]	2190	2400
Titanium alloy:		
Volume [mm ³]	2003	2003
Weight [g]	3.04	3.04
Contact surface [mm ²]	745	745
Silicone (membrane):		
Volume [mm ³]	198	---
Weight [g]	0.20	---
Contact surface [mm ²]	248	---
Silicone (catheter 60 cm optional):		
Volume [mm ³]	2400	2400
Weight [g]	2.40	2.40
Contact surface [mm ²]	7000	7000

3.09 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxx indicates the identifier of the responsible notified body.
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Serial number
	Unique Device Identifier
UDI-DI	UDI-DI number
	Sterilised using ethylene oxide
	Do not resterilise
	Do not reuse
	Double sterile barrier system
	Do not use if package is damaged and consult instructions for use.
	Keep dry
	Temperature limit
	Atmospheric pressure limitation
	Consult instructions for use / electronic instructions for use

Symbol	Explanation
	Caution
	Non-pyrogenic
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians.
	MR Conditional
	Electric and electronic waste
	Patient identification
	Date of implantation
	Health care centre or doctor
	Patient information website
	Model number / European Medical Device Nomenclature Code

4.00 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions. You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany