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TABLE OF CONTENTS

READER UNIT FOR SENSOR RESERVOIR	31
INDICATIONS	33
INTENDED PURPOSE	33
DEVICE DESCRIPTION	33
SAFE HANDLING	33
INTENDED USE	33
CONTRAINDICATIONS	33
OPERATING PRINCIPLE	34
SYMBOLS ON THE READER UNIT DISPLAY	34
USING THE READER UNIT	35
SYSTEM START-UP	36
FUNCTION CHECKS	36
OPERATION OF THE READER UNIT	37
GENERAL OPERATING INSTRUCTIONS	37
INFO MENU	38
TAKING MEASUREMENTS	38
DATA MANAGEMENT/DISPLAY OF MEASURED DATA	40
SETTINGS	44
CLEANING/DISINFECTION	48
MANUAL CLEANING/DISINFECTION	48
CALIBRATION OF THE BARD SENSOR	48
TECHNICAL SERVICE	49
TROUBLESHOOTING	49
MALFUNCTIONS DISPLAYED WITH AN ERROR MESSAGE	50
OTHER MALFUNCTIONS/POSSIBLE FAULTS	53
TECHNICAL DATA	54
AMBIENT CONDITIONS	55
CE MARKING	55
STORAGE AND TRANSPORT	55
DISPOSAL	55

Reader Unit for SENSOR RESERVOIR

- 1 ON/OFF switch
- Display screen
 Antenna socket
 Antenna
 SD card slot

- 6 Function keys
- Particular Reys
 Connectorsocket for power supply unit
 Plug
 Control light
 Power supply unit

Symbol	Explanation
\triangle	Caution, general warning sign Caution, read the documentation supplied with the product
(LE	Follow the instructions for use
((<u>``</u>))	Non-ionising radiation
	Labelling of electric and electronic devices in accordance with Directive 2002/96/EG (WEEE), see chapter on "Disposal"
	ON/OFF switch
	Function key: The 4 function keys are used to perform the options shown on the display screen
IP44	Reader unit 7502 0000
IP44	Antenna 7503 0000
IP40	Power supply unit 7504 0000
	Protection class II (protective insulation)
¥	Antenna symbol Antenna socket
	SD card SD card slot
DCIN	DCIN connector for power supply unit

INDICATIONS

- Functional test of the shunt system implanted into the patient
- Technical functional diagnosis of the shunt system, e.g. detection and localisation of occlusions

INTENDED PURPOSE

The Reader Unit set for SENSOR RESERVOIR is used for the pressure-based functional control of shunt systems.

DEVICE DESCRIPTION

Description	Article number
Reader Unit Set for	75100000
SENSOR RESERVOIR	(AAG ref. FV905X)

SAFE HANDLING

Risk of injury due to incorrect operation of the product!

Appropriate product training prior to using the product is recommended. Please contact Christoph Miethke GmbH & Co. KG for information about product training.

- Clean factory-fresh product thoroughly after removing transport packaging (see chapter on "Manual Cleaning/Disinfection").
- Ensure that the product is in proper condition and fully functional prior to use (see chapter on "Inspection").
- To prevent nosocomial infection and multiple resistances, the device should be disinfected after each use. Disinfection is preferably carried out by wiping with the disinfectant described in the section on "Cleaning/Disinfection".
- To prevent damage as the result of improper assembly or operation and to avoid compromising warranty and liability, please note the following:
- Use the product only as specified in these instructions for use.
- Observe safety information and maintenance instructions.

- Use only products by Christoph Miethke GmbH & Co. KG in combination, see chapter on "Reader Unit System Components".
- Device and accessories must only be operated and used by persons with the requisite training, knowledge or experience.
- Keep the instructions for use accessible for the user.

INTENDED USE

In conjunction with the SENSOR RESERVOIR / SENSOR PRECHAMBER, the Reader Unit Set for the SENSOR RESERVOIR / SENSOR PRECHAMBER (henceforth: the reader unit) is used for the functional diagnostic of shunts and for the reliable non-invasive detection of shunt systems (without intervention or radiological diagnostics) as well as for the identification of mechanical damage to shunt valves in the treatment of hydrocephalus.

CONTRAINDICATIONS

- · Intracranial pressure measurement
- Intracranial pressure measurement without a shunt system
- Intracranial pressure management based on measured values
- Functional diagnosis of extracorporeal fluid drainage system (traumatic brain injury)

Fig. 1: The working distance for telemetric data transmission, i.e. the distance between the antenna and the SENSOR RESERVOIR / SENSOR PRECHAMBER is 10-30 mm

OPERATING PRINCIPLE

The telemetric system serves to measure CSF pressure readings whose relative, time-dependent changes and characteristics may allow conclusions as to the function of an implanted ventricular drainage system. Diagnostic methods can be used to initiate a pressure change in the shunt system, which can then be interpreted. With the reader unit set, measured data from the measuring cell in the SENSOR RESERVOIR / SENSOR PRECHAMBER can be read and displayed. Measured data are automatically stored on the SD card and can thus be evaluated at a later stage. The relative behaviour of the CSF pressure in the shunt system can provide information about its function. This permits the non-invasive detection and localisation of an occlusion within the shunt system and also of a mechanical loss of function

SYMBOLS ON THE READER UNIT DISPLAY

Info	For displaying additional information or starting with the information menu
Menü	For starting with the selection menu
●]])	Fast measurement - Fast measurement can be started without a further selection
	Upward navigation of cursor
	Downward navigation of cursor
	For quitting the current menu
ОК	Activates or confirms the selected function
Start	Starts measurement
Stop	Stops the current measurement
6	For deleting a saved measurement in the <data management=""> menu</data>
*	For placing a marker during continuous measurement
*~>	For opening the diagram processing menu
	For activating the zoom function
♦ 	For activating the cursor function
	For adjusting the display of min. and max. values of the pressure axis

∽⊨	For returning to the diagram processing menu
⊨ ••	For returning to the diagram processing menu
	Navigate forward
-	Navigate backward
on	On
off	Off
\checkmark	Confirm
*	For maximizing the display of a sequence in the <data management=""> menu</data>
	For minimizing the display of a sequence in the <data management=""> menu</data>
×	For deactivating the acoustic signal

USING THE READER UNIT

Caution: Limited accuracy and no storage of measured data for measurements without SD card!

• Insert the SD card designated for use with the SENSOR RESERVOIR / SENSOR PRE-CHAMBER into the reader unit. The implant (see patient data card) is associated with the respective SD card via the identification number (ID). Measured data can only be save on the designated SD card.

Caution: Unexpected termination of measurement process if SD card is full!

 Check memory status prior to starting measurement and delete data no longer required.

Caution: Danger from use of MRI! Only operate reader unit outside of MRI field of use!

Caution: The reader unit cannot be operated with low battery status!

 Check battery status. Connect to power supply if battery status is low. Caution: No or bad communication! The telemetric link between the reader unit and the implant may be disrupted if a further reader unit is used.

• Increase the distance between reader units.

Caution: The telemetric link between the antenna and the implant may be disrupted by a further implant within the vicinity of the antenna.

Increase the distance between implants.

Caution: The telemetric link between the antenna and the implant may be disrupted by metal components within the vicinity of the implant.

• Increase the distance from metal components!

Note: If the temperature of the SENSOR RE-SERVOIR / SENSOR PRECHAMBER falls below 20°C, use procedure for "Measurement without SD card". The device will be operated with limited accuracy and measured data cannot be stored. The following error message can occur during measurement with SD card:

- During continuous measurement: pressure values out of range
- During a single or fast measurement: temperature values out of range

For a more detailed description of errors, see chapter on "Troubleshooting".

SYSTEM START-UP

Caution: Allow reader unit to acclimatize for approx. three hours at room temperature.

The reader unit is equipped with a battery which on a 100% charge permits independent operation without external power supply for up to five hours. To that end, the battery needs to be charged when first using the device. The power supply unit is used to charge the reader unit battery. It takes approx. 6 hours for the battery to be completely charged.

The reader unit can also be operated during the charging process (with connected power supply unit).

Note: Charging is not possible at an ambient temperature of more than 35°C!

Connecting the power supply

Caution: Incorrect accessories can cause damage to the reader unit and pose a safety risk to user and patient!

· Only use the original power supply unit for the reader unit.

Mains voltage must correspond to the voltage range stated on the label of the reader unit power supply.

- Insert the connector of the power supply unit into the socket of the reader unit, see Fig. 2.
- Plug power supply unit into a wall power socket.

Fig. 2

Connecting Antenna

Connect the antenna to the reader unit, see Fig. 3

Inserting the SD card

Caution: Improper use can damage the SD card!

Do not touch the contact strips on the SD card

Insert the SD card designated for use with the respective SENSOR RESERVOIR / SENSOR PRECHAMBER (see patient data card) into the reader unit until it locks into place (see Fig. 4).

Note: Press SD card slightly to remove

Fig. 4

FUNCTION CHECKS

- Check that the entire unit is in good condition and fully functional prior to each use.
- To check battery status, perform a function check without power supply.
- Ensure that neither the device nor the antenna used for the function test show any visible signs of damage.
- Check the function of the following components in the given order:
 - 1. Activate the On/Off switch 1.
 - 2. Automatic self-test after activation inclu-
ding display and speaker volume

 Check battery status; if necessary, plug in power supply unit and charge battery.

The following is displayed on the screen: Fade-in "Selftest ..." Fade-in "booting ..."

Fig. 5a

The following is displayed on the screen of the reader unit during operation with the mains power supply unit:

Fig. 5b

Note: The control lamp of the reader unit power supply must be on!

- Check system time, correct if necessary, see chapter on "Settings".
- Activate the On/Off switch ①.

OPERATION OF THE READER UNIT

Key

- A Date
- B Time
- C SD card D Assignment of memory space (on SD card)
- E Battery status
- F Info/OK key
- G Downward navigation
- H Menu key

Activate reader unit with On/Off switch O Start screen content is displayed (see Fig. 7).

GENERAL OPERATING INSTRUCTIONS

The menu-controlled operation of the device uses the four function keys. The context-related function of these soft keys is indicated by the symbols shown on the display for each of these keys. The exact meaning of these symbols is explained in the list of "Reader unit display symbols" on p. 34/35.

The preselected submenu is highlighted by a frame. For improved navigation, the currently selected position in the submenu is highlighted by a scroll bar in addition to the frame. The keys

for standard operation

Note: When not in use, the reader unit switches to stand-by mode within 1-5 minutes, depending on setting.

INFO MENU

The following information can be called up from the info menus:

- · Hotline & Service
- · System Details
- · Manufacturer Info
- · Language SD-Card

Press Info to open the info menu.

The following information is displayed on the screen:

Fig. 7

Preselect the desired submenu with the keys and , and confirm with OK. Return to the previous menu by pressing .

Submenu content

Hotline & Service

- Hotline & Service
- · Tel: +49 331 620 83-0

System Details

- Product name: Reader Unit for Sensor Reservoir
- · Article number: 7502 0000
- Serial number: XXXXX
- · Software version: 2.0X
- · Service date: dd/mm/yyyy

Manufacturer Info

Manufacturer: Christoph Miethke GmbH & Co. KG Ulanenweg 2 | 14469 Potsdam Germany

TAKING MEASUREMENTS

Press Menü to open the menu.

The following information is displayed on the screen:

Fig. 8

Press **OK** to open the submenu <measurement>.

The following information is displayed on the screen:

Fig. 9

There are three types of measurement:

1. Single measurement: The mean of 8-10 measurements is displayed and saved to the SD card.

 Continuous measurement: The sequence of single measurements within the selected time interval is shown as a time diagram and save to the SD card.

3. Fast measurement: The sequence of unaveraged single measurements at the maximum available measurement rate (approx. 40 measurements per second) is shown as a time diagram and save to the SD card.

Single measurement

Press **OK** to select the <single measurement> menu.

The following information is displayed on the screen:

Fig. 10

Press Start to start the single measurement. Press to return to the selection menu.

The following information is displayed on the screen for a single measurement:

Fig. 11

Continuous measurement

In the <measurement> menu, press and to select the <continuous measurement> menu and confirm with OK .

The following information is displayed on the screen:

Fig. 12

Press Start to start the continuous measurement. Press to return to the selection menu.

The following information is displayed on the screen: Symbols for reception quality

Press to place a marker. Multiple markers can be placed during a measurement process. Press Stop to stop the measurement.

Note: Markers permit the circumstantial analysis of the measured data

Symbols for reception quality

Symbol	Explanation
	Communication started
*	Distance from antenna to measuring cell: - good
\leftrightarrow	Distance from antenna to measuring cell too close: - increase distance
→←	Distance from antenna to measuring cell too large: - decrease distance

Fast measurement

There are two ways to start a fast measurement:

On the start screen, press ())) and then start measurement by pressing **Start**, or press **Menü** to open the following submenu:

Fig. 14

Press **OK** to open the menu <measurement>.

In the <measurement> menu, press and and to select the <fast measurement> menu The following is displayed on the screen:

Abb. 15

and confirm with OK

The following is displayed on the screen:

Fig. 16

During fast measurement, the following is displayed on the screen:

Symbols for reception quality

Press Stop to stop the measurement.

For an explanation of the symbols for reception quality: see chapter on "Continuous measurement".

DATA MANAGEMENT/DISPLAY OF MEASURED DATA

On the start screen, press Menü to open the following submenu:

Fig. 18

Press and to select the <data management> menu and confirm with **OK**.

Fig. 19

Measured data are listed chronologically (from the start of measurement) and shown as follows:

Fig. 20

Press and to select the desired file and confirm with **OK**.

Measured data files are labelled as follows:

Symbol	Explanation
Ť	Single measurement
	Continuous measurement
*****	Fast measurement

Note: An hour glass is shown when loading a continuous or fast measurement and the loading time displayed in seconds.

Single measurement

Fig. 22

Deletion has to be confirmed by pressing . Press to return to the previous view. Press again to return to the selection menu.

Continuous measurement

Fig. 23

Press to delete the data file. Deletion has to be confirmed by pressing . Press to return to the previous view. Press again to return to the selection menu. Press to open the <diagram processing> menu.

The following information is displayed on the screen:

Press to open the <zoom> menu. The following information is displayed on the screen:

Fig. 25

Press or to display an expanded or compressed time diagram (zoom) of the measured sequence. Press to return to the <diagram processing> menu. Press to return to the main <continuous measurement> menu.

Press in the <diagram processing> to open the following screen:

Fig. 26

Press return to the <diagram processing> menu. Press return to the <diagram proto move the time axis to the left or right. Press to return to the main <continuous measurement> menu.

Press **Press** in the <diagram processing> to open the following screen:

Fig. 27

Cursor

Use the and keys for manual setting of maximum or minimum scale values. Make sure that the "automatic" function is switched off. If the "automatic" function is selected, the setting of maximum or minimum values is irrelevant. In this function, axis scaling is automatic according to the measured data. Press to toggle cursor between maximum, minimum and automatic.

Press to return to the display of measured values.

Fast measurement

Fig. 28

Press to delete the measurement. Deletion has to be confirmed by pressing Press to return to the previous view. Press to return to the selection menu. Press to open the <diagram processing> menu.

The following information is displayed on the screen:

Fig. 29

Press to open the <zoom> menu.

Fig. 30

Press or to display an expanded or compressed time diagram (zoom) of the measured sequence. Press to return to the <diagram processing> menu. Press to return to the main <fast measurement> menu.

Note: Switching between displays can take a few seconds, depending on file size.

Press in the <diagram processing> to open the following screen:

Fig. 31

Press Press to return to the <diagram processing> menu. Press or to move the time axis to the left or right. Press to return to the main <fast measurement> menu. Press in the <diagram processing> to open the following screen:

Use the and keys for manual setting of maximum or minimum scale values. Make sure that the "automatic" function is switched off. If the "automatic" function is selected, the setting of maximum or minimum values is irrelevant. In this function, axis scaling is automatic according to the measured data. Press

to toggle cursor between maximum, minimum and automatic.

Press to return to the display of measured values.

SETTINGS

On the start screen, press Menü to open the following submenu:

Fig. 33

Press and to select the <settings> menu.

Fig. 34

Press **OK** to open the menu <settings>. The following is displayed on the screen:

Fig. 35 44

Date/Time

Press **OK** to open the menu <date/time>. The following information is displayed on the screen:

Press to change the cursor position. Use the and keys to change the values at the cursor position. Press to return to the display of measured values.

Note: Modified values are saved immediately.

Units

Press and to select the <units> menu.

Fig. 37

Press **OK** to open the menu <units>. The following is displayed on the screen:

Fig. 38

Press and to select the desired unit and confirm with **OK**. Press **H** to return to the <settings> menu.

Measurement interval

These settings function only in the continuous measurement mode.

In the <settings> menu, press and and to select the <measurement interval> menu.

Fig. 39

Press **OK** to open the menu <measurement interval>. The following information is displayed on the screen:

Fig. 40

Press and to select the desired measurement interval and confirm with **OK**. Available settings range from 1-300 s. The available memory space on the SD card is also shown. Press **C** to return to the <settings> menu.

Measurement display

In the <settings> menu, press and to select the <display data> menu.

Fig. 41

Press **OK** to open the menu <display value>. The following is displayed on the screen:

Fig. 42

Press or to select <display value>. Press to return to the <settings> menu.

Language

Fia. 43

Press **OK** to open the menu <language>. The following information is displayed on the screen.

language and confirm with OK. Press to return to the <settings> menu.

Volume

Fig. 45

Press **OK** to open the menu <volume>. The following information is displayed on the screen:

Press and to select the desired volume and confirm with **OK**. Available settings are 1-5. When setting the volume, a corresponding acoustic signal is emitted. Press to return to the <settings> menu.

Note: The setting determines the volume of acoustic alerts.

Exception: Error alerts are generally emitted at volume level 5.

Standby mode

In the <settings> menu, press and to select the <standby-mode> menu.

Fig. 47

Press **OK** to open the menu <standbymode>. The following information is displayed on the screen:

Fig. 48

Press and to select the desired time after which the device switches to standby mode and confirm with **OK**. Available settings are 1-5 minutes. Press **I C** to return to the <settings> menu.

Brightness

In the <settings> menu, press and to select the
brightness> menu.

Fig. 49

Press **OK** to open the menu
-brightness>. The following information is displayed on the screen:

Press and to select the desired background brightness and confirm with **OK**. Available settings are 000-100 (in steps of 5). Press **Here** to return to the <settings> menu.

Contrast

In the <settings> menu, press and to select the <contrast> menu.

Fig. 51

Press **OK** to open the menu <contrast>. The following information is displayed on the screen:

Fig. 52

Press and to select the desired contrast and confirm with **ok**. Available settings are 000-100 (in steps of 5). Press to return to the <settings> menu.

CLEANING/DISINFECTION

Caution: Risk of electric shock and fire!

- · Disconnect from mains before cleaning.
- Do not use flammable or explosive detergents or disinfectants.
- Ensure that no fluids can infiltrate into the product.

Caution: Mechanical cleaning or disinfection as well as inappropriate detergents or disinfectants can damage or destroy the product!

- · Only disinfect the product manually
- · Never sterilize the product

Use detergents/disinfectants approved for surface cleaning only as specified by the manufacturer. Observe specifications regarding concentration, temperature and contact time.

MANUAL CLEANING/DISINFECTION

Wipe disinfection for electrical devices without sterilization

Phase	1
Step	Wipe disinfection
T (°C/°F)	RT (room temperature)
t (min)	≥1
Conc. (%)	-
Water quality	-
Chemical	Meliseptol HBV Tissues 50% propan-1-ol

Phase I

- Remove any visible residue with a disposable disinfectant wipe.
- Wipe all surfaces of the visually clean product with a fresh disposable disinfectant wipe.
- Observe the specified contact time (1 minute minimum).

Inspection

- Check product for damages after each cleaning/disinfection.
- · Immediately set aside product if damaged.

Storage

 Store reprocessed products protected from dust in a dry, dark and cool room with as little bacterial contamination as possible.

CALIBRATION OF THE BARD SENSOR

The reader unit contains a barometric pressure sensor (BARD sensor). An annual calibration is required to ensure that predefined tolerance limits are maintained.

Note: As part of the calibration, a comprehensive functional and safety test of the device is also performed. Omitting to perform the annual calibration can result in a drift of the BARD sensor outside of the tolerance limits. To this end, the device has to be sent to the technical service section on an annual basis. The due date for the next calibration can be found in the menu under System Details > Service date.

TECHNICAL SERVICE

Caution: Risk on injury and/or malfunction! Do not modify the product.

For service and maintenance, please contact Christoph Miethke GmbH & Co. KG.

Modifications of medical devices can result in a loss of guarantee or warranty claims. Christoph Miethke GmbH & Co. KG only accepts responsibility for the safety, reliability and performance of the device if

- the device is used in accordance with the instructions for use.
- new settings, modifications and repairs have only been carried out by persons authorized by us.
- the electrical installation of the respective room conforms to national standards (IEC definition).

For service, maintenance, repair and disposal, please contact Christoph Miethke GmbH & Co. KG.

Technical service:

Christoph Miethke GmbH & Co. KG Ulanenweg 2 14469 Potsdam, Germany Tel.: +49 331 62083-0 Fax: +49 331 62083-40 E-mail: info@miethke.com

TROUBLESHOOTING

If an error occurs, it is shown on the screen.

Example of an error message:

Error message

Fig. 53

MALFUNCTIONS DISPLAYED WITH AN ERROR MESSAGE

Message displayed	Cause	Fault detection/rectification
Battery flat - auto off	Battery capacity exhausted (0 %)	All data are saved after 2 min. Reader unit switches off automatically. Connect to original power supply.
Battery voltage incompatible - use original power supply	Battery voltage of the reader unit is too low	Reader unit switches off automatically after 20 seconds. Connect to original power supply.
Low Battery Voltage	Battery voltage is too low	Backlight switches off after 3 seconds. Con- nect to power supply. No interruption to on- going measurements.
Flashing of crossed-out antenna symbol Info : Antenna error	Antenna faulty	Deactivate and reactivate device. If the error persists, contact Technical Service
Flashing of crossed-out antenna symbol	Antenna not connected at the start of measurement	Connect antenna. Measurement starts again.
Info : Antenna not connected	- or - Antenna was disconnected duri- ng measurement	-or- Connect antenna. Measurement continues.
Flashing of crossed-out antenna symbol Info: No communication	Data recording interrupted du- ring continuous measurement (telemetric link disrupted)	Measurement resumes automatically once communication is restored.
Flashing of crossed-out anten- na symbol Info: SD card has been removed. Measurement possible	SD card was removed during continuous measurement	Insert SD card. Start measurement again.
Data error! Info: File opening not possible	Validation of file failed	File cannot be opened. Try again.
Continuous key activation Keypad error	Persistent activation of a key > 60 seconds	Release the key.
Pressure values out of range	Pressure data measured from implant not plausible - physiologically meaningless data	Measurement continues. If the error persists, contact Technical Service.
Input voltage incompatible.	Power supply voltage is too high	Reader unit switches off automatically after 20 seconds. Use original power supply.

Message displayed	Cause	Fault detection/rectification
Wrong SD card inserted! Remove SD card! Info : Measurement without data sto- rage possible - or - Insert SD card with correct ID	Measurement starts without SD card. During measurement, an SD card not designated for use with that implant is used during measurement	Insert the SD card designated for use with the SENSOR RESERVOIR / SENSOR PRECHAMBER into the reader unit. The implant (see patient card) is associated with the respective SD card via the identifica- tion number (ID).
Wrong SD card inserted! Remove SD card! Info: Measurement without data sto- rage possible - or - Insert SD card with correct ID XXXXXXXXXX!	SD card with different ID to the implant has been used	Insert the SD card designated for use with the SENSOR RESERVOIR / SENSOR PRECHAMBER into the reader unit. The implant (see patient card is associated with the respective SD card via the identifica- tion number (ID).
Wrong implant - restart measurement! Info Switching to another implant during continuous measure- ment not possible!	After the start of a continuous measurement, data from another implant have been received	Increase the distance between the two implants.
Unit-temperature incorrect Info: Unit temperature from 10°C to 40°C valid!	Temperature in the reader unit exceeds calibrated interval	The reader unit can only be used at device temperatures from 10°C to 40°C. Current measurement is stopped.
Internal voltage incompatible	Internal voltage of device too high/ too low	Reader unit switches off automatically after 20 seconds. Contact Technical Service.
SD-card faulty! Info: Measurement w/o storage of data possible	SD card not writeable or readable (dirt, corrosion, faulty contacts)	Use of implant's internal calibration data. Data are not saved.
SD-card faulty!	SD card not readable (dirt, corro- sion, faulty contacts)	Check SD card for damage or dirt.
SD-card inserted! Restart measurement. Info: Storage of reading after restart of measurement possible!	Measurement starts without SD card. During measurement, the designated SD card for use with that implant is inserted	Start measurement again.
SD-card missing! Info: Insert SD card!	No SD card inserted in data ma- nagement mode	Insert SD card.

Message displayed	Cause	Detection/remedy
SD-card missing! Info : Measurement w/o storage of data possible - or - Insert SD-card with correct ID XXXXXXXXX!	SD card not inserted	Insert SD card designated for use with the SENSOR RESERVOIR / SENSOR PRECHAMBER into the reader unit. The implant (see patient card) is associated with the respective SD card via the identification number (ID).
SD-card not readable! Info: Measurement w/o data storage possible!	Insertion of wrongly formatted or unformatted SD card	Insert SD card designated for use with the SENSOR RESERVOIR / SENSOR PRECHAMBER into the reader unit. The implant (see patient card) is associated with the respective SD card via the identification number (ID).
SD-card not readable! Info: Measurement w/o storage of data possible!	SD card missing - or - SD card ID not readable - or - SD card contains no calibration data	Insert SD card designated for use with the SENSOR RESERVOIR / SENSOR PRECHAMBER into the reader unit. Measurement with data storage possible, but with limited accuracy.
SD-card memory full. Option to measure w/o saving	SD card memory full (100%) during continuous measurement	Delete any data no longer required, measure- ment without data storage possible.
SD-card memory full	SD card memory full (100%) during continuous measurement	Delete any data no longer required.
SD-card memory almost full	SD card memory almost full (99%) during continuous or fast measurement	Delete any data no longer required.
System error Info: Ambient pressure not readable	Barometric (BARD) pressure not readable	Device function locked. Deactivate and reac- tivate device. If the error persists, contact Technical Service.
System error - Incompatibility	Device hardware and software version not compatible	Device function locked. Deactivate and re- activate device. If the error persists, contact Technical Service.
System error - Antenna incompatible	Hardware version of antenna and reader unit not compatible	Device function locked. Deactivate and re- activate device. If the error persists, replace antenna and/or contact Technical Service.
System error - ID data incorrect	Implant identification data corrupted	Measurement is stopped. Deactivate and re- activate device. If the error persists, contact Technical Service.
System error - Implant voltage out of range	Implant voltage outside of permissible range	Measurement is stopped. Deactivate and re- activate device. If the error persists, contact Technical Service.

Message displayed	Cause	Fault detection/rectification
System error - Calibration data incorrect	Calibration data in implant cor- rupted or not readable (only applies if no SD card has been inserted)	Measurement is stopped. Deactivate and re- activate device. If the error persists, contact Technical Service. Measurement with inserted SD card possible.
System error - contact Service	System test detected an error	Device function locked. Deactivate and re- activate device. If the error persists, contact Technical Service.
Keypad faulty	A pressed key is registered upon activation of the reader unit	Release the key. Deactivate and reactivate device.
Temperature increase out of range	A spontaneous temperature increase in implant above 2 K or beyond 39.0 °C	Measurement is stopped. Leave to rest for 10 minutes. If the error persists, contact Technical Service.
Temperature value out of range	Temperature data measured from implant not plausible - physiologically meaningless data	Measurement is stopped. Deactivate and re- activate device. If the error persists, contact Technical Service.
Timeout! Restart measurement!	The time between starting measurement and successfully establishing communication was exceeded (60 seconds).	Restart measurement. Optimize distance between antenna and im- plant.
Ambient pressure incorrect Info: Ambient pressure from 800 to 1.100 mbar valid	The permissible barometric pres- sure during measurement is too low or too high.	The reader unit can only be used at barome- tric pressures from 800 to 1100 hPa. Current measurement is stopped.

OTHER MALFUNCTIONS/POSSIBLE FAULTS

Message displayed	Cause	Fault detection/rectification
Device cannot be activated	Battery is completely empty	Connect to power supply. It takes approx. 6 hours to charge the battery fully. The reader unit can also be operated during the charging process (with connected power supply unit). Note: Charging is not possible at an ambient temperature of more than 35°C.
Device switches off	Unfavourable operating condi- tions (e.g. low humidity or unsu- itable floor coverings)	Switch device back on. If the error persists, contact Technical Service.

TECHNICAL DATA

Description	Specifications and standards
Voltage range Reader unit Power supply unit	6 V (DC) 100-240 V (50–60 Hz)
Power consumption Reader unit Power supply unit	1.4 A (DC) 0.25 A (50–60 Hz)
Operating frequency	133 kHz
Pressure measuring range	-66.67 hPa bis 133.33 hPa±5.32 hPa
Temperature range for pressure measurement	10°C to 40°C (Reader unit) 20°C to 39°C (SENSOR RESERVOIR / SENSOR PRECHAMBER)
Functional distance between antenna and SENSOR RESERVOIR / SENSOR PRECHAMBER	10-30 mm
Battery life	5 years minimum
Battery charging cycle	at least 250 charging cycles
Battery self-discharge	Remaining charge after 3 months (storage) $> 70\%$
Flammability category of housing	UL 94 HB
Moisture protection or protection against ingress of liquids Reader unit Antenna Power supply unit	IP44 IP44 IP40
Durability - impact test	in accordance with IEC 60601-1 (3rd edition) 15.3.2
Durability - drop test	in accordance with IEC 60601-1 (3rd edition) 15.3.4.1
Weight: Reader unit Antenna Power supply unit	0.600 kg 0.215 kg 0.115 kg
Dimensions (W x H x D) Reader unit Antenna (without cable) Power supply unit	144 x 270 x 65 mm 100 x 250 x 25 mm 55 x 65 x 70 mm
Protection class (power supply unit)	11
Conformity to standards	IEC/DIN EN60601-1/ EN 45502-1
EMC	EN 60601-1-2:2007

AMBIENT CONDITIONS

Operating conditions	
Relative humidity	30 % to 75 %
Ambient temperature	10°C to 40°C, no condensation
Atmospheric pressure	800 hPa to 1100 hPa
Transport conditions	
Relative humidity	15 % to 95 %
Ambient temperature	0°C to 50°C
Atmospheric pressure	500 hPa to 1100 hPa
Storage conditions	
Relative humidity	15 % to 95 %
Ambient temperature	10°C to 40°C
Atmospheric pressure	500 hPa to 1100 hPa

CE MARKING

Approval for marking with the CE mark for active implantable medical devices (in accordance with Directive 90/385/EEC) was granted for the first time in 2011.

STORAGE AND TRANSPORT

- Store reprocessed products protected from dust in a dry, dark and cool room with as little bacterial contamination as possible.
- Ship product only in original packaging.
- For storage and transport conditions, see chapter on "Ambient conditions".

Caution: The product may be damaged if used too soon after storage/transport at temperatures below 10°C.

· Leave reader unit to acclimatize for approx. three hours at room temperature.

DISPOSAL

 Adhere to national regulations when disposing of or recycling the product, its components or its packaging!

Products carrying this symbol are have to be disposed of in the separate collection of electrical and electronic devices. Within the European Union, disposal is arranged free of charge by the manufacturer.

For recycling, please return the reader unit to the manufacturer.

If you have any questions regarding the disposal of the product, please contact Christoph Miethke GmbH & Co. KG; see also chapter on "Technical Service".

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- 1. Actionner la touche On/Off (©).