THOMALE GUIDE

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Ⓜ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
CONTENT

INDICATION	14
TECHNICAL DESCRIPTION	14
ASSEMBLY AND COMPONENTS	16
INSTRUCTION FOR MEASUREMENT OF THE INSERTION ANGLE	16
ANGLE ADJUSTMENT	18
APPLICATION	18
STERILE FIRST DELIVERY	19
RE-USE	19
RE-STERILIZATION	20
CONTRAINDICATION	21
ACCESSORIES	21
REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE RL 93/42/EEC	21
MEDICAL PRODUCT CONSULTANT	21
GENERAL INFORMATION	22

INDICATION

The *THOMALE GUIDE* (fig. 2) is a surgical instrument to control the placement of a ventricular catheter (VC) for the treatment of hydrocephalus. It is only used as navigation aid.

Therefore, the following conditions should be adhered to:

The *THOMALE GUIDE* is designated for the application of a ventricular catheter via a frontal pre-coronary burrhole. The usage of navigation devices like the *THOMALE GUIDE* is advisable for the implantation of a ventricular catheter.

It is possible that the application of a VC is possible without the *THOMALE GUIDE* in case of enlarged ventricles. It is for the surgeon to decide. In case of narrow ventricles (FOHR-value smaller than 0.55) the usage of the *THOMALE GUIDE* is appropriate. Especially for slit ventricles with a FOHWR-value smaller than 0.05 a navigation device is recommended.

TECHNICAL DESCRIPTION

The purpose of the VC is to allow the drainage of cerebrospinal fluid post surgical placement into the ventricles. A relevant risk factor for obstruction of the VC is the non-optimal placement and contact with adjacent tissue. In the past, catheters were placed mostly free-hand. In the literature it shows, that in 40-50% of treatments, an incorrect positioning of the VC can occur. Hence, in recent years X-rays or neuro navigation was needed in order to control the implantation procedure. This means a highly technical effort compared to a relatively small surgical intervention.

The *THOMALE GUIDE* was designed for an exact placement of the VC, using a simple technique, to avoid complications (fig.1 and fig. 2).

The *THOMALE GUIDE* allows a guided puncture of the ventricles with a VC, via a frontal pre-coronary burrhole, with a relatively low-tech technique.

Abb. 1: *THOMALE GUIDE* with ventricular catheter


Abb. 2: *THOMALE GUIDE*

It is based on the principle, that in the sagittal plane a right angle to the surface (tangent) always reaches the ventricles correctly (Fig. 3), whereas in the coronal plane an individual angle has to be adjusted in order to allow a correct placement (fig. 4 and fig.5). Therefore the *THOMALE GUIDE* is placed directly over the

burrhole with its three feet on the skull bones and with its line markings parallel to the centre line, so that a right angle in the sagittal plane is guaranteed. In the coronal plane an individual angle can be set with an angle adjustment on the *THOMALE GUIDE*. The rotation of the catheter guide, which is directed via guide rails, results from a very close fulcrum over the burr-hole.

The *THOMALE GUIDE* was designed for the Miethke-VC. For the Miethke-VC only the catheter guide with an inner diameter of 2,6 mm (ID2.6) should be used. Optionally, the *THOMALE GUIDE* is supplied with catheter guides with an inner diameter of 3,1 mm (ID3.1) and 3,6 mm (ID3.6) so that a range of ventricular catheters with different diameters can be used. It is very important to take into account that the catheters should have an external diameter 0.1 mm less than the catheter guide's inner diameter. Christoph Miethke GmbH & Co. KG does not assume any liability for the improper use of the *THOMALE GUIDE*.


Abb. 4: *THOMALE GUIDE* in coronal plane (flat skull)


Abb. 3: *THOMALE GUIDE* in sagittal plane



Abb. 5: *THOMALE GUIDE* in coronal plane (pointed skull)

ASSEMBLY AND COMPONENTS


The *THOMALE GUIDE* is initially supplied sterile and preassembled with the catheter guide ID2.6 (see chapter first delivery). After use the *THOMALE GUIDE* should be disassembled, cleaned and reprocessed (see chapter Re-use).

Step 1.


Install the catheter guide (1) from below into the module consisting of a base ring (3) and guide rails (2), so that the sliding surfaces of the catheter guide rest against the inner radii of the guide rails.

Step 2.

Slip the slider (4) on the top of the catheter guide so that the inner radii corresponds and rests directly on the outer radii.

Step 3.

Screw the locking nut (5) hand-tight onto the screw thread of the catheter guide with the small ledge of the locking nut pointing toward the slider.

The steps above are performed in reverse order for disassembly (for example to change or to clean the catheter guide).

INSTRUCTION FOR MEASUREMENT OF THE INSERTION ANGLE

Basically, for the measurement of the insertion angle, determine the selected angle of the ventricular catheter to the tangent in the coronal plane above the burrhole. Therefore we propose two possibilities. It is in the responsibility of the user to decide for the most practical method.

1. The most accurate method to determine the individual insertion angle of the ventricular catheter guide to the tangent is done by importing biparietal data of the skull into a planning software (for example, Osirix or a similar navigation software). An optimal trajectory for the ventricular catheter can thereby be defined. The coordinates of the entry point, e.g. into the nasion and the centre line can be measured. Draw the tangent of the entry point on the projected coronal plane and measure the angle of the trajectory to the tangent (Fig. 6 and 7). Now the insertion

parameters can be applied on to the patient data and intra-operatively the angle can be applied to the *THOMALE GUIDE* for catheter insertion. This method is especially recommended in case of very narrow ventricles.

Caution: The accuracy of the determination of the angle depends on the used image / the underlying image.

get. Now, draw the tangent. Proceed by setting two points, both with the same distance medial and laterally to the entry point. The relation between those two points correspond to a tangent. Accordingly, the angle between trajectory and tangent can be determined and set on the *THOMALE GUIDE*.

Fig. 6: CCT in a reconstructed layer along the catheter


Fig. 8: CCT in coronary layer


Fig. 7: MRI in a reconstructed layer along the catheter


Fig. 9: MRI in coronary layer

2. To allow an easier unique measurement of the trajectory angle, import a coronary cross-section. Preferably, the image should show the level of the anterior horn in front of the foramen Monro (Fig. 8 and 9). Mark on this image the entry point, the projected target location in the ventricle and thereby the resulting trajectory. The distance to the centre line and the depth of the ventricle need to be considered when marking the projected tar-

get. Additionally to these methods and this instruction for use, an application for smart phones is available.

For example: by importing an image, the angle of the trajectory to the tangent can now easily be determined with the help of the application.

ANGLE ADJUSTMENT

After loosening the lock nut slightly (a quarter turn in counterclockwise direction) the catheter guide can slide along the guide rails and the required angle can be set.


Fig. 10: Angle adjustment on the THOMALE GUIDE

The scale (Fig. 11) of the guide rail is marked by thick solid lines in increments of 10° and the thin half-lines at increments of 5° . It is recommended to use a lens or a magnifying glass for the adjustment of the angle. The setting of an angle is done by aligning the arrow on the slider to the desired marking on one of the two rails. Our recommendation is to start the adjustment always at the 0° marking. The 0° mark is characterized by a higher lineweight (Fig. 11).


Fig. 11: Angle adjustment at the THOMALE GUIDE

The sequent 5° - or 10° markings are added until the required angle is reached. Angles which are not multiples of 5° can approximately be estimated between the markings. Angles between 0° and 60° can be set on both sides. The overall THOMALE GUIDE, except for the base ring, is completely symmetrical. This allows the angle to be set on either side of the

scale. Once the correct angle has been set, it can be fixed by lightly tightening the locking nut clockwise.

APPLICATION

Prior to implantation of ventricular catheters using the THOMALE GUIDE, clear the skull from the periosteum in a circular area of about 20 to 25 mm around the burrhole.

Caution: We therefore recommend to remove the periosteum completely from the skull bone. The correct placement of the catheter in the ventricle can only be guaranteed if the feet are placed directly on the skull bone. Close attention has to be taken that the skull does not show any unevenness resulting from previous operations.

Place the THOMALE GUIDE on the skull bone and aligned in such a way, that its rotation axis with the open side of the base ring is aligned parallel to the median plane (Fig. 12).


Fig. 12: THOMALE GUIDE aligned parallel to the median plane

To help find and maintain the correct position, especially with a patient positioned on their side, marks can be applied to the skin. Particular attention should be taken during alignment that the coronary insertion angle setting is correct and that the lateral opening of the base ring is aligned towards the median plane (Fig.12). The ventricular catheter can now be inserted and pushed through the catheter guide up to the dura. The positioning of the catheter

tip can be controlled by the lateral opening in the base ring.

Please note that while cautiously pushing the ventricular catheter, the catheter tip points directly to the previously created dura opening and that the dura boarders will not be affected. Each contact with the dura bears the risk of deflecting the catheter. Control the insertion depth whilst moving the catheter in the direction of the ventricle. Through the lateral openings of the base ring, the markings on the VC can be identified and read. The standard insertion length in adult patients is about 60 to 65 mm. Ideally, the length can be determined individually by pre-operative measurement of a coronary image. When positioned correctly, the catheter is then fixed at the burrhole level with the aid of atraumatic forceps through the lateral opening of the base ring. Remove the stylet and check for correct CSF flow. Now the *THOMALE GUIDE* can be removed from the site without changing the position of the catheter.

STERILE FIRST DELIVERY

The *THOMALE GUIDE* is supplied sterile and it is re-sterilizable.

The products are sterilized with steam under closely monitored conditions. The double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

Note: in any case the position of the VC should be monitored post-operatively. The Christoph Miethke GmbH & Co.KG assumes no liability in case of misplacement the ventricular catheter.

RE-USE

Cleaning and disinfection of the *THOMALE GUIDE*

To ensure validated cleaning of the *THOMALE GUIDE* it has to be disassembled into its four components (see chapter Assembly and Components).

Avoid damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for surgical steels according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55°C.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning for products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect the *THOMALE GUIDE* mechanically, provided it can be securely fixed in machines or storage devices in such a way that the *THOMALE GUIDE* will be thoroughly cleaned.

Manual cleaning/disinfecting

- Check visible surfaces for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Mechanical cleaning/disinfecting

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).

Control, care and inspection

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged and does not have any loose, bent, broken, cracked, worn,

or fractured components.

- Set aside the product if it is damaged. We recommend putting the *THOMALE GUIDE* in a double sterile bag after cleaning. So the next use in the theatre (sterile area) will be prepared as its best.

Mechanical cleaning/disinfecting with manual pre-cleaning

Manual pre-cleaning with ultrasound

Stage	Step	T (°C/°F)	t (min)	Conc. (%)	Waterquality	Chemical
i	Ultrasonic disinfection cleaning	RT (cold)	15	2	T-W	BBraun Stabimed; aldehydphenol- and QAV-free; pH = 9
ii	Intermediate rinse	RT (cold)	1		T-W	

T-W: Drinkingwater, RT: Roomtemperature

Mechanical alkaline cleaning and thermal disinfection

Machinetype: Single-chamber-cleaning/ disinfection machine without ultrasound

- Place product on a tray that is suitable for cleaning (avoid rinsing blind spots)

Stage I

- Clean the product in a ultrasonic cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage	Step	T (°C/°F)	t (min)	Waterquality	Chemical
I	Pre-rinse	<25/77	3	-	
II	Cleaning	55/131	10	VE-W	- Concentrate alkaline: pH = 10,9 < 5 % anionic tenside - 1% solution pH = 10,5
III	Neutralization	20/68	2	VE-W	- Concentrate, acid: pH = 2,6 Basis: Citric acid - 1% solution pH = 3,0
IV	Intermediate rinse	70/158	1	VE-W	
V	Thermodisinfected	94/201	10	VE-W	
VI	Drying	90/194	40	-	

T-W: Drinkingwater, VE-W: Deminiralized water

RE-STERILIZATION

The *THOMALE GUIDE* should be sterilized in its sterile double wrapping in an autoclave (steam sterilization, fractionated vacuum method) at 134° C and 5 minutes cycling time. Demount all parts prior to sterilization. Since different catheter guides can be applied during operation, corresponding to the used ventricular catheter, we recommend to pack each catheter guide in a separate sterile bag and to mark them accordingly.

CONTRAINDICATION

The usage of the *THOMALE GUIDE* is advised against if a frontal way to the ventricles is not possible or not indicated, a parietal or occipital entry to the ventricle is not possible, an irregular anatomy of the skull surface or a FOHWR smaller than 0,05 exist.

Explanation: At a smaller FOHWR-value particularly close slit ventricles exist. In this case use a navigation guide for the positioning of the VC.

ACCESSORIES

The *THOMALE GUIDE Application* was designed as an accessory for the *THOMALE GUIDE*.

It serves for determination and calculation of the insertion angle, of the catheter length and of the position of the burrhole when puncturing the pre-coronary ventricle with a VC.

REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE RL 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied on human beings. The requirements on product identification and traceability has to be maintained also at the user's area.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available up to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

MEDICAL PRODUCT CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

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
GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product	<i>THOMALE GUIDE</i>
Intended use	Controlled placement of a ventricular catheter (VC) for the treatment of hydrocephalus
Re-sterilizable	
Store in a dry, clean place	
<p>Schematic representation of the <i>THOMALE GUIDE</i> with its external dimensions:</p> 	

CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
CE marking according to directive 93/42/EEC
Label CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE

Technische Änderungen vorbehalten
Technical alterations reserved
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifichie tecniche

Manufacturer acc. MDD 93/42/EEC:

■ **CHRISTOPH MIETHKE GMBH & CO. KG**

Christoph Miethke GmbH & Co. KG | Ulanerweg 2 | 14469 Potsdam | Germany
Phone +49 (0) 331 62 083-0 | Fax +49 (0) 331 62 083-40 | www.miethke.com

Distributed by:

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 74 61 95-26 00 | www.aesculap.com

Aesculap - a B. Braun company

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AESCLAP AG