Tunneller

**A NEW FLEXIBILITY IN SHUNT SURGERY.
YOU DEFINE THE PATH.**

DE Gebrauchsanweisung | EN Instructions for use | FR Mode d'emploi | ES Instrucciones de manejo
IT Istruzioni per l'uso | PT Instruções de utilização | NL Gebruiksaanwijzing

 www.miethke.com

US This Instruction for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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info@miethke.com

CONTENTS

0.00 PREFACE AND IMPORTANT INFORMATION	15
1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	15
1.01 Explanation of the warnings	15
1.02 Display conventions	15
1.03 Other accompanying documents and additional information	16
1.04 Feedback on the instructions for use	16
1.05 Copyright, disclaimer, warranty and other information	16
2.00 DESCRIPTION OF THE <i>Tunneller</i>	16
2.01 Medical purpose	16
2.02 Clinical use	16
2.03 Indications	16
2.04 Contraindications	16
2.05 Intended patient groups	16
2.06 Intended users	16
2.07 Intended use environment	17
2.08 Technical description	17
2.09 System components	18
3.00 <i>Tunneller</i> PROPERTIES	18
3.01 Product description	18
3.02 Important safety information	19
3.03 Transport and storage	20
3.04 Use of the product	20
3.05 Disposal	21
3.06 Troubleshooting	22
3.07 Technical information	23
3.08 Symbols used for labelling	24
4.00 MEDICAL DEVICE CONSULTANTS	24

0.00 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing the *Tunneller*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product.

Your team at Christoph Miethke GmbH & Co. KG.

Relevance of the instructions for use

WARNING

Improper handling and use of this product can cause risks and damage. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

The *Tunneller* is comprised of the following components:

- ▶ Tunnelling rod
- ▶ Handle

Length	Standard variant	Rigid variant
30 cm	FX005SU	FX001SU
45 cm	FX006SU	FX002SU
60 cm	FX007SU	FX003SU
70 cm	FX008SU	FX004SU

1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

1.01 EXPLANATION OF THE WARNINGS

DANGER

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

WARNING

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

CAUTION

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

NOTE

Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for risk, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

1.02 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates product names

1.03 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

1.04 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

1.05 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

2.00 DESCRIPTION OF THE *Tunneller*

2.01 MEDICAL PURPOSE

The *Tunneller* is a surgical instrument that is used to create a subcutaneous tunnel in which a Shunt System or shunt components are placed to drain cerebrospinal fluid (CSF).

2.02 CLINICAL USE

Treatment options using the variants

- ▶ 8 variants of the *Tunneller*
- ▶ Ergonomic handpiece with different handle options
- ▶ 2 different tips (sharp and blunt) for each *Tunneller*
- ▶ Quick and easy placement of the catheter into the subcutaneous tunnel
- ▶ Two-way tunnelling possible from retroauricular or paraumbilical

2.03 INDICATIONS

The *Tunneller* is used for the following indications:

Indications

- ▶ Treatment of hydrocephalus

2.04 CONTRAINDICATIONS

The following contraindications apply to the *Tunneller*:

Contraindications

- ▶ Coagulation disorders (risk of postoperative bleeding)
- ▶ Infections or a suspected infection in the regions of the body affected by the implantation (e.g. skin infection, meningitis, ventriculitis, bacteraemia, septicaemia, also peritonitis if the catheter is integrated into a Shunt System)
- ▶ Allergy to the materials or procedures used

2.05 INTENDED PATIENT GROUPS

Patients who are treated with a CSF-draining Shunt System on account of their clinical picture.

2.06 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product

must only be used by users with the following qualifications:

- ▶ Neurosurgeons
- ▶ Trained medical professionals with knowledge of the function and intended use of the product

2.07 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Use under sterile operating theatre conditions

2.08 TECHNICAL DESCRIPTION


Fig. 1: Tunneller, comprising a handle (1) and tunnelling rod (2) with a transport protection cover (3)

The Tunneller is made up of two parts (Fig. 1): A handle (1) and a tunnelling rod (2) fitted with a transport protection cover (3) on the accessible end. The handle (Fig. 2) is made of plastic and is comprised of a recessed handle (1a) and a locking rocker (1b).


Fig. 2: Handle, comprised of a recessed handle (1a) and locking rocker (1b)

Pressing on the back of the rocker ('OPEN') unlocks the handle and releases the tunnelling rod. The tunnelling rod can be removed and assembled lengthwise and crosswise (Fig. 3). The tunnelling rod must be inserted into the handle all the way to the stop. The tunnelling rod is correctly secured in place when the user feels it click into place upon pressing onto the front of the locking rocker.


Fig. 3: Tunneller handle lengthwise (left) and crosswise (right)

The tunnelling rod (Fig. 4) is made of surgical stainless steel and is available in the effective lengths of 30, 45, 60 and 70 cm. The tunnelling rod is available in two strengths; the more flexible version is referred to as "standard" and the less flexible version as "rigid".


Fig. 4: Tunnelling rod, left, with the sharp (2a) and, right, with the blunt (2b) end, catheter hole (2c) and marking (2d) only on standard rods

The standard tunnelling rods feature a three-part ring marking (2d) near the tunnelling tips. Each tunnelling rod has a sharp (2a) and a blunt (2b) end. Both ends have a drill hole (2c) that is suitable for holding catheters with an outer diameter of up to 2.5 +/- 0.1 mm.

CAUTION

Catheters with a diameter of > 2.6 mm cannot be used with the Tunneller as this may damage the catheter material.

When delivered, the Tunneller from Christoph Miethke GmbH & Co. KG is configured with a lengthwise handle and a sharp-ended tunnelling rod. The sharp end has a transport protection cover (Fig. 1, Pos. 3) that is to be removed directly before use. The user can configure the Tunneller based on the individual patient's anatomy and the preferred tunnelling technique within the scope of the options above.

**WARNING**

The transport protection cover is used to ensure the integrity of the packaging. It must be removed immediately before tunnelling and is not intended to come into contact with patient tissue.

2.09 SYSTEM COMPONENTS

The *Tunneller* can be safely used to implant shunt components from Christoph Miethke GmbH & Co. KG. We recommend only using products from Christoph Miethke GmbH & Co. KG in combination with the *Tunneller*.

The *Tunneller* is compatible with all catheter products from Christoph Miethke GmbH & Co. KG, including all shunt configurations with a catheter.

3.00 Tunneller PROPERTIES**3.01 PRODUCT DESCRIPTION****3.01.01 Tunneller VARIANTS**

The *Tunneller* is available in the following variants:

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3.01.02 SCOPE OF DELIVERY

Box content	Number
<i>Tunneller</i> in sterile packaging	10
Instructions for use of the <i>Tunneller</i>	1

3.01.03 STERILITY

The Christoph Miethke GmbH & Co. KG *Tunneller* is sterilised with steam under controlled conditions and is packed in a double sterile barrier. The sterile barriers have an indicator that turns brown after sterilisation.

**WARNING**

The products must not be used if the packaging or the product is damaged or after the expiry date.

3.01.04 REPEATED USE AND RESTERILISATION

Tunnellers that have already been used to create a subcutaneous tunnel must not be reused either on the same or a different patient in order to minimise the risk of infection. The *Tunneller* is intended for single use and must not be re-sterilised.

WARNING

The *Tunneller* must not be re-sterilised or reprocessed in any other way as the safe functioning and sterility of the product cannot be guaranteed.

3.01.05 OPERATING CONDITIONS

Sterile operating environment

3.01.06 PRODUCT LIFE

These medical devices have been constructed in such a way as to ensure their precise and reliable operation. The expected service life of the *Tunneller* is limited to the duration of use during a single surgical procedure.

3.01.07 DISPOSABLE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the *Tunneller*. No guarantee can be assumed for the functional safety of re-sterilised products.

3.01.08 PRODUCT CONFORMITY

The product meets current regulatory requirements.

The Medical Device Directive requires the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the medical device should therefore be recorded in the patient's medical records to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

3.02 IMPORTANT SAFETY INFORMATION

3.02.01 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ If the sterile packaging is damaged, the products must not be used under any circumstances.
- ▶ The product must not be implanted after its expiry date.
- ▶ Products that have already been used must not subsequently be reused on the same or another patient.

CAUTION

- ▶ Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be read carefully and understood before the product is used for the first time.
- ▶ Prior to use, it is essential to check that the *Tunneller* is complete and intact.

3.02.02 COMPLICATIONS, SIDE EFFECTS AND RESIDUAL RISKS

Various side effects and complications can occur associated with the *Tunneller* and the implantation of a Shunt System, particularly along the tunnel channel and/or at the entry and exit points:

- ▶ Redness/irritation of the skin and tightness around the subcutaneous tunnel as an indication of a possible infection
- ▶ Haematoma
- ▶ Infection
- ▶ Wound healing disorders
- ▶ Injuring or perforating the surrounding tissue or organs
- ▶ Allergic reaction
- ▶ Scarring

A physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the *Tunneller*:

- ▶ Serious infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Damaging/puncturing the tissue
- ▶ Skin irritation
- ▶ Allergic reaction

3.02.03 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible authority of the EU member state in which you work.

3.02.04 PATIENT INFORMATION

The attending doctor is responsible for informing the patient and/or his/her proxy in advance. This includes a comprehensive description of the course of the operation, the surgical technique and the medical devices to be used. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product.

3.03 TRANSPORT AND STORAGE

3.03.01 TRANSPORT

Transport conditions

Ambient temperature	-18 °C to +50 °C
---------------------	---------------------

3.03.02 STORAGE

The medical devices must always be stored in a clean and dry place.

3.04 USE OF THE PRODUCT

3.04.01 INTRODUCTION

The *Tunneller* is a surgical instrument that is used to create a subcutaneous tunnel into which a Shunt System or shunt components are placed to drain CSF.

The *Tunneller* is suitable for tunnelling from a proximal to distal direction and vice versa.

3.04.02 SAFETY AND WARNING NOTICES

CAUTION

- ▶ Prior to use, it is essential to check that the *Tunneller* is complete and intact.
- ▶ In the event of improper use of the *Tunneller* (e.g. too much force or incorrect length of the rod), the handle or parts of the handle may break or become damaged, the tunnelling rod may bend or the connection between the handle and the rod may fail. In this case, the proper functioning of the *Tunneller* can no longer be guaranteed, and the instrument must be replaced with a new one. Improper use may lead to injuries to the user and/or the patient and impair the functionality of the Shunt System to be implanted.
- ▶ Selecting a layer of skin that is unsuitable for tunnelling may injure the surrounding tissue or organs and increase the risk of infection.
- ▶ Do not reshape the tunnelling rod during its subcutaneous insertion in order to avoid injuries to the skin.
- ▶ Correct tunnelling is essential when passing through particularly critical areas (collar bone, thorax, fascia) in order to avoid injuring the surrounding tissue structures or organs and accidentally piercing the skin.

**CAUTION**

- ▶ While tunnelling, the user should always grip the handle firmly in order to ensure a safe guidance of the tunnelling rod through the tissue.
- ▶ The use of catheters with an outer diameter of > 2.6 mm can cause damage to the catheter material, which can result in the impairment or loss of function of the Shunt System.
- ▶ When pulling through the catheter, care must be taken that the catheter is not under tension in order to prevent damage to the catheter material and to avoid breaking the connection between the tunnelling rod and the catheter. When pulling back the *Tunneller*, the user should always grip the handle firmly in order to ensure a safe guidance of the tunnelling rod through the tissue.

3.04.03 REQUIRED MATERIALS

The *Tunneller* is designed for safe use with the shunt components described in chapter "2.09 System components".

3.04.04 PREPARING FOR SURGERY**Checking the sterile packaging**

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Checking the integrity of the *Tunneller*

Prior to surgery, the integrity of the *Tunneller* must be checked. To do so, check the connection and handling of the tunnelling rod and the handle.

Remove the protection cover of the tunnelling rod before use.

Position of the patient

The patient is to be placed in a suitable position to ensure that the tunnel is as straight as possible.

3.04.05 TUNNELLING**Surgical technique**

The surgeon is responsible for selecting the correct surgical technique and the right instrument configuration, which depends on the patient. The patient is to be placed in a suitable position to ensure that the tunnel is as straight as possible. An additional relief incision may be necessary. The tunnelling rod is malleable and can be adjusted by the user to the patient's anatomy and the procedure.

The tunnelling rod is inserted into the subcutaneous tissue through an incision in the skin in accordance with the indications of the surgical procedure. The tunnelling rod is pushed forwards until it reaches the desired exit point.

Placement of the catheter

The catheter is pulled through the hole at the end of the tunnelling rod and secured with a ligature. The *Tunneller* is carefully pulled back by the handle, which places the catheter inside the subcutaneous tunnel.

Depending on the selected tunnelling technique, the position and the direction of the handle to pull through the catheter can be changed. Asepsis must be ensured at all times when handling the product. Only use the product for the intended purpose.

3.05 DISPOSAL

Product components that have been used and have been in contact with patients must be disposed of correctly as potentially infectious material in accordance with medical practice as well as respective regional laws and regulations.

Tunnellers that have been used once must not be reused.

3.06 TROUBLESHOOTING

The table below lists possible errors or malfunctions, their cause and how to remedy them:

Error/malfunction	Cause	Remedy
The handle cannot be separated from the tunnelling rod	The tunnelling rod is locked in the handle	Press on the back of the locking rocker (OPEN) to release the handle
The tunnelling rod cannot be inserted into the handle	The handle is not or not fully open	Press down on the back of the locking rocker (OPEN) as far as it will go
The tunnelling rod cannot be locked into the handle	The tunnelling rod is not fully inserted into the handle	Insert the tunnelling rod into the handle as far as it will go and press the locking rocker down
The tip is not suitable for the intended patient	The sharp or blunt side of the tunnelling rod is in the handle	Remove the tunnelling rod and turn it around
The direction of the handle is not suitable for the planned procedure	The handle can be positioned lengthwise and crosswise	Release the handle, remove the tunnelling rod and reinsert it in the required direction
The tunnelling rod is not suitable for the intended tunnel	The tunnelling rod does not correspond to the anatomical shape of the intended tunnel	The tunnelling rod can be bent to fit the anatomical shape


3.07 TECHNICAL INFORMATION

3.07.01 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>Tunneller</i>
Intended purpose	Subcutaneous tunnelling to place a catheter or shunt system
Sterilisability	Cannot be resterilised
Storage	Store in a clean and dry place
Material	Tunnelling rod: 1.4301 surgical stainless steel Handle: Polyoxymethylene (POM) Transport protection cover: TPE

For single use only

Sketch with outer dimensions:

3.07.02 MATERIALS THAT COME INTO CONTACT WITH BODILY TISSUE/FLUIDS

- ▶ Tunnelling rod: 1.4301 surgical stainless steel
- ▶ Handle: Polyoxymethylene (POM)

Quantitative material properties

Tunnelling rod	
1.4301 surgical stainless steel:	
Volume [mm ³]	4.2 – 9.2
Weight [g]	32.8 – 72.5
Contact surface [mm ²]	4,310 – 9,340

Handle	
Polyoxymethylene (POM):	
Volume [mm ³]	17
Weight [g]	22
Contact surface [mm ²]	7,650

3.08 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxx indicates the number of the responsible Notified Body
	Medical device
	Manufacturer
	Date of manufacture
	Use by
	Production lot number, batch
	Article number
	UDI number (Unique Device Identifier)
	Steam-sterilised
	Do not resterilise
	Do not reuse
	Do not use if package is damaged
	Store in a dry place
	Temperature limitation
	Consult instructions for use / electronic instructions for use
	Caution, read the documentation supplied with the product
	Pyrogen-free
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians

Symbol	Explanation
	MR unsafe
	Date
	Model number / European Medical Device Nomenclature Code

4.00 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany