M.flow[®]

EXPAND FLEXIBILITY AND MANAGE THE FLOW

 Gebrauchsanweisung |  Instructions for use |  Mode d'emploi
 Instrucciones de manejo

 www.miethke.com

 This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0
info@miethke.com

CONTENTS

1	PREFACE AND IMPORTANT INFORMATION	22
2	INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	22
2.1	EXPLANATION OF THE WARNINGS	22
2.2	DISPLAY CONVENTIONS	22
2.3	OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION	22
2.4	FEEDBACK ON THE INSTRUCTIONS FOR USE	22
2.5	COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION	23
3	DESCRIPTION OF THE <i>M.flow</i>	23
3.1	MEDICAL PURPOSE	23
3.2	CLINICAL BENEFITS	23
3.3	INDICATIONS	23
3.4	CONTRAINDICATIONS	23
3.5	INTENDED PATIENT GROUPS	23
3.6	INTENDED USERS	23
3.7	INTENDED USE ENVIRONMENT	23
3.8	TECHNICAL DESCRIPTION	24
3.9	OPERATING PRINCIPLE OF THE VALVE	25
3.10	PRESSURE-FLOW CHARACTERISTICS	26
3.11	FLOW-REDUCTION LEVEL DETERMINATION IN X-RAY IMAGES	27
3.12	APPLICATION OF <i>M.flow Instruments</i>	28
3.13	SYSTEM COMPONENTS	31
3.14	FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	32
4	PROPERTIES OF THE <i>M.flow</i>	32
4.1	PRODUCT DESCRIPTION	32
4.2	IMPORTANT SAFETY INFORMATION	33
4.3	PATIENT EDUCATION	34
4.4	TRANSPORT AND STORAGE	34
4.5	USE OF THE PRODUCT	34
4.6	EXPLANATION AND DISPOSAL	37
4.7	TECHNICAL INFORMATION	37
4.8	SYMBOLS USED FOR LABELLING	38
5	MEDICAL DEVICE CONSULTANTS	38

1 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing the medical device *M.flow*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product. Your team at Christoph Miethke GmbH & Co. KG

Relevance of the instructions for use

WARNING

Improper handling and non-intended use of this product can cause risks and damages. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

The *M.flow* includes the following components:

- ▶ *M.flow*

Additional options:

- ▶ MIETHKE valve
- ▶ *Burrhole Reservoirs*
(including paediatric versions)
- ▶ *Prechambers*
(including paediatric versions)
- ▶ *Burrhole Deflectors*
(including paediatric version)
- ▶ *Ventricular Catheters* with stylet
- ▶ *Peritoneal Catheters*

2 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

2.1 EXPLANATION OF THE WARNINGS

DANGER

Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE

Indicates information considered important, but not hazard-related (e.g. messages relating to property damage).

The symbols for danger, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

2.2 DISPLAY CONVENTIONS

Display	Description
Italics	Indicates <i>product names</i>

2.3 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

The latest version of the instructions for use as well as translations into additional languages can be found on our website:

<https://miethke.com/downloads/>

Each delivery is accompanied by a patient card and patient labels with information about the product. The patient labels are intended to provide the treating physician with all the product information in a compact form for the patient record.

If you still need additional information despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

2.4 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests or criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

2.5 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products from another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

3 DESCRIPTION OF THE *M.flow*

3.1 MEDICAL PURPOSE

The product *M.flow* is intended for the shunting of cerebrospinal fluid (CSF).

3.2 CLINICAL BENEFITS

The clinical benefit of the MIETHKE shunt component *M.flow* as part of a MIETHKE shunt system is improvement of symptoms in hydrocephalus therapy.

3.3 INDICATIONS

The following indications apply to the *M.flow*:

- ▶ Hydrocephalus

3.4 CONTRAINDICATIONS

The following contraindications apply to the *M.flow*:

- ▶ Infections in the implantation area
- ▶ Pathological concentration levels (e.g. of blood components and/or protein) within the cerebrospinal fluid
- ▶ Intolerance to materials of the shunt system

3.5 INTENDED PATIENT GROUPS

- ▶ Patients who are treated with a cerebrospinal fluid shunt system on account of their clinical picture

3.6 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical professionals, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product. These can be acquired, e.g., in training offers organised by Christoph Miethke GmbH & Co. KG (see Ch. 5).

3.7 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions

3.8 TECHNICAL DESCRIPTION


Fig. 1: Cross section of the M.flow

1. Adjustable flow-reducing unit

- | | | |
|------------------|---------------------------|--------------------------------------|
| 2. Rotor | 5. Shunt switch-off range | 7. Inlet connector with check valve |
| 3. Channel inlet | 6. Outlet | 8. Outlet connector with check valve |
| 4. Flow channel | | |

The *M.flow* is a valve made of titanium with an adjustable flow-reducing unit and an additional shunt switch-off function (Fig. 1).

Flow resistance can be set by adjusting the flow-reducing unit. Increasing the flow resistance makes it possible to reduce the CSF diversion volume per time unit in the shunt system and thus control the CSF quantity to be shunted. The *M.flow*'s shunt switch-off function allows cerebrospinal fluid shunting to be fully halted.

The flow-reducing unit (1) consists of a rotor (2) with an aperture, the channel inlet (3) and a circular flow channel (4) with a channel interruption, the shunt switch-off range (5) and an outlet (6). Depending on the position of the rotatable mounted rotor, a connection between the channel inlet and the flow channel is established or prevented in the area of the channel interruption (shunt switch-off function).

The position of the channel inlet in the rotor to the flow channel determines its effective length and thus the set flow resistance (see Ch. 3.9).

The flow resistance, i.e. the flow-reduction level, can be adjusted with the *M.flow Instruments* postoperatively and non-invasively through the skin using the clickable valve membrane.

In addition to the adjustable flow-reducing unit, check valves are integrated in the inlet and outlet connectors (7,8), which ensure a directed flow in the distal direction (Fig. 2).


Fig. 2: Inlet connector with check valve
a) open, b) closed

In addition, the volume inside the *M.flow* can be reduced by pushing in the clickable valve membrane. Due to the arrangement of the check valves in the connectors combined with the clickable valve membrane, the *M.flow* has a distal pumping property. The volume per pump cycle is approx. 0.04 ml.

WARNING

When pumping the M.flow, an adjustment of the flow-reduction level by pressing the valve membrane several times cannot be ruled out. Checking the flow-reduction level after the pump cycle is therefore required! The patient should be properly informed about this risk.

3.9 OPERATING PRINCIPLE OF THE VALVE

The M.flow is a posture-independent flow-reducing valve. The flow-reducing unit allows the flow resistance to be adjusted and thus the shunted cerebrospinal fluid volume per unit of time to be controlled. This makes the M.flow significantly different from conventional valves, where the pressure value is adjusted.

The flow resistance, and thus the cerebrospinal fluid shunt volume per unit of time, is influenced by the effective length of the flow channel. This effective length is determined by the length of the flow channel from the channel inlet to the outlet and can be changed by the rotor position.

In the rotor position in which the channel inlet is congruent with the outlet (Fig. 3), the length of the flow channel has no effect. Cerebrospinal fluid flows directly through the flow-reducing unit without the detour via the flow channel. The flow resistance, and thus the flow-reducing effect, is minimal in this setting (flow-reduction level 0).


Fig. 3: Rotor position for minimal flow resistance (flow-reduction level 0)

In the rotor position in which the channel inlet at the other end of the flow channel (Fig. 4), the flow resistance, and thus the flow-reducing effect, is at a maximum (flow-reduction level 10). In this setting, the liquid must first flow through the entire flow channel before it reaches the flow-reducing unit outlet. This allows the shunted cerebrospinal fluid volume per time unit to be significantly reduced.


Fig. 4: Rotor position for maximum flow resistance (flow-reduction level 10)

Any desired setting can be chosen between the two extremes (Fig. 5).


Fig. 5: Rotor position for flow-reduction level 7

To temporarily interrupt the connection of the channel inlet and outlet of the M.flow, the channel inlet in the rotor can be completely closed. To do this, place the rotor in the position where the channel inlet of the rotor is in the range of the flow channel interruption (shunt switch-off range, Fig. 6). The channel inlet is now closed by the flat surface. No cerebrospinal fluid shunting is possible with this setting 'M.flow off' (⊗).


Fig. 6: Rotor position of the shunt switch-off range, no flow ('M.flow off' ⊗)

Flow-reduction level	Flow resistance	Cerebrospinal fluid shunt volume per time
0	minimal	almost unaffected
10	maximal	significantly reduced
⊗	-	no shunting, 'M.flow off', shunt impermeable

3.10 PRESSURE-FLOW CHARACTERISTICS

When applying a flow rate of 5 to 50 ml/h, the following pressure-flow characteristics are obtained. The pressure at the reference flow rate of 20 ml/h is highlighted in bold. The reference flow rate is highlighted by a dotted line (...). Note that testing the product may provide different results depending on the test setup.

The valve design allows continuously variable adjustment (shown by double arrows: ←→) between the illustrated pressure values to find the individual pressure setting for each patient.

The combination with *Catheters* (inner diameter > 1 mm) and shunt components (*Burrhole Deflectors*, *Reservoirs*, *Titanium Connectors*) does not fundamentally affect the pressure-flow characteristics.

The pressure-flow characteristics of the continuously adjustable *M.flow* are shown below exemplarily for the flow-reduction levels 0, 2, 6 and 10, irrespective of body position:


Fig. 7: Pressure-flow characteristics of selected flow-reduction levels of the *M.flow*

Flow-reduction level	Flow [ml/h]	Pressure [cmH ₂ O]	Tolerance [cmH ₂ O]
0	5	0–2	+ 1 / – 1
	50	1–11	+ 5.5 / – 4.5
2	5	0–8	+ 4.5 / – 3.5
	50	10–30	+ 10 / – 10
6	5	2–18	+ 9 / – 7
	50	30–65	+ 17 / – 18
10	5	3–25	+ 11 / – 11
	50	55–95	+ 18 / – 22

Influence on the shunt system

For integration of the *M.flow* into a shunt system, the pressure values of the individual components are added for physical reasons. Fig. 8 shows an example of the combination of a valve with a gravity and differential pressure unit (opening pressure = 25 cmH₂O) with an *M.flow* (flow-reduction level 6) in a vertical position.


Fig. 8: Pressure-flow characteristics of a shunt system, consisting of M.flow and M.blue 5
 a) M.blue 5 (gravitational unit = 20 cmH₂O), vertical
 b) M.flow (flow-reduction level 6)
 c) Shunt system (combination of both valves a + b)

WARNING

The combination of the M.flow with another valve leads to significantly higher pressure values in the shunt system, depending on the selected flow-reduction level.

3.11 FLOW-REDUCTION LEVEL DETERMINATION IN X-RAY IMAGES

The selected flow-reduction level of the M.flow should always be monitored using the M.flow Compass, but it can also be checked using an X-ray image (Fig. 9).


Fig. 9: X-ray image for M.flow with flow-reduction level 4 set

The position of the channel inlet (1) is crucial (Fig. 10), which is determined by the position of the rotor.


Fig. 10: Schematic representation of rotor in an X-ray image

1. Channel inlet of the flow-reducing unit
2. Flow-reduction level, 3. Rotor stop
4. Shunt switch-off range 'M.flow off' (⊗)
5. Non-adjustable range
6. Outlet side of the flow-reducing unit
7. Proximal, 8. Distal

For determining the flow-reduction level, the X-ray template supplied with the M.flow Instruments is placed over the X-ray image. When aligning the X-ray template, the position of the rotor stop (3) and the outlet side (6) must match the X-ray image.

The position of the channel inlet (1) shows the selected flow-reduction level (2) (Fig. 10).

If the position of the channel inlet is in the dark grey range (4), the flow-reducing unit, and thus the shunt system, is closed (shunt switch-off range, 'M.flow off' ⊗). To prevent reading the flow-reduction level in reverse, the rotor stop (3) serves as orientation - in a top view onto the implanted valve as in Fig. 10, the rotor stop is visible on the lower right-hand side.

3.12 APPLICATION OF M.flow Instruments

WARNING

Because of the magnets inside the *M.flow Instruments*, *M.flow Instruments* must not be used in the vicinity of active implants that are influenced by magnetism, such as cardiac pacemakers. When adjusting with the *M.flow Instruments*, the safety distance specified by the manufacturer of the active implant should be observed.

CAUTION

- ▶ The *M.flow Instruments* must not be taken into an MRI facility as this may pose a safety risk to the patient and/or the user.
- ▶ Only *M.flow Instruments* approved for that purpose must be used to determine, change and monitor the flow-reduction level of the *M.flow*.

NOTICE

The *M.blue plus Adjustment Ring* emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.

The *M.flow Instruments* may only be used by trained specialists.

The selected flow-reduction level of the *M.flow* can be determined, adjusted and monitored using the *M.flow Instruments*.

The *M.flow Compass* (Fig. 11) is used for localising and reading the adjustable flow-reducing unit of the *M.flow*.


Fig. 11: *M.flow Compass*

a) open, 1. Scale ring

b) closed, 2. Float gauge

The *M.blue plus Adjustment Ring* (Fig. 12) can be used to set the flow-reduction level of the *M.flow* from 0 to 10 as well as the shunt switch-off range 'M.flow off' (⊗).

Fig. 12: *M.blue plus Adjustment Ring*

The flow-reduction level for the adjustable flow-reducing unit of the *M.flow* can be changed before or after implantation. It is pre-set to 0 (minimum flow-reducing effect) by the manufacturer.

The following steps must be performed in order to set the flow-reduction level of the *M.flow*:

1. Localisation

CAUTION

- ▶ The *M.flow Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus Ad-*

justment Ring should not be placed in the immediate vicinity of the **M.flow Compass** when determining the flow-reduction level. We recommend a minimum distance of 30 cm.

CAUTION

- ▶ Swelling of the skin may make adjustment difficult for a few days after surgery. If the valve setting cannot be checked conclusively using the **M.flow Compass**, we recommend checking it with an imaging method.
- ▶ The **M.flow Compass** should be placed as centrally as possible over the valve, as the determined flow-reduction level may otherwise be incorrect.
- ▶ Possible air pockets in the **M.flow Compass** do not affect the compass function.

When the **M.flow Compass** is opened, a circular cut-out becomes visible, which is used to localise the valve in the patient's head as centrally as possible using the index finger (Fig. 13).

Fig. 13: Localising the valve

The directional markings (proximal and distal) must be aligned according to the inlet and outlet connector of the valve. The directional markings show the flow direction.

2. Verification procedure

In order to determine the selected flow-reduction level, the **M.flow Compass** is then closed again. The float gauge should now be centred within the designated circular marking by moving the **M.flow Compass** (Fig. 14). Once the float gauge is centred, the currently set flow-

reduction level of the **M.flow** can be read off via the line mark on the float gauge (Fig. 15).


Fig. 14: Determining the flow-reduction level with the **M.flow Compass**

The scale ring is marked with two scales (Fig. 15).

The grey marked adjustment range from 0 (lowest flow-reduction level) to 10 (highest flow-reduction level) on the outer scale applies to the flow-reduction level of the **M.flow**. The setting '**M.flow off**' (⊗) (shunt switch-off range) causes the flow to be blocked. No cerebrospinal fluid shunting is possible with this setting.


Fig. 15: **M.flow Compass** scale ring

1. Outer: Scale for the **M.flow** from 0 to 10 with additional shunt switch-off range '**M.flow off**' (⊗) (flow-reduction level in the illustrated example: 0).
2. Inner: Scale for the **M.blue** gravitational unit from 0 to 40 cmH_2O (the opening pressure of the **M.blue** gravitational unit in the illustrated example: 28 cmH_2O)

3. Adjustment process

WARNING

The position 'M.flow off' (⊗), represented by a circle with a cross, completely stops the flow through the valve. Danger of underdrainage!

CAUTION

Between flow-reduction level 0 (minimal flow-reducing effect) and the setting 'M.flow off' (⊗) (valve impermeable/cerebrospinal fluid shunting prevented) there is a stop that prevents the extreme positions from being accidentally set incorrectly.

When adjusting, care must be taken to change the valve by a maximum of 4 levels per adjustment process; otherwise, errors may occur.

Example: The flow-reduction level is to be changed from 0 to 8. The correct method is an adjustment in two stages: an initial adjustment from 0 to 4 and subsequently from 4 to 8.

Example: The flow-reduction level is to be changed from 'M.flow off' (⊗) to 0. The correct method is an adjustment in three stages: an initial adjustment from 'M.flow off' (⊗) to 8, from 8 to 4 and subsequently from 4 to 0.

3a. Adjusting with the M.blue plus Adjustment Ring

In order to adjust the flow-reduction level of the M.flow, the M.flow Compass is opened, but without changing the position of the scale ring. The M.blue plus Adjustment Ring is now inserted into the scale ring in such a manner that the line marking points to the desired value on the scale of the scale ring. (Fig. 16)


Fig. 16: Inserting the M.blue plus Adjustment Ring:
1. M.blue plus Adjustment Ring, 2. Scale ring


Fig. 17: Setting the flow-reduction level of the M.flow in the example illustrated to 3.

By applying slight pressure with the index finger to the valve membrane located in the centre of the M.blue plus Adjustment Ring under the skin, the rotor brake is released and the flow-reduction level of the M.flow is changed to the desired value (Fig. 18).


Fig. 18: Adjustment using the M.blue plus Adjustment Ring

The flow-reducing unit of the M.flow is equipped with a feedback mechanism.

Due to the valve housing design, targeted pressure exerted on the valve membrane produces an audible acoustic signal (a clicking sound) and/or palpable resistance as soon as the rotor brake has been released. The valve thus shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released, the rotor is once again adjustment-proof. Although the click caused by releasing the rotor brake is easily audible before implantation, it may be considerably reduced

after implantation and the filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patient or by using a stethoscope.

3b. Adjustment using the *M.blue plus Adjustment Assistant*

Alternatively, the *M.blue plus Adjustment Assistant* can be used to adjust the flow-reduction level. To do this, insert the *M.blue plus Adjustment Assistant* into the *M.blue plus Adjustment Ring* aligned to the desired value and press it with your index finger (Fig. 19).

Fig. 19: *M.blue plus Adjustment Assistant* inserted in the *M.blue plus Adjustment Ring*

4. Checking after adjustment

After adjusting the flow-reduction level, it is recommended to check the setting. To do this, proceed as in Points 1 and 2. Should the measured value not correspond with the required flow-reduction level, the adjustment process should be repeated. To do this, start again at Point 3.

5. Checking and adjusting using the *M.flow Checkmate*

The *M.flow Checkmate* (Fig. 20) is supplied sterile and can be re-sterilised. The *M.flow Checkmate* can be used to change and check the flow-reduction level adjustment before and during valve implantation directly on the *M.flow*. To determine the flow-reduction level, the *M.flow Checkmate* is centred over the *M.flow*. The *M.flow Checkmate* automatically aligns itself over the valve. The flow-reduction level can be read from the direction of the proximal catheter (at the inlet connector). If the flow-reduction level is to be adjusted, the *M.flow Checkmate* is placed centrally onto the *M.flow*. When do-

ing so, the desired flow-reduction level must point in the direction of the proximal catheter (leading towards the valve). By slightly pressing the *M.flow Checkmate* onto the valve, the rotor brake in the *M.flow* is released and the flow-reduction level set.

When adjusting, make sure that the valve is changed by a maximum of 4 flow-reduction levels per adjustment; otherwise, errors may occur (see Point '3. Adjustment process').

Fig. 20: *M.flow Checkmate*

3.13 SYSTEM COMPONENTS

Combination with shunt components

The product *M.flow* can be safely combined with our range of implantable shunt components. We recommend using the Christoph Miethke GmbH & Co. KG products in combination with *M.flow*.

Reservoirs

The use of shunt systems with a *Reservoir* makes it possible to puncture the *Reservoir* and check the shunt function.

Any puncture of the *Reservoir* should be performed as perpendicular as possible to the reservoir surface with a sharp cannula of max. \varnothing 0.9 mm. The silicone membrane can be punctured 30 times without any restrictions. A stable titanium base prevents the reservoir base from being pierced.

Thanks to an integrated check valve in the *SPRUNG RESERVOIR* and the *CONTROL RESERVOIR*, cerebrospinal fluid can be pumped in the shunt direction, thus making it possible to perform a shunt function check of the distal part of the shunting system as well as the ventricular catheter. If the reservoir can be easily squeezed out, then the valve and distal catheter are patent. If the reservoir fills quickly after being squeezed out, then the ventricular catheter is patent. The volume per pump cycle is approx. 0.2 ml for the normal reservoirs and

approx. 0.1 ml for the paediatric *Reservoirs* with check valve. Access to the ventricular catheter is closed off during pumping.

WARNING

The integration of the *M.flow* into the shunt increases the flow resistance depending on the setting and therefore also the behaviour of the reservoir pump test. Depending on the flow-reduction level, increased resistance can be observed when the silicone membrane is pressed in.

CAUTION

Frequent pumping of the *Reservoir* can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.

Burrhole Deflectors

Because of the tight fit on the *Ventricular Catheter*, the *Burrhole Deflector* makes it possible to choose the length of catheter penetrating into the skull prior to implantation. The *Ventricular Catheter* is deflected at a right angle in the burr hole (Ch. 4.5.5).

Tube systems

The *M.flow* can be ordered as an individual valve unit and can be combined with MIETHKE catheters (internal diameter 1.2 mm, external diameter 2.5 mm). These catheters are designed for MIETHKE valves or MIETHKE shunt components with a connector outer diameter of 1.9 mm.

The catheters are made from silicone and contain barium sulphate. This ensures that the X-ray image is clearly visible. For a new connection of catheters and connectors, the catheters must be carefully secured with a ligature to the valve's connectors or shunt components.

3.14 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

- ▶ The *M.flow* and the entire shunt system are safely able to resist positive and negative pressures up to 100 cmH₂O occurring during and after surgery.
- ▶ Magnetic resonance imaging (MRI) examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations

can be performed without danger or impairment to the valve function. The *M.flow* is MR Conditional. The catheters are MR Safe. *Burrhole Reservoirs*, *Prechambers*, *Burrhole Deflectors* and *Titanium Connectors* are MR Conditional. The stylet is MR Unsafe.

CAUTION

- ▶ In MRI imaging *M.flow* creates artefacts that are larger than the valve itself.
- ▶ If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, an adjustment of the valve cannot be ruled out.

The conditions for MRI safety of the products are found on our website:

<https://miethke.com/downloads/>

4 PROPERTIES OF THE *M.flow*

4.1 PRODUCT DESCRIPTION

4.1.1 *M.flow* VARIANTS

The *M.flow* has an adjustable flow-reducing unit that can be used to adjust the flow resistance, as well as a shunt switch-off function.

Adjustable flow-reducing unit	
Flow-reduction levels	0 - 10
Shunt switch-off function	⊗ (' <i>M.flow</i> off')

4.1.2 SCOPE OF DELIVERY

Box contents	Number
<i>M.flow</i> in sterile packaging	1
Instructions for use for <i>M.flow</i>	1
Patient card	1
Patient label	10

4.1.3 STERILITY

WARNING

The product must not be used if the sterile packaging or the product is damaged or after the expiry date.

The products are sterilised with steam under strictly controlled conditions. The respective expiry date is printed on the packaging.

4.1.4 REPEATED USE AND RESTERILISATION

WARNING

The product must not be resterilised or reprocessed in any other way as the safe functioning and sterility of the product cannot be guaranteed.

In order to minimise the risk of infection, products that have already been implanted in a patient must not be reimplanted either on the same or a different patient.

4.1.5 PRODUCT LIFETIME

The technical lifetime of the medical device *M.flow* is 15 years. Patient-specific circumstances or known side effects such as biological deposits or infection can shorten the lifetime and make revision necessary. In a favourable case, the real lifetime can also exceed the expected time span.

4.1.6 SINGLE-USE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the product *M.flow*. No guarantee can be assumed for the functional safety of resterilised products.

4.1.7 PRODUCT CONFORMITY

The product meets current regulatory requirements.

The requirements stipulate the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the medical device should therefore be recorded in the patient's medical records to ensure complete traceability.

4.2 IMPORTANT SAFETY INFORMATION

4.2.1 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety

instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ The product must not be used if the sterile packaging or the product is damaged or after the expiry date.
- ▶ Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be carefully read and understood before the product is used for the first time.
- ▶ Prior to use, it is essential to check the product for completeness and integrity.

4.2.2 COMPLICATIONS, SIDE EFFECTS, PRECAUTIONS AND RESIDUAL RISKS

The following complications and side effects can occur in conjunction with the *M.flow*:

- ▶ Headaches, dizzy spells, mental confusion, vomiting in cases of possible leakage from the shunt system and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood components in the cerebrospinal fluid
- ▶ Allergic reaction/intolerance to product materials
- ▶ Overdrainage/underdrainage
- ▶ Noise development

Violent external shocks (accident, fall, etc.) may put the integrity of the shunt system at risk.

As a precaution, a physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the product *M.flow*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis)/allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture, skin irritation
- ▶ Local shunt irritation/allergic reaction

4.2.3 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the competent state authority.

4.3 PATIENT EDUCATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product. (see Ch. 4.1.5, 4.2, 4.5, 4.7).

4.4 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place.

4.4.1 TRANSPORT

Transport conditions

Temperature range for transport	-30 °C ... +40 °C
---------------------------------	-------------------

4.4.2 STORAGE

Storage conditions

Temperature range for storage	room temperature
-------------------------------	------------------

4.5 USE OF THE PRODUCT

4.5.1 INTRODUCTION

The *M.flow* is a valve with an adjustable flow-reducing unit and an additional shunt switch-off function. The flow-reducing unit allows the adjustment of the flow resistance. Its shunt switch-off function allows cerebrospinal fluid shunting to be fully halted.

The *M.flow* is used for shunting cerebrospinal fluid into the peritoneum in the treatment of hydrocephalus. Valves and Reservoirs are placed in suitable positions along the course of the shunt.

4.5.2 SAFETY AND WARNING NOTICES

WARNING

- ▶ The correct flow direction must be strictly observed as per the arrows. Incorrect implantation will result in flow prevention and underdrainage.
- ▶ The *M.flow* may only be used in combination with another MIETHKE valve in a shunt system, as the *M.flow* alone does not completely prevent backflow.
- ▶ The combination of the *M.flow* with another valve leads to significantly higher pressure values in the shunt system, depending on the selected flow-reduction level.
- ▶ The position '*M.flow* off' (⊗), represented by a circle with a cross, completely stops the flow through the valve. Danger of underdrainage!
- ▶ When pumping the *M.flow*, an adjustment of the flow-reduction level by pressing the valve membrane several times cannot be ruled out. Checking the flow-reduction level after the pump cycle is therefore required! The patient should be properly informed about this risk.
- ▶ The integration of the *M.flow* into the shunt increases the flow resistance increases depending on the setting and therefore also the behaviour of the reservoir pump tests. Depending on the flow-reduction level, increased resistance can be observed when the silicone membrane is pressed in.
- ▶ The adjustable MIETHKE valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e. g. underneath heavily scarred tissue).
- ▶ During placement, it must be ensured that the adjustable flow-reducing valve is positioned a maximum of 10 mm below the skin surface to make it easier to locate the valve as well as to read and adjust the flow-reduction level. If the implantation site is unfavourably chosen or the skin above the valve is too thick (skin and tissue thickness greater than 10 mm), it may no longer be possible to locate the valve as well as read and adjust the flow-reduction level. The valve then works with a fixed flow-reduction level.

**WARNING**

- ▶ **Silicone is extremely electrostatic. Care must be taken to avoid the catheters coming into contact with dry cloths, talcum powder or rough surfaces. Clinging particles could lead to tissue reactions.**
- ▶ **When using sharp instruments, care should be taken to avoid cuts and scratches in the silicone elastomer.**
- ▶ **It must be ensured that the ligature is not tightened excessively. Damage may result in a loss of integrity of the shunt and thus necessitate a revision.**
- ▶ **The catheters should only be blocked with an atraumatic clamp and not directly behind the shunt component (e.g. valve, reservoir or connector) as they might be damaged otherwise.**

**CAUTION**

- ▶ **Frequent pumping of the Reservoir can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.**
- ▶ **If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, an adjustment of the valve cannot be ruled out.**
- ▶ **In MRI imaging M.flow creates artefacts that are larger than the valve itself.**

4.5.3 REQUIRED MATERIALS

The *M.flow* is designed so that it can be safely used with the shunt components described in Chapter 3.13.

Catheters with an inner diameter of 1.2 mm and an outer diameter of 2.5 mm should be used for connection. In any case, catheters have to be carefully secured to the connectors of the shunt components with a ligature. Any kinks in the catheter must be avoided.

4.5.4 PREPARING FOR IMPLANTATION**Checking the sterile packaging**

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Preoperative valve test

The *M.flow* should be vented and checked for permeability before implantation. The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological saline solution. The valve is patent if saline solution can be extracted (Fig. 21).


Fig. 21: Patency test

**WARNING**

Contaminants in the solution used for testing can impair the product's performance.

**WARNING**

Pressurisation with a single-use syringe should be avoided both at the proximal and the distal end (Fig. 22).

Fig. 22: Avoidance of pressurisation

4.5.5 PERFORMING THE IMPLANTATION

Positioning of the *Ventricular Catheter*

Several surgical techniques are available for placement of the *Ventricular Catheter*. The required skin incision should be semicircular, with the closed side towards the shunting catheter or by a straight skin cut. If a *Burrhole Reservoir* – or a *SPRUNG RESERVOIR* – is used, the skin incision should not be located right above the reservoir. To avoid cerebrospinal fluid leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole.

The *M.flow* is available in a range of different configurations: If a *Burrhole Reservoir* – or a *SPRUNG RESERVOIR* – is used, then the *Ventricular Catheter* is implanted first. Once the stylet has been removed, the patency of the *Ventricular Catheter* can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the *Burrhole Reservoir* or *SPRUNG RESERVOIR* is connected and the connection is secured with a ligature.

When using a shunt system with a *CONTROL RESERVOIR* or a *Prechamber*, a *Burrhole Deflector* is included. The *Burrhole Deflector* is used to adjust the length of the catheter to be implanted and to position it inside the ventricle. The *Ventricular Catheter* is deflected by 90° and the *CONTROL RESERVOIR* or *Prechamber* is put into place. The position of the *Ventricular Catheter* should be inspected after the procedure by imaging (such as CT or MRI).

Valve positioning

For ventriculoperitoneal shunting, a location behind the ear is suitable as an implantation position. The implantation height has no influence on the valve function.

The adjustable valves should be touching the bone or the periosteum since pressure must be exerted on the valve during any later adjustment. A large arch-shaped or a small straight skin cut with a pocket for the valve should be made. The catheter is pushed forward from the burr hole to the selected valve implantation location, shortened if necessary, and secured to the *M.flow* by ligature. The valve should not be located directly under the skin incision. The valve unit has an arrow in the flow direction (arrow towards distal or downwards). The surface of the valve with the arrow markings points to the outside.

WARNING

- ▶ The correct flow direction must be strictly observed as per the arrows. Incorrect implantation will result in flow prevention and underdrainage.
- ▶ During placement, it must be ensured that the adjustable flow-reducing valve is positioned a maximum of 10 mm below the skin surface to make it easier to locate the valve as well as to read and adjust the flow-reduction level. If the implantation site is unfavourably chosen or the skin above the valve is too thick (skin and tissue thickness greater than 10 mm), it may no longer be possible to locate the valve as well as read and adjust the flow-reduction level. The valve then works with a fixed flow-reduction level.

Positioning of the *Peritoneal Catheter*

The access site for the *Peritoneal Catheter* is left to the surgeon's discretion. For example, it can be used in a paraumbilical application or applied at the level of the epigastrium. Likewise, various surgical techniques are available for placing the *Peritoneal Catheter*.

The recommendation is to pull the *Peritoneal Catheter* from the valve to the intended position using a subcutaneous Tunneller, if necessary with the aid of an auxiliary incision. The *Peritoneal Catheter*, usually securely attached to the valve, has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the *Peritoneal Catheter* (shortened, if necessary) is pushed forward into the open space of the abdominal cavity.

4.5.6 POSTOPERATIVE TEST

Postoperative valve test

The *M.flow* has been constructed as a reliably functioning unit. The valve test can be performed by flushing, pressure measurement or pumping via a reservoir.

4.6 EXPLANTATION AND DISPOSAL

4.6.1 EXPLANTATION

The explantation of the *M.flow* should be performed according to the state of the art and in compliance with medical practice.

4.6.2 DISPOSAL

***M.flow* and shunt components**

Products and product parts not used in the implantation or that have been surgically removed must be disposed of correctly as potentially infectious material in accordance with medical practice as well as respective regional laws and regulations.

4.7 TECHNICAL INFORMATION


4.7.1 MATERIALS THAT COME INTO CONTACT WITH BODILY TISSUE/FLUIDS

Quantitative material properties

All implantable parts of a shunt system with direct contact to body tissues or fluids are biocompatible and are composed as follows:

- ▶ Titanium alloy TiAl6V4 (ASTM F136)
- ▶ Titanium grade 1 (ASTM F67)
- ▶ Silicone elastomer (long-term implantable)

4.7.2 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>M.flow</i>
Medical purpose	Shunting of cerebrospinal fluid (CSF)
Sterilizability	Cannot be resterilised
Storage	Store in a clean and dry place at room temperature
For single use only	
Sketch with outer dimensions:	
	

4.8 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxxx indicates the identifier of the responsible notified body.
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Serial number
	Quantity
	Sterilised using steam
	Do not resterilise
	Do not reuse
	Double sterile barrier system
	Do not use if package is damaged and consult instructions for use
	Keep dry
	Consult instructions for use / electronic instructions for use
	Caution

Symbol	Explanation
	Non-pyrogenic
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians.
	MR Conditional

5 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions. You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany