proGAV®

Ⓛ Gebrauchsanweisung | Ⓜ Instructions for use | Ⓛ Mode d'emploi
Ⓛ Instrucciones de manejo | Ⓛ Istruzioni per l'uso
Ⓛ 사 용 방 법 | Ⓛ 使用适应症

Ⓛ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000.

A paper copy will be provided to you upon request at no additional cost.

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CONTENT

INDICATION	21
TECHNICAL DESCRIPTION OF THE VALVE	21
FUNCTION OF THE VALVE	22
PHYSICS BACKGROUND	23
SELECTING THE APPROPRIATE SHUNT	24
ACCESSORIES	25
<i>proGAV Verification Tool</i>	25
<i>proGAV Masterdisc</i>	25
<i>proGAV Compass</i>	25
<i>proGAV Adjustment Tool</i>	26
<i>proGAV Adjustment Disc</i>	26
<i>proGAV Check-mate</i> 	27
ADJUSTING THE ADJUSTABLE DP-UNIT	27
READING THE PRESSURE SETTING FROM AN X-RAY IMAGE	30
POSSIBLE SHUNT COMPONENTS	30
TUBE SYSTEMS	31
SURGICAL PROCEDURE	31
PREOPERATIVE VALVE TEST	32
PRESSURE-FLOW CHARACTERISTICS	32
CONTRAINDICATIONS	33
INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS	33
RE-IMPLANTATION	33
SAFETY MEASURES	33
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	34
POSTOPERATIVE VALVE TEST	34
FUNCTIONAL SAFETY	34
ADVERSE REACTION	34
STERILIZATION	34
RESTERILISATION	34
NOTE ON THE INSTRUCTIONS FOR USE	34
REQUIREMENTS OF THE MDD 93/42/EEC	34
MEDICAL PRODUCTS CONSULTANT	35
GENERAL INFORMATION	35
VARIATIONS	36

INDICATION

The Miethke *proGAV* is intended to shunt cerebrospinal fluid (CSF) from the lateral ventricles of the brain into the peritoneum.

TECHNICAL DESCRIPTION OF THE VALVE

The *proGAV* is a posture dependent hydrocephalus valve. It comprises a ball-cone valve with adjustable opening pressure (adjustable DP-unit) and a fixed pressure gravitational unit. In this way, optimal CSF drainage is ensured for each individual patient in any body position.

Fig. 1: The *proGAV* is a combination of an adjustable DP-unit and a gravitational unit.

The function of both valve units is described as follows. Fig. 2 shows a schematic cross section of the *proGAV*:

The adjustable DP-unit is composed of a solid titanium body with a well-ried ball-cone valve (1) integrated in its proximal part. A bow spring (2) defines the opening pressure of the ball-cone valve. The pretensioning of the spring, and thus the valve opening pressure, can be adjusted by turning a rotor (3), with the valve implanted under the patient's skin.

The gravitational unit contains a tantalum ball (4), which defines the opening pressure of this valve, and a sapphire ball (5), which ensures the precise closure of the valve. A connector under the silicone catheter (6) always allows a later disconnection of the catheter if necessary.


Fig. 2: Schematic cross section of the *proGAV*

FUNCTION OF THE VALVE

The opening pressure of the *proGAV* is composed of the opening pressure of the adjustable DP-unit and the opening pressure of the gravitational unit.

Horizontal position

When the patient is lying down, the gravitational unit is always open and therefore does not present any resistance to the fluid flow. Hence, the opening pressure of the *proGAV* is defined by the adjustable DP-unit. The operational principle of the adjustable DP-unit is illustrated in figures 3 and 4.


Fig. 3: Gravitational unit in horizontal position

In fig. 4a, the ball-cone valve is closed. The drainage is blocked. In fig. 4b, the adjustable DP-unit is shown in the open condition. The patient's IVP is increased and the spring force, which otherwise keeps closed the ball-cone valve, is overcome. The closing ball moves out of the cone and a gap opens to allow drainage.


Fig. 4: Adjustable DP-unit
a) closed b) open

Vertical position

As soon as the patient moves into an upright position, the gravitational unit closes, the opening pressure of the *proGAV* is significantly increased and CSF drainage is blocked (fig. 5a). Only when the sum of the IVP and the hydrostatic pressure exceeds the opening pressure of the *proGAV*, drainage will be possible again (fig. 5b). The opening pressure of the *proGAV* in the vertical position is the sum of the opening pressures of both the adjustable DP-unit and the gravitational unit (weight of the tantalum ball).


Fig. 5: Gravitational unit in vertical position
a) closed, b) open

PHYSICS BACKGROUND

The intraventricular pressure is positive in a healthy human in a horizontal position. To adjust this pressure through shunt drainage, one has to choose the appropriate pressure range, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig. 6).

The ventricular pressure in a healthy human in a vertical position becomes slightly negative. To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure.


Fig. 6: Calculating the intraventricular pressure for horizontal and vertical body position

SELECTING THE APPROPRIATE SHUNT

The *proGAV* is a position-dependent shunt, meaning the opening pressure changes depending on the position of the patient. To choose the suitable *proGAV* for an individual patient, one opening pressure is set for the horizontal position (patient lying down), and one for the vertical position (patient standing upright).

Horizontal position

The opening pressure for the horizontal position is defined by the adjustable DP-unit. The pressure level should be chosen according to the clinical situation and indication. The unit can be adjusted to a pressure setting between 0 cmH₂O and 20 cmH₂O.

Vertical position

The opening pressure of the *proGAV* for the vertical position is calculated by the sum of the opening pressure of the adjustment and of the gravitational unit. The selection of the gravitational unit depends on the activity and the abdominal pressure (adiposity).

The following opening pressure ranges for the gravitational unit are possible, the pressure range selected can be checked postoperatively on X-ray image:

Opening pressure for vertical posture	Coding of gravitational unit
10 cmH ₂ O	small, no ring
15 cmH ₂ O	large, no ring
20 cmH ₂ O	large, 1 ring
25 cmH ₂ O	large, 2 rings
30 cmH ₂ O	large, 3 rings
35 cmH ₂ O	large, 4 rings


Fig. 7: X-ray image of the gravitational unit
a) large, 1 ring = 20 cmH₂O,
b) small = 10 cmH₂O

ACCESSORIES

Caution: Due to magnets inside the *proGAV Tools*, do not use the *proGAV Tools* nearby pacemakers.

Further more do not use the *proGAV Tools* nearby MRI scanner, since there is a danger of damaging the MRI-scanner.

proGAV Verification Tool

The *proGAV Verification Tool* is used for reading the valve opening pressure setting. First, the *proGAV Verification Tool* should be placed on the valve. The marking (3) on the *proGAV Verification Tool* must be in line with the proximal (ventricular) catheter. The tool contains two magnets. As soon as the button on the *proGAV Verification Tool* is pushed (2) the magnets in the tool align with the magnets in the valve. The opening pressure is shown on the scale (1).

Fig. 8a: *proGAV Verification Tool*

proGAV Masterdisc

The *proGAV Verification Tool* can be easily checked by using the *proGAV Masterdisc* before measuring the opening pressure of the valve. On the disc the positions 0, 5, 10, 15 and 20 cmH₂O are indicated. If the *proGAV Verification Tool* is placed on the disc the opening pressure shown by the instrument should be aligned to the value of the *proGAV Masterdisc*.

Example: The *proGAV Verification Tool* is put on the *proGAV Masterdisc* so that the marking on the instrument is in line with the value 10 cmH₂O on the *proGAV Masterdisc*.

The *proGAV Verification Tool* should indicate the value of 10 cmH₂O.

Fig.8b: *proGAV Masterdisc*

proGAV Compass

Alongside the *proGAV Compass* there is an additional device for measuring the adjusted opening pressure. The *proGAV Compass* can be used to locate the valve when palpation is not possible. The *proGAV Compass* is set on the skin over the implanted valve and moved in a circling motion until the internal disc is fixed over the valve. The opening pressure corresponds to the value indicated towards the direction of the ventricular catheter.


Fig. 9: *proGAV Compass*

Airbubbles inside the *proGAV Compass* do not affect its functionality!

proGAV Adjustment Tool

The *proGAV Adjustment Tool* is used for adjusting the valve opening pressure. First the intended pressure setting is selected at the knurled dial (1), the opening pressure is shown on a scale (2). Then the *proGAV Adjustment Tool* is placed central on the valve. By pushing the button (3), the adjustment tip (4) appears, the brake is decoupled, the rotor turns and the adjusted pressure is achieved.

The marking (5) on the *proGAV Adjustment Tool* has to point towards the proximal catheter (leading to the ventricle).

Warning note: If the implantation site is poorly selected or if the skin over the valve is too thick, adjustment of the *proGAV* can be difficult or sometimes impossible. The valve works then as gravitational valve with constant opening pressure. The patient should be informed about the risk.


Fig. 10a: *proGAV Adjustment Tool*

proGAV Adjustment Disc

The *proGAV Adjustment Disc* offers an other option to adjust the pressure setting.

The *proGAV Adjustment Disc* is placed centrally over the valve. The desired pressure setting should be aligned with the proximal catheter (leading to the ventricle). By pressing down the *proGAV Adjustment Disc* on the valve, the brake is decoupled and the opening pressure of the *proGAV* is changed.

Ensure that the pressure setting is changed no more than 8 cmH₂O. (see Chapter „Adjusting the adjustable DP-unit“)

Use a longer tip (size L) for thicker skin (see Fig. 10b).


Fig. 10b: *proGAV Adjustment Disc*

a) Size S b) Size M c) Size L

proGAV Check-mate  0297

The *proGAV Check-mate* is delivered sterile and is intended to be re-sterilised. It is possible to change and to verify an applied pressure setting on the valve directly. To verify the actual pressure setting the *proGAV Check-mate* has to be put centrally over the valve. The *proGAV Check-mate* will immediately start to move. If it remains stable, the pressure setting can be read in alignment to the inlet connector.

To adjust a new pressure setting, the *proGAV Check-mate* has to be placed centrally over the valve. The new pressure setting has to point towards the proximal catheter (leading to the ventricle). By pressing down slightly the *proGAV Check-mate*, the brake of the valve is decoupled, the rotor turns and the opening pressure of the *proGAV* is changed.

Please be aware that the steps for changing the pressure setting should not be more than 8 cmH₂O per step.

Fig. 11: *proGAV Check-mate*

ADJUSTING THE ADJUSTABLE DP-UNIT

Please verify specifically before using any tool for verifying or adjusting the opening pressure:

for the adjustable DP-unit, use only *proGAV Tools*

and

for the adjustable gravitational-unit, use only *proSA Tools*

Each *proGAV* is calibrated under strict quality control procedures. The presetting of the adjustable DP-unit is 5 cmH₂O, but it must be checked before implantation. The setting is changed in the following steps:

1. Locating the valve

The valve is located under the skin (fig. 12).


Fig. 12: Locating the valve

2. Verifying the opening pressure

The *proGAV Verification Tool* is positioned centrally on the valve. The marking on the instrument has to point towards the proximal catheter (leading to the ventricle). The arrow on the bottom side of the *proGAV Verification Tool* indicates the direction of CSF-flow. The button is pushed and the pressure setting is read (fig. 13).


Fig. 13: Verifying the pressure setting
a) correctly, b) incorrectly

Caution: Placing the *proGAV Verification Tool* in a non-central position on the valve can lead to erroneous readings!

3. Adjusting the opening pressure

3a. Adjustment with the *proGAV Adjustment Tool*.
The *proGAV Adjustment Tool* is set to the required opening pressure by turning the knurled dial of the unit (fig. 14a).

Fig. 14a: Adjustment of the *proGAV Adjustment Tool*

Caution: The new opening pressure setting of the valve must not differ from the measured opening pressure by more than 8 cmH₂O in any one setting (see Nr.4).

The *proGAV Adjustment Tool* is positioned centrally on the valve. The marking on the instrument has to point towards the proximal catheter (leading to the ventricle). As soon as the *proGAV Adjustment Tool* has been positioned centrally on the valve, the button is pushed onto the valve and the adjustment tip appears. This triggers the mechanical decoupling and the valve setting is adjusted to the required opening pressure setting (fig. 14b).


Fig. 14b: Adjusting the pressure setting

Caution: Ensure that the instrument remains close to the valve during the adjustment procedure. When treating patients who have a low tolerance to pain, local anaesthesia (e.g. applied through a plaster) should be considered, in cases where contraindication can be excluded.

If the site of implantation is poorly selected or if the skin over the valve is too thick, an adjustment of the *proGAV* can be difficult or sometimes impossible.

The adjustable DP-unit then behaves like a valve with one opening pressure. Patients should discuss the risks involved with their surgeon.

When adjusting the *proGAV* preoperatively through the packaging, moderate force with the *proGAV Adjustment Tool* is sufficient. **DO NOT USE THE BUTTON!** Strong pressure can cause damage to the housing, which might affect valve function.

Caution: Due to postoperative swelling of the skin the adjustment of the valve setting may be difficult within the first few days.

3b. Adjustment with the proGAV Adjustment Disc

Center the *proGAV Adjustment Disc* over the adjustable DP-unit of the *proGAV* and align the desired pressure setting (b) on top of the *proGAV Adjustment Disc* in direction of the ventricular catheter (c), see Fig. 15a.


Fig. 15a: Adjustment with the *proGAV Adjustment Disc*

For changing the opening pressure, press down the *proGAV Adjustment Disc* and release (Fig. 15b). Do not press and turn. Finally, remove the *proGAV Adjustment Disc* and confirm the setting with the *proGAV Verification Tool*.


Fig. 15b: Positioning the *proGAV Adjustment Disc*

4. Verifying the adjustment

After the adjustment, the valve opening pressure is measured again, as described in step 2. If the pressure measured now differs from the intended pressure level, the adjustment procedure has to be repeated from step 3.

During the adjustment the opening pressure of the adjustable DP-unit should not be changed more than 8 cmH₂O per adjustment procedure. Example: Opening pressure is to be changed from 3 to 18 cmH₂O. With only one adjustment procedure the rotor would turn in a counter clockwise direction (short way) and would stop at the position 0 cmH₂O (fig. 16a).


Fig. 16: Rotor rotation during adjustment
a) false b) correct

The correct adjustment is in 2 steps: Adjustment from 3 to 11, and from 11 to 18 cmH₂O. The rotor turns now in a clockwise direction (fig. 17b).

Caution: If the pressure configuration of the valve cannot be determined with complete certainty by the *proGAV Verification Tool*, the use of imaging techniques is recommended (excluding MRI: danger of artefacts).

MRI examinations must be performed at field strengths no greater than 3.0 tesla.

READING THE PRESSURE SETTING FROM AN X-RAY IMAGE

The pressure setting of the *proGAV* should be measured with the *proGAV Verification Tool*, but if there is any discrepancy between the desired adjustment setting and the setting that is read by the *proGAV Verification Tool*, then a radiographic confirmation (x-ray) can be performed to confirm the actual valve setting. The magnets are visible in the image as white dots. The direction of the rotor tip indicates the pressure setting. The rotor tip can take any position outside region (Fig. 17). Thus, the opening pressure of the *proGAV* can be adjusted in increments of 1 cmH₂O between 0 and 20 cmH₂O. In order to avoid misidentification of the adjusted opening pressure in the X-ray image, the valve is flattened on one side.


Fig. 17: Schematic X-ray image


Fig. 18: Radiograph pressure setting 0 cmH₂O

POSSIBLE SHUNT COMPONENTS

The *proGAV* is available with different shunt accessories. These variations are comprised of a variety of components, which are described briefly introduced below:

The *Burrhole Reservoir* is positioned in the cranial burrhole. It allows measurement of intraventricular pressure, injection of drugs and extraction of CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. A special *Burrhole Reservoir* is the *SPRUNG RESERVOIR*. An additional new feature of this reservoir is that CSF can be flushed towards the *proGAV* because of a non-return valve in the bottom of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The *prechamber* is positioned on the cranium. It allows measurement of intraventricular pressure, injection of drugs, extraction of CSF and palpatory ventricle inspection. Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of maximum diameter 0,9 mm. 30 punctures are possible without any restrictions. A special prechamber is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the *proGAV* because of a non-return valve in the proximal inlet of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

Warning note: Frequent pumping can lead to overdrainage and thus to pressure conditions outside the normal physiological range. The patient should discuss the risks (involved) with their surgeon.

Tight tolerancing of the deflector ensures a good fit with the ventricular catheter. By adjusting the deflector (prior to implantation) the length of catheter penetrating into the skull can be optimised. The ventricular catheter is "deflected" at a right angle in the burrhole (see chapter "Variations").

TUBE SYSTEMS

The *proGAV* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. Regardless, the catheters must be carefully fixed, with a ligature, to the valve connectors. It is essential that kinks in the catheter are avoided. The included catheters have virtually no effect on the Pressure-flow characteristics.

SURGICAL PROCEDURE

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burrhole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *proGAV* is available in different shunt variations:

When using a *proGAV SHUNTSYSTEM with Burrhole Reservoir* or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the *Burrhole Reservoir* is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir.

The *proGAV SHUNTSYSTEM with prechamber* or *CONTROL RESERVOIR* comes with a *deflector*. This deflector is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the valve

The adjustable DP-unit of the *proGAV* is supplied with a factory setting of 5 cmH₂O. This opening pressure can be set to a different value prior to implantation (see chapter "Adjusting the *proGAV*"). The gravitational unit of the *proGAV* is a posture-dependent valve. Therefore, care must be taken that the unit is implanted parallel to the body axis. A suitable implantation site is behind the ear.

After the skin incision and tunneling under the skin, the catheter is pushed forward, from the burrhole to the intended shunt implantation site. The catheter is shortened, if necessary, and secured at the *proGAV* with a ligature. The shunt should not be located directly under the skin incision. The valve is marked with an arrow pointing in the direction of flow (arrow pointing to distal or downward).

Warning note: The adjustable DP-unit must be placed over a hard boney surface and should not be implanted within an area, which makes finding and feeling the valve more difficult (e. g. under a scar).

Frequent pumping can lead to overdrainage and thus to unphysiological pressure conditions. The patient should be informed about the risk.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction or transrectally at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to the intended position of the catheter. The peritoneal catheter, which is usually securely attached to the *proGAV*, has an open distal end, but no wall slits. Following the exposure of, and the entry into, the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

PREOPERATIVE VALVE TEST


Fig. 19: Patency test

The *proGAV* can be filled by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see Fig. 19).

Caution: Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end.

Contaminations in the solution used for the test can impair the product's performance.

PRESSURE-FLOW CHARACTERISTICS

Horizontal position

The following diagrams show the pressure-flow characteristics of the adjustable DP-unit of the *proGAV* for the pressure settings 0, 10 and 20 cmH₂O.








Fig. 20: Pressure-flow characteristics for some pressure settings of the adjustable DP-unit

Vertical position

The opening pressure of the *proGAV* in the vertical position is the sum of the opening pressure of the adjustable DP-unit and the gravitational unit.

The following diagrams show the pressure-flow characteristics for some pressure settings in the vertical body position.

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.

CONTRAINDICATIONS

The adjustable DP-unit should not be implanted within an area which makes locating and sensing the valve more difficult (e. g. under a scar). The valve should lie on the periosteum or the bone to make an adjustment after implantation possible.

INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS

The *proGAV* with gravitational unit should not be used under any circumstances in conjunction with hydrostatic valves, as this can bring about abnormally high ventricular pressure outside of the normal physiological range. Hydrostatic valves allow for changes in hydrostatic pressure in the drainage system caused by changes in position. If in doubt, please contact the medical product consultants at Christoph Miethke GmbH & CO. KG.

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried out without endangering or impairing the functionality of the shunt. The *proGAV* is MR Conditional (ASTM-F2503-13). All components are visible via X-ray. The provided catheters are MRI Safe. Reservoirs, deflectors and connectors are MR Conditional.

Warning note: When using a magnetic field and simultaneous pressing on the valve an adjustment of the valve cannot be excluded. The *proGAV* will produce artifacts or signal-intensity voids in MR images larger than the physical size of the device.

POSTOPERATIVE VALVE TEST

The *proGAV* has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the unit if a shunt system with a prechamber or a *Burrhole Reservoir* is used. Valve tests can be carried out by flushing or pressure measurements.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. Still, it cannot be excluded that the shunt system needs to be replaced for technical or medical reasons. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

Warning note for carriers of pacemakers: Due to the implantation of a *proGAV* the function of a pacemaker can be affected.

ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

STERILIZATION

The products are sterilized with steam under closely monitored conditions. The double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILISATION

The functional safety and reliability of resterilized products cannot be guaranteed, therefore re-sterilisation is not recommended.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co.KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz
Michaela Funk-Neubarth


Christoph Miethke GmbH & Co. KG
Ulanenweg 2
D-14469 Potsdam
Phone: +49(0) 7000 6438453 or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap Platz
D-78532 Tuttlingen
Phone: +49 (0) 7461 95-0
Fax: +49 (0) 7461 95-26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042
AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>proGAV</i>
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the valve with its external dimensions: 	

VARIATIONS

The *proGAV* is available as a single valve or as a shunt system comprising various components.

proGAV

proGAV SHUNTSYSTEM

proGAV SHUNTSYSTEM with SPRUNG RESERVOIR or with (pediatric) Burrhole Reservoir

proGAV SHUNTSYSTEM with CONTROL RESERVOIR or with (pediatric) prechamber

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany