proGAV[®]

 Gebrauchsanweisung |  Instructions for use |  Mode d'emploi

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 www.miethke.com

 This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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CONTENT

INDICATION	21
TECHNICAL DESCRIPTION OF THE VALVE	21
FUNCTION OF THE VALVE	22
PHYSICS BACKGROUND	23
SELECTING THE APPROPRIATE SHUNT	24
READING THE PRESSURE SETTING FROM AN X-RAY IMAGE	24
ACCESSORIES	25
APPLICATION OF INSTRUMENTS	27
POSSIBLE SHUNT COMPONENTS	30
TUBE SYSTEMS	30
SURGICAL PROCEDURE	30
PREOPERATIVE VALVE TEST	31
POSTOPERATIVE VALVE TEST	31
PRESSURE-FLOW CHARACTERISTICS	32
CONTRAINDICATIONS	33
INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS	33
RE-IMPLANTATION	33
SAFETY MEASURES	33
FUNCTIONAL SAFETY AND COMPATIBILITY	
WITH DIAGNOSTIC PROCEDURES	33
ADVERSE REACTION	34
STERILISATION	34
RESTERILISATION	34
NOTE ON THE INSTRUCTIONS FOR USE	34
REQUIREMENTS OF THE MDD 93/42/EEC	34
MEDICAL PRODUCTS CONSULTANT	34
GENERAL INFORMATION	34
VARIATIONS	35

INDICATION

The Miethke proGAV is intended to shunt cerebrospinal fluid (CSF) from the lateral ventricles of the brain into the peritoneum.

TECHNICAL DESCRIPTION OF THE VALVE

The proGAV is a posture dependent hydrocephalus valve. It comprises a ball-cone valve with adjustable opening pressure (adjustable DP-unit) and a fixed pressure gravitational unit. In this way, optimal CSF drainage is ensured for each individual patient in any body position.

Fig. 1: The proGAV is a combination of an adjustable DP-unit and a gravitational unit.

The function of both valve units is described as follows. Fig. 2 shows a schematic cross section of the proGAV:

The adjustable DP-unit is composed of a solid titanium body with a well-tried ball-cone valve (1) integrated in its proximal part. A bow spring (2) defines the opening pressure of the ball-cone valve. The pretensioning of the spring, and thus the valve opening pressure, can be adjusted by turning a rotor (3), with the valve implanted under the patient's skin.

The gravitational unit contains a tantalum ball (4), which defines the opening pressure of this valve, and a sapphire ball (5), which ensures the precise closure of the valve. A connector under the silicone catheter (6) always allows a later disconnection of the catheter if necessary.


Fig. 2: Schematic cross section of the proGAV

FUNCTION OF THE VALVE

The opening pressure of the proGAV is composed of the opening pressure of the adjustable DP-unit and the opening pressure of the gravitational unit.

Horizontal position

When the patient is lying down, the gravitational unit is always open and therefore does not present any resistance to the fluid flow.

Hence, the opening pressure of the proGAV is defined by the adjustable DP-unit. The operational principle of the adjustable DP-unit is illustrated in figures 3 and 4.


Fig. 3: Gravitational unit in horizontal position

In fig. 4a, the ball-cone valve is closed. The drainage is blocked. In fig. 4b, the adjustable DP-unit is shown in the open condition. The patient's IVP is increased and the spring force, which otherwise keeps closed the ball-cone valve, is overcome. The closing ball moves out of the cone and a gap opens to allow drainage.


Fig. 4: Adjustable DP-unit

a) closed b) open

Vertical position

As soon as the patient moves into an upright position, the gravitational unit closes, the opening pressure of the proGAV is significantly increased and CSF drainage is blocked (fig. 5a). Only when the sum of the IVP and the hydrostatic pressure exceeds the opening pressure of the proGAV, drainage will be possible again (fig. 5b). The opening pressure of the proGAV in the vertical position is the sum of the opening pressures of both the adjustable DP-unit and the gravitational unit (weight of the tantalum ball).


Fig. 5: Gravitational unit in vertical position

a) closed, b) open

PHYSICS BACKGROUND

The intraventricular pressure is positive in a healthy human in a horizontal position. To adjust this pressure through shunt drainage, one has to choose the appropriate pressure range, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig. 6).

The ventricular pressure in a healthy human in a vertical position becomes slightly negative. To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure.


Fig. 6: Calculating the intraventricular pressure for horizontal and vertical body position

SELECTING THE APPROPRIATE SHUNT

The proGAV is a position-dependent shunt, meaning the opening pressure changes depending on the position of the patient.

To choose the suitable proGAV for an individual patient, one opening pressure is set for the horizontal position (patient lying down), and one for the vertical position (patient standing upright).


Fig. 7: Schematic X-ray image

Horizontal position

The opening pressure for the horizontal position is defined by the adjustable DP-unit. The pressure level should be chosen according to the clinical situation and indication. The unit can be adjusted to a pressure setting between 0 cmH₂O and 20 cmH₂O.

Vertical position

The proGAV opening pressure for the vertical body position is calculated from the sum of the opening pressure of the adjustable differential pressure unit and the gravitational unit. Patient size, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure for this position (see pressure level recommendations at <https://www.miethke.com/downloads/>).


Fig. 8: Radiograph pressure setting 0 cmH₂O

The following opening pressure ranges for the gravitational unit are possible, the pressure range selected can be checked postoperatively on X-ray image:

READING THE PRESSURE SETTING FROM AN X-RAY IMAGE

The pressure setting of the proGAV (*adjustable DP-unit*) should be measured with the proGAV Verification Tool, but if there is any discrepancy between the desired adjustment setting and the setting that is read by the proGAV Verification Tool, then a radiographic confirmation (x-ray) can be performed to confirm the actual valve setting. The magnets are visible in the image as white dots. The direction of the rotor tip indicates the pressure setting. The rotor tip can take any position outside region (Fig. 7). Thus, the opening pressure of the proGAV can be adjusted in increments of 1 cmH₂O between 0 and 20 cmH₂O.

In order to avoid misidentification of the adjusted opening pressure in the X-ray image, the valve is flattened on one side.

Opening pressure for vertical posture	Coding of gravitational unit
10 cmH ₂ O	
15 cmH ₂ O	
20 cmH ₂ O	
25 cmH ₂ O	
30 cmH ₂ O	
35 cmH ₂ O	


Fig. 9: X-ray image of the gravitational unit

a) large, 1 ring = 20 cmH₂O

b) small = 10 cmH₂O

ACCESSORIES

Instruments used for the adjustable differential pressure unit:

proGAV Verification Tool

proGAV Masterdisc

proGAV Compass

proGAV Adjustment Tool

proGAV Adjustment Disc

proGAV Checkmate

NOTE

There are different testing and adjusting instruments for locating, testing and adjusting the adjustable differential pressure unit.

proGAV Verification Tool

The proGAV Verification Tool is used for reading the valve opening pressure setting. First, the proGAV Verification Tool should be placed on the valve. The marking (3) on the proGAV Verification Tool must be in line with the proximal (ventricular) catheter. The tool contains two magnets. As soon as the button on the proGAV Verification Tool is pushed (2) the magnets in the tool align with the magnets in the valve. The opening pressure is shown on the scale (1).

CAUTION

The proGAV Tools emit a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.

CAUTION

Because of the magnets inside the proGAV Tools, the instruments must not be used in the vicinity of cardiac pacemakers. Furthermore, there is risk of damage to MRI instruments in their vicinity. Therefore, the proGAV Tools must not be used in such locations!

Fig. 10a: proGAV Verification Tool

proGAV Masterdisc

The *proGAV Verification Tool* can be easily checked by using the *proGAV Masterdisc* before measuring the opening pressure of the valve. On the disc the positions 0, 5, 10, 15 and 20 cmH₂O are indicated. If the *proGAV Verification Tool* is placed on the disc the opening pressure shown by the instrument should be aligned to the value of the *proGAV Masterdisc*.

Example: The *proGAV Verification Tool* is put on the *proGAV Masterdisc* so that the marking on the instrument is in line with the value 10 cmH₂O on the *proGAV Masterdisc*.

The *proGAV Verification Tool* should indicate the value of 10 cmH₂O.

Fig. 10b: *proGAV Masterdisc*

proGAV Compass

Alongside the *proGAV Compass* there is an additional device for measuring the adjusted opening pressure. The *proGAV Compass* can be used to locate the valve when palpation is not possible. The *proGAV Compass* is set on the skin over the implanted valve and moved in a circling motion until the internal disc is fixed over the valve. The opening pressure corresponds to the value indicated towards the direction of the ventricular catheter.

Fig. 11: *proGAV Compass*

NOTE

Airbubbles inside the *proGAV Compass* do not affect its functionality!

proGAV Adjustment Tool

The *proGAV Adjustment Tool* is used for adjusting the valve opening pressure. First the intended pressure setting is selected at the knurled dial (1), the opening pressure is shown on a scale (2). Then the *proGAV Adjustment Tool* is placed central on the valve. By pushing the button (3), the adjustment tip (4) appears, the brake is decoupled, the rotor turns and the adjusted pressure is achieved.

The marking (5) on the *proGAV Adjustment Tool* has to point towards the proximal catheter (leading to the ventricle).


Fig. 12a: *proGAV Adjustment Tool*

proGAV Adjustment Disc

The *proGAV Adjustment Disc* offers an other option to adjust the pressure setting.

The *proGAV Adjustment Disc* is placed centrally over the valve. The desired pressure setting should be aligned with the proximal catheter (leading to the ventricle). By pressing down the *proGAV Adjustment Disc* on the valve, the brake is decoupled and the opening pressure of the *proGAV* is changed.

Ensure that the pressure setting is changed no more than 8 cmH₂O. (see Chapter „Adjusting the adjustable DP-unit“)

Use a longer tip (size L) for thicker skin (see Fig. 12b).


Fig. 12b: proGAV Adjustment Disc
a) Size S b) Size M c) Size L

proGAV Checkmate

The proGAV Checkmate is delivered sterile and is intended to be re-sterilised. It is possible to change and to verify an applied pressure setting on the valve directly. To verify the actual pressure setting the proGAV Checkmate has to be put centrally over the valve. The proGAV Checkmate will immediately start to move. If it remains stable, the pressure setting can be read in alignment to the inlet connector.

To adjust a new pressure setting, the proGAV Checkmate has to be placed centrally over the valve. The new pressure setting has to point towards the proximal catheter (leading to the ventricle). By pressing down slightly the proGAV Checkmate, the brake of the valve is decoupled, the rotor turns and the opening pressure of the proGAV is changed.

Please be aware that the steps for changing the pressure setting should not be more than 8 cmH₂O per step.

Fig. 13: proGAV Checkmate

APPLICATION OF INSTRUMENTS

The proGAV Tools may only be used by trained specialist personnel. The selected pressure level of proGAV can be determined, adjusted and monitored using the proGAV Tools.

CAUTION

It is absolutely essential that only proGAV Tools are used to determine, change and monitor the opening pressure of the differential pressure unit of the proGAV.

proGAV

If the proGAV is combined with an adjustable gravitational unit (proSA), proSA Tools must be used to adjust and check the adjustable pressure level of the proSA valve.

proSA

NOTE

If the implantation site is unfavourable or the skin above the valve is too thick, the adjustable differential pressure unit may no longer be able to be adjusted. The valve then works with unchanged pressure levels. The patient should be made aware of this risk.

**NOTE**

Swelling of the skin can make adjustment difficult for a few days after the operation. During preoperative adjustment through the sterile packaging, the *Adjustment Tool* is only placed on the valve. Do not press the push-button. Excessive pressure with the *Adjustment Tool* can damage

**NOTE**

If the valve setting cannot be checked conclusively using the *Verification Tool* or the *Verification Compass*, we recommend carrying out a check using an imaging method.


Fig. 15: Verifying the pressure setting
a) correctly, b) incorrectly

Each proGAV is calibrated under strict quality control procedures. The presetting of the adjustable DP-unit is 5 cmH₂O, but it must be checked before implantation. The setting is changed in the following steps:

1. Locating the valve

The valve is located under the skin (fig. 14).

Fig. 14: Locating the valve

2. Verifying the opening pressure

The proGAV Verification Tool is positioned centrally on the valve. The marking on the instrument has to point towards the proximal catheter (leading to the ventricle). The arrow on the bottom side of the proGAV Verification Tool indicates the direction of CSF-flow. The button is pushed and the pressure setting is read (fig. 15).

**WARNING**

Placing the proGAV Verification Tool in a non-central position on the valve can lead to erroneous readings!

3. Adjusting the opening pressure

**WARNING**

The new opening pressure setting of the valve must not differ from the measured opening pressure by more than 8 cmH₂O in any one setting (see Nr.4).

3a. Adjustment with the proGAV Adjustment Tool. The proGAV Adjustment Tool is set to the required opening pressure by turning the knurled dial of the unit (fig. 16a).

Fig. 16a: Adjustment of the proGAV Adjustment Tool

The proGAV Adjustment Tool is positioned centrally on the valve. The marking on the instrument has to point towards the proximal catheter (leading to the ventricle). As soon as the proGAV Adjustment Tool has been positioned centrally on the valve, the button is pushed onto the valve and the adjustment tip appears. This

triggers the mechanical decoupling and the valve setting is adjusted to the required opening pressure setting (fig. 16b).


Fig. 16b: Adjusting the pressure setting

WARNING

Ensure that the instrument remains close to the valve during the adjustment procedure.

NOTE

When treating patients who have a low tolerance to pain, local anaesthesia (e.g. applied through a plaster) should be considered, in cases where contraindication can be excluded.

3b. Adjustment with the proGAV Adjustment Disc

Center the *proGAV Adjustment Disc* over the adjustable DP-unit of the *proGAV* and align the desired pressure setting (b) on top of the *proGAV Adjustment Disc* in direction of the ventricular catheter (c), see Fig. 17a.


Fig. 17a: Adjustment with the *proGAV Adjustment Disc*

For changing the opening pressure, press down the *proGAV Adjustment Disc* and release (Fig. 17b). Do not press and turn. Finally, remove the *proGAV Adjustment Disc* and confirm the setting with the *proGAV Verification Tool*.


Fig. 17b: Positioning the *proGAV Adjustment Disc*

4. Verifying the adjustment

After the adjustment, the valve opening pressure is measured again, as described in step 2. If the pressure measured now differs from the intended pressure level, the adjustment procedure has to be repeated from step 3.

During the adjustment the opening pressure of the adjustable DP-unit should not be changed more than 8 cmH₂O per adjustment procedure. Example: Opening pressure is to be changed from 3 to 18 cmH₂O. With only one adjustment procedure the rotor would turn in a counter clockwise direction (short way) and would stop at the position 0 cmH₂O (fig. 18a).


Fig. 18: Rotor rotation during adjustment
a) false b) correct

The correct adjustment is in 2 steps: Adjustment from 3 to 11, and from 11 to 18 cmH₂O. The rotor turns now in a clockwise direction (fig. 18b).

If the pressure configuration of the valve cannot be determined with complete certainty by the *proGAV Verification Tool*, the use of imaging techniques is recommended (excluding MRI: danger of artefacts).

POSSIBLE SHUNT COMPONENTS

The *proGAV* is available with different shunt accessories. These variations are comprised of a variety of components, which are described briefly introduced below:

The *Burhole Reservoir* is positioned in the cranial burrhole. It allows measurement of intraventricular pressure, injection of drugs and extraction of CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. A special *Burhole Reservoir* is the *SPRUNG RESERVOIR*. An additional new feature of this reservoir is that CSF can be flushed towards the *proGAV* because of a non-return valve in the bottom of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The *prechamber* is positioned on the cranium. It allows measurement of intraventricular pressure, injection of drugs, extraction of CSF and palpatory ventricle inspection. Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of maximum diameter 0.9 mm. 30 punctures are possible without any restrictions. A special prechamber is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the *proGAV* because of a non-return valve in the proximal inlet of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

WARNING

Frequent pumping can lead to overdrainage and thus to pressure conditions outside the normal physiological range. The patient should discuss the risks (involved) with their surgeon.

Tight tolerancing of the deflector ensures a good fit with the ventricular catheter. By adjusting the deflector (prior to implantation) the length of catheter penetrating into the skull can be optimised. The ventricular catheter is "deflected" at a right angle in the burrhole (see chapter "Variations").

TUBE SYSTEMS

The *proGAV* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. Regardless, the catheters must be carefully fixed, with a ligature, to the valve connectors. It is essential that kinks in the catheter are avoided. The included catheters have virtually no effect on the Pressure-flow characteristics.

SURGICAL PROCEDURE

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burrhole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *proGAV* is available in different shunt variations:

When using a *proGAV SHUNTSYSTEM with Burhole Reservoir* or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the *Burhole Reservoir* is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir.

The proGAV SHUNTSYSTEM with prechamber or CONTROL RESERVOIR comes with a deflector. This deflector is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the valve

The adjustable DP-unit of the proGAV is supplied with a factory setting of 5 cmH₂O. This opening pressure can be set to a different value prior to implantation (see chapter "APPLICATION OF INSTRUMENTS"). The gravitational unit of the proGAV is a posture-dependent valve. Therefore, care must be taken that the unit is implanted parallel to the body axis. A suitable implantation site is behind the ear.

After the skin incision and tunneling under the skin, the catheter is pushed forward, from the burrhole to the intended shunt implantation site. The catheter is shortened, if necessary, and secured at the proGAV with a ligature. The shunt should not be located directly under the skin incision. The valve is marked with an arrow pointing in the direction of flow (arrow pointing to distal or downward).

NOTE

The gravitational unit of proGAV is position-dependent. You must therefore ensure that the gravitational unit is implanted parallel to the body axis.

Therefore, if a shunt system in which the valve has been pre-fitted with a *Burrhole Reservoir* is being used, only the occipital access should be used.

WARNING

The adjustable DP-unit must be placed over a hard boney surface and should not be implanted within an area, which makes finding and feeling the valve more difficult (e. g. under a scar).

CAUTION

The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction or at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to the intended position of the catheter. The peritoneal catheter, which is usually securely attached to the proGAV, has an open distal end, but no wall slits. Following the exposure of, and the entry into, the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

PREOPERATIVE VALVE TEST


Fig. 19: Patency test

The proGAV can be filled by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see Fig. 19).

WARNING

Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end. Contaminations in the solution used for the test can impair the product's performance.

POSTOPERATIVE VALVE TEST

The proGAV has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the unit if a shunt system with a prech-

amber or a *Burhole Reservoir* is used. Valve tests can be carried out by flushing or pressure measurements.

PRESSURE-FLOW CHARACTERISTICS

Horizontal position

The following diagrams show the pressure-flow characteristics of the adjustable DP-unit of the proGAV for the pressure settings 0, 10 and 20 cmH₂O.








Fig. 20: Pressure-flow characteristics for some pressure settings of the adjustable DP-unit

Vertical position

The opening pressure of the proGAV in the vertical position is the sum of the opening pressure of the adjustable DP-unit and the gravitational unit.

The following diagrams show the pressure-flow-characteristics for some pressure settings in the vertical body position.

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.

CONTRAINDICATIONS

The adjustable DP-unit should not be implanted within an area which makes locating and sensing the valve more difficult (e. g. under a scar). The valve should lie on the periost or the bone to make an adjustment after implantation possible.

INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS

The proGAV with gravitational unit should not be used under any circumstances in conjunction with hydrostatic valves, as this can bring about abnormally high ventricular pressure outside of the normal physiological range. Hydrostatic valves allow for changes in hydrostatic pressure in the drainage system caused by changes in position. If in doubt, please contact the medical product consultants at Christoph Miethke GmbH & CO. KG.

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successfull decontamination of the device cannot be reached without functional degradation.

SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system.

FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons. These medical devices are able to resist positive and negative pressures up to 100 cmH₂O during and after implantation. These medical devices have to be stored in a clean and dry environment at all times.

Nuclear magnetic resonance examinations up to a field strength of 3 Tesla or computed tomography examinations can be performed without risk or impairment to the valve function. The valve is MR Conditional. Supplied catheters are MR Safe. Reservoirs, deflectors and connectors are MR Conditional. The conditions for MR compatibility of the products can be found on our website:
<https://www.miethke.com/downloads>

WARNING

When using a magnetic field and simultaneous pressing on the valve an adjustment of the valve cannot be excluded. The proGAV will produce artifacts or signal-intensity voids in MR images larger than the physical size of the device.

CAUTION

For carriers of pacemakers: Due to the implantation of a proGAV the function of a pacemaker can be affected.

ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

STERILISATION

The products are sterilised with steam under closely monitored conditions. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILISATION

The functional safety and reliability of resterilised products cannot be guaranteed, therefore resterilisation is not recommended.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the implanted valve should therefore be recorded in the patient's medical records and patient data card to ensure complete traceability.

Translations of these instructions for use into additional languages can be found on our website: <https://www.miethke.com/downloads>


MEDICAL PRODUCTS CONSULTANTS

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical products consultants as contacts for all product-related questions:

You can contact our medical devices consultants via:

Tel. +49 331 62083-0
info@miethke.com

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	proGAV
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the valve with its external dimensions:	
	

VARIATIONS

The proGAV is available as a single valve or as a shunt system comprising various components.

proGAV

proGAV SHUNTSYSTEM

proGAV SHUNTSYSTEM with
SPRUNG RESERVOIR or with
(pediatric) Burrhole Reservoir

proGAV SHUNTSYSTEM with
CONTROL RESERVOIR or with
(pediatric) prechamber

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany