proSA[®]

DE Gebrauchsanweisung | **EN** Instructions for use | **FR** Mode d'emploi | **ES** Instrucciones de manejo
 www.miethke.com

US This Instructions for Use is NOT intended for United States users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
TABLE OF CONTENTS

INDICATION	23
TECHNICAL DESCRIPTION	23
PHYSICS BACKGROUND	24
FUNCTION OF THE VALVE	25
SELECTION OF APPROPRIATE PRESSURE LEVEL	27
READING THE PRESSURE SETTING FROM AN X-RAY IMAGE	27
ACCESSORIES	29
APPLICATION OF INSTRUMENTS	32
POSSIBLE SHUNT COMPONENTS	35
TUBE SYSTEMS	36
SURGICAL PROCEDURE	36
VALVE TEST	37
PRESSURE-FLOW CHARACTERISTICS	38
PRECAUTIONARY MEASURES AND CONTRAINDICATIONS	38
INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS	38
RE-IMPLANTATION	39
FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	39
ADVERSE REACTIONS AND INTERACTIONS	39
STERILISATION	39
NOTE ON THE INSTRUCTIONS FOR USE	39
REQUIREMENTS OF THE MDD 93/42/EEC	39
MEDICAL PRODUCTS CONSULTANT	39
GENERAL INFORMATION	40
VARIANTS	41

INDICATION

The proSA is intended to divert cerebrospinal fluid (CSF) from the lateral ventricles of the brain into the peritoneum. It is the first hydrocephalus valve that is adjustable for the vertical body position, respectively for inclined body position.

TECHNICAL DESCRIPTION

The proSA combines a fixed differential pressure (DP) unit or an adjustable DP unit and an adjustable gravitational unit (fig. 1-3). Each valve unit consists of a robust titanium housing and uses the reliable ball-cone principle.

A spiral spring (1) defines the fixed opening pressure of the fixed DP unit. The sapphire ball (2) ensures the precise closure of the valve.


Fig. 1: Schematic cross-section of the proSA with DP unit

A bow spring (3) defines the opening pressure of the adjustable DP unit. The pretensioning of the spring and thus the opening pressure can be adjusted by turning the rotor (4) with the valve implanted under the patient's skin.

Fig. 2: Schematic cross-section of the adjustable DP unit

The adjustable gravitational unit contains a tantalum weight (5), which is connected to the bow spring (3). The weight keeps the sapphire ball in its position. Depending on the body position, the influence of the tantalum weight changes and therefore the opening pressure. By turning the rotor (4) the pretension of the bow spring can be adjusted and the opening pressure can be changed with the valve implanted under the patient's skin.

Fig. 3: Schematic cross-section of the adjustable gravitational unit

PHYSICS BACKGROUND

The intraventricular pressure (IVP) is positive in a healthy human in the horizontal position. To adjust this pressure through shunt drainage, the appropriate pressure range must be chosen, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig. 4).

The ventricular pressure in a healthy human when moving to the vertical position becomes slightly negative.

To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate for the change in pressure differential. This is equal to the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure. Conventional shunts open immediately as soon as the patient stands up, which can lead to critical overdrainage.


Fig. 4: Pressure ratios for the lying and the upright body position

IVP : Intraventricular pressure

P_L : Opening pressure in horizontal position (DP unit only)

P_s : Opening pressure in vertical position (DP unit + gravitational unit)

P_B : Pressure in the abdominal cavity

P_H : Hydrostatic pressure

Horizontal: $IVP = P_L + P_B$

Vertical: $IVP = P_H - P_s - P_B$

FUNCTION OF THE VALVE

proSA is a posture-dependent hydrocephalus valve. The opening pressure for proSA consists of the opening pressures for the adjustable gravitational unit and the differential pressure unit combined.

Horizontal body position:

proSA with DP unit

The operating principle of the DP unit is illustrated in fig. 5a and 5b. Fig. 5 shows the fixed DP unit in the horizontal position. The ballcone valve is closed and drainage is prevented. If the patient's IVP increases and continues to rise, the spring pressure of the ball-cone unit will be overcome. The sapphire ball will move away from the cone and a gap will open for fluid drainage.


Fig. 5: DP unit
a) closed
b) open

proSA with adjustable DP unit

In fig. 6a the ball-cone unit of the adjustable DP unit is closed, hence the drainage is blocked. Fig. 6b shows the adjustable DP unit in an open condition. The patient's IVP is increased and the spring force, which otherwise keeps the ball-cone unit closed, is overcome. The sapphire ball moves out of the cone and a gap opens up to allow drainage.


Fig. 6: Adjustable DP unit

- a) closed
- b) open

With the proSA in the horizontal position (fig. 7), the weight (1) does not affect the sapphire ball (2). Hence in this configuration, the proSA is open and the force brought to bear by the weight does not counteract the fluid pressure. The ventricular pressure is now controlled solely by the differential pressure valve installed.


Fig. 7: Adjustable gravitational unit in horizontal position

Vertical body position:

In the vertical position, the DP unit and the adjustable gravitational unit of the proSA work together. The total opening pressure is the sum of the opening pressures of the DP unit and the adjustable gravitational unit. If the IVP (intraventricular pressure) of the patient and the hydrostatic pressure exceed this opening pressure, the closing ball of the gravitational unit is pushed away from the cone seal and a gap opens for fluid drainage (fig. 8).

The valve system prevents increases in intraventricular pressure above the physiological range.

To ensure optimal CSF drainage for each individual patient in any body position, the opening pressure of the adjustable gravitational unit can be adjusted from 0 to 40 cmH₂O. The adjustment can be achieved non-invasively using the proSA Tools, eliminating the need for further surgery.


Fig. 8: Gravitational unit in vertical body position

a) closed

b) open

SELECTION OF APPROPRIATE PRESSURE LEVEL

The proSA is a position-dependent valve, that means the opening pressure changes depending on the body position of the patient. Hence, one opening pressure for the horizontal position and one for the vertical position is set to fit the proSA individually to the patient.

Horizontal position

In the lying position, the gravitational unit has no influence on the valve opening pressure. In this position, the valve opening pressure is thus determined exclusively by the differential pressure unit. In this case, the pressure level should be set in accordance with the clinical picture and indications. By default, a differential pressure unit with an opening pressure of 5 cmH₂O is recommended.

Vertical position

The proSA opening pressure for the vertical body position is calculated from the sum of the opening pressure in the differential pressure unit and the adjustable gravitational unit. Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure for the gravitational unit (see pressure level recommendations at <https://www.miethke.com/downloads/>). This is a non-binding recommendation for the attending physician. According to his diagnosis, the physician decides each case independently, without instructions and individually.

READING THE PRESSURE SETTING FROM AN X-RAY IMAGE

Pressure levels for the differential pressure unit

A DP unit or an adjustable DP unit is available for the lying position.

The coding of the DP unit can be identified according to the shape of the valve's housing.

Pressure setting (cmH ₂ O)	Coding
0	
5	
10	
15	

Example: The DP unit with an opening pressure of 5 cmH₂O has a concave proximal part (curved inwards) and a convex distal part (curved outwards).


Fig. 9: Radiographic image of the fixed DP unit: pressure rating 5 cmH₂O

The pressure setting of the adjustable DP unit should be checked with the *proGAV Verification Tool*. If there is any discrepancy between the desired adjustment setting and the setting that is read by the *proGAV Verification Tool*, then a radiographic confirmation can be performed in addition to confirm the actual valve setting. The position of the rotor tip indicates the opening pressure. The rotor can take any position outside the region indicated shown below (fig. 10).


Fig. 10: Schematic X-ray image of adjustable DP unit

a) non-adjustable range

b) valve flattened on one side

The opening pressure of the adjustable DP unit can be adjusted from 0 to 20 cmH₂O, in increments of 1 cmH₂O. To further assist correct identification of the adjusted opening pressure, the housing of the valve has a flat profile on one side (fig. 11).


Fig. 11: X-ray image adjustable DP unit: setting 0 cmH₂O

b) valve flattened on one side

Pressure settings for the adjustable gravitational unit

The pressure setting of the adjustable gravitational unit should be checked with the proSA Verification Tool. If there is any discrepancy between the desired adjustment setting and the setting that is read by the proSA Verification Tool, then a radiographic confirmation can be performed in addition to confirm the actual valve setting. The position of the rotor tip indicates the opening pressure. The rotor can take

any position outside the indicated region shown below (fig. 12).


Fig. 12: Schematic X-ray image of adjustable gravitational unit

a) non-adjustable range

b) cut-out

c) rotor tip

The opening pressure of the adjustable gravitational unit can be adjusted from 0 to 40 cmH₂O, in increments of 1 cmH₂O. To further assist correct identification of the adjusted opening pressure, a cut-out in the casing ring is visible under X-ray imaging (fig. 13).


Fig. 13: X-ray image adjustable gravitational unit: setting

31 cmH₂O

b) cut-out

ACCESSORIES

NOTE

In order to locate, to verify and to adjust either the adjustable DP unit or the adjustable gravitational unit, different verification and adjustment tools are used:

Tools used for the adjustable DP unit:

proGAV Verification Tool
proGAV Masterdisc
proGAV Compass
proGAV Adjustment Tool
proGAV Adjustment Disc
proGAV Checkmate

Tools used for the adjustable gravitational unit:

proSA Verification Tool
proSA Masterdisc
proSA Compass
proSA Adjustment Tool
proSA Adjustment Disc
proSA Checkmate

CAUTION

The proGAV Tools and proSA Tools emit a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.

CAUTION

The magnets inside the proGAV and proSA Tools can alter the function of a pacemaker.

CAUTION

Do not use the proGAV and proSA Tools nearby an MRI scanner, since there is a danger of damaging the MRI scanner due to the magnets inside proGAV and proSA Tools.

The modes of operation of proGAV Tools and proSA Tools are identical. Hence, the instructions on the following pages are classified for both proGAV Tools and proSA Tools.

VERIFICATION TOOL

The Verification Tool is used for reading the valve opening pressure setting (fig. 14). Firstly, it is essential that the Verification Tool is placed centrally over the valve. The notch (3) on the Verification Tool must be in line with the proximal (ventricular) catheter. The tool contains two magnets. As soon as the button (2) on the instrument is pushed the magnets in the tool align with the magnets in the valve. The opening pressure is shown on the scale (1). The arrow on the underside of the Verification Tool should be aligned with the direction of flow.

Fig. 14: proSA Verification Tool

MASTERDISC

The *Verification Tool* can be easily checked by using the *Masterdisc* before measuring the opening pressure of the valve. On the *proGAV Masterdisc* the positions 0, 5, 10, 15, 20 cmH₂O (fig. 15) and on the *proSA Masterdisc* the positions 0, 10, 20, 30 and 40 cmH₂O (fig. 16) are indicated. If the *Verification Tool* is placed on the disc the opening pressure shown by the instrument should be aligned to the value of the *Masterdisc*.

Example: The *Verification Tool* is put on the *Masterdisc* so that the marking on the instrument is in line with the value 10 cmH₂O on the *Masterdisc*. The *Verification Tool* should indicate the value 10 cmH₂O.

Fig. 15: *proGAV Masterdisc*

Fig. 16: *proSA Masterdisc*

COMPASS

Alongside the *Verification Tool* there is an additional device for checking the adjusted opening pressure. The *Compass* can be used to locate the valve when palpation is not possible. The *Compass* is set to the skin above the implanted valve. The opening pressure corresponds to the value indicated towards the direction of the ventricular catheter.

Fig. 17: *proGAV Compass*

Fig. 18: *proSA Compass*

NOTE

Airbubbles inside the compass do not affect its functionality.

ADJUSTMENT TOOL

The *Adjustment Tool* (fig. 19) is used for adjusting the valve opening pressure. First the intended pressure setting is selected at the knurled dial (1), the opening pressure is shown on a scale (2). Then the *Adjustment Tool* is placed centrally on the valve. By pushing the button (3), the adjustment tip (4) appears, the brake is decoupled, the rotor turns and the adjusted pressure is set. The marking (5) on the *Adjustment Tool* has to point towards the proximal catheter (leading to the ventricle).


Fig. 19: proSA Adjustment Tool

ADJUSTMENT DISC

The *Adjustment Disc* offers another option to adjust the pressure setting (fig. 20-21). The *Adjustment Disc* is placed centrally over the valve. The desired pressure setting should be aligned with the proximal catheter (leading to the ventricle). By pressing down the *Adjustment Disc* on the valve, the brake is decoupled and the opening pressure of the valve is changed.

The *proGAV Adjustment Disc* has a blue casing and is available in three different sizes (fig. 20). The scale shows the selectable pressure settings from 0 to 20 cmH₂O.


Fig. 20: proGAV Adjustment Disc

- a) size S
 - b) size M
 - c) size L
- Colour: blue

The *proSA Adjustment Disc* has a red casing and is only available in one size (fig. 21). The scale shows the selectable pressure settings from 0 to 40 cmH₂O.

Fig. 21: proSA Adjustment Disc

- One size only
Colour: red

CHECKMATE

The *Checkmate* is delivered sterile and is intended to be re-sterilised. It is possible to change and to verify an applied pressure setting on the adjustable DP unit and on the adjustable gravitational unit directly in the OR.

To verify the actual pressure setting the *Checkmate* has to be put centrally over the adjustable DP unit or the adjustable gravitational unit. The *Checkmate* will immediately start to move. When it remains stable, the pressure setting can be read in alignment to the inlet connector.

To adjust a new pressure setting, the *Checkmate* has to be placed centrally over the adjustable DP unit or the adjustable gravitational unit. The new pressure setting has to point towards the proximal catheter (leading to the ventricle). By pressing the *Checkmate* down slightly, the brake of the adjustable DP unit or adjustable gravitational unit is decoupled, the rotor turns and the opening pressure of the valve is changed.

*Fig. 22: proGAV Checkmate
Pressure range 0 to 20 cmH₂O
Colour: titanium*

*Fig. 23: proSA Checkmate
Pressure range 0 to 40 cmH₂O
Colour: gold*

APPLICATION OF INSTRUMENTS

The *proSA Tools* may only be used by trained specialist personnel.

NOTE

The modes of operation of *proGAV Tools* and *proSA Tools* are identical. Hence, the instructions on the following pages are shown using the example of *proSA Tools* but the instructions are classified for both *proGAV Tools* and *proSA Tools*.

CAUTION


- **It is absolutely essential that only the *proSA Tools* are used to determine, change and monitor the opening pressure of the gravitational unit of the *proSA*.**

proSA

- **If the *proSA* is combined with an adjustable differential pressure unit (*proGAV*), *proGAV Tools* must be used to adjust and check the adjustable pressure level of the *proGAV* valve.**

proGAV

The pressure setting of the adjustable DP unit and the adjustable gravitational unit can be adjusted pre- and postoperatively.

Each *proSA* is calibrated under strict quality control procedures. The adjustable DP unit is preset to 5 cmH₂O and the adjustable gravitational unit is preset to 20 cmH₂O. Each unit must be checked before implantation.

The adjustment procedure for the adjustable DP unit and the adjustable gravitational unit (in the following named adjustable unit) is outlined in the steps below:

1. Localisation

The adjustable unit is located under the skin (fig. 24).


Fig. 24: Locating the adjustable unit

2. Verifying

The proSA Verification Tool is positioned centrally over the valve (fig. 25). The notch on the instrument must point towards the proximal catheter (leading to the ventricle).

Fig. 25: Place centrally over the valve

The arrow on the bottom side of the Verification Tool indicates the direction of CSF-flow. The button is pushed and the pressure setting is read (fig. 26).

WARNING

Failing to centre the Verification Tool over the valve can lead to erroneous readings!


Fig. 26: Measuring the pressure setting

- a) correctly
- b) incorrectly

3. Adjustment of the opening pressure

NOTE

When adjusting the proSA preoperatively through the packaging, moderate force with the *Adjustment Tool* is sufficient. DO NOT USE THE BUTTON. Strong pressure can cause damage to the housing, which might affect valve function.

CAUTION

Due to postoperative swelling of the skin, adjusting the valve may be difficult within the first few days!

WARNING

Please ensure that the opening pressure is changed by no more than 8 cmH₂O per each adjustment for the adjustable DP unit and no more than 16 cmH₂O per each adjustment for the adjustable gravitational unit.


Fig. 27: Rotor rotation during adjustment

- a) false 1. non-adjustable range

b) correct 1. non-adjustable range

CAUTION

Ensure that the instrument remains close to the valve during the adjustment procedure.

NOTE

When treating patients who have a low pain tolerance, local anaesthesia (e. g. applied through a plaster) should be considered, in cases where contraindication can be excluded.

3a. Adjustment with the Adjustment Tool

The *Adjustment Tool* is set to the required opening pressure by turning the knurled dial of the unit (fig. 28).


Fig. 28: Setting the Adjustment Tool

The *Adjustment Tool* is positioned centrally on the valve. The notch on the instrument and the scale (a) must point towards the proximal catheter (b) (leading to the ventricle), see fig. 29.


Fig. 29: Positioning the Adjustment Tool

As soon as the *Adjustment Tool* has been positioned centrally over the valve, the button is pushed, the adjustment tip emerges and applies pressure to the valve. This triggers the mechanical decoupling of the rotor and the valve setting is adjusted to the required opening pressure setting (fig. 30).


Fig. 30: Adjusting the pressure setting

3b. Adjustment with the Adjustment Disc

Center the *Adjustment Disc* over the valve and align the desired pressure setting (b) on top of the disc in direction of the *Ventricular Catheter* (c), see fig. 31.


Fig. 31: Adjustment with the Adjustment Disc

For changing the opening pressure, press down the *Adjustment Disc* and release (fig. 32). Do not press and turn.


Fig. 32: Press the Adjustment Disc down slightly and release

Finally, remove the *Adjustment Disc* and confirm the setting with the *Verification Tool*.

4. Verifying after adjustment

After the adjustment, the valve opening pressure has to be measured again, as described in step 2. If the pressure measured now differs from the intended pressure level, the adjustment procedure has to be repeated from step 3.

NOTE

If the pressure configuration of the valve cannot be determined with complete certainty by the *Verification Tool*, the use of imaging techniques is recommended (excluding MRI: danger of artefacts).

WARNING

MRI examinations must be performed at field strengths no greater than 3.0 tesla.

CAUTION

If the site of implantation is poorly selected or if the skin over the valve is too thick, an adjustment of the adjustable unit can be difficult or sometimes impossible. The adjustable gravitational unit then behaves like a gravitational unit with a fixed opening pressure for a given position.

POSSIBLE SHUNT COMPONENTS

The proSA is available in combination with different accessory components for the treatment of adult and pediatric hydrocephalus which are briefly described below:

The DP unit is available in pressure configurations of 0, 5, 10 and 15 cmH₂O. The adjustable differential pressure unit is adjustable from 0 to 20 cmH₂O, in increments of 1 cmH₂O.

The burr hole reservoir is positioned in the cranial burr hole. It allows measurement of the intraventricular pressure and a site for the injection of drugs and extraction of CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors.

A special burr hole reservoir is the *SPRUNG RESERVOIR*. As an additional new feature of this reservoir CSF can be flushed towards the valve because of a one-way valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The prechamber or the *CONTROL RESERVOIR* are positioned on the cranium. It allows measurement of the intraventricular pressure and a site for the injection of drugs and extraction of CSF, as well as palpative inspection of the ventricle. Its solid titanium base is highly puncture-resistant. If required, a puncture of the prechamber or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. diameter 0.9 mm. Up to 30 punctures are possible without any restrictions.

A special prechamber is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the valve because of a one-way valve in the proximal inlet of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is prevented during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

**WARNING**

Frequent pumping can lead to overdrainage and thus pressure conditions outside of the normal physiological range. The risks involved should be discussed with the patient.

Tight tolerancing of the *Deflector* ensures a good fit with the ventricular catheter. By adjusting the *Deflector* (prior to implantation) the length of catheter penetrating into the skull can be optimised. The ventricular catheter is "deflected" at a right angle in the burr hole.

TUBE SYSTEMS

The *proSA* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as an individual valve unit with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows the use of catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. Regardless, the catheters must be fixed carefully, with a ligature, to the valve connectors. It is essential that kinks in the catheter are avoided. The included catheters have virtually no effect on the pressure-flow characteristics.

SURGICAL PROCEDURE

If a *proSA* without integrated catheters is being used, the catheters must be secured (proximally and distally) with a ligature. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

proSA is available in a range of different configurations: When using a *proSA* shunt system with burr hole reservoir or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the burr hole reservoir is connected, with the connection secured by a ligature. The skin incision should not be located directly above the reservoir. The *proSA* shunt system with prechamber or *CONTROL RESERVOIR* comes with a *Deflector*. This *Deflector* is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the valve

The *proSA* is supplied with a factory setting of

- ▶ 5 cmH₂O for the adjustable DP unit
- ▶ 20 cmH₂O for the adjustable gravitational unit.

This opening pressure can be set to a different value prior to implantation (see chapter "Application of instruments").

The *proSA* is a posture-dependent valve. Therefore, care must be taken that the unit is implanted parallel to the body axis. A suitable implantation site is behind the ear, whereby the implantation height has no influence on the valve function. The adjustable valve should be touching the bone or the periosteum since pressure must be exerted on the valve during any later adjustment. A large arch-shaped or a small straight skin incision with a pocket for the valve should be made.

After the skin incision and tunneling under the skin, the catheter is pushed forward, from the burr hole to the intended shunt implantation site. The catheter is shortened, if necessary, and secured at the *proSA* with a ligature. The shunt should not be located directly under the skin incision. The valve is marked with an arrow pointing in the direction of flow (arrow pointing to distal or downward).

**NOTE**

proSA is position-dependent. For that reason, care must be taken to implant the valve parallel to the body axis.

Therefore, if a shunt system in which the valve has been pre-fitted with a *Burrhole Reservoir* is being used, only the occipital access should be used.

**WARNING**

The adjustable gravitational unit must be placed over a hard bony surface and should not be implanted within an area that makes locating the valve more difficult (e. g. under a scar).

**NOTE**

To prevent damage to the catheter, clamps with protective coverslips should be used. The catheter should not be tied off directly behind the valve.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction or at the height of the epigastrum. Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to the intended position of the catheter. The peritoneal catheter, which is usually attached securely to the proSA, has an open distal end, but no wall slits. Following the exposure of and the entry into the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

VALVE TEST

Preoperative valve test


Fig. 33: Patency test
a) isotonic sterile saline solution

proSA should be vented before implantation and checked for permeability. The proSA can be filled by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (fig. 33).

**CAUTION**

Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end. Contaminations in the solution used for the test can impair the product's performance.

Postoperative valve test

The proSA has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the unit if a shunt system with a prechamber or a burr hole reservoir is used. Valve tests can be carried out by flushing or pressure measurements.

PRESSURE-FLOW CHARACTERISTICS

Horizontal position:

DP unit

The following diagram shows the pressure-flow characteristics of the DP unit for the opening pressures of 0, 5, 10 and 15 cmH₂O.


Fig. 34: Pressure-flow characteristics for the pressure settings of the fixed DP unit: pressure (cmH₂O), flow (ml/h)

Adjustable DP unit

The following diagrams show the pressure-flow characteristics of the adjustable DP unit for the pressure settings of 0, 5, 10 and 15 cmH₂O.


Fig. 35: Pressure-flow characteristics for the pressure settings of the adjustable DP unit: pressure (cmH₂O), flow (ml/h)

Vertical position:

Adjustable gravitational unit

The total opening pressure in the vertical position is the sum of the opening pressure of the DP unit and the adjustable gravitational unit. The following diagram shows the pressure-flow characteristics of the adjustable gravitational unit as an example for the pressure settings 0, 20 and 40 cmH₂O.


Fig. 36: Pressure-flow characteristics for selected pressure settings of the adjustable gravitational unit: pressure (cmH₂O), flow (ml/h)

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrate reaches 20 ml/h, the opening pressures are approximately 1 to 2 cmH₂O higher.

PRECAUTIONARY MEASURES AND CONTRAINDICATIONS

Patients must be carefully monitored after implantation. Reddening of skin or tightness in the area of the drained tissue may be indications of infections at the shunt system. Symptoms such as headache, dizziness, confusion or vomiting often occur in conjunction with shunt dysfunction. These symptoms and a leakage within the shunt system require the immediate replacement of the affected shunt component or the entire shunt system.

The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e. g. meningitis, ventriculitis, peritonitis, bacteraemia, septicaemia) in the region affected by the implantation.

INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS

The proSA should not be used under any circumstances in conjunction with hydrostatic valves, as this can bring about abnormally high ventricular pressure outside of the normal physiological range. Hydrostatic valves allow for changes in hydrostatic pressure in the drainage system caused by changes in position. If in doubt, please contact the medical product consultants at Christoph Miethke GmbH & Co. KG.

RE-IMPLANTATION

Under no circumstances should products that had previously been implanted in a patient be subsequently reimplanted in the same or another patient, because a successful decontamination of the device cannot be reached without functional degradation.

FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons. These medical devices are able to resist positive and negative pressures up to 100 cmH₂O during and after implantation. These medical devices have to be stored in a clean and dry environment at all times.

Nuclear magnetic resonance (MRI) examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. The valve is MR Conditional. The supplied catheters are MR Safe. Reservoirs, deflectors and connectors are MR Conditional. The conditions for MR compatibility of the products can be found on our website: <https://www.miethke.com/downloads/>

WARNING

If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, it is not possible to rule out valve adjustment. In MRI imaging proSA creates artefacts which are larger than the valve itself.

CAUTION

For people using cardiac pacemakers: It is possible that the function of the pacemaker is influenced by the implantation of a proSA.

ADVERSE REACTIONS AND INTERACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in

very rare cases noise development. Violent shocks from the outside (accident, fall etc.) may put the integrity of the shunt system at risk.

STERILISATION

The products are sterilised with steam under strictly controlled conditions. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used under any circumstances. No guarantee can be given for the functional safety and reliability of resterilised products.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

REQUIREMENTS OF THE MDD 93/42/EEC

The Medical Device Directive requires comprehensive documentation of the fate of medical devices used in humans, especially for implants. The individual implant identification number should therefore be recorded in the patient's medical records and patient card to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website: <https://www.miethke.com/downloads/>

MEDICAL PRODUCTS CONSULTANT


In compliance with the European Directive concerning medical devices (Directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical devices consultants as contacts for all product-related questions.

You can contact our medical devices consultants via:

Tel. +49 331 62083-0

info@miethke.com

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	proSA
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the proSA with its external dimensions:	
proSA	
proSA with adjustable DP unit	

VARIANTS

The proSA is available as a single valve or as a proSA shunt system with DP unit or adjustable DP unit. All figures in this chapter are to scale (1:1).

Fig. 37: proSA

Fig. 39: proSA with distal catheter


Fig. 38: proSA with differential pressure unit

Fig. 40: proSA with adjustable differential pressure unit

proSA shunt system with (pediatric) prechamber

proSA shunt system with (pediatric) burrhole reservoir

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany