proSA

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GB Instructions for use

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Indication	24
Technical description	24
Physics background	25
Function of the <i>proSA</i>	26
Selecting the appropriate <i>proSA</i>	28
Accessories	28
Adjusting the <i>proSA</i>	29
Reading the pressure setting from an X-ray image	32
Possible shunt components	32
Tube systems	33
Interactions with products from other manufacturers	33
Testing valve patency	34
Valve test prior to implantation	34
Pressure-flow-characteristics	37
Test on backflow safety	38
Surgical procedure	38
Contraindications	39
Re-implantations	39
Safety measures	39
Compatibility with diagnostic procedures	40
Postoperative valve tests	40
Functional Safety	40
Sterilisation	40
Resterilisation	40
Medical products consultant	41
Requirements of the MDD 93/42/EEC	41
Note on the instructions for use	41
General information	41
Variants	42

Indication

The Miethke *proSA* is intended to shunt cerebrospinal fluid (CSF) from the lateral ventricles of the brain into the peritoneum.

Technical description

The *proSA* is the first hydrocephalus valve that enables the adjustment of the opening pressure in the upright position, as well as for all angles of elevation between the upright and the supine positions. This opening pressure depends on the body position of the patient. The *proSA* combines an adjustable gravitational unit and a differential pressure unit (Fig. 1).


The *proSA* comprises a robust titanium casing (1) whose proximal end contains a ball-cone valve. The sapphire ball (2) ensures the precise closure of the valve and a bow spring (4) determines the opening pressure of the ball-in-cone valve. The weight (3) is connected to the spring and keeps the sapphire ball in its position. In this way, optimal CSF drainage is ensured for each individual patient in any body position.

The tension of the spring, and thus the valve opening pressure, can be adjusted by turning the rotor (5), with the valve implanted under the patient's skin. The cerebrospinal fluid flows in through the inlet port (6) and flows out by way of the outlet port (7).

- 1 titanium casing
- 2 sapphire ball
- 3 weight
- 4 bow spring
- 5 rotor
- 6 inlet connector
- 7 outlet connector
- 8 spiral spring
- 9 ligature
- 10 catheter

differential pressure unit

adjustable gravitational unit

proSA

Fig. 1: Schematic cross section of the *proSA*

- IVP Intraventricular pressure
- P_{Vli} Opening pressure in horizontal position
(DP-unit only)
- P_{Vst} Opening pressure in vertical position
(DP-unit + gravitational unit)
- P_B Pressure in the abdominal cavity
- P_{Hyd} Hydrostatic pressure

horizontal: $IVP = P_{Vli} + P_B$
 vertical: $IVP = P_{Hyd} - P_{Vst} - P_B$


Fig. 2: Calculating the intraventricular pressure

Physics background

The intraventricular pressure is positive in a healthy human in the horizontal position. To adjust this pressure through shunt drainage, the appropriate pressure range must be chosen, taking into account the abdominal cavity pressure.

The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (Fig. 2).

The ventricular pressure in a healthy human when moving to the vertical position becomes slightly negative.

To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate for the change in pressure differential.


This is equal to the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure.

Conventional shunts open immediately as soon as the patient stands up, which can lead to critical overdrainage.

Function of the proSA

Horizontal body position

With the *proSA* in the horizontal position (Fig. 3), the weight (2) does not affect the sapphire ball. Hence, in this configuration the *proSA* is open, and the force brought to bear by the weight does not counteract the fluid pressure. The ventricular pressure is now controlled solely by the differential-pressure valve installed.

b)
Fig. 3: DP-unit
a) closed and b) open

The operating principle of the DP-unit is illustrated in Fig. 3. a) and b)

Fig. 3 a) shows the DP-unit in the horizontal position. The ball-cone valve is closed and drainage is prevented.

If the patient's IVP increases and continues to rise, the spring pressure of the ball-cone unit will be overcome. The sealing ball will move away from the cone and a gap will open for fluid drainage.

Vertical body position

In the vertical position, the gravitational unit of the *proSA* and the DP-unit work together. The total opening pressure is the sum of the opening pressures of the DP-unit and the *proSA*. If the IVP (intraventricular pressure) of the patient and the hydrostatic pressure exceed this opening pressure, the closing ball of the *proSA* is pushed away from the cone seal and a gap opens for fluid drainage (Fig. 4).

The valve system prevents increases in intraventricular pressure above the physiological range.

To ensure optimal CSF drainage for each individual patient in any body position, the opening pressure of the adjustable gravitational unit can be adjusted from 0 to 40 cmH₂O.

The adjustment can be achieved non-invasively using the *proSA TOOLSET*, preventing the need for further surgery.


Fig. 4: gravitational unit in vertical body position a) closed b) open

Selecting the appropriate proSA

The *proSA* is a position-dependent shunt, meaning the opening pressure changes depending on the position of the patient. We recommend implantation of the gravitational unit in combination with either a simple DP-unit or with the *proGAV* (adjustable DP-unit).

Horizontal body position

In the supine position the opening pressure is determined solely by the DP-valve. The *proSA* does not provide any resistance to flow in this position. The opening pressure of the valve should be chosen according to the clinical needs of the patient.

Vertical body position

The choice of an appropriate *proSA* is determined by the ventricular pressure required for the vertical position. The sum total of the opening pressures of the differential pressure unit and the gravitational unit used determines the ventricular pressure for the patient in the upright position.

Accessories

proSA VERIFICATION TOOL

The *proSA* VERIFICATION TOOL is used for reading the valve opening pressure setting (Fig. 5). Firstly, it is essential that the *proSA* VERIFICATION TOOL is placed centrally over the valve. The notch (3) on the *proSA* VERIFICATION TOOL must be in line with the proximal (ventricular) catheter. The tool contains two magnets. As soon as the button (2) on the instrument is pushed the magnets in the tool align with the magnets in the valve. The opening pressure is shown on the scale (1). The arrow on the underside of the *proSA* VERIFICATION TOOL should be aligned with the direction of flow.

Fig.5: *proSA* VERIFICATION TOOL

proSA ADJUSTMENT TOOL

The *proSA ADJUSTMENT TOOL* is used for adjusting the valve opening pressure (Fig. 6). First the intended pressure setting is selected at the knurled dial (1). The opening pressure is shown on a scale (2). The *proSA ADJUSTMENT TOOL* is then placed centrally over the valve. By pushing the button (3), the adjustment tip (4) appears, the brake is decoupled, the rotor turns and the opening pressure is adjusted accordingly.

The notch (5) on the *proSA ADJUSTMENT TOOL* must point towards the proximal catheter (leading to the ventricle). The arrow on the underside of the *proSA ADJUSTMENT TOOL* should be aligned with the direction of flow.


Fig. 6: *proSA ADJUSTMENT TOOL*

Adjusting the proSA

Each *proSA* is calibrated under strict quality control procedures.

The adjustable gravitational unit is preset to 20 cmH₂O, but must be checked before implantation. The adjustment procedure is outlined in the steps below:

1. The adjustable gravitational unit is located under the skin (Fig. 7).


Fig. 7: Locating the valve

2. The *proSA VERIFICATION TOOL* is positioned centrally over the valve. (Fig. 8a). The notch on the instrument must point towards the proximal catheter (leading to the ventricle).


Fig. 8a: Place centrally over the valve

The button is pushed and the pressure setting is read (Fig. 8b).


Fig. 8b: Measuring the pressure setting
a) correctly, b) incorrectly

Caution:

Failing to centre the proSA VERIFICATION TOOL over the valve can lead to erroneous readings!

3. The proSA ADJUSTMENT TOOL is set to the required opening pressure by turning the knurled dial of the unit (Fig. 9).

Fig. 9:
Setting the proSA ADJUSTMENT TOOL

Please ensure that the opening pressure is changed by no more than 16 cmH₂O per each adjustment (see page 31).

4. The proSA ADJUSTMENT TOOL is positioned centrally on the valve. The notch on the instrument and the scale (a) must point towards the proximal catheter (b) (leading to the ventricle), see Fig. 10.


Fig. 10:
Placing the proSA ADJUSTMENT TOOL

5. As soon as the proSA ADJUSTMENT TOOL has been positioned centrally over the valve, the button is pushed and the adjustment tip appears to apply pressure to the valve. This triggers the mechanical decoupling of the rotor and the valve setting is adjusted to the required opening pressure setting (Fig. 11).

Caution:

Ensure that the instrument remains close to the valve during the adjustment procedure.
When treating patients who have a low tolerance to pain, local anaesthesia (e.g. applied through a plaster) should be considered, in cases where contraindication can be excluded.

Abb. 11: Adjusting the pressure setting

6. After the adjustment, the valve opening pressure is measured again, as described in step 2. If the pressure measured now differs from the intended pressure level, the adjustment procedure has to be repeated from step 3.

Caution: Due to postoperative swelling of the skin the adjustment of the valve setting may be difficult within the first few days.

Warning note:

If the site of implantation is poorly selected or if the skin over the valve is too thick, an adjustment of the *proSA* can be difficult or sometimes impossible. The valve then behaves like a **SHUNTASSISTANT** with one opening pressure for a given position. Patients should discuss the risks involved with their surgeon.

During the adjustment, the opening pressure should not be changed more than 16 cmH₂O at a time.

Example:

Opening pressure is to be changed from 6 to 36 cmH₂O. With only one adjustment step the rotor would turn in a counter clockwise direction (shortest path) and would stop at the position 0 cmH₂O (Fig. 12a).

The correct adjustment is performed in 2 steps: Adjustment from 6 to 22 cmH₂O, and then from 22 to 36 cmH₂O. The rotor turns now in a clockwise direction (Fig. 12b).


Fig. 12: Rotor rotation during adjustment
a) false b) correct

Reading the pressure setting from an X-ray image

The pressure setting of the *proSA* should be measured with the *proSA VERIFICATION TOOL*, but if there is any discrepancy between the desired adjustment setting and the setting that is read by the *proSA VERIFICATION TOOL*, then a radiographic confirmation (X-ray) can be performed to confirm the actual valve setting. The magnets of the rotor are visible in the X-ray image as white dots (Fig. 13). The orientation of the rotor tip indicates the pressure setting. The rotor tip can take any position outside the region indicated shown below (Fig. 13). Thus, the opening pressure of the adjustable gravitational unit can be adjusted from 0 to 40 cmH₂O, in increments of 1 cmH₂O.

To further assist correct identification of the adjusted opening pressure, a cut-out in the casing ring is visible under X-ray imaging (Fig. 14).


Fig. 13: Schematic X-ray image


Fig. 14: X-ray image position 31 cmH₂O

MRI examinations must be performed at field strengths no greater than 3.0 tesla.

Note:

If the pressure configuration of the valve cannot be determined with complete certainty by the *proSA VERIFICATION TOOL*, the use of imaging techniques is recommended (excluding MRI: danger of artefacts).

Possible shunt components

The *proSA* is available in combination with different accessory components, which are briefly described below:

A simple differential pressure unit is available in pressure configurations of 0, 5, 10 and 15 cmH₂O.

The adjustable differential pressure unit is adjustable from 0 to 20 cmH₂O, in increments of 1 cmH₂O.

The *borehole reservoir* is positioned in the cranial borehole. It allows measurement of the intraventricular pressure and a site for the injection of drugs and extraction of CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. A special *borehole reservoir* is the *SPRUNG RESERVOIR*. As an additional new feature of this reservoir CSF can be flushed towards the valve because of a one-way valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The *FLUSHING RESERVOIR* is a pre-chamber and is positioned on the cranium. It allows measurement of the intraventricular pressure and a site for the injection of drugs and extraction of CSF, as well as palpatory inspection of the ventricle. Its solid titanium base is highly puncture-resistant. If required, a puncture of the *FLUSHING RESERVOIR* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. diameter 0,9 mm. Up to 30 punctures are possible without any restrictions.

A special prechamber is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the valve because of a one-way valve in the proximal inlet of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is prevented during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

Warning note:

Frequent pumping can lead to over-drainage and thus pressure conditions outside of the normal physiological range. Patients should discuss the risks involved with their surgeon.

Tight tolerancing of the *deflector* ensures a good fit with the ventricular catheter. By adjusting the *deflector* (prior to implantation) the length of catheter penetrating into the skull can be optimised. The *ventricular catheter* is “deflected” at a right angle in the borehole (see chapter “Variants”).


Tube systems

The *proSA* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as an individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows the use of catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. Regardless, the catheters must be carefully fixed, with a ligature, to the valve connectors. It is essential that kinks in the catheter are avoided.

The included catheters have virtually no effect on the pressure-flow characteristics.

Interactions with products from other manufacturers

The *proSA* should not be used under any circumstances in conjunction with hydrostatic valves, as this can bring about abnormally high ventricular pressure outside of the normal physiological range. Hydrostatic valves allow for changes in hydrostatic pressure in the drainage system caused by changes in position. If in doubt, please contact the medical product consultants at Christoph Miethke GmbH & CO. KG.

isotonic sterile sodium chloride solution

Fig. 15: Patency test

Testing valve patency

The *proSA* can be filled most gently by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see Fig. 15).

Caution: Applying pressure through the single-use syringe should be avoided, both at the proximal and the distal end.

Contamination in the solution used for the test can impair the product's performance.

Valve test prior to implantation

Each *proSA* valve has been tested to ensure that the performance specifications given on the label are always met. The dynamic performance characteristics of the shunt cannot be tested in a static test performed in the operating room.

If the surgeon wishes to verify, prior to implantation, that the shunt meets the specifications given by the manufacturer, the following test can be performed in the operating room:

Caution: Always take care that sterility is maintained and particle contamination is avoided.

Test method

Equipment required for this test:

- sterile fluid reservoir or water bath
- sterile fluid 60-cm water manometer with millimeter grading and three-branch faucet at the base
- sterile syringe, 30 cc to 50 cc
- sterile 5- μ tip filter
- sterile tube adapter
- sterile silicone tube

Setting up the equipment

- Position the manometer and the water bath in such a way that the zero point of the manometer and the fluid level of the water bath are at the same height (see Fig. 16).
- Fill the syringe, with the 5- μ tip filter attached, with sterile water (always use the 5- μ tip filter when topping up the syringe). Remove the tip filter when the syringe is full.
- Connect the syringe, the manometer and the silicone tube together. Use the tube adapter if necessary (see Fig. 16).
- To release all air from the test assembly, turn the three-way faucet as shown in Fig. 17.
- Immerse the silicone tube in the sterile water bath and rinse it with the sterile water from the syringe.

Calibrating the equipment

- a) Turn the three-way faucet as shown in Fig. 18 and fill the manometer to at least 5 cmH₂O.
- b) With the silicone tube immersed in the water bath, turn the three-way faucet so that the syringe is isolated from the manometer (see Fig. 19).
- c) Allow the water column in the manometer to drop.
- d) The water column should stop dropping at the zero point. Adjust the zero point of the manometer to fluid level of the water bath, if necessary.
- e) The manometer has now been calibrated to the zero-level of the water bath. Fixate the manometer to maintain its position in relation to the water bath.


Fig. 16: Test setup

- 1 proSA a) horizontal, b) vertical
- 2 water bath
- 3 constant water level
- 4 silicone tube
- 5 three-way tap
- 6 single-use syringe with syringe filter
- 7 manometer


Fig. 17


Fig. 18


Fig. 19

Test procedure

Please note: During the test the shunt must be submerged in the water bath. The zero point of the manometer must be aligned with the water level of the water bath in order to obtain correct results.

a) Connect the sterile shunt to be tested to the already assembled, sterile test equipment.

b) Turn the three-way faucet as shown in Fig. 18 and fill the manometer to 10 cmH₂O above the expected opening pressure. (Example: For testing a *proSA* in combination with a *miniNAV* having an opening pressure setting of 5 cmH₂O in the horizontal and 25 cmH₂O in the vertical position, the manometer is filled to 15 cmH₂O with the shunt in the horizontal position and to 40 cmH₂O with the shunt in the vertical position.)

c) Turn the three-way faucet as shown in Fig. 17 so that the manometer is isolated

d) Remove all air from the shunt and the test setup by carefully rinsing it through with sterile water from the syringe.

e) Immerse the sterile shunt in the sterile water bath. The distal part of the shunt must be under water to obtain valid test results.

f) Carefully maintain a flow through the shunt and turn the three-way faucet as shown in Fig. 19 to isolate the syringe. As soon as the three-way faucet is in the correct position, the water column should begin to drop.

The syringe is now isolated from the valve and it is not necessary anymore to maintain its flow. Repeat steps b) to f) if the water column fails to drop.

g) Allow the water level in the manometer to drop for 2 to 2.5 minutes. Read the resulting pressure at the manometer.

Test results of preimplantation test

The following table shows results, which should be achieved by this method, for some selected pressure levels:

DP-unit (*proSA* adjusted to „0“)

Opening pressure cmH ₂ O	Acceptable pressure range
0	0-5 cmH ₂ O
5	2-9 cmH ₂ O
10	7-15 cmH ₂ O
15	12-20 cmH ₂ O

proSA (without DP-unit)

Opening pressure cmH ₂ O	Acceptable pressure range
0	0-5 cmH ₂ O
10	2-14 cmH ₂ O
20	8-24 cmH ₂ O
30	15-34 cmH ₂ O
40	22-44 cmH ₂ O

Pressure-flow characteristics

Horizontal position

The following diagrams show the pressure-flow characteristics of the *proSA* in combination with the *miniNAV* for the pressure settings 0, 5, 10 and 15 cmH₂O


Fig. 20: Pressure-flow characteristics for the pressure settings of the DP-unit

Vertical position

The total opening pressure in the vertical position is the sum of the opening pressure of the DP-valve and the *proSA*. The following diagram shows the pressure-flow-characteristics in the vertical body position for several pressure settings.

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrate reaches 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.




Fig. 21: Pressure-flow characteristics for adjustments of the gravitational unit

Testing backflow safety

This test is carried out with the same equipment used for the pre-implantation test. The shunt is carefully filled with sterile saline solution from the syringe before the air is removed from it (Fig. 22). The shunt is connected against the direction of flow (see arrow on the shunt). The outlet of the shunt has to be at the zero level of the manometer. The manometer is filled up to 14 cmH₂O (Fig. 23).

The three-way faucet is used for unblocking the flow to the shunt and blocking the flow to syringe. In this setup, no more than two drops (0.1 cc) per minute should emerge from the proximal part of the shunt (Fig. 24).


Fig. 22 isotonic sterile sodium chloride solution


Fig. 23


Fig. 24

Caution: Be careful to maintain sterility and to avoid particle contamination.

Surgical procedure

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the borehole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *proSA* is available in different shunt variants:

When using a *proSA-SHUNTSYSTEM with borehole reservoir* or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the *borehole reservoir* is connected, with the connection secured by a ligature. The skin incision should not be located directly above the reservoir.

The *proSA-SHUNTSYSTEM with FLUSHING* or *CONTROL RESERVOIR* comes with a *deflector*. This *deflector* is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the *proSA*

The *proSA* is supplied with a factory setting of 20 cmH₂O. This opening pressure can be set to a different value prior to implantation (see chapter “Adjusting the *proSA*”). The *proSA* is a posture-dependent valve. Therefore, care must be taken that the unit is implanted parallel to the body axis. A suitable implantation site is behind the ear.

After the skin incision and tunneling under the skin, the catheter is pushed forward, from the borehole to the intended shunt implantation site. The catheter is shortened, if necessary, and secured at the *proSA* with a ligature. The shunt should not be located directly under the skin incision. The valve is marked with an arrow pointing in the direction of flow (arrow pointing to distal or downward).

Warning note:

The adjustable gravitational unit must be placed over a hard boney surface and should not be implanted within an area that makes locating the valve more difficult (e. g. under a scar).

Frequent pumping can lead to over-drainage and thus pressure conditions outside of the normal physiological range. Patients should discuss the risks (involved) with their surgeon. To prevent damage to the catheter, clamps with protective coverslips should be used. The catheter should not be tied off directly behind the valve.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction or transrectally at the height of the epigastrium.

Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to the

intended position of the catheter. The peritoneal catheter, which is usually securely attached to the *proSA*, has an open distal end, but no wall slits. Following the exposure of, and the entry into the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

Contraindications

The adjustable gravitational unit should not be implanted within an area which makes locating and sensing the valve more difficult (e. g. under a scar). The valve should lie on the peristome or the bone to make an adjustment after implantation possible.

Re-implantation

Shunt components that have already been implanted must not be implanted again into another patient.

Safety measures

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system.

Compatibility with diagnostic procedures

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried out without altering or impairing the functionality of the shunt.

The *proSA* is MRI safe (ASTM-F2503-05). All components are visible via X-ray. The provided catheters are MRI compatible. Reservoirs, deflectors and connectors are MRI safe.

Warning note:

When using a magnetic field and simultaneous pressing on the valve, the possibility of valve adjustment can not be excluded.

Postoperative valve test

The *proSA* has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the unit if a shunt system with a prechamber or a *BOREHOLE RESERVOIR* is used. Valve tests can be carried out by flushing or pressure measurements.

Functional safety

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded.

The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

Warning note for carriers of pacemakers:

Implantation of a *proSA* unit can alter the function of a pacemaker.

Sterilisation

The products are sterilised with steam under closely monitored conditions. Double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from damaged wrapping must not be used under any circumstances.

Resterilisation

The functional safety and reliability of resterilised products cannot be guaranteed, therefore resterilisation is not recommended.

Medical products consultant

In compliance with the requirements of the European law MDD 93/42/EEC, the following individuals are listed by Christoph Miethke GmbH&Co. KG as medical product consultants, to whom all queries concerning the products should be addressed:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
D-14469 Potsdam
Phone: +49(0) 7000 MIETHKE or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
E-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap Platz
D-78532 Tuttlingen
Phone: +49 (0) 7461 95-9
Fax: +49 (0) 7461 95-26 00
E-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1 (800) 282-9000
www.aesculapusa.com

Requirements of the MDD 93/42/EEC

Medical product laws call for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is out-fitted with a sticker for this purpose.

Note on the instructions for use

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.


General information

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	proSA
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the proSA with its external dimensions:	
scale 1:1	

Variants

The *proSA* is available as a single valve or as a *SHUNTSYSTEM* comprising various components.

proSA


proSA with adjustable differential pressure unit

proSA with distal catheter

proSA-SHUNTSYSTEM with BOREHOLE or SPRUNG RESERVOIR (adult and pediatric)

proSA-SHUNTSYSTEM with FLUSHING RESERVOIR or CONTROL RESERVOIR (adult and pediatric)

Scale 1:1

CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
CE marking according to directive 93/42/EEC
Label CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE

Technische Änderungen vorbehalten
Technical alterations reserved
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche

Distributed by

B | BRAUN

AESULAP AG
Am Aesculap-Platz
D-78532 Tuttlingen
Tel.: +49 (0) 7461 95-0
Fax: +49 (0) 7461 95-26 00
E-mail: information@aesculap.de
www.aesculap.de

Manufacturer acc. MDD 93/42/EEC

CHRISTOPH MIETHKE GMBH & CO. KG

Ulanenweg 2
D-14469 Potsdam
Tel.: +49 (0) 7000 MIETHKE or
Tel.: +49 (0) 331 62083-0
Fax: +49 (0) 331 62083-40
E-mail: info@miethke.com
www.miethke.com