# *Ventricular Catheter*

# *Peritoneal Catheter*

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 [www.miethke.com](http://www.miethke.com)

**US** This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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**TABLE OF CONTENTS**

|  |   |
|--|---|
| INDICATION   | 6 |
| TECHNICAL DESCRIPTION                              | 6 |
| IMPLANTATION                                       | 7 |
| SAFETY MEASURES AND CONTRAINDICATIONS              | 7 |
| INTERACTION WITH PRODUCTS FROM OTHER MANUFACTURERS | 7 |
| COMPATIBILITY TO DIAGNOSTIC PROCEDURES             | 7 |
| FUNCTIONAL SAFETY                                  | 7 |
| ADVERSE REACTIONS AND INTERACTIONS                 | 7 |
| STERILIZATION   PACKAGING                          | 8 |
| MEDICAL DEVICES CONSULTANTS                        | 8 |
| REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)      | 8 |

## INDICATION

The catheters are used for cerebrospinal fluid (CSF) drainage in the treatment of hydrocephalus.

## TECHNICAL DESCRIPTION

The catheters from the company Christoph Miethke GmbH & Co. KG have an inner diameter of 1.2 mm and an outer diameter of 2.5 mm. With these dimensions, the pressure characteristics of a shunt system are only slightly affected. The catheters are dimensioned to work safely with Miethke valves suitable for ventricular peritoneal CSF drainage.

### Ventricular Catheter

The Ventricular Catheters are silicone catheters and contain barium sulfate. This ensures that the X-ray image is clearly visible. In order to achieve a safe placement of the Ventricular Catheter it has number markings after 3, 5, 7, 10 cm and 13 cm distance from the tip.

The Ventricular Catheters are produced in such a way that the mandrin can be pushed forward into the tip. This prevents buckling during implantation (fig. 1). The drill holes begin 3 mm from the tip and end after 14 mm. This prevents drainage from the ventricles to the brain during implantation in very narrow ventricles.

1. Silicone catheter
2. Introducing stylet / mandrin
3. Burrhole


Fig. 1: Ventricular Catheter in cross section  
Length: 180 mm or 250 mm

### THOMALE CATHETER

The THOMALE CATHETER is a ventricular catheter with 6 burrholes next to the tip of the catheter. It is used preferably in the treatment in case of narrow ventricles.

#### 1. Burrhole


Fig. 2: THOMALE CATHETER in a schematic cross section

Length: 180 mm or 250 mm

### Peritoneal Catheter

The peritoneal catheters are silicone catheters with 9-20 % BaSO<sub>4</sub>. Therefore a good identification on X-ray image is warranted.


Fig. 3: Peritoneal catheter in a schematic cross section

Length: 50 mm, 600 mm, 900 mm or 1200 mm

### Peritoneal Catheter, barium striped

The peritoneal catheter, barium striped is a silicone catheter with a stripe of BaSO<sub>4</sub> which is completely enclosed by silicone (fig. 4).


Fig. 4: Peritoneal catheter, barium striped in a schematic cross section

1. Silicone
2. BaSO<sub>4</sub>

## IMPLANTATION

When implanting a catheter, the physician have to ensure that the catheter is secured with ligatures. The ventricular catheter can be implanted with the aid of an introducing stylet. After removing the stylet, the patency of the ventricular catheter can be checked by dripping out CSF. The length of the catheter that extends into the cranium can be pre-determined with the help of the deflector, after which it can be brought through the 90° angle. The deflector (*Burrhole Deflector*) can be ordered separately.

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

## SAFETY MEASURES AND CONTRAINDICATIONS

To avoid cuts and scratches in the silicone elastomer close attention should be paid while using sharp instruments. Likewise attention should be paid not to tighten the ligature too firmly. Damage may result in a loss of integrity of the shunt and therefore may require a revision.

Silicone is extremely electrostatic. You shall ensure that the catheter does not come into contact with dry towels, talcum or rough surfaces. Adherent particles may lead to tissue reactions.

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflamed skin and tension in the area of the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation, or vomiting often occur in cases of shunt dysfunction. In the event such symptoms occur, or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

Do not implant if an infection (f.e. meningitis, ventriculitis, peritonitis, bacteremia, septicemia) or the suspicion of an infection is diagnosed in a region of the patient's body where the shunt system will be implanted.

## INTERACTION WITH PRODUCTS FROM OTHER MANUFACTURERS

When using accessory products of the Christoph Miethke GmbH & Co. KG in combination with valves from other manufacturers it should be ensured that the connection between products is secure. The combination of products from different manufacturers in a shunt system is generally not recommended.

## COMPATIBILITY TO DIAGNOSTIC PROCEDURES

All catheters are made of silicone. MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt (ASTM F2503-2020). The catheters are MR Safe. The provided Burrhole Deflectors are MR Conditional.

## FUNCTIONAL SAFETY

The medical products have been designed for long-term reliable and precise operation. Still, it cannot be excluded that the medical products need to be replaced for technical or medical reasons. The medical products are able to resist positive and negative pressure up to 200 cmH<sub>2</sub>O during and after implantation.

The medical products should be stored dry and clean. The products should not be removed from the packaging until they are required.

## ADVERSE REACTIONS AND INTERACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): Infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in very rare cases noise development. Violent shocks for the outside (accident, fall) may put the integrity of the shunt system at risk.

## STERILIZATION | PACKAGING

All products are carefully and thoroughly steam sterilized. Owing to the fact that the product is packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

The functional safety and reliability of resterilized products cannot be guaranteed.

## MEDICAL DEVICES CONSULTANTS

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical products consultants as contacts for all product-related questions.

You can reach our medical devices consultants at:

Phone +49 331 62083-0

info@miethke.com

## REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the whereabouts of medical devices used in humans, especially for implants. The individual identification number of the implanted valve should therefore be recorded in the patient's medical records and patient passport to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website (<https://www.miethke.com/en/products/downloads>).

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany