# *Ventricular Catheter*

# *Peritoneal Catheter*

Ⓛ Gebrauchsanweisung | Ⓜ Instructions for use | Ⓧ Mode d'emploi  
Ⓝ Instrucciones de manejo | ⓘ Istruzioni per l'uso

Ⓤ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000.

A paper copy will be provided to you upon request at no additional cost.

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## DESCRIPTION

The catheter is used for draining cerebrospinal fluid (CSF) from the ventricles into an appropriate cavity. Due to the dimensions (inner diameter 1.2mm and outer diameter 2.5mm) the pressure-flow-characteristic of the shunt system is only an insignificant influence.

### Ventricular Catheter

The ventricular catheters of the Christoph Miethke GmbH & Co. KG are silicone catheters added with 9-14% BaSO<sub>4</sub>. Therefore a good identification on X-ray image is guaranteed. For a safe placement the ventricular catheter has got markings of 3, 5, 7, 10 and 13 cm seen from the catheter tip (see variations). By using the ventricular catheter of the Christoph Miethke GmbH & Co. KG the introducing stylet can be pushed into the tip of the catheter (Fig. 1). Therefore kinking of the catheter tip is avoided. The holes start at 3 mm from the tip and end after 14 mm. In treatment of slitlike ventricles drainage from the ventricles into the brain is therefore avoided.


Fig. 1: Ventricular catheter

### THOMALE CATHETER

The THOMALE CATHETER is a ventricular catheter with 6 holes next to the tip of the catheter. It is used preferably in the treatment in case of narrow ventricles.


Fig. 2: THOMALE CATHETER

### Peritoneal Catheter

The peritoneal catheter of the Christoph Miethke GmbH & Co. KG are silicone catheters with 9-14% BaSO<sub>4</sub>. Therefore a good identification on X-ray image is warranted.

### Peritoneal Catheter, barium striped

The peritoneal catheter, barium striped is a silicone catheter with a stripe of BaSO<sub>4</sub> which is completely enclosed by silicone (Fig. 3).


Fig. 3: Schematic cross section of the peritoneal catheter, barium striped

## IMPLANTATION AND REIMPLANTATION

When implanting a catheter, the physician have to ensure that the catheter is secured with ligatures. The ventricular catheter can be implanted with the aid of an introducing stylet. The length of the catheter that extends into the cranium can be predetermined with the help of the deflector, after which it can be brought through the 90° angle. (The deflector can be ordered separately also).

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

## TOLERABILITY TO DIAGNOSTIC PROCEDURES

All catheters manufactured by Christoph Miethke GmbH & Co. KG are made of silicone. MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt (ASTM F2503-13).The catheters are MR Safe.

## **SAFETY MEASURES**

To avoid cuts and scratches in the silicone elastomer close attention should be paid while using sharp instruments. Likewise attention should be paid not to tighten the ligature too firmly. Damage may result in a loss of integrity of the shunt and therefore may require a revision.

Silicone is extremely electrostatic. You shall ensure that the catheter does not come into contact with dry towels, talcum or rough surfaces. Adherent particles may lead to tissue reactions.

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflamed skin and tension in the area of the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation, or vomiting often occur in cases of shunt dysfunction. In the event such symptoms occur, or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

## **ADVERSE REACTION**

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage as well as intercranial and interabdominal injuries.

## **STERILISATION**

All products are carefully and thoroughly steam sterilized. Owing to the fact that the product is packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

## **RESTERILISATION**

The functional safety and reliability of resterilized products cannot be guaranteed.

## **REQUIREMENTS OF THE MDD 93/42/EEC**

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

## **NOTE ON THE INSTRUCTIONS FOR USE**

The descriptions and explanations given in this document are based on the clinical experience available up to date. It is for the surgeon's to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

**MEDICAL PRODUCTS CONSULTANT**

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke  
Dipl.-Ing. Roland Schulz  
Michaela Funk-Neubarth

Christoph Miethke GmbH & Co. KG  
Ulanenweg 2  
D-14469 Potsdam  
Phone: +49(0) 7000 6438453 or  
Phone: +49(0) 331 620 83 0  
Fax: +49(0) 331 620 83 40  
e-mail: info@miethke.com

Please address any enquiries to:  
AESCULAP AG  
Am Aesculap Platz  
D-78532 Tuttlingen  
Phone: +49 (0) 7461 95-0  
Fax: +49 (0) 7461 95-26 00  
e-mail: information@aesculap.de  
Service address in the US  
AESCULAP Inc.  
Attn. AESCULAP Technical Services  
615 Lambert Pointe Road  
Hazelwood, MO, 63042

AESCULAP Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada  
AESCULAP Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
Phone: +1-800-282-9000  
www.aesculapusa.com

**GENERAL INFORMATION**

<b>Manufacturer</b>	<b>Christoph Miethke GmbH &amp; Co. KG</b>
Product name	Catheters
Intended application	Treatment of hydrocephalus
Intended for one-time use only (disposable)	
Store in a clean and dry place	

## VARIATIONS

The catheters are available as ventricular catheters and peritoneal catheters.

Inside measurement: 1.2 mm  
Outside measurement: 2.5 mm

*Ventricular Catheter,*  
length 180 mm or 250 mm (adult and pediatric)


*THOMALE CATHETER,*  
length 180 mm or 250 mm

*Ventricular Catheter with deflector,*  
length 180 mm or 250 mm  
(adult and pediatric)

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany