# XABO<sup>®</sup> Catheters

## **XTRA PROTECTION** AGAINST INFECTION

- DE** Gebrauchsanweisung | **EN** Instructions for use | **FR** Mode d'emploi |  
**ES** Instrucciones de manejo | **PT** Instruções de utilização | **IT** Istruzioni per l'uso |  
**NL** Gebruiksaanwijzing |

 [www.miethke.com](http://www.miethke.com)

**US** This Instructions for Use is NOT intended for United States users. Please discard.  
The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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**Tel. +49 331 62083-0**  
**info@miethke.com**

**CONTENTS**

1	PREFACE AND IMPORTANT INFORMATION	13
2	INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	13
2.1	EXPLANATION OF THE WARNINGS	13
2.2	DISPLAY CONVENTIONS	13
2.3	OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION	13
2.4	FEEDBACK ON THE INSTRUCTIONS FOR USE	13
2.5	COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION	13
3	DESCRIPTION OF THE <i>XABO Catheters</i>	14
3.1	MEDICAL PURPOSE	14
3.2	CLINICAL BENEFITS	14
3.3	INDICATIONS	14
3.4	CONTRAINDICATIONS	14
3.5	INTENDED PATIENT GROUPS	14
3.6	INTENDED USERS	14
3.7	INTENDED USE ENVIRONMENT	14
3.8	TECHNICAL DESCRIPTION	14
3.9	SYSTEM COMPONENTS	16
3.10	FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	16
4	PROPERTIES OF THE <i>XABO Catheters</i>	16
4.1	PRODUCT DESCRIPTION	16
4.2	IMPORTANT SAFETY INFORMATION	17
4.3	PATIENT EDUCATION	17
4.4	TRANSPORT AND STORAGE	17
4.5	USE OF THE PRODUCT	18
4.6	EXPLANTATION AND DISPOSAL	19
4.7	TECHNICAL INFORMATION	19
4.8	SYMBOLS USED FOR LABELLING	20
5	MEDICAL DEVICE CONSULTANTS	20

## 1 PREFACE AND IMPORTANT INFORMATION

### Preface

Thank you for purchasing the medical device *XABO Catheters*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product. Your team at Christoph Miethke GmbH & Co. KG

### Relevance of the instructions for use

#### WARNING

**Improper handling and use of this product can cause hazards and damages. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.**

### Scope

*XABO Catheters* includes the following variants:

- ▶ *XABO Ventricular Catheter* with mandrin
- ▶ *XABO Ventricular Catheter* with mandrin and *Burrhole Deflector* or *Pediatric Burrhole Deflector*
- ▶ *XABO Peritoneal Catheter*
- ▶ *XABO Catheter Set*

## 2 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

### 2.1 EXPLANATION OF THE WARNINGS

#### DANGER

**Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.**

#### WARNING

**Indicates a possible risk. If not avoided, it could lead to death or serious injuries.**

#### CAUTION

**Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.**

#### NOTE

**Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.**

The symbols for danger, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

### 2.2 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates <i>product names</i>

### 2.3 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

### 2.4 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

### 2.5 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies

solely to jurisdictions in which it has trademark rights.

### 3 DESCRIPTION OF THE XABO Catheters

#### 3.1 MEDICAL PURPOSE

*XABO Catheters* are intended for the shunting of cerebrospinal fluid (CSF).

#### 3.2 CLINICAL BENEFITS

##### Clinical benefits of the *XABO Catheters* :

- ▶ Implantation of a long-term implant for the extraction of cerebrospinal fluid from the CSF circulation and its shunting into the peritoneum
- ▶ Treatment of hydrocephalus, e.g. by alleviation of the clinical symptoms
- ▶ Reducing the risk of infection with gram-positive bacteria thanks to antibiotic-impregnated catheters

#### 3.3 INDICATIONS

The following indications apply to *XABO Catheters*:

- ▶ Treatment of hydrocephalus

#### 3.4 CONTRAINDICATIONS

The following contraindications apply to *XABO Catheters*:

- ▶ Infections in the implantation area
- ▶ Pathological concentrations levels (e.g. of blood components and/or protein) within the cerebrospinal fluid
- ▶ Intolerance to materials of the shunt system
- ▶ Hypersensitivity to rifampicin and/or clindamycin hydrochloride

#### 3.5 INTENDED PATIENT GROUPS

- ▶ Patients who are treated with a CSF shunt system on account of their clinical picture

#### 3.6 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical professionals, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product
- ▶ Successful participation in product training

#### 3.7 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions

#### 3.8 TECHNICAL DESCRIPTION

*XABO Catheters* from Christoph Miethke GmbH & Co. KG have an inner diameter of 1.2 mm and an outer diameter of 2.5 mm. With these dimensions, the Pressure-Flow Characteristics of a shunt system are only slightly affected. *XABO Catheters* are dimensioned to work safely with Miethke valves suitable for ventricular peritoneal cerebrospinal fluid shunting. *XABO Catheters* are made from silicone, they are impregnated with antibiotics in a special treatment process and contain 0.054% rifampicin and 0.15% clindamycin hydrochloride. Laboratory test have shown *XABO Catheters* to reduce colonisation with gram-positive bacteria of the silicone catheter surface. The laboratory tests were performed with *Staphylococcus aureus*, *Staphylococcus epidermidis* and *Bacillus subtilis*. Systemic therapeutic effects are highly unlikely as the amounts of rifampicin and clindamycin hydrochloride contained in the catheter represent only a fraction of the therapeutic dose of these antibiotics.

##### ***XABO Ventricular Catheter***

*XABO Ventricular Catheters* are silicone catheters and contain barium sulphate. This ensures that the X-ray image is clearly visible. To ensure safe placement of the *XABO Ventricular Catheter*, it has numerical markings after 3, 5, 7, 10 and 13 cm from the tip.

*XABO Ventricular Catheters* are produced in such a way that the mandrin can be pushed-forward into the tip. This prevents buckling during implantation (Fig. 1). The drill holes begin 3 mm from the tip and end after 14 mm. This prevents drainage from the ventricles to the brain during implantation in very narrow ventricles.


Fig. 1: XABO Ventricular Catheter, cross-section of catheter tip

(lengths: 180 mm or 250 mm):

1. Silicone catheter, 2. Mandrin, 3. Drainage burr hole


Fig. 2: Position of the XABO Ventricular Catheter:

1: Burrhole Deflector (optional accessory)  
 2: Proximal side of the catheter  
 3: Distal side of the catheter

**XABO Peritoneal Catheter**

XABO Peritoneal Catheters are silicone catheters and contain barium sulphate. This ensures that the X-ray image is clearly visible.


Fig. 3: XABO Peritoneal Catheter in cross section  
 Lengths: 600 mm, 900 mm or 1200 mm


Fig. 4: Position of the XABO Peritoneal Catheter (1)

**XABO Catheter Set**

The XABO Catheter Set (Fig. 5) is a combination of the XABO Ventricular Catheter (Fig. 1) and the XABO Peritoneal Catheter (Fig. 3) as described above.


Fig. 5: Position of the XABO Catheter Set:

1 Burrhole Deflector (optional accessory)  
 2 XABO Ventricular Catheter  
 3 Valve (not included in the XABO Catheter Set)  
 4 XABO Peritoneal Catheter

### 3.9 SYSTEM COMPONENTS

#### Combination with shunt components

The *XABO Catheters* can be safely combined with our range of implantable shunt components. We recommend using the Christoph Miethke GmbH & Co. KG products in combination with the *XABO Catheters*.

#### 3.10 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

- ▶ The medical device has been designed to work accurately and reliably over long periods of time. We guarantee that our medical devices are free from defects and fully functional at the time of dispatch. Excluded from this guarantee are cases in which the medical device has been exchanged or rather explanted for technical or medical reasons outside of our sphere of responsibility.
- ▶ *XABO Catheters* together with the entire shunt system can safely withstand pressure of up to 100 cm H<sub>2</sub>O occurring during and after surgery.
- ▶ Nuclear magnetic resonance examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment. *XABO Catheters* are MR Safe. *Burhole Deflectors* are MR Conditional. The conditions for MRI security of the products are found on our website: <https://miethke.com/downloads/>

### 4 PROPERTIES OF THE XABO Catheters

#### 4.1 PRODUCT DESCRIPTION

##### 4.1.1 PRODUCT VARIANTS

*XABO Catheters* are available in the following variants:

*XABO Ventricular Catheter*:

- ▶ Length 180 mm
- ▶ Length 250 mm

*XABO Peritoneal Catheter*:

- ▶ Length 600 mm
- ▶ Length 900 mm
- ▶ Length 1200 mm

#### 4.1.2 SCOPE OF DELIVERY

Box content	Number
Sterile packaging with <i>XABO Catheters</i>	1
Instructions for use for <i>XABO Catheters</i>	1

#### 4.1.3 STERILITY

#### WARNING

**The product must not be used if the sterile packaging or the product is damaged or after the expiry date.**

*XABO Catheters* are sterilised under strictly controlled conditions using irradiation. The respective expiry date is printed on the packaging. The inner sterile packaging consists of a gas- and liquid-tight barrier to protect the antibiotics.

#### 4.1.4 REPEATED USE AND RESTERILISATION

#### WARNING

**The product must not be resterilised or reprocessed in any other way as the safe functioning and sterility of the product cannot be guaranteed.**

Products that have already been implanted in a patient must not be reused either on the same or a different patient in order to minimise the risk of infection.

#### 4.1.5 SINGLE-USE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the *XABO Catheters*. No guarantee can be assumed for the functional safety of resterilised products.

#### 4.1.6 PRODUCT CONFORMITY

The product meets current regulatory requirements.

The requirements stipulate the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the medical device should therefore be recorded in the patient's medical records to ensure complete traceability.

**4.2 IMPORTANT SAFETY INFORMATION**

**4.2.1 SAFETY INSTRUCTIONS**

**Important!** Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

- ▶ **The product must not be used if the sterile packaging or the product is damaged or after the expiry date.**
- ▶ **Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be carefully read and understood before the product is used for the first time.**
- ▶ **Prior to use, it is essential to check the product for completeness and integrity.**

**4.2.2 COMPLICATIONS, SIDE EFFECTS, PRECAUTIONS AND RESIDUAL RISKS**

The following complications can occur in conjunction with the *XABO Catheters*:

- ▶ Headaches, dizzy spells, mental confusion, vomiting in cases of possible leakage from the shunt system and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid
- ▶ Allergic reaction / intolerance to product materials
- ▶ Overdrainage/underdrainage

Violent external shocks (accident, fall, etc.) may put the integrity of the shunt system at risk. As a precaution, a physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the *XABO Catheters*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation
- ▶ Local shunt irritation
- ▶ Allergic reactions to catheter components, particularly to the active antibiotic substances rifampicin and clindamycin hydrochloride

**4.2.3 REPORTING OBLIGATION**

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible state authority.

**4.3 PATIENT EDUCATION**

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product (Ch. 4.2, 4.5).

**4.4 TRANSPORT AND STORAGE**

The medical devices must always be transported and stored in a clean and dry place. *XABO Catheters* are to be protected from direct sunlight. The products should not be removed from the packaging until they are required.

**4.4.1 TRANSPORT**

**Transport conditions**

<b>Ambient temperature</b>	≤ 40 °C
----------------------------	---------

**4.4.2 STORAGE**

**Storage conditions**

<b>Temperature range for storage</b>	≤ 30 °C
--------------------------------------	---------

## 4.5 USE OF THE PRODUCT

### 4.5.1 INTRODUCTION

XABO Catheters are used for cerebrospinal fluid (CSF) shunting in the treatment of hydrocephalus.

### 4.5.2 SAFETY NOTICES AND WARNINGS

#### WARNING

- ▶ XABO Catheters must not be immersed in antibiotic solutions. The contact of the catheter with sterile water or sterile saline solution should be limited to a minimal period of time and should be done immediately prior to implantation. The solution may take on a faint orange colour.

#### CAUTION

- ▶ Silicone is extremely electrostatic. Care must be taken to avoid the catheter coming into contact with dry cloths, talcum powder or rough surfaces. Clinging particles could lead to tissue reactions.
- ▶ When using sharp instruments, care should be taken to avoid cuts and scratches in the silicone elastomer.
- ▶ It must be ensured that the ligature is not tightened excessively. Damage may result in a loss of integrity of the shunt and thus necessitate a revision.
- ▶ The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

### 4.5.3 REQUIRED MATERIALS

The XABO Catheters is designed for safe use with the described shunt components in chapter 3.9. In any case, catheters have to be carefully fixed with a ligature to the connectors of the shunt components. Any kinks in the catheter must be avoided.

### 4.5.4 PREPARING FOR IMPLANTATION

#### Checking the sterile packaging

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier

system. The products should only be removed from the packaging immediately prior to use.

#### Checking the integrity of the XABO Catheters

The integrity of the XABO Catheters must be checked prior to implantation. To this end, XABO Catheters must be visually inspected.

### 4.5.5 PERFORMING THE IMPLANTATION

#### CAUTION

The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

When implanting the XABO Catheters, care must be taken to ensure that the catheter is secured with a ligature (Fig. 6). To that end, insert the catheter into the connector up to the limit stop. At the centre of the connector, a bulge of the catheter becomes visible. A ligature is placed in this area, secured with a surgical knot.


Fig. 6: Ligature-fixed catheter

1: Catheter, 2: Ligature, 3: Limit stop

#### Positioning of the XABO Ventricular Catheter

Several surgical techniques are available for the positioning of the XABO Ventricular Catheter. The required skin incision should be made in form of a lobule pedicled towards the shunting catheter or by a straight skin incision. If a *Burrhole Reservoir* - or a *SPRUNG RESERVOIR* - is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole.

XABO Catheters is available in a range of different configurations: If using a *Burrhole Reser-*

voir - or a *SPRUNG RESERVOIR* -, the *XABO Ventricular Catheter* is implanted first. Once the mandrin has been removed, the patency of the *XABO Ventricular Catheter* can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the *Burrhole Reservoir* - or the *SPRUNG RESERVOIR* - connected and the connection secured with a ligature.

When using a shunt system with a *CONTROL RESERVOIR*, a *Burrhole Deflector* is included. The *Burrhole Deflector* is used to adjust the length of the catheter to be implanted and to position it inside the ventricle. The *Ventricular Catheter* is deflected by 90° and the *CONTROL RESERVOIR* put into place. The position of the *XABO Ventricular Catheter* should be inspected after the procedure by imaging (such as CT or MRI).

#### **Positioning of the *XABO Peritoneal Catheter***

The access site for the *XABO Peritoneal Catheter* is left to the surgeon's discretion. For example, it can be used in a paraumbilical application or applied at the level of the epigastrium. Likewise, various surgical techniques are available for placing the *XABO Peritoneal Catheter*. The recommendation is to pull the

*XABO Peritoneal Catheter* from the valve to the intended position using a subcutaneous *Tunneller*, if necessary with the aid of an auxiliary incision. The *XABO Peritoneal Catheter*, usually securely attached to the valve, has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the *XABO Peritoneal Catheter* (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

## **4.6 EXPLANTATION AND DISPOSAL**

### **4.6.1 EXPLANTATION**

The explantation of the *XABO Catheters* should be performed according to the state of the art and in compliance with medical practice.

### **4.6.2 DISPOSAL**

#### ***XABO Catheters* and shunt components**

Products and product parts not used in the implantation or surgically removed must be disposed of correctly as potentially infectious material in accordance with medical practice as well as respective regional laws and regulations.

Explanted medical devices must not be reused.

## **4.7 TECHNICAL INFORMATION**

### **4.7.1 TECHNICAL DATA**

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>XABO Catheters</i>
Medical purpose	Shunting of cerebrospinal fluid (CSF)
Sterilisability	Cannot be resterilised
Storage	Store in a clean and dry place at ≤ 30 °C
For single use only	

## 4.8 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxx indicates the identifier of the responsible notified body
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Sterilised using irradiation
	Do not resterilise
	Do not reuse
	Do not use if package is damaged and consult instructions for use
	Keep dry
	Upper limit of temperature
	Keep away from sunlight
	Consult instructions for use / electronic instructions for use
	Caution
	Contains a medical substance
	Non-pyrogenic

Symbol	Explanation
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians.
	MR Safe
	MR Conditional
	Patient identification
	Date
	Healthcare centre or doctor
	Patient information website
	Model number / European Medical Device Nomenclature Code

## 5 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany