# *M.blue*<sup>®</sup> *XABO*<sup>®</sup>

## **THE BALANCED WAY OF LIFE INSPIRED BY YOU**

- Ⓛ DE Gebrauchsanweisung | Ⓛ EN Instructions for use | Ⓛ FR Mode d'emploi |  
Ⓛ ES Instrucciones de manejo | Ⓛ PT Instruções de utilização | Ⓛ IT Istruzioni per l'uso |  
Ⓛ NL Gebruiksaanwijzing |

 [www.miethke.com](http://www.miethke.com)

Ⓛ US This Instructions for Use is NOT intended for United States users. Please discard.  
The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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**info@miethke.com**

**CONTENTS**

1	PREFACE AND IMPORTANT INFORMATION	22
2	INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	22
2.1	EXPLANATION OF THE WARNINGS	22
2.2	DISPLAY CONVENTIONS	22
2.3	OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION	22
2.4	FEEDBACK ON THE INSTRUCTIONS FOR USE	22
2.5	COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION	22
3	DESCRIPTION OF THE <i>M.blue XABO</i>	23
3.1	MEDICAL PURPOSE	23
3.2	CLINICAL BENEFITS	23
3.3	INDICATIONS	23
3.4	CONTRAINDICATIONS	23
3.5	INTENDED PATIENT GROUPS	23
3.6	INTENDED USERS	23
3.7	INTENDED USE ENVIRONMENT	23
3.8	TECHNICAL DESCRIPTION	24
3.9	OPERATING PRINCIPLE OF THE VALVE	25
3.10	SELECTION OF THE APPROPRIATE PRESSURE SETTING	25
3.11	PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES	26
3.12	PRESSURE-FLOW CHARACTERISTICS	26
3.13	APPLICATION OF THE M.BLUE PLUS INSTRUMENTS	28
3.14	SYSTEM COMPONENT	32
3.15	FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	32
4	PROPERTIES OF THE <i>M.blue XABO</i>	32
4.1	PRODUCT DESCRIPTION	32
4.2	IMPORTANT SAFETY INFORMATION	33
4.3	PATIENT EDUCATION	34
4.4	TRANSPORT AND STORAGE	34
4.5	USE OF THE PRODUCT	34
4.6	EXPLANTATION AND DISPOSAL	36
4.7	TECHNICAL INFORMATION	37
4.8	SYMBOLS USED FOR LABELLING	38
5	MEDICAL DEVICE CONSULTANTS	38

## 1 PREFACE AND IMPORTANT INFORMATION

### Preface

Thank you for purchasing the medical device *M.blue XABO*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product. Your team at Christoph Miethke GmbH & Co. KG

### Relevance of the instructions for use

#### WARNING

**Improper handling and use of this product can cause hazards and damages. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.**

### Scope

*M.blue XABO* includes the following components:

- ▶ *M.blue XABO*

Additional options:

- ▶ *Reservoirs* (including paediatric version)
- ▶ *Prechamber* (including paediatric version)
- ▶ *Burrhole Deflector* (including paediatric version)
- ▶ *Ventricular Catheter*
- ▶ *Titanium Connectors*

## 2 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

### 2.1 EXPLANATION OF THE WARNINGS

#### DANGER

**Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.**

#### WARNING

**Indicates a possible risk. If not avoided, it could lead to death or serious injuries.**

#### CAUTION

**Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.**

#### NOTE

**Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.**

The symbols for danger, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

### 2.2 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates <i>product names</i>

### 2.3 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

Each delivery is accompanied by a patient card with information about the product. The patient card is to provide the treating physician with all the product information in a compact form for the patient record.

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

### 2.4 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

### 2.5 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

### 3 DESCRIPTION OF THE *M.blue XABO*

#### 3.1 MEDICAL PURPOSE

*M.blue XABO* is intended for shunting of cerebrospinal fluid (CSF).

#### 3.2 CLINICAL BENEFITS

##### Clinical benefits of the *M.blue XABO* :

- ▶ Implantation of a long-term implant for the extraction of cerebrospinal fluid from the CSF circulation and its shunting into the peritoneum
- ▶ Treatment of hydrocephalus, e.g. by alleviation of the clinical symptoms
- ▶ Reducing the risk of infection with gram-positive bacteria thanks to antibiotic-impregnated catheters

#### 3.3 INDICATIONS

The following indications apply to *M.blue XABO*:

- ▶ Treatment of hydrocephalus

#### 3.4 CONTRAINDICATIONS

The following contraindications apply to *M.blue XABO*:

- ▶ Infections in the implantation area
- ▶ Pathological concentrations levels (e.g. of blood components and/or protein) within the cerebrospinal fluid
- ▶ Intolerance to materials of the shunt system
- ▶ Hypersensitivity to rifampicin and/or clindamycin hydrochloride

#### 3.5 INTENDED PATIENT GROUPS

- ▶ Patients who are treated with a CSF shunt system on account of their clinical picture

#### 3.6 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical professionals, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product
- ▶ Successful participation in product training

#### 3.7 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions

### 3.8 TECHNICAL DESCRIPTION


Fig. 1: M.blue cross section

#### 1. Adjustable gravitational unit

2. Sapphire ball
3. Tantalum weight
4. Rotor

#### 5. Differential pressure unit

6. Sapphire ball
7. Micro-coil spring

The *M.blue* is a valve made of titanium. It consists of an adjustable gravitational unit and a differential pressure unit (Fig. 1).

The adjustable gravitational unit (1) in the proximal part of the valve contains a tantalum weight (3), which holds a sapphire ball in the ball seat via a lever (2). Depending on the patient's body position, the influence of the tantalum weight on the sapphire ball changes, which also affects the valve opening pressure. Via a rotor (4), the pretension of the torsion spring connected to the lever can be adjusted through the skin after surgery. This way, the impact of the tantalum weight on the sapphire ball can be influenced and, consequently, the valve opening pressure can be adjusted.

In the distal part of the valve, a micro-coil spring (7) controls the opening pressure of the differential pressure unit (5). A sapphire ball (6)

ensures precise opening and closing of the ball-cone unit.

The *M.blue XABO* consists of an *M.blue* valve in combination with *XABO Catheters*.

*XABO Catheters* are made from silicone, they are impregnated with antibiotics in a special treatment process and contain 0.054% rifampicin and 0.15% clindamycin hydrochloride. Laboratory test have shown *XABO Catheters* to reduce colonisation with gram-positive bacteria of the silicone catheter surface. The laboratory tests were performed with *Staphylococcus aureus*, *Staphylococcus epidermidis* and *Bacillus subtilis*. Systemic therapeutic effects are highly unlikely as the amounts of rifampicin and clindamycin hydrochloride contained in the catheter represent only a fraction of the therapeutic dose of these antibiotics.

### 3.9 OPERATING PRINCIPLE OF THE VALVE

The *M.blue XABO* is a posture-dependent hydrocephalus valve. The opening pressure of the *M.blue XABO* is composed of the opening pressures of the adjustable gravitational unit and the differential pressure unit.

#### Horizontal body position

In the horizontal position, the gravitational unit is always open and does not present any resistance.

The opening pressure of the *M.blue XABO* in the horizontal position is thus determined by the adjustable differential pressure unit. The basic operating principle of the differential pressure unit is shown in Fig. 2 a) and b).


Fig. 2: Functional principle of the differential pressure unit

a) closed b) open

In Fig. 2 a), the valve is closed; hence, shunting is not possible.

If the patient's intraventricular pressure (IVP) exceeds the spring force of the micro-coil spring, which otherwise keeps the differential pressure unit closed, the sealing ball moves out of the cone, leaving a gap for shunting of CSF (Fig. 2 b).

#### Vertical body position

As the patient sits or stands up, the gravitational unit closes the discharge channel in the proximal part of the valve (Fig. 3 a). Thus, the opening pressure of the *M.blue XABO* is increased in the upright position, because now the weight of the tantalum weight (opening pressure of the gravitational unit) must be overcome in addition to the opening pressure of the differential pressure unit. Shunting is only possible again when the sum of IVP and hydrostatic suction is greater than the opening pressure of both units (Fig. 3 b). For individual adaptation of the

opening pressure to the patient, a valve opening pressure between 0 and 40 cmH<sub>2</sub>O can be selected for the adjustable gravitational unit.


Fig. 3: Gravitational unit in the vertical body position  
a) closed b) open

#### NOTE

**During physical activity associated with shock (e.g. jogging), the opening pressure of the *M.blue XABO* may decrease temporarily according to laboratory results. Fundamentally, though, functionality remains the same. At the end of physical activity, the opening pressure returns to its original level and remains stable.**

### 3.10 SELECTION OF THE APPROPRIATE PRESSURE SETTING

Recommended pressure settings for the *M.blue XABO* are available at:

<https://www.miethke.com/downloads/>

This is a non-binding recommendation for the treating physician. The physician decides on each case independently, individually and autonomously in accordance with his/her diagnosis.

The adjustable gravitational unit of the *M.blue XABO* is set to an opening pressure of 20 cmH<sub>2</sub>O upon delivery.

This pre-set opening pressure can be changed to a different pressure prior to implantation.

#### Horizontal body position

The opening pressure in the horizontal body position is determined by the differential pressure unit.

In this case, the pressure level should be set in accordance with the patient's clinical picture, indication and age. The standard settings for the horizontal position are 5 to 10 cmH<sub>2</sub>O.

#### Vertical body position

The opening pressure of the *M.blue XABO* for the vertical body position is calculated from the

sum of the opening pressure of the differential pressure unit and the gravitational unit.

Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure level for the gravitational unit.

### 3.11 PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES

#### **M.blue** gravitational unit

The selected pressure level of the *M.blue* gravitational unit should always be monitored using the *M.blue plus Compass*, but it can also be checked using X-ray (Fig. 4).


Fig. 4: X-ray image (adjustable gravitational unit pre-set to 20 cmH<sub>2</sub>O; differential pressure unit: 0 cmH<sub>2</sub>O)

The rotor setting is decisive in this case. The four magnets in the rotor appear in the X-ray image as white dots and are located opposite each other in pairs. On one side of the rotor, two additional burr holes (right and left next to the magnet pairs) serve as orientation. They appear as black dots in the X-ray image. This side can be described as the rear side of the rotor. The two front magnets are on the opposite side.

The space between these two magnets can be considered as the apex of the triangle. The pressure level can be read off using the orientation of this intermediate space (Fig. 5). The apex of the triangle can take up any position except the space marked as the non-adjustable range in Fig. 5. This means that the opening pressure of the *M.blue* can be infinitely variably adjusted from 0 up to 40 cmH<sub>2</sub>O. To prevent reading the pressure level in reverse, the top view of the valve shows an obvious recess in the housing ring with the tantalum weight to the right of the inlet connector (Fig. 4).


Fig. 5: Schematic representation of rotor in an X-ray image

1. Coding holes for the adjustable gravitational unit
2. Coding holes for the differential pressure unit
3. Triangle apex, 4. Non-adjustable range
5. proximal, 6. distal

#### **Differential pressure unit of the M.blue**

The pressure level of the pre-set *M.blue* differential pressure unit can be recognised in the X-ray image by the encoding (Fig. 6). The following pressure levels are possible for the differential pressure unit:


Fig. 6: Pressure level encodings of the differential pressure unit (1.) with coding holes (2.)

### 3.12 PRESSURE-FLOW CHARACTERISTICS

The opening pressure refers to a reference flow of 20 ml/h. The supplied XABO Catheters do

not fundamentally change the pressure-flow characteristics.

**Horizontal body position**

The pressure flow characteristics of the M.blue differential pressure unit are shown below for pressure settings of 0, 5, 10 and 15 cmH<sub>2</sub>O in the horizontal valve position.









Fig. 7: Pressure-flow characteristics of selected M.blue pressure levels in the horizontal body position; pressure (cmH<sub>2</sub>O), flow rate (ml/h):

a) M.blue 0, b) M.blue 5, c) M.blue 10, d) M.blue 15

**Vertical body position**

In the vertical body position, the M.blue opening pressure is composed of the setting of the differential pressure unit and the adjustable gravitational unit. The pressure flow characteristics for various pressure levels in the vertical body position are shown below.

### 3.13 APPLICATION OF THE M.blue plus Instruments

#### CAUTION

Only M.blue plus Instruments approved for that purpose must be used to determine, change and monitor the opening pressure of the gravitational unit of the (M.blue).

The M.blue plus Instruments may only be used by trained specialists.

The selected pressure level of the M.blue can be determined, adjusted and monitored using the M.blue plus Instruments.

The M.blue plus Compass (Fig. 9) is used for localising and reading the adjustable gravitational unit of the M.blue.


Fig. 9: M.blue plus Compass

a) open, 1. scale ring

b) closed, 2. float gauge compass

The M.blue plus Adjustment Ring (Fig. 10) is used to set the opening pressure of the gravitational unit of the M.blue from 0 to 40 cmH<sub>2</sub>O.

Fig. 10: M.blue plus Adjustment Ring

The opening pressure of the adjustable gravitational unit of the M.blue can be changed before

Fig. 8: Pressure-flow characteristics of selected M.blue pressure levels in the vertical body position; pressure (cmH<sub>2</sub>O), flow rate (ml/h):

a) M.blue 0, b) M.blue 5, c) M.blue 10, d) M.blue 15

or after implantation. It is pre-set by the manufacturer to 20 cmH<sub>2</sub>O.

In order to set the opening pressures, the following steps must be performed:

## 1. Localisation

### WARNING

The *M.blue plus Compass* should be placed as centrally as possible over the valve, as the determined opening pressure may otherwise be incorrect.

### NOTE

- ▶ Possible air pockets in the *M.blue plus Compass* do not affect its function.
- ▶ The *M.blue plus Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *M.blue plus Adjustment Ring* should not be placed in the immediate vicinity of the *M.blue plus Compass* when determining the opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin may make adjustment difficult for a few days after surgery. If the valve setting cannot be checked conclusively using the *M.blue plus Compass*, we recommend checking it with an imaging method.

When the *M.blue plus Compass* is opened, a circular cut-out becomes visible, which is used to localise the valve in the patient's head as centrally as possible using the index finger (Fig. 11).

Fig. 11: Localising the valve

The direction markings show the flow direction.

## 2. Test procedure

In order to determine the selected pressure level, the *M.blue plus Compass* is then closed again. The float gauge should now be centred by moving the instrument within the designated circular marking (Fig. 12). Once the float gauge is centred, the current setting of the opening pressure the gravitational unit (*M.blue*) can be read from the line marking on the float gauge (Fig. 12).


Fig. 12: Determining the pressure setting with the *M.blue plus Compass*

There are two scales on the scale ring (Fig. 13). The blue adjustment range from 0 to 40 cmH<sub>2</sub>O on the inner scale applies to the opening pressure of the gravitational unit of the *M.blue*.


Fig. 13: Scale ring of the M.blue plus Compass

1. Inner: Scale for the M.blue gravitational unit from 0 to 40 cmH<sub>2</sub>O (the opening pressure of the M.blue gravitational unit in the illustrated example is 16 cmH<sub>2</sub>O)
2. Outer: Scale for the proGAV 2.0 differential pressure unit from 0 to 20 cmH<sub>2</sub>O (the opening pressure of the proGAV 2.0 differential pressure unit in the illustrated example is 17 cmH<sub>2</sub>O)

### 3. Adjustment process

#### WARNING

When adjusting the gravitational unit of the M.blue, care must be taken to change the opening pressure by a maximum of 16 cmH<sub>2</sub>O per adjustment process; otherwise, errors can result.

**Example: The opening pressure is to be changed from 6 to 36 cmH<sub>2</sub>O. The correct method is an adjustment in two stages: an initial adjustment from 6 to 22 cmH<sub>2</sub>O and subsequently from 22 to 36 cmH<sub>2</sub>O.**

In order to adjust the opening pressure, the M.blue plus Compass is opened, but without changing the position of the scale ring. The M.blue plus Adjustment Ring is now inserted into the scale ring in such a manner that the line marking points to the desired value on the scale of the scale ring (Fig. 14).


Fig. 14: Inserting the M.blue plus Adjustment Ring

1. M.blue plus Adjustment Ring, 2. Scale ring


Fig. 15: setting the M.blue gravitational unit to 32 cmH<sub>2</sub>O.

By applying slight pressure with the index finger to the valve diaphragm located in the centre of the M.blue plus Adjustment Ring and under the skin, the rotor brake is released and the opening pressure of the differential pressure or gravitational unit is changed to the desired value (Fig. 16).


Fig. 16: Adjustment using the M.blue plus Adjustment Ring

The gravitational unit of the M.blue is equipped with a feedback mechanism.

Due to the valve housing design, targeted pressure on the valve produces an audible acoustic signal (a clicking sound) and/or palpable resistance as soon as the rotor brake has been released. The valve thus shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released, the rotor is once again adjustment-proof. Although the click caused by releasing the rotor brake is easily audible before implantation, it may be considerably reduced after implantation and the filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patient or by using a stethoscope.

### Adjustment using the *M.blue plus Adjustment Assistant*

Alternatively, the *M.blue plus Adjustment Assistant* can be used to adjust the opening pressure. To do this, insert the *M.blue plus Adjustment Assistant* into the *M.blue plus Adjustment Ring* aligned to the desired value and press it with your index finger (Fig. 17).


Fig. 17: *M.blue plus Adjustment Assistant*

### Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the set pressure level. To do this, proceed as in Points 1 and 2. Should the measured value not correspond with the required pressure level, the adjustment process should be repeated. To do this, start again at Point 3.

### *M.blue Checkmate*

The *M.blue Checkmate* (Fig. 18) is supplied sterile and can be re-sterilised. The *M.blue Checkmate* can be used to change and check the pressure level before and during valve implantation directly at the *M.blue*. To determine the pressure level, the *M.blue Checkmate*

is placed centrally onto the *M.blue*. The *M.blue Checkmate* automatically aligns itself over the valve. The pressure level can be read from the direction of the proximal catheter (leading towards the valve). If the pressure level is to be adjusted, the *M.blue Checkmate* is placed centrally onto the *M.blue*. When doing so, the required pressure level must point towards the proximal catheter (leading towards the valve). By slightly pressing the *M.blue Checkmate* onto the valve, the rotor brake in the *M.blue* is released and the pressure level set.

Fig. 18: *M.blue Checkmate*, colour: blue  
Pressure levels: 0 to 40 cmH<sub>2</sub>O

### Checking and adjusting without unpacking

To protect the impregnation, the *M.blue XABO* is wrapped in opaque packaging. However, it is still possible to check and adjust the setting of the *M.blue XABO*.

The position of the *M.blue XABO* is indicated by an arrow mark. The arrow indicates the flow direction.

For localisation, the *M.blue plus Compass* is placed centrally on the marking of the valve that needs to be checked. Because of the packaging, the *M.blue plus Compass* is slightly tilted on the packaging.

Checking and/or adjusting without unpacking retains the sterility of the product and the integrity of the packaging.


Fig. 19: Packaging of the *M.blue XABO*: marking of the adjustable valve with an arrow


Fig. 20: Position of the M.blue plus Compass on the packaging of the M.blue XABO

### 3.14 SYSTEM COMPONENT

#### Combination with shunt components

The M.blue XABO can be safely combined with our range of implantable shunt components. We recommend using the Christoph Miethke GmbH & Co. KG products in combination with the M.blue XABO.

#### Reservoirs

The use of shunt systems with a reservoir provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to an integrated check valve in the CONTROL RESERVOIR and the SPRUNG RESERVOIR, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the shunt system as well as the Ventricular Catheter.

During the pump action, access to the Ventricular Catheter is closed. The use of this Reservoir does not increase the opening pressure of the shunt system. Puncturing the Reservoir should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. A firm titanium floor prevents the bottom surface from being pierced. 30 punctures are possible without any restrictions.

**WARNING**

**Frequent pumping of the Reservoir can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.**

#### Burrhole Deflector

Because of the tight fit on the Ventricular Catheter, the Burrhole Deflector makes it possible to choose the length of catheter penetrat-

ing into the skull prior to implantation. The Ventricular Catheter is deflected at a right angle in the burr hole (see chapter 4.5.5).

#### Tube systems

The M.blue XABO is supplied as a shunt system with integrated catheters (internal diameter 1.2 mm, external diameter 2.5 mm).

For a new connection of catheter and connector, the catheters have to be carefully fixed with a ligature to the valve's Titanium Connectors.

### 3.15 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

- ▶ The medical device has been designed to work accurately and reliably over long periods of time. We guarantee that our medical devices are free from defects and fully functional at the time of dispatch. Excluded from this guarantee are cases in which the medical device has been exchanged or rather explained for technical or medical reasons outside of our sphere of responsibility.
- ▶ M.blue XABO together with the entire shunt system can safely withstand pressure of up to 100 cm H<sub>2</sub>O occurring during and after surgery.
- ▶ Nuclear magnetic resonance examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. M.blue XABO is MR Conditional. The supplied catheters are MR Safe. Reservoirs, Burrhole Deflectors and Titanium Connectors are MR Conditional.

The conditions for MRI security of the products are found on our website: <https://miethke.com/downloads/>

## 4 PROPERTIES OF THE M.blue XABO

### 4.1 PRODUCT DESCRIPTION

#### 4.1.1 M.blue XABO VARIANTS

The M.blue XABO is available in a range of different variants. They differ in the pre-set pressure of the differential pressure unit.

Differential pressure unit	Adjustable gravitational unit
0 cmH <sub>2</sub> O	0–40 cmH <sub>2</sub> O
5 cmH <sub>2</sub> O	0–40 cmH <sub>2</sub> O
10 cmH <sub>2</sub> O	0–40 cmH <sub>2</sub> O
15 cmH <sub>2</sub> O	0–40 cmH <sub>2</sub> O

These variants are also available as shunt systems and can come with the following components: *XABO Ventricular Catheter*, *Prechambers*, *Reservoirs*.

#### 4.1.2 SCOPE OF DELIVERY

Box content	Number
Sterile packaging with <i>M.blue XABO</i> shunt system	1
Instructions for use for <i>M.blue XABO</i>	1
Patient card	1
Recommended pressure level	1

#### 4.1.3 STERILITY

#### WARNING

**The product must not be used if the sterile packaging or the product is damaged or after the expiry date.**

*M.blue XABO* is sterilised under strictly controlled conditions using irradiation. The respective expiry date is printed on the packaging. The inner sterile packaging consists of a gas- and liquid-tight barrier to protect the antibiotics.

#### 4.1.4 REPEATED USE AND RESTERILISATION

#### WARNING

**The product must not be resterilised or reprocessed in any other way as the safe functioning and sterility of the product cannot be guaranteed.**

Products that have already been implanted in a patient must not be reused either on the same or a different patient in order to minimise the risk of infection.

#### 4.1.5 SINGLE-USE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the *M.blue XABO*. No guarantee can be assumed for the functional safety of resterilised products.

#### 4.1.6 PRODUCT CONFORMITY

The product meets current regulatory requirements.

The requirements stipulate the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the medical device should therefore be recorded in the patient's medical records to ensure complete traceability.

#### 4.2 IMPORTANT SAFETY INFORMATION

##### 4.2.1 SAFETY INSTRUCTIONS

**Important!** Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

#### WARNING

- ▶ **The product must not be used if the sterile packaging or the product is damaged or after the expiry date.**
- ▶ **Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be carefully read and understood before the product is used for the first time.**
- ▶ **Prior to use, it is essential to check the product for completeness and integrity.**

##### 4.2.2 COMPLICATIONS, SIDE EFFECTS, PRECAUTIONS AND RESIDUAL RISKS

The following complications can occur in conjunction with the *M.blue XABO*:

- ▶ Headaches, dizzy spells, mental confusion, vomiting in cases of possible leakage from the shunt system and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid

- ▶ Allergic reaction / intolerance to product materials
- ▶ Overdrainage/underdrainage
- ▶ Noise development

Violent external shocks (accident, fall, etc.) may put the integrity of the shunt system at risk.

As a precaution, a physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the *M.blue XABO*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation
- ▶ Local shunt irritation
- ▶ Allergic reactions to catheter components, particularly to the active antibiotic substances rifampicin and clindamycin hydrochloride

#### 4.2.3 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible state authority.

#### 4.3 PATIENT EDUCATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product (Ch. 4.2, 4.5).

#### 4.4 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place. *M.blue XABO* is to be protected from direct sunlight. The products should not be removed from the packaging until they are required.

#### 4.4.1 TRANSPORT

##### Transport conditions

Ambient temperature	≤ 40 °C
---------------------	---------

#### 4.4.2 STORAGE

##### Storage conditions

Temperature range for storage	≤ 30 °C
-------------------------------	---------

#### 4.5 USE OF THE PRODUCT

##### 4.5.1 INTRODUCTION

The *M.blue XABO* is a posture-dependent valve with an adjustable gravitational unit and a pre-set differential pressure unit combined with antibiotic impregnated *XABO Catheters*.

*M.blue XABO* is used for shunting of cerebrospinal fluid (CSF) in the treatment of hydrocephalus. Valves and Reservoirs are placed in suitable positions along the course of the shunt.

##### 4.5.2 SAFETY NOTICES AND WARNINGS

##### WARNING

- ▶ *M.blue XABO* must not be immersed in antibiotic solutions. The contact of the catheter with sterile water or sterile saline solution should be limited to a minimal period of time and should be done immediately prior to implantation. The solution may take on a faint orange colour.
- ▶ Frequent pumping of the *Reservoir* can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.
- ▶ The gravitational unit of the *M.blue XABO* is posture-dependent. For that reason, care must be taken to implant the gravitational unit is implanted parallel to the body axis.

**WARNING**

- ▶ The adjustable valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e. g. underneath heavily scarred tissue). If the implantation site is unfavourably chosen or the skin above the valve is too thick, it may no longer be possible to adjust the adjustable unit. The valve then works with unchangeable pressure levels.
- ▶ If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, and adjustment of the valve cannot be ruled out.

**CAUTION**

- ▶ Silicone is extremely electrostatic. Care must be taken to avoid the catheters coming into contact with dry cloths, talcum powder or rough surfaces. Clinging particles could lead to tissue reactions.
- ▶ When using sharp instruments, care should be taken to avoid cuts and scratches in the silicone elastomer.
- ▶ It must be ensured that the ligature is not tightened excessively. Damage may result in a loss of integrity of the shunt and thus necessitate a revision.
- ▶ The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.
- ▶ For people using cardiac pacemakers: It is possible that the function of the pacemaker is influenced by the implantation of a *M.blue XABO*.

**NOTE**

- ▶ In MRI imaging *M.blue XABO* creates artefacts that are larger than the valve itself.

**4.5.3 REQUIRED MATERIALS**

The *M.blue XABO* is designed so that it can be safely used with the shunt components described in chapter 3.14. Catheters with an inner diameter of 1.2 mm and an outer diameter of approx. 2.5 mm should be used for connection. In any case, catheters have to be care-

fully fixed with a ligature to the connectors of the shunt components. Any kinks in the catheter must be avoided.

**4.5.4 PREPARING FOR IMPLANTATION****Checking the sterile packaging**

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

**Preoperative valve test**

The *M.blue XABO* should be vented before implantation and checked for permeability. The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological salt solution. The valve is patent if saline solution can be extracted (Fig. 21).


Fig. 21: Patency test

**WARNING**

- ▶ No antibiotic solution must be used for valve testing and/or venting, as this may result in a reaction with the active substances in the impregnation.
- ▶ Contaminants in the solution used for testing can impair the product's performance.
- ▶ Pressurisation with a single-use syringe should be avoided both at the proximal and the distal end. (Fig. 22).

Fig. 22: Avoidance of pressurisation

## 4.5.5 PERFORMING THE IMPLANTATION

### Positioning of the XABO Ventricular Catheter

Several surgical techniques are available for the positioning of the XABO Ventricular Catheter. The required skin incision should be made in form of a lobule pedicled towards the shunting catheter or by a straight skin incision. If a Burrhole Reservoir - or a SPRUNG RESERVOIR - is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole.

M.blue XABO is available in a range of different configurations: If using a Burrhole Reservoir - or a SPRUNG RESERVOIR -, the XABO Ventricular Catheter is implanted first. Once the mandrin has been removed, the patency of the XABO Ventricular Catheter can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the Burrhole Reservoir - or the SPRUNG RESERVOIR - connected and the connection secured with a ligature.

When using a shunt system with a CONTROL RESERVOIR, a Burrhole Deflector is included. The Burrhole Deflector is used to adjust the length of the catheter to be implanted and to position it inside the ventricle. The Ventricular Catheter is deflected by 90° and the CONTROL RESERVOIR put into place. The position of the XABO Ventricular Catheter should be inspected after the procedure by imaging (such as CT or MRI).

### Positioning of the valve

A location behind the ear is suitable as an implantation position, whereby the implantation height has no influence on the function of the valve.

The adjustable valve should be touching the bone or the periosteum since pressure must be exerted on the valve during any later adjustment.

A large arch-shaped or a small straight skin cut with a pocket for the valve should be made. The catheter is pushed forward from the burr hole to the selected valve implantation location, shortened if necessary, and secured to the M.blue XABO by ligation. The valve should not be located directly under the skin incision.

The valve unit has an arrow in the flow direction (arrow towards distal or downwards). The surface of the valve with the arrow markings points to the outside.

### WARNING

**The gravitational unit of the M.blue XABO is posture-dependent. For that reason, care must be taken to implant the gravitational unit is implanted parallel to the body axis.**

### Positioning of the XABO Peritoneal Catheter

The access site for the XABO Peritoneal Catheter is left to the surgeon's discretion. For example, it can be used in a paraumbilical application or applied at the level of the epigastrium. Likewise, various surgical techniques are available for placing the XABO Peritoneal Catheter. The recommendation is to pull the XABO Peritoneal Catheter from the valve to the intended position using a subcutaneous Tunneler, if necessary with the aid of an auxiliary incision. The XABO Peritoneal Catheter, usually securely attached to the valve, has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the XABO Peritoneal Catheter (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

## 4.5.6 POSTOPERATIVE VALVE TEST

### Postoperative valve test

The M.blue XABO has been constructed as a reliably functioning unit without pump or test function. The valve test can be performed by flushing, pressure measurement or pumping via a Reservoir or a Prechamber.

## 4.6 EXPLANTATION AND DISPOSAL

### 4.6.1 EXPLANTATION

The explantation of the M.blue XABO should be performed according to the state of the art and in compliance with medical practice.

### 4.6.2 DISPOSAL

#### M.blue XABO and shunt components


Products and product parts not used in the implantation or surgically removed must be disposed of correctly as potentially infectious material in accordance with medical practice

as well as respective regional laws and regulations.

Explanted medical devices must not be reused.


## 4.7 TECHNICAL INFORMATION

### 4.7.1 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>M.blue XABO</i>
Medical Purpose	Shunting of cerebrospinal fluid (CSF)
Sterilisability	Cannot be resterilised
Storage	Store in a clean and dry place at $\leq 30$ °C
For single use only	
Sketch with outer dimensions:	
 <p>The sketch shows two views of the device. The top view is a circular blue disc with a diameter of 25 mm and a central circular area with a diameter of 16.6 mm. The side view shows a blue cylindrical component with a height of 4.2 mm. The device has two small protrusions on each side, one on the top view and one on the side view.</p>	

## 4.8 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxxx indicates the identifier of the responsible notified body
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Serial number
	Sterilised using irradiation
	Do not re-sterilise
	Do not reuse
	Do not use if package is damaged and consult instructions for use
	Keep dry
	Upper limit of temperature
	Keep away from sunlight
	Consult instructions for use / electronic instructions for use
	Caution
	Contains a medical substance

Symbol	Explanation
	Non-pyrogenic
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians.
	MR Conditional
	Patient identification
	Date
	Healthcare centre or doctor
	Patient information website
	Model number / European Medical Device Nomenclature Code

## 5 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany