M.blue plus[®] *XABO*[®]

THE BALANCED WAY OF LIFE
INSPIRED BY YOU

- DE Gebrauchsanweisung | EN Instructions for use | FR Mode d'emploi |
ES Instrucciones de manejo | PT Instruções de utilização | IT Istruzioni per l'uso |
NL Gebruiksaanwijzing |

 www.miethke.com

US This Instructions for Use is NOT intended for United States users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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info@miethke.com

CONTENTS

1	PREFACE AND IMPORTANT INFORMATION	22
2	INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	22
2.1	EXPLANATION OF THE WARNINGS	22
2.2	DISPLAY CONVENTIONS	22
2.3	OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION	22
2.4	FEEDBACK ON THE INSTRUCTIONS FOR USE	22
2.5	COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION	22
3	DESCRIPTION OF THE <i>M.blue plus XABO</i>	23
3.1	MEDICAL PURPOSE	23
3.2	CLINICAL BENEFITS	23
3.3	INDICATIONS	23
3.4	CONTRAINDICATIONS	23
3.5	INTENDED PATIENT GROUPS	23
3.6	INTENDED USERS	23
3.7	INTENDED USE ENVIRONMENT	23
3.8	TECHNICAL DESCRIPTION	24
3.9	FUNCTIONING OF THE VALVE SYSTEM	25
3.10	SELECTION OF THE APPROPRIATE PRESSURE SETTING	26
3.11	PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES	26
3.12	PRESSURE-FLOW CHARACTERISTICS	27
3.13	APPLICATION OF THE M.BLUE PLUS INSTRUMENTS	28
3.14	SYSTEM COMPONENT	33
3.15	FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	33
4	PROPERTIES OF THE <i>M.blue plus XABO</i>	33
4.1	PRODUCT DESCRIPTION	33
4.2	IMPORTANT SAFETY INFORMATION	34
4.3	PATIENT EDUCATION	35
4.4	TRANSPORT AND STORAGE	35
4.5	USE OF THE PRODUCT	35
4.6	EXPLANTATION AND DISPOSAL	37
4.7	TECHNICAL INFORMATION	38
4.8	SYMBOLS USED FOR LABELLING	39
5	MEDICAL DEVICE CONSULTANTS	39

1 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing the medical device *M.blue plus XABO*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product. Your team at Christoph Miethke GmbH & Co. KG

Relevance of the instructions for use

WARNING

Improper handling and use of this product can cause hazards and damages. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

M.blue plus XABO includes the following components:

- ▶ *M.blue plus XABO*

Additional options:

- ▶ *Reservoirs* (including paediatric version)
- ▶ *Prechamber* (including paediatric version)
- ▶ *Burrhole Deflector* (including paediatric version)
- ▶ *Ventricular Catheter*
- ▶ *Titanium Connectors*

2 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

2.1 EXPLANATION OF THE WARNINGS

DANGER

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

WARNING

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

CAUTION

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

NOTE

Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for danger, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

2.2 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates <i>product names</i>

2.3 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

Each delivery is accompanied by a patient card with information about the product. The patient card is to provide the treating physician with all the product information in a compact form for the patient record.

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

2.4 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

2.5 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

3 DESCRIPTION OF THE *M.blue plus XABO*

3.1 MEDICAL PURPOSE

M.blue plus XABO is intended for shunting of cerebrospinal fluid (CSF).

3.2 CLINICAL BENEFITS

Clinical benefits of the *M.blue plus XABO*
:

- ▶ Implantation of a long-term implant for the extraction of cerebrospinal fluid from the CSF circulation and its shunting into the peritoneum
- ▶ Treatment of hydrocephalus, e.g. by alleviation of the clinical symptoms
- ▶ Reducing the risk of infection with gram-positive bacteria thanks to antibiotic-impregnated catheters

3.3 INDICATIONS

The following indications apply to *M.blue plus XABO*:

- ▶ Treatment of hydrocephalus

3.4 CONTRAINDICATIONS

The following contraindications apply to *M.blue plus XABO*:

- ▶ Infections in the implantation area
- ▶ Pathological concentrations levels (e.g. of blood components and/or protein) within the cerebrospinal fluid
- ▶ Intolerance to materials of the shunt system
- ▶ Hypersensitivity to rifampicin and/or clindamycin hydrochloride

3.5 INTENDED PATIENT GROUPS

- ▶ Patients who are treated with a CSF shunt system on account of their clinical picture

3.6 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical professionals, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product
- ▶ Successful participation in product training

3.7 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions

3.8 TECHNICAL DESCRIPTION


Fig. 1: M.blue plus cross section

1. Adjustable differential pressure unit

2. Sapphire ball
3. Torsion spring
4. Rotor

5. Adjustable gravitational unit

6. Sapphire ball
7. Tantalum weight
8. Rotor

The *M.blue plus* is a valve system made of titanium. It consists of an adjustable differential pressure unit (hereinafter also *proGAV 2.0*) and an adjustable gravitational unit (hereinafter also *M.blue*) (Fig. 1).

The adjustable differential pressure unit (1) in the proximal part of the valve system consists of a stable titanium housing with an integrated ball-cone valve (2) in the front section. A torsion spring (3) determines the opening pressure of this unit. The pretension of the spring, and thus the valve opening pressure, can be adjusted post-surgically through the skin using a pivoted rotor (4).

The adjustable gravitational unit (5) contains a tantalum weight (7), which holds a sapphire ball in the ball seat via a lever (6). Depending on the patient's body position, the influence of the tantalum weight on the sapphire ball changes, which also affects the valve opening pressure.

Via a rotor (8), the pretension of the torsion spring connected to the lever can be adjusted through the skin after surgery. This way, the impact of the tantalum weight on the sapphire ball can be influenced and, consequently, the valve opening pressure can be adjusted.

The *M.blue plus XABO* consists of an *M.blue plus* valve system in combination with *XABO Catheters*.

XABO Catheters are made from silicone, they are impregnated with antibiotics in a special treatment process and contain 0.054% rifampicin and 0.15% clindamycin hydrochloride. Laboratory tests have shown *XABO Catheters* to reduce colonisation with gram-positive bacteria of the silicone catheter surface. The laboratory tests were performed with *Staphylococcus aureus*, *Staphylococcus epidermidis* and *Bacillus subtilis*. Systemic therapeutic effects are highly unlikely as the amounts

of rifampicin and clindamycin hydrochloride contained in the catheter represent only a fraction of the therapeutic dose of these antibiotics.

3.9 FUNCTIONING OF THE VALVE SYSTEM

The *M.blue plus XABO* is a posture-dependent hydrocephalus valve system. The opening pressure of the *M.blue plus XABO* is composed of the opening pressures of the adjustable differential pressure unit and the adjustable gravitational unit.

Horizontal body position

In the horizontal position, the gravitational unit is always open and does not present any resistance.

The opening pressure of the *M.blue plus XABO* in the horizontal position is thus determined by the adjustable differential pressure unit. The basic operating principle of the differential pressure unit is shown in Fig. 2 a) and b).


Fig. 2: Functional principle of the differential pressure unit

a) closed b) open

In Fig. 2 a), the valve is closed; hence, shunting is not possible.

If the patient's intraventricular pressure (IVP) exceeds the spring force of the torsion spring, which otherwise keeps the differential pressure unit closed, the sealing ball moves out of the cone, leaving a gap for shunting of CSF (Fig. 2 b).

Vertical body position

As the patient sits or stands up, the gravitational unit closes the discharge channel in the proximal part of the valve (Fig. 3 a).

Thus, the opening pressure of the *M.blue plus XABO* is increased in the upright position, because now the weight of the tantalum weight (opening pressure of the gravitational unit) must be overcome in addition to the opening pressure of the differential pressure unit. Shunting is only possible again when the sum of IVP and hydrostatic suction is greater than the opening pressure of both units (Fig. 3 b).

For individual adaptation of the opening pressure to the patient, an opening pressure between 0 and 20 cmH₂O can be selected for the adjustable differential pressure unit, and additionally a valve opening pressure between 0 and 40 cmH₂O can be selected for the adjustable gravitational unit.


Fig. 3: Gravitational unit in the vertical body position
a) closed b) open

NOTE

During physical activity associated with shock (e.g. jogging), the opening pressure of the *M.blue plus XABO* may decrease temporarily according to laboratory results. Fundamentally, though, functionality remains the same. At the end of physical activity, the opening pressure returns to its original level and remains stable.

3.10 SELECTION OF THE APPROPRIATE PRESSURE SETTING

Recommended pressure settings for the *M.blue plus XABO* are available at:

<https://www.miethke.com/downloads/>

This is a non-binding recommendation for the treating physician. The physician decides on each case independently, individually and autonomously in accordance with his/her diagnosis.

The adjustable differential pressure unit of the *M.blue plus XABO* is set to an opening pressure of 5 cmH₂O upon delivery.

The adjustable gravitational unit of the *M.blue plus XABO* is set to an opening pressure of 20 cmH₂O upon delivery.

This pre-set opening pressure can be changed to a different pressure prior to implantation.

Horizontal body position

The opening pressure in the horizontal body position is determined by the differential pressure unit.

In this case, the pressure level should be set in accordance with the patient's clinical picture, indication and age. The standard settings for the horizontal position are 5 to 10 cmH₂O.

Vertical body position

The opening pressure of the *M.blue plus XABO* for the vertical body position is calculated from the sum of the opening pressure of the differential pressure unit and the gravitational unit. Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure level for the gravitational unit.

3.11 PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES

proGAV 2.0 differential pressure unit

The selected pressure level of the *M.blue plus* differential pressure unit should always be monitored using the *proGAV 2.0 Compass* or the *M.blue plus Compass*, but it can also be checked with the aid of an X-ray image (Fig. 4).


Fig. 4: X-ray image (*proGAV 2.0* adjustable differential pressure unit, set to 14 cmH₂O)

The rotor setting is decisive in this case. The four magnets in the rotor appear in the X-ray image as white dots and are located opposite each other in pairs. On one side of the rotor, two additional burr holes (right and left next to the magnet pairs) serve as orientation. They appear as black dots in the X-ray image. This side can be described as the rear side of the rotor. The two front magnets are on the opposite side.


Fig. 5: Schematic representation of rotor in an X-ray image

1. Admittance markings, 2. Valve marking
3. Triangle apex, 4. Non-adjustable range
5. proximal, 6. distal

The space between these two magnets can be considered as the apex of the triangle. The pressure level can be read off using the orientation of this intermediate space (Fig. 5). The apex of the triangle can take up any position except the space marked as a the non-adjustable range in Fig. 5. The opening pres-

sure of the *proGAV 2.0* can thus be infinitely variably adjusted from 0 up to 20 cmH₂O. To prevent reading the pressure level in reverse, the valve is fitted with a valve marking on one side, which appears black in the X-ray image – in a top view onto the implanted valve as in Fig. 4, the recess is visible on the right-hand side.

M.blue gravitational unit

The selected pressure level of the *M.blue* gravitational unit should always be monitored using the *M.blue plus Compass*, but it can also be checked using X-ray (Fig. 6).


Fig. 6: X-ray image (adjustable gravitational unit pre-set to 20 cmH₂O; differential pressure unit: 0 cmH₂O)

The rotor setting is decisive in this case. The four magnets in the rotor appear in the X-ray image as white dots and are located opposite each other in pairs. On one side of the rotor, two additional burr holes (right and left next to the magnet pairs) serve as orientation. They appear as black dots in the X-ray image. This side can be described as the rear side of the rotor. The two front magnets are on the opposite side. The space between these two magnets can be considered as the apex of the triangle. The pressure level can be read off using the orientation of this intermediate space (Fig. 7). The apex of the triangle can take up any position except the space marked as a the non-adjustable range in Fig. 7. This means that the opening pressure of the *M.blue* can be infinitely variably adjusted from 0 up to 40 cmH₂O. To prevent reading the pressure level in reverse, the top view of the valve shows an obvious recess in the housing ring with the tantalum weight to the right of the inlet connector (Fig. 6).


Fig. 7: Schematic representation of rotor in an X-ray image

- 1. Coding holes for the adjustable gravitational unit
- 2. Coding holes for the differential pressure unit
- 3. Triangle apex, 4. Non-adjustable range
- 5. proximal, 6. distal

3.12 PRESSURE-FLOW CHARACTERISTICS

The opening pressure refers to a reference flow of 20 ml/h. The supplied *XABO Catheters* do not fundamentally change the pressure-flow characteristics.

Horizontal body position

Fixed M.blue 0 differential pressure unit

The pressure flow characteristic for the fixed *M.blue 0* differential pressure unit is presented below.


Fig. 8: Pressure-flow characteristic of the M.blue 0 in the horizontal body position; pressure (cmH₂O), flow rate (ml/h)


Fig. 10: Pressure-flow characteristics of the M.blue in the vertical body position; pressure (cmH₂O), flow rate (ml/h)

Adjustable differential pressure unit of the proGAV 2.0

The pressure flow characteristics of the adjustable differential pressure unit are shown below as an example for pressure settings 0, 10 and 20 cmH₂O.


Fig. 9: Pressure-flow characteristics for selected pressure levels of the adjustable differential pressure unit; pressure (cmH₂O), flow rate (ml/h)

Vertical body position

M.blue adjustable gravitational unit

The pressure flow characteristics for various pressure level settings of the M.blue in the vertical body position are shown below. Only values of the adjustable gravitational unit are shown without consideration of the adjustable differential pressure unit.

3.13 APPLICATION OF THE M.blue plus Instruments

CAUTION

Only M.blue plus Instruments approved for that purpose must be used to determine, change and monitor the opening pressure of the gravitational unit of the (M.blue).

The opening pressure of the differential pressure unit (proGAV 2.0) can be determined, changed and monitored with the M.blue plus Instruments approved for that purpose or with the proGAV 2.0 Tools.

The M.blue plus Instruments may only be used by trained specialists.

The selected pressure level of the M.blue plus can be determined, adjusted and monitored using the M.blue plus Instruments.

The M.blue plus Compass (Fig. 11) is used for localising and reading the adjustable gravitational unit of (M.blue) as well as the adjustable differential pressure unit (proGAV 2.0) of the M.blue plus.


Fig. 11: M.blue plus Compass

a) open, 1. scale ring

b) closed, 2. float gauge compass

The M.blue plus Adjustment Ring (Fig. 12) is used to set the opening pressure of the gravitational unit of the M.blue from 0 to 40 cmH₂O. The M.blue plus Adjustment Ring can also be used to set the opening pressure of the differential pressure unit of the proGAV 2.0 from 0 to 20 cmH₂O.

Fig. 12: M.blue plus Adjustment Ring

The opening pressure of the adjustable gravitational unit of the M.blue can be changed before or after implantation. It is pre-set by the manufacturer to 20 cmH₂O.

The opening pressure of the adjustable differential pressure unit of the proGAV 2.0 can be changed before or after implantation. It is pre-set by the manufacturer to 5 cmH₂O.

In order to set the opening pressures, the following steps must be performed:

1. Localisation

WARNING

The M.blue plus Compass should be placed as centrally as possible over the valve, as the determined opening pressure may otherwise be incorrect.

NOTE

- ▶ Possible air pockets in the M.blue plus Compass do not affect its function.
- ▶ The M.blue plus Compass reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the M.blue plus Adjustment Ring should not be placed in the immediate vicinity of the M.blue plus Compass when determining the opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin may make adjustment difficult for a few days after surgery. If the valve setting cannot be checked conclusively using the M.blue plus Compass, we recommend checking it with an imaging method.

When the M.blue plus Compass is opened, a circular cut-out becomes visible, which is used to localise the valve in the patient's head as centrally as possible using the index finger (Fig. 13).

Fig. 13: Localising the valve

The direction markings show the flow direction.

2. Test procedure

In order to determine the selected pressure level, the M.blue plus Compass is then closed

again. The float gauge should now be centred by moving the instrument within the designated circular marking (Fig. 14). Once the float gauge is centred, the current setting of the opening pressure of the differential pressure unit (*proGAV 2.0*) or the gravitational unit (*M.blue*) can be read from the line marking on the float gauge (Fig. 14).


Fig. 14: Determining the pressure setting with the *M.blue plus Compass*

There are two scales on the scale ring (Fig. 15). The grey adjustment range from 0 to 20 cmH₂O on the outer scale applies to the opening pressure of the adjustable differential pressure unit of the *proGAV 2.0*.

The blue adjustment range from 0 to 40 cmH₂O on the inner scale applies to the opening pressure of the gravitational unit of the *M.blue*.


Fig. 15: Scale ring of the *M.blue plus Compass*

1. Inner: Scale for the *M.blue* gravitational unit from 0 to 40 cmH₂O (the opening pressure of the *M.blue* gravitational unit in the illustrated example is 16 cmH₂O)
2. Outer: Scale for the *proGAV 2.0* differential pressure unit from 0 to 20 cmH₂O (the opening pressure of the *proGAV 2.0* differential pressure unit in the illustrated example is 17 cmH₂O)

3. Adjustment process

WARNING

When adjusting the differential pressure unit of the *proGAV 2.0*, care must be taken to change the opening pressure by a maximum of 8 cmH₂O per adjustment process; otherwise, errors can result. Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in two stages: an initial adjustment from 3 to 11 cmH₂O and subsequently from 11 to 18 cmH₂O.

WARNING

When adjusting the gravitational unit of the *M.blue*, care must be taken to change the opening pressure by a maximum of 16 cmH₂O per adjustment process; otherwise, errors can result. Example: The opening pressure is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: an initial adjustment from 6 to 22 cmH₂O and subsequently from 22 to 36 cmH₂O.

In order to adjust the opening pressure, the *M.blue plus Compass* is opened, but without changing the position of the scale ring. The *M.blue plus Adjustment Ring* is now inserted into the scale ring in such a manner that the line marking points to the desired value on the scale of the scale ring (Fig. 16).


Fig. 16: Inserting the M.blue plus Adjustment Ring
1. M.blue plus Adjustment Ring, 2. Scale ring


Fig. 17: The example illustrates setting the proGAV 2.0 differential pressure unit to 1 cmH₂O or setting the M.blue gravitational unit to 32 cmH₂O.

By applying slight pressure with the index finger to the valve diaphragm located in the centre of the M.blue plus Adjustment Ring and under the skin, the rotor brake is released and the opening pressure of the differential pressure or gravitational unit is changed to the desired value (Fig. 18).

Fig. 18: Adjustment using the M.blue plus Adjustment Ring

Both the differential pressure unit of the proGAV 2.0 and the gravitational unit of the M.blue are equipped with a feedback mechanism.

Due to the valve housing design, targeted pressure on the valve produces an audible acoustic signal (a clicking sound) and/or palpable resistance as soon as the rotor brake has been released. The valve thus shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released, the rotor is once again adjustment-proof. Although the click caused by releasing the rotor brake is easily audible before implantation, it may be considerably reduced after implantation and the filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patient or by using a stethoscope.

Adjustment using the M.blue plus Adjustment Assistant

Alternatively, the M.blue plus Adjustment Assistant can be used to adjust the opening pressure. To do this, insert the M.blue plus Adjustment Assistant into the M.blue plus Adjustment Ring aligned to the desired value and press it with your index finger (Fig. 19).


Fig. 19: M.blue plus Adjustment Assistant

Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the set pressure level. To do this, proceed as in Points 1 and 2. Should the measured value not correspond with the required pressure level, the adjustment process should be repeated. To do this, start again at Point 3.

M.blue Checkmate

The M.blue Checkmate (Fig. 20) is supplied sterile and can be re-sterilised. The M.blue Checkmate can be used to change and check

the pressure level before and during valve implantation directly at the *M.blue*. To determine the pressure level, the *M.blue Checkmate* is placed centrally onto the *M.blue*. The *M.blue Checkmate* automatically aligns itself over the valve. The pressure level can be read from the direction of the proximal catheter (leading towards the valve). If the pressure level is to be adjusted, the *M.blue Checkmate* is placed centrally onto the *M.blue*. When doing so, the required pressure level must point towards the proximal catheter (leading towards the valve). By slightly pressing the *M.blue Checkmate* onto the valve, the rotor brake in the *M.blue* is released and the pressure level set.

Fig. 20: *M.blue Checkmate*, colour: blue
Pressure levels: 0 to 40 cmH₂O

proGAV Checkmate

The *proGAV Checkmate* (Fig. 21) is supplied sterile and can be re-sterilised. The *proGAV Checkmate* can be used to change and check the pressure level before and during valve implantation directly on the *proGAV 2.0*. To determine the pressure level, the *proGAV Checkmate* is placed centrally on the *proGAV 2.0*. The *proGAV Checkmate* automatically aligns itself over the valve. The pressure level can be read from the direction of the proximal catheter (leading towards the valve). For pressure level adjustment, the *proGAV Checkmate* is placed centrally on the *proGAV 2.0*. When doing so, the required pressure level must point towards the proximal catheter (leading towards the valve). Lightly pressing the *proGAV Checkmate* onto the valve releases the rotor brake in the *proGAV 2.0* and sets the pressure level.

Fig. 21: *proGAV Checkmate*,
Pressure levels 0-20 cmH₂O

Checking and adjusting without unpacking

To protect the impregnation, the *M.blue plus XABO* is wrapped in opaque packaging. However, it is still possible to check and adjust the setting of the *M.blue plus XABO*.

The positions of the valves of the *M.blue plus XABO* (*M.blue* and *proGAV 2.0*) are indicated by two arrow markings. The arrows indicate the flow direction.

For localisation, the *M.blue plus Compass* is placed centrally on the marking of the valve that needs to be checked. Because of the packaging, the *M.blue plus Compass* is slightly tilted on the packaging.

Checking and/or adjusting without unpacking retains the sterility of the product and the integrity of the packaging.


Fig. 22: Packaging of the *M.blue plus XABO*: marking of the adjustable valves with an arrow indicating the flow direction


Fig. 23: Position of the *M.blue plus Compass* on the packaging of the *M.blue plus XABO*

3.14 SYSTEM COMPONENT

Combination with shunt components

The *M.blue plus XABO* can be safely combined with our range of implantable shunt components. We recommend using the Christoph Miethke GmbH & Co. KG products in combination with the *M.blue plus XABO*.

Reservoirs

The use of shunt systems with a reservoir provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to an integrated check valve in the *CONTROL RESERVOIR* and the *SPRUNG RESERVOIR*, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the shunt system as well as the *Ventricular Catheter*.

During the pump action, access to the *Ventricular Catheter* is closed. The use of this Reservoir does not increase the opening pressure of the shunt system. Puncturing the Reservoir should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. A firm titanium floor prevents the bottom surface from being pierced. 30 punctures are possible without any restrictions.

WARNING

Frequent pumping of the Reservoir can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.

Burrhole Deflector

Because of the tight fit on the *Ventricular Catheter*, the *Burrhole Deflector* makes it possible to choose the length of catheter penetrating into the skull prior to implantation. The *Ventricular Catheter* is deflected at a right angle in the burr hole (see chapter 4.5.5).

Tube systems

The *M.blue plus XABO* is supplied as a shunt system with integrated catheters (internal diameter 1.2 mm, external diameter 2.5 mm).

For a new connection of catheter and connector, the catheters have to be carefully fixed with a ligature to the valve's *Titanium Connectors*.

3.15 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

- ▶ The medical device has been designed to work accurately and reliably over long periods of time. We guarantee that our medical devices are free from defects and fully functional at the time of dispatch. Excluded from this guarantee are cases in which the medical device has been exchanged or rather explanted for technical or medical reasons outside of our sphere of responsibility.
- ▶ *M.blue plus XABO* together with the entire shunt system can safely withstand pressure of up to 100 cm H₂O occurring during and after surgery.
- ▶ Nuclear magnetic resonance examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. *M.blue plus XABO* is MR Conditional. The supplied catheters are MR Safe. *Reservoirs*, *Burrhole Deflectors* and *Titanium Connectors* are MR Conditional.

The conditions for MRI security of the products are found on our website: <https://miethke.com/downloads/>

4 PROPERTIES OF THE *M.blue plus XABO*

4.1 PRODUCT DESCRIPTION

4.1.1 *M.blue plus XABO* VARIANTS

The *M.blue plus XABO* has an adjustable differential pressure unit and an adjustable gravitational unit.

Adjustable differential pressure unit	Adjustable gravitational unit
0–20 cmH ₂ O	0–40 cmH ₂ O

The *M.blue plus XABO* is also available as a shunt system and can come with the follow-

ing components: *XABO Ventricular Catheter*, *Prechambers*, *Reservoirs*.

4.1.2 SCOPE OF DELIVERY

Box content	Number
Sterile packaging with <i>M.blue plus XABO</i> shunt system	1
Instructions for use for <i>M.blue plus XABO</i>	1
Patient card	1
Recommended pressure level	1

4.1.3 STERILITY

WARNING

The product must not be used if the sterile packaging or the product is damaged or after the expiry date.

M.blue plus XABO is sterilised under strictly controlled conditions using irradiation. The respective expiry date is printed on the packaging. The inner sterile packaging consists of a gas- and liquid-tight barrier to protect the antibiotics.

4.1.4 REPEATED USE AND RESTERILISATION

WARNING

The product must not be resterilised or reprocessed in any other way as the safe functioning and sterility of the product cannot be guaranteed.

Products that have already been implanted in a patient must not be reused either on the same or a different patient in order to minimise the risk of infection.

4.1.5 SINGLE-USE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the *M.blue plus XABO*. No guarantee can be assumed for the functional safety of resterilised products.

4.1.6 PRODUCT CONFORMITY

The product meets current regulatory requirements.

The requirements stipulate the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the medical device should therefore be recorded in the patient's medical records to ensure complete traceability.

4.2 IMPORTANT SAFETY INFORMATION

4.2.1 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ **The product must not be used if the sterile packaging or the product is damaged or after the expiry date.**
- ▶ **Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be carefully read and understood before the product is used for the first time.**
- ▶ **Prior to use, it is essential to check the product for completeness and integrity.**

4.2.2 COMPLICATIONS, SIDE EFFECTS, PRECAUTIONS AND RESIDUAL RISKS

The following complications can occur in conjunction with the *M.blue plus XABO*:

- ▶ Headaches, dizzy spells, mental confusion, vomiting in cases of possible leakage from the shunt system and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid
- ▶ Allergic reaction / intolerance to product materials
- ▶ Overdrainage/underdrainage
- ▶ Noise development

Violent external shocks (accident, fall, etc.) may put the integrity of the shunt system at risk. As a precaution, a physician must be consulted immediately if the patient suffers from

skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the *M.blue plus XABO*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation
- ▶ Local shunt irritation
- ▶ Allergic reactions to catheter components, particularly to the active antibiotic substances rifampicin and clindamycin hydrochloride

4.2.3 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible state authority.

4.3 PATIENT EDUCATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product (Ch. 4.2, 4.5).

4.4 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place. *M.blue plus XABO* is to be protected from direct sunlight. The products should not be removed from the packaging until they are required.

4.4.1 TRANSPORT

Transport conditions

Ambient temperature	≤ 40 °C
----------------------------	---------

4.4.2 STORAGE

Storage conditions

Temperature range for storage	≤ 30 °C
--------------------------------------	---------

4.5 USE OF THE PRODUCT

4.5.1 INTRODUCTION

The *M.blue plus XABO* is a posture-dependent valve system with an adjustable gravitational unit and an adjustable differential pressure unit combined with antibiotic impregnated *XABO Catheters*.

M.blue plus XABO is used for shunting of cerebrospinal fluid (CSF) in the treatment of hydrocephalus. Valves and Reservoirs are placed in suitable positions along the course of the shunt.

4.5.2 SAFETY NOTICES AND WARNINGS

WARNING

- ▶ *M.blue plus XABO* must not be immersed in antibiotic solutions. The contact of the catheter with sterile water or sterile saline solution should be limited to a minimal period of time and should be done immediately prior to implantation. The solution may take on a faint orange colour.
- ▶ Frequent pumping of the *Reservoir* can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.
- ▶ The gravitational unit of the *M.blue plus XABO* is posture-dependent. For that reason, care must be taken to implant the gravitational unit (*M.blue*) is implanted parallel to the body axis.
- ▶ The adjustable valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e. g. underneath heavily scarred tissue). If the implantation site is unfavourably chosen or the skin above the valve is too thick, it may no longer be possible to adjust the adjustable unit. The valve then works with unchangeable pressure levels.

⚠ WARNING

- ▶ If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, and adjustment of the valve cannot be ruled out.

⚠ CAUTION

- ▶ Silicone is extremely electrostatic. Care must be taken to avoid the catheters coming into contact with dry cloths, talcum powder or rough surfaces. Clinging particles could lead to tissue reactions.
- ▶ When using sharp instruments, care should be taken to avoid cuts and scratches in the silicone elastomer.
- ▶ It must be ensured that the ligature is not tightened excessively. Damage may result in a loss of integrity of the shunt and thus necessitate a revision.
- ▶ The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.
- ▶ For people using cardiac pacemakers: It is possible that the function of the pacemaker is influenced by the implantation of a *M.blue plus XABO*.

! NOTE

- ▶ In MRI imaging *M.blue plus XABO* creates artefacts that are larger than the valve itself.

4.5.3 REQUIRED MATERIALS

The *M.blue plus XABO* is designed so that it can be safely used with the shunt components described in chapter 3.14. Catheters with an inner diameter of 1.2 mm and an outer diameter of approx. 2.5 mm should be used for connection. In any case, catheters have to be carefully fixed with a ligature to the connectors of the shunt components. Any kinks in the catheter must be avoided.

4.5.4 PREPARING FOR IMPLANTATION

Checking the sterile packaging

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier

system. The products should only be removed from the packaging immediately prior to use.

Preoperative valve test

The *M.blue plus XABO* should be vented before implantation and checked for permeability. The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological salt solution. The valve is patent if saline solution can be extracted (Fig. 24).


Fig. 24: Patency test

⚠ WARNING

- ▶ No antibiotic solution must be used for valve testing and/or venting, as this may result in a reaction with the active substances in the impregnation.
- ▶ Contaminants in the solution used for testing can impair the product's performance.
- ▶ Pressurisation with a single-use syringe should be avoided both at the proximal and the distal end. (Fig. 25).


Fig. 25: Avoidance of pressurisation

4.5.5 PERFORMING THE IMPLANTATION

Positioning of the *XABO Ventricular Catheter*

Several surgical techniques are available for the positioning of the *XABO Ventricular Catheter*. The required skin incision should be made in form of a lobule pedicled towards the shunting catheter or by a straight skin incision. If a *Burhole Reservoir* - or a *SPRUNG RESERVOIR* - is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage,

care should be taken that the dura opening is kept as small as possible after applying the burr hole.

M.blue plus XABO is available in a range of different configurations: If using a *Burrhole Reservoir* - or a *SPRUNG RESERVOIR* -, the *XABO Ventricular Catheter* is implanted first. Once the mandrin has been removed, the patency of the *XABO Ventricular Catheter* can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the *Burrhole Reservoir* - or the *SPRUNG RESERVOIR* - connected and the connection secured with a ligature.

When using a shunt system with a *CONTROL RESERVOIR*, a *Burrhole Deflector* is included. The *Burrhole Deflector* is used to adjust the length of the catheter to be implanted and to position it inside the ventricle. The *Ventricular Catheter* is deflected by 90° and the *CONTROL RESERVOIR* put into place. The position of the *XABO Ventricular Catheter* should be inspected after the procedure by imaging (such as CT or MRI).

Positioning of the valve system

A location behind the ear is suitable as an implantation position, whereby the implantation height has no influence on the function of the valve system.

The adjustable valve should be touching the bone or the periosteum since pressure must be exerted on the valve during any later adjustment.

A large arch-shaped or a small straight skin cut with a pocket for the valve system should be made. The catheter is pushed forward from the burr hole to the selected valve implantation location, shortened if necessary, and secured to the *M.blue plus XABO* by ligation. The valve system should not be located directly under the skin incision. The valve unit has an arrow in the flow direction (arrow towards distal or downwards). The surface of the valve with the arrow markings points to the outside.

WARNING

The gravitational unit of the *M.blue plus XABO* is posture-dependent. For that reason, care must be taken to implant the gravitational unit (*M.blue*) is implanted parallel to the body axis.

Positioning of the *XABO Peritoneal Catheter*

The access site for the *XABO Peritoneal Catheter* is left to the surgeon's discretion. For example, it can be used in a paraumbilical application or applied at the level of the epigastrium. Likewise, various surgical techniques are available for placing the *XABO Peritoneal Catheter*. The recommendation is to pull the *XABO Peritoneal Catheter* from the valve to the intended position using a subcutaneous *Tunneller*, if necessary with the aid of an auxiliary incision. The *XABO Peritoneal Catheter*, usually securely attached to the valve, has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the *XABO Peritoneal Catheter* (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

4.5.6 POSTOPERATIVE VALVE TEST

Postoperative valve test

The *M.blue plus XABO* has been constructed as a reliably functioning unit without pump or test function. The valve test can be performed by flushing, pressure measurement or pumping via a *Reservoir* or a *Prechamber*.

4.6 EXPLANTATION AND DISPOSAL

4.6.1 EXPLANTATION

The explantation of the *M.blue plus XABO* should be performed according to the state of the art and in compliance with medical practice.

4.6.2 DISPOSAL


***M.blue plus XABO* and shunt components**

Products and product parts not used in the implantation or surgically removed must be disposed of correctly as potentially infectious material in accordance with medical practice as well as respective regional laws and regulations.

Explanted medical devices must not be reused.


4.7 TECHNICAL INFORMATION

4.7.1 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>M.blue plus XABO</i>
Medical Purpose	Shunting of cerebrospinal fluid (CSF)
Sterilisability	Cannot be resterilised
Storage	Store in a clean and dry place at $\leq 30\text{ }^{\circ}\text{C}$
For single use only	
Sketch with outer dimensions:	
 <p>The sketch illustrates the components and dimensions of the M.blue plus XABO shunt. It includes two views of the valve heads: a side view of the grey valve with a diameter of 17 mm and a top view of the blue valve with a diameter of 16.6 mm. Below these, a side view of the complete shunt assembly is shown, featuring a grey 'proGAV 2.0' valve with a diameter of 4.5 mm and a blue 'M.blue plus' valve with a diameter of 4.2 mm.</p>	

4.8 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxx indicates the identifier of the responsible notified body
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Serial number
	Sterilised using irradiation
	Do not resterilise
	Do not reuse
	Do not use if package is damaged and consult instructions for use
	Keep dry
	Upper limit of temperature
	Keep away from sunlight
	Consult instructions for use / electronic instructions for use
	Caution
	Contains a medical substance

Symbol	Explanation
	Non-pyrogenic
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians.
	MR Conditional
	Patient identification
	Date
	Healthcare centre or doctor
	Patient information website
	Model number / European Medical Device Nomenclature Code

5 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany