GAV 2.0 XABO[®]

- ⒹE Gebrauchsanweisung | ⒺN Instructions for use | ⒻR Mode d'emploi |
ⒺS Instrucciones de manejo | ⒻT Instruções de utilização | ⒿT Istruzioni per l'uso |
ⒹL Gebruiksaanwijzing |

 www.miethke.com

ⒺS This Instructions for Use is NOT intended for United States users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0
info@miethke.com

CONTENTS

1	PREFACE AND IMPORTANT INFORMATION	16
2	INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	16
2.1	EXPLANATION OF THE WARNINGS	16
2.2	DISPLAY CONVENTIONS	16
2.3	OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION	16
2.4	FEEDBACK ON THE INSTRUCTIONS FOR USE	16
2.5	COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION	16
3	DESCRIPTION OF THE GAV 2.0 XABO	17
3.1	MEDICAL PURPOSE	17
3.2	CLINICAL BENEFITS	17
3.3	INDICATIONS	17
3.4	CONTRAINDICATIONS	17
3.5	INTENDED PATIENT GROUPS	17
3.6	INTENDED USERS	17
3.7	INTENDED USE ENVIRONMENT	17
3.8	TECHNICAL DESCRIPTION	18
3.9	OPERATING PRINCIPLE OF THE VALVE	19
3.10	SELECTION OF THE APPROPRIATE PRESSURE SETTING	19
3.11	PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES	20
3.12	PRESSURE-FLOW CHARACTERISTICS	20
3.13	SYSTEM COMPONENT	21
3.14	FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	22
4	PROPERTIES OF THE GAV 2.0 XABO	22
4.1	PRODUCT DESCRIPTION	22
4.2	IMPORTANT SAFETY INFORMATION	23
4.3	PATIENT EDUCATION	23
4.4	TRANSPORT AND STORAGE	23
4.5	USE OF THE PRODUCT	24
4.6	EXPLANATION AND DISPOSAL	26
4.7	TECHNICAL INFORMATION	26
4.8	SYMBOLS USED FOR LABELLING	27
5	MEDICAL DEVICE CONSULTANTS	27

1 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing the medical device GAV 2.0 XABO. Please contact us if you have any questions about the contents of these instructions for use or the use of the product. Your team at Christoph Miethke GmbH & Co. KG

Relevance of the instructions for use

WARNING

Improper handling and use of this product can cause hazards and damages. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

GAV 2.0 XABO includes the following components:

- ▶ GAV 2.0 XABO

Additional options:

- ▶ Reservoirs (including paediatric version)
- ▶ Prechamber (including paediatric version)
- ▶ Burrhole Deflector (including paediatric version)
- ▶ Ventricular Catheter
- ▶ Titanium Connectors

2 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

2.1 EXPLANATION OF THE WARNINGS

DANGER

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

WARNING

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

CAUTION

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

NOTE

Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for danger, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

2.2 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates <i>product names</i>

2.3 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

Each delivery is accompanied by a patient card with information about the product. The patient card is to provide the treating physician with all the product information in a compact form for the patient record.

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

2.4 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

2.5 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

3 DESCRIPTION OF THE GAV 2.0 XABO

3.1 MEDICAL PURPOSE

GAV 2.0 XABO is intended for shunting of cerebrospinal fluid (CSF).

3.2 CLINICAL BENEFITS

Clinical benefits of the GAV 2.0 XABO :

- ▶ Implantation of a long-term implant for the extraction of cerebrospinal fluid from the CSF circulation and its shunting into the peritoneum
- ▶ Treatment of hydrocephalus, e.g. by alleviation of the clinical symptoms
- ▶ Reducing the risk of infection with gram-positive bacteria thanks to antibiotic-impregnated catheters

3.3 INDICATIONS

The following indications apply to GAV 2.0 XABO:

- ▶ Treatment of hydrocephalus

3.4 CONTRAINDICATIONS

The following contraindications apply to GAV 2.0 XABO:

- ▶ Infections in the implantation area
- ▶ Pathological concentrations levels (e.g. of blood components and/or protein) within the cerebrospinal fluid
- ▶ Intolerance to materials of the shunt system
- ▶ Hypersensitivity to rifampicin and/or clindamycin hydrochloride

3.5 INTENDED PATIENT GROUPS

- ▶ Patients who are treated with a CSF shunt system on account of their clinical picture

3.6 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical professionals, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product
- ▶ Successful participation in product training

3.7 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions

3.8 TECHNICAL DESCRIPTION


Fig. 1: GAV 2.0 cross section

1. Differential pressure unit

- 2. Sapphire ball
- 3. Micro-coil spring

The GAV 2.0 is a valve made of titanium. It consists of a differential pressure unit (1) and a gravitational unit (4). This ensures physiological intraventricular pressure (IVP) in any body position.

In the proximal part of the GAV 2.0, a micro-coil spring (3) ensures the opening pressure for the differential pressure unit. The gravitational unit in the distal part of the valve consists of a tantalum ball (5), which determines the opening pressure of this unit, and a sapphire ball (6), which ensures its precise closure. Coding (7)

4. Gravitational unit

- 5. Tantalum ball
- 6. Sapphire ball
- 7. X-ray coding

facilitates identification of pressure levels in an X-ray image.

The GAV 2.0 XABO consists of an GAV 2.0 valve in combination with XABO Catheters.

XABO Catheters are made from silicone, they are impregnated with antibiotics in a special treatment process and contain 0.054% rifampicin and 0.15% clindamycin hydrochloride. Laboratory test have shown XABO Catheters to reduce colonisation with gram-positive bacteria of the silicone catheter surface. The laboratory tests were performed with *Staphylococcus aureus*, *Staphylococcus epidermidis* and *Bacillus subtilis*. Systemic therapeutic effects are highly unlikely as the amounts

of rifampicin and clindamycin hydrochloride contained in the catheter represent only a fraction of the therapeutic dose of these antibiotics.

3.9 OPERATING PRINCIPLE OF THE VALVE

The GAV 2.0 XABO is a posture-dependent valve.

Horizontal body position

In the horizontal position, the gravitational unit is always open and does not present any resistance.

The opening pressure of the GAV 2.0 XABO in this body position is thus determined by the adjustable differential pressure unit. The basic operating principle of the differential pressure unit is shown in Fig. 2 a) and b).


Fig. 2: GAV 2.0 in the horizontal position
a) closed b) open

In Fig. 2 a), the differential pressure unit is closed; hence, shunting is not possible. If the patient's intraventricular pressure (IVP) exceeds the spring force of the micro-coil spring, the sealing ball moves out of the cone, leaving a gap for shunting of CSF (Fig. 2 b).

Vertical body position

As the patient sits or stands up, the gravitational unit is activated and the opening pressure of the GAV 2.0 XABO is significantly increased (Fig. 3 a). In addition to the opening pressure of the differential pressure unit, the weight of the tantalum ball (opening pressure of the gravitational unit) has to be exceeded. Shunting is only possible again when the sum of intraventricular pressure (IVP) and hydrostatic suction is greater than the opening pressure of both units (Fig. 3 b).


Fig. 3: GAV 2.0 in the vertical position
a) closed b) open

NOTE

During physical activity associated with shock (e.g. jogging), the opening pressure of the GAV 2.0 XABO may decrease temporarily according to laboratory results. Fundamentally, though, functionality remains the same. At the end of physical activity, the opening pressure returns to its original level and remains stable.

3.10 SELECTION OF THE APPROPRIATE PRESSURE SETTING

Recommended pressure settings for the GAV 2.0 XABO are available at:

<https://www.miethke.com/downloads/>

This is a non-binding recommendation for the treating physician. The physician decides on each case independently, individually and autonomously in accordance with his/her diagnosis.

Horizontal body position

The opening pressure in the horizontal body position is determined by the differential pressure unit.

Depending on the patient's clinical picture, indication and age, the opening pressure for this position can be selected between pressure levels from 5 to 10 cmH₂O.

Vertical body position

The opening pressure of the GAV 2.0 XABO for the vertical body position is calculated from the sum of the opening pressure of the differential pressure unit and the gravitational unit.

Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure level for the gravitational unit.




3.11 PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES

The pressure levels of the GAV 2.0 can be identified post-surgically by the coding in the X-ray image:

Pressure level		Coding
horizontal	vertical	
5 cmH ₂ O	20 cmH ₂ O	
5 cmH ₂ O	25 cmH ₂ O	
5 cmH ₂ O	30 cmH ₂ O	
5 cmH ₂ O	35 cmH ₂ O	
10 cmH ₂ O	25 cmH ₂ O	
10 cmH ₂ O	30 cmH ₂ O	

3.12 PRESSURE-FLOW CHARACTERISTICS

The following diagrams show the pressure-flow characteristics for the available pressure levels of the GAV 2.0 XABO. The opening pressure refers to a reference flow of 5 ml/h. For flow rates of 20 ml/h, stated pressures are approx. 1–2 cmH₂O higher. The supplied XABO catheters do not fundamentally change the pressure-flow characteristics.


Fig. 4: Pressure-flow characteristics of selected pressure levels of the GAV 2.0 in the horizontal and vertical body position; pressure (cmH₂O), flow rate (ml/h): a) GAV 2.0 5/20, b) GAV 2.0 5/25, c) GAV 2.0 5/30, d) GAV 2.0 5/35, e) GAV 2.0 10/25, f) GAV 2.0 10/30

3.13 SYSTEM COMPONENT

Combination with shunt components

The GAV 2.0 XABO can be safely combined with our range of implantable shunt components. We recommend using the Christoph Miethke GmbH & Co. KG products in combination with the GAV 2.0 XABO.

Reservoirs

The use of shunt systems with a reservoir provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to an integrated check valve in the CONTROL RESERVOIR and the SPRUNG RESERVOIR, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the shunt system as well as the Ventricular Catheter.

During the pump action, access to the Ventricular Catheter is closed. The use of this Reservoir does not increase the opening pressure of the shunt system. Puncturing the Reservoir should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. A firm titanium floor prevents the bottom surface from being pierced. 30 punctures are possible without any restrictions.

 **WARNING**

Frequent pumping of the Reservoir can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.

Burrhole Deflector

Because of the tight fit on the Ventricular Catheter, the Burrhole Deflector makes it possible to choose the length of catheter penetrating into the skull prior to implantation. The Ventricular Catheter is deflected at a right angle in the burr hole (see chapter 4.5.5).

Tube systems

The GAV 2.0 XABO is supplied as a shunt system with integrated catheters (internal diameter 1.2 mm, external diameter 2.5 mm).

For a new connection of catheter and connector, the catheters have to be carefully fixed with a ligature to the valve's *Titanium Connectors*.

3.14 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

- ▶ The medical device has been designed to work accurately and reliably over long periods of time. We guarantee that our medical devices are free from defects and fully functional at the time of dispatch. Excluded from this guarantee are cases in which the medical device has been exchanged or rather explained for technical or medical reasons outside of our sphere of responsibility.
- ▶ GAV 2.0 XABO together with the entire shunt system can safely withstand pressure of up to 100 cm H₂O occurring during and after surgery.
- ▶ Nuclear magnetic resonance examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. GAV 2.0 XABO is MR Conditional. The supplied catheters are MR Safe. *Reservoirs*, *Burhole Deflectors* and *Titanium Connectors* are MR Conditional.

The conditions for MRI security of the products are found on our website: <https://miethke.com/downloads/>

4 PROPERTIES OF THE GAV 2.0 XABO

4.1 PRODUCT DESCRIPTION

4.1.1 GAV 2.0 XABO VARIANTS

The GAV 2.0 XABO is available in a range of different variants. They differ in the pre-set pressure level depending on body position.

Horizontal	Vertical
5 cmH ₂ O	20 cmH ₂ O
5 cmH ₂ O	25 cmH ₂ O
5 cmH ₂ O	30 cmH ₂ O
5 cmH ₂ O	35 cmH ₂ O
10 cmH ₂ O	25 cmH ₂ O

Horizontal	Vertical
10 cmH ₂ O	30 cmH ₂ O

These variants are also available as shunt systems and can come with the following components: *XABO Ventricular Catheter*, *Prechambers*, *Reservoirs*.

4.1.2 SCOPE OF DELIVERY

Box content	Number
Sterile packaging with GAV 2.0 XABO shunt system	1
Instructions for use for GAV 2.0 XABO	1
Patient card	1
Recommended pressure level	1

4.1.3 STERILITY

WARNING

The product must not be used if the sterile packaging or the product is damaged or after the expiry date.

GAV 2.0 XABO is sterilised under strictly controlled conditions using irradiation. The respective expiry date is printed on the packaging. The inner sterile packaging consists of a gas- and liquid-tight barrier to protect the antibiotics.

4.1.4 REPEATED USE AND RESTERILISATION

WARNING

The product must not be resterilised or reprocessed in any other way as the safe functioning and sterility of the product cannot be guaranteed.

Products that have already been implanted in a patient must not be reused either on the same or a different patient in order to minimise the risk of infection.

4.1.5 SINGLE-USE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the GAV 2.0 XABO. No guarantee can be assumed for the functional safety of resterilised products.

4.1.6 PRODUCT CONFORMITY

The product meets current regulatory requirements.

The requirements stipulate the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the medical device should therefore be recorded in the patient's medical records to ensure complete traceability.

4.2 IMPORTANT SAFETY INFORMATION

4.2.1 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ **The product must not be used if the sterile packaging or the product is damaged or after the expiry date.**
- ▶ **Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be carefully read and understood before the product is used for the first time.**
- ▶ **Prior to use, it is essential to check the product for completeness and integrity.**

4.2.2 COMPLICATIONS, SIDE EFFECTS, PRECAUTIONS AND RESIDUAL RISKS

The following complications can occur in conjunction with the GAV 2.0 XABO:

- ▶ Headaches, dizzy spells, mental confusion, vomiting in cases of possible leakage from the shunt system and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid
- ▶ Allergic reaction / intolerance to product materials
- ▶ Overdrainage/underdrainage
- ▶ Noise development

Violent external shocks (accident, fall, etc.) may put the integrity of the shunt system at risk.

As a precaution, a physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the GAV 2.0 XABO:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation
- ▶ Local shunt irritation
- ▶ Allergic reactions to catheter components, particularly to the active antibiotic substances rifampicin and clindamycin hydrochloride

4.2.3 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible state authority.

4.3 PATIENT EDUCATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product (Ch. 4.2, 4.5).

4.4 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place. GAV 2.0 XABO is to be protected from direct sunlight. The products should not be removed from the packaging until they are required.

4.4.1 TRANSPORT

Transport conditions

Ambient temperature	≤ 40 °C
----------------------------	---------

4.4.2 STORAGE

Storage conditions

Temperature range for storage	≤ 30 °C
-------------------------------	---------

4.5 USE OF THE PRODUCT

4.5.1 INTRODUCTION

The GAV 2.0 XABO is a posture-dependent valve with a pre-set differential pressure unit and a pre-set gravitational unit combined with antibiotic impregnated XABO Catheters.

GAV 2.0 XABO is used for shunting of cerebrospinal fluid (CSF) in the treatment of hydrocephalus. Valves and Reservoirs are placed in suitable positions along the course of the shunt.

4.5.2 SAFETY NOTICES AND WARNINGS

WARNING

- ▶ GAV 2.0 XABO must not be immersed in antibiotic solutions. The contact of the catheter with sterile water or sterile saline solution should be limited to a minimal period of time and should be done immediately prior to implantation. The solution may take on a faint orange colour.
- ▶ Frequent pumping of the Reservoir can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.
- ▶ The gravitational unit of the GAV 2.0 XABO is posture-dependent. For that reason, care must be taken to implant the gravitational unit is implanted parallel to the body axis.

CAUTION

- ▶ Silicone is extremely electrostatic. Care must be taken to avoid the catheters coming into contact with dry cloths, talcum powder or rough surfaces. Clinging particles could lead to tissue reactions.
- ▶ When using sharp instruments, care should be taken to avoid cuts and scratches in the silicone elastomer.

CAUTION

- ▶ It must be ensured that the ligature is not tightened excessively. Damage may result in a loss of integrity of the shunt and thus necessitate a revision.
- ▶ The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

4.5.3 REQUIRED MATERIALS

The GAV 2.0 XABO is designed so that it can be safely used with the shunt components described in chapter 3.13. Catheters with an inner diameter of 1.2 mm and an outer diameter of approx. 2.5 mm should be used for connection. In any case, catheters have to be carefully fixed with a ligature to the connectors of the shunt components. Any kinks in the catheter must be avoided.

4.5.4 PREPARING FOR IMPLANTATION

Checking the sterile packaging

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Preoperative valve test

The GAV 2.0 XABO should be vented before implantation and checked for permeability. The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological salt solution. The valve is patent if saline solution can be extracted (Fig. 5).


Fig. 5: Patency test

**WARNING**

- ▶ **No antibiotic solution must be used for valve testing and/or venting, as this may result in a reaction with the active substances in the impregnation.**
- ▶ **Contaminants in the solution used for testing can impair the product's performance.**
- ▶ **Pressurisation with a single-use syringe should be avoided both at the proximal and the distal end. (Fig. 6).**


Fig. 6: Avoidance of pressurisation

4.5.5 PERFORMING THE IMPLANTATION

Positioning of the XABO Ventricular Catheter

Several surgical techniques are available for the positioning of the XABO Ventricular Catheter. The required skin incision should be made in form of a lobule pedicled towards the shunting catheter or by a straight skin incision. If a *Burrhole Reservoir* - or a *SPRUNG RESERVOIR* - is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole.

GAV 2.0 XABO is available in a range of different configurations: If using a *Burrhole Reservoir* - or a *SPRUNG RESERVOIR* -, the XABO Ventricular Catheter is implanted first. Once the mandrin has been removed, the patency of the XABO Ventricular Catheter can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the *Burrhole Reservoir* - or the *SPRUNG RESERVOIR* - connected and the connection secured with a ligature.

When using a shunt system with a *CONTROL RESERVOIR*, a *Burrhole Deflector* is included. The *Burrhole Deflector* is used to adjust the length of the catheter to be implanted and to position it inside the ventricle. The *Ventricular Catheter* is deflected by 90° and the *CONTROL RESERVOIR* put into place. The position of the XABO Ventricular Catheter should be inspected

after the procedure by imaging (such as CT or MRI).

Positioning of the valve

A location behind the ear is suitable as an implantation position, whereby the implantation height has no influence on the function of the valve.

A large arch-shaped or a small straight skin cut with a pocket for the valve should be made. The catheter is pushed forward from the burr hole to the selected valve implantation location, shortened if necessary, and secured to the GAV 2.0 XABO by ligation. The valve should not be located directly under the skin incision. The valve unit has an arrow in the flow direction (arrow towards distal or downwards). The surface of the valve with the arrow markings points to the outside.

**WARNING**

The gravitational unit of the GAV 2.0 XABO is posture-dependent. For that reason, care must be taken to implant the gravitational unit is implanted parallel to the body axis.

Positioning of the XABO Peritoneal Catheter

The access site for the XABO Peritoneal Catheter is left to the surgeon's discretion. For example, it can be used in a paraumbilical application or applied at the level of the epigastrium. Likewise, various surgical techniques are available for placing the XABO Peritoneal Catheter. The recommendation is to pull the XABO Peritoneal Catheter from the valve to the intended position using a subcutaneous *Tunneller*, if necessary with the aid of an auxiliary incision. The XABO Peritoneal Catheter, usually securely attached to the valve, has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the XABO Peritoneal Catheter (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

4.5.6 POSTOPERATIVE VALVE TEST

Postoperative valve test

The GAV 2.0 XABO has been constructed as a reliably functioning unit without pump or test function. The valve test can be performed by

flushing, pressure measurement or pumping via a *Reservoir* or a *Prechamber*.

4.6 EXPLANTATION AND DISPOSAL

4.6.1 EXPLANTATION

The explantation of the GAV 2.0 XABO should be performed according to the state of the art and in compliance with medical practice.

4.6.2 DISPOSAL

GAV 2.0 XABO and shunt components

Products and product parts not used in the implantation or surgically removed must be disposed of correctly as potentially infectious material in accordance with medical practice as well as respective regional laws and regulations.

Explanted medical devices must not be reused.

4.7 TECHNICAL INFORMATION

4.7.1 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	GAV 2.0 XABO
Medical Purpose	Shunting of cerebrospinal fluid (CSF)
Sterilisability	Cannot be resterilised
Storage	Store in a clean and dry place at ≤ 30 °C
For single use only	
Sketch with outer dimensions:	
<p>GAV 2.0 (VP shuntina)</p> <p>Technical drawing of the GAV 2.0 VP shuntina showing outer dimensions:</p> <ul style="list-style-type: none"> Left port diameter: $\varnothing 1.9$ mm Right port diameter: $\varnothing 4.2$ mm Distance between port centers: 13.4 mm Total length: 22.2 mm 	

4.8 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxx indicates the identifier of the responsible notified body
	Medical device
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Serial number
	Sterilised using irradiation
	Do not resterilise
	Do not reuse
	Do not use if package is damaged and consult instructions for use
	Keep dry
	Upper limit of temperature
	Keep away from sunlight
	Consult instructions for use / electronic instructions for use
	Caution
	Contains a medical substance

Symbol	Explanation
	Non-pyrogenic
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians.
	MR Conditional
	Patient identification
	Date
	Healthcare centre or doctor
	Patient information website
	Model number / European Medical Device Nomenclature Code

5 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany