proGAV[®] 2.0 XABO[®]

IN TOUCH WITH YOU

- DE** Gebrauchsanweisung | **EN** Instructions for use | **FR** Mode d'emploi |
ES Instrucciones de manejo | **PT** Instruções de utilização | **IT** Istruzioni per l'uso |
NL Gebruiksaanwijzing |

 www.miethke.com

US This Instructions for Use is NOT intended for United States users. Please discard.
The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0
info@miethke.com

CONTENTS

| | | |
|------|----------------------------------------------------------------|----|
| 1 | PREFACE AND IMPORTANT INFORMATION | 21 |
| 2 | INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE | 21 |
| 2.1 | EXPLANATION OF THE WARNINGS | 21 |
| 2.2 | DISPLAY CONVENTIONS | 21 |
| 2.3 | OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION | 21 |
| 2.4 | FEEDBACK ON THE INSTRUCTIONS FOR USE | 21 |
| 2.5 | COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION | 21 |
| 3 | DESCRIPTION OF THE <i>proGAV 2.0 XABO</i> | 22 |
| 3.1 | MEDICAL PURPOSE | 22 |
| 3.2 | CLINICAL BENEFITS | 22 |
| 3.3 | INDICATIONS | 22 |
| 3.4 | CONTRAINDICATIONS | 22 |
| 3.5 | INTENDED PATIENT GROUPS | 22 |
| 3.6 | INTENDED USERS | 22 |
| 3.7 | INTENDED USE ENVIRONMENT | 22 |
| 3.8 | TECHNICAL DESCRIPTION | 23 |
| 3.9 | FUNCTIONING OF THE VALVE SYSTEM | 24 |
| 3.10 | SELECTION OF THE APPROPRIATE PRESSURE SETTING | 25 |
| 3.11 | PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES | 25 |
| 3.12 | PRESSURE-FLOW CHARACTERISTICS | 26 |
| 3.13 | APPLICATION OF THE PROGAV 2.0 TOOLS | 27 |
| 3.14 | SYSTEM COMPONENT | 30 |
| 3.15 | FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES | 30 |
| 4 | PROPERTIES OF THE <i>proGAV 2.0 XABO</i> | 31 |
| 4.1 | PRODUCT DESCRIPTION | 31 |
| 4.2 | IMPORTANT SAFETY INFORMATION | 31 |
| 4.3 | PATIENT EDUCATION | 32 |
| 4.4 | TRANSPORT AND STORAGE | 32 |
| 4.5 | USE OF THE PRODUCT | 32 |
| 4.6 | EXPLANTATION AND DISPOSAL | 35 |
| 4.7 | TECHNICAL INFORMATION | 35 |
| 4.8 | SYMBOLS USED FOR LABELLING | 36 |
| 5 | MEDICAL DEVICE CONSULTANTS | 36 |

1 PREFACE AND IMPORTANT INFORMATION

Preface

Thank you for purchasing the medical device *proGAV 2.0 XABO*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product. Your team at Christoph Miethke GmbH & Co. KG

Relevance of the instructions for use

WARNING

Improper handling and use of this product can cause hazards and damages. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

proGAV 2.0 XABO includes the following components:

- ▶ *proGAV 2.0 XABO*

Additional options:

- ▶ *Reservoirs* (including paediatric version)
- ▶ *Prechamber* (including paediatric version)
- ▶ *Burrhole Deflector* (including paediatric version)
- ▶ *Ventricular Catheter*
- ▶ *Titanium Connectors*

2 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE

2.1 EXPLANATION OF THE WARNINGS

DANGER

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

WARNING

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

CAUTION

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

NOTE

Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for danger, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

2.2 DISPLAY CONVENTIONS

| Display | Description |
|----------------|--------------------------------|
| <i>Italics</i> | Indicates <i>product names</i> |

2.3 OTHER ACCOMPANYING DOCUMENTS AND ADDITIONAL INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

Each delivery is accompanied by a patient card with information about the product. The patient card is to provide the treating physician with all the product information in a compact form for the patient record.

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

2.4 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

2.5 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

3 DESCRIPTION OF THE *proGAV 2.0 XABO*

3.1 MEDICAL PURPOSE

proGAV 2.0 XABO is intended for shunting of cerebrospinal fluid (CSF).

3.2 CLINICAL BENEFITS

Clinical benefits of the *proGAV 2.0 XABO* :

- ▶ Implantation of a long-term implant for the extraction of cerebrospinal fluid from the CSF circulation and its shunting into the peritoneum
- ▶ Treatment of hydrocephalus, e.g. by alleviation of the clinical symptoms
- ▶ Reducing the risk of infection with gram-positive bacteria thanks to antibiotic-impregnated catheters

3.3 INDICATIONS

The following indications apply to *proGAV 2.0 XABO*:

- ▶ Treatment of hydrocephalus

3.4 CONTRAINDICATIONS

The following contraindications apply to *proGAV 2.0 XABO*:

- ▶ Infections in the implantation area
- ▶ Pathological concentrations levels (e.g. of blood components and/or protein) within the cerebrospinal fluid
- ▶ Intolerance to materials of the shunt system
- ▶ Hypersensitivity to rifampicin and/or clindamycin hydrochloride

3.5 INTENDED PATIENT GROUPS

- ▶ Patients who are treated with a CSF shunt system on account of their clinical picture

3.6 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical professionals, e. g. neurosurgeons
- ▶ Knowledge of mode of operation and intended use of the product
- ▶ Successful participation in product training

3.7 INTENDED USE ENVIRONMENT

Professional healthcare facilities

- ▶ Implantation under sterile operating theatre conditions

3.8 TECHNICAL DESCRIPTION


Fig. 1: proGAV 2.0 cross section

- 1 Adjustable differential pressure unit**
- 2 Torsion spring
- 3 Sapphire ball
- 4 Rotor

- 5 Gravitational unit (SHUNTASSISTANT 2.0)**
- 6 Tantalum ball
- 7 Sapphire ball

The *proGAV 2.0* is a valve system made of titanium. It consists of an adjustable differential pressure unit (1) and a gravitational unit (5) (Abb. 1).

The adjustable differential pressure unit (1) in the proximal part of the valve system consists of a stable titanium housing with an integrated ball-cone valve (3) in the front section. A torsion spring (2) determines the opening pressure of this unit. The pretension of the spring, and thus the valve opening pressure, can be adjusted post-surgically through the skin using a pivoted rotor (4).

Major components of the gravitational unit (5) are a tantalum ball (6), which defines the opening pressure for this valve depending on body

position, and a sapphire ball (7), which ensures precise closure.

The *proGAV 2.0 XABO* consists of an *proGAV 2.0* valve system in combination with *XABO Catheters*.

XABO Catheters are made from silicone, they are impregnated with antibiotics in a special treatment process and contain 0.054% rifampicin and 0.15% clindamycin hydrochloride. Laboratory test have shown *XABO Catheters* to reduce colonisation with gram-positive bacteria of the silicone catheter surface. The laboratory tests were performed with *Staphylococcus aureus*, *Staphylococcus epidermidis* and *Bacillus subtilis*. Systemic therapeutic effects are highly unlikely as the amounts of rifampicin and clindamycin hydrochloride

contained in the catheter represent only a fraction of the therapeutic dose of these antibiotics.

3.9 FUNCTIONING OF THE VALVE SYSTEM

The proGAV 2.0 XABO is a posture-dependent hydrocephalus valve system. The opening pressure of the proGAV 2.0 XABO is composed of the opening pressures of the adjustable differential pressure unit and the gravitational unit.

Horizontal body position

In the horizontal body position, the gravitational unit is always open and does not present any resistance (Fig. 2).


Fig. 2: The gravitational unit in the horizontal body position

The opening pressure of the proGAV 2.0 XABO in the horizontal body position is thus determined by the adjustable differential pressure unit. The basic operating principle of the differential pressure unit is shown in Fig. 3 a) and b).


Fig. 3: The adjustable differential pressure unit in the horizontal body position
a) closed b) open

In Fig. 3 a), the differential pressure unit is closed; hence, shunting is not possible.

If the patient's intraventricular pressure (IVP) exceeds the spring force of the torsion spring, which otherwise keeps the differential pressure unit closed, the sealing ball moves out of the cone, leaving a gap for shunting of CSF (Fig. 3 b).

Vertical body position

As soon as the patient sits or stands up, the gravitational unit immediately closes (Fig. 4 a). The opening pressure of the proGAV 2.0 XABO is thus increased significantly, because now the weight of the tantalum ball (opening pressure of the gravitational unit) must be overcome in addition to the opening pressure of the differential pressure unit. Shunting is only possible again when the sum of IVP and hydrostatic suction is greater than the opening pressure of both units (Fig. 4 b).

For individual adaptation of the opening pressure to the patient, an opening pressure between 0 and 20 cmH₂O can be selected for the adjustable differential pressure unit.


Fig. 4: Gravitational unit in the vertical body position
a) closed b) open

NOTE

During physical activity associated with shock (e.g. jogging), the opening pressure of the proGAV 2.0 XABO may decrease temporarily according to laboratory results. Fundamentally, though, functionality remains the same. At the end of physical activity, the opening pressure returns to its original level and remains stable.

3.10 SELECTION OF THE APPROPRIATE PRESSURE SETTING

Recommended pressure settings for the proGAV 2.0 XABO are available at:
<https://www.miethke.com/downloads/>

This is a non-binding recommendation for the treating physician. The physician decides on each case independently, individually and autonomously in accordance with his/her diagnosis.

The adjustable differential pressure unit of the proGAV 2.0 XABO is set to an opening pressure of 5 cmH₂O upon delivery.

This pre-set opening pressure can be changed to a different pressure prior to implantation.

Horizontal body position

The opening pressure in the horizontal body position is determined by the differential pressure unit.

Depending on the patient's clinical picture, indication and age, the opening pressure for this position can be selected between pressure levels from 0 to 20 cmH₂O.

Vertical body position

The opening pressure of the proGAV 2.0 XABO for the vertical body position is calculated from the sum of the opening pressure of the differential pressure unit and the gravitational unit. Patient height, activity level and potentially increased abdominal pressure (obesity) should be taken into account in selecting the opening pressure level for the gravitational unit.

3.11 PRESSURE LEVEL IDENTIFICATION IN X-RAY IMAGES

proGAV 2.0 differential pressure unit

The selected pressure level of differential pressure unit should always be monitored using the proGAV 2.0 Compass or the M.blue plus Com-

pass, but it can also be checked with the aid of an X-ray image (Fig. 5).


Fig. 5: X-ray image (proGAV 2.0 adjustable differential pressure unit, set to 14 cmH₂O)

The rotor setting is decisive in this case. The four magnets in the rotor appear in the X-ray image as white dots and are located opposite each other in pairs. On one side of the rotor, two additional burr holes (right and left next to the magnet pairs) serve as orientation. They appear as black dots in the X-ray image. This side can be described as the rear side of the rotor. The two front magnets are on the opposite side.


Fig. 6: Schematic representation of rotor in an X-ray image

- 1. Admittance markings, 2. Valve marking
- 3. Triangle apex, 4. Non-adjustable range
- 5. proximal, 6. distal

The space between these two magnets can be considered as the apex of the triangle. The pressure level can be read off using the orientation of this intermediate space (Fig. 6). The apex of the triangle can take up any position except the space marked as a the non-

adjustable range in Fig. 6. The opening pressure of the *proGAV 2.0* can thus be infinitely variably adjusted from 0 up to 20 cmH₂O.

To prevent reading the pressure level in reverse, the valve is fitted with a valve marking on one side, which appears black in the X-ray image – in a top view onto the implanted valve as in Fig. 5, the recess is visible on the right-hand side.

SA 2.0 gravitational unit

The pressure levels of the gravitational unit can be identified post-surgically by the coding:

| Pressure level | Coding |
|-----------------------|--------|
| 10 cmH ₂ O | |
| 15 cmH ₂ O | |
| 20 cmH ₂ O | |
| 25 cmH ₂ O | |
| 30 cmH ₂ O | |
| 35 cmH ₂ O | |

3.12 PRESSURE-FLOW CHARACTERISTICS

The following diagrams show the pressure-flow characteristics for the *proGAV 2.0 XABO*. The opening pressure refers to a reference flow of 5 ml/h. For flow rates of 20 ml/h, stated pressures are approx. 1–2 cmH₂O higher. The supplied *XABO* catheters do not fundamentally change the pressure-flow characteristics.

Horizontal valve position

The pressure flow characteristics for the adjustable differential pressure unit of the *proGAV 2.0 XABO* are shown below as an example for the pressure settings of 0, 10 and 20 cmH₂O in the horizontal valve position.




Fig. 7: Pressure-flow characteristics for selected pressure levels of the adjustable differential pressure unit; pressure (cmH₂O), flow rate (ml/h)

Vertical valve position

In the vertical body position, the opening pressure of the *proGAV 2.0* is composed of the setting of the adjustable differential pressure unit and the gravitational unit. The pressure flow characteristics for various pressure levels in the vertical body position are shown below:





Fig. 8: Pressure-flow characteristics of the available pressure levels of the proGAV 2.0 in the vertical body position; pressure (cmH₂O), flow rate (ml/h): a) 10 cmH₂O, b) 15 cmH₂O, c) 20 cmH₂O, d) 25 cmH₂O, e) 30 cmH₂O, f) 35 cmH₂O

3.13 APPLICATION OF THE proGAV 2.0 Tools

Only proGAV 2.0 Tools or M.blue plus Instruments approved for that purpose must be used to determine, change and monitor the opening pressure of the differential pressure unit of the proGAV 2.0.

The proGAV 2.0 Tools may only be used by trained specialists.

The selected pressure level of proGAV 2.0 can be determined, adjusted and monitored using the proGAV 2.0 Tools.

The proGAV 2.0 Compass (Fig. 9) is used for localising and reading out the adjustable unit of the proGAV 2.0.

Fig. 9: proGAV 2.0 Compass

The proGAV 2.0 Adjustment Instrument (Fig. 10) is used to set the opening pressure of the adjustable unit of the proGAV 2.0 from 0 to 20 cmH₂O.

Fig. 10: proGAV 2.0 Adjustment Instrument

The opening pressure of the adjustable differential pressure unit of the *proGAV 2.0* can be changed before or after implantation. It is pre-set by the manufacturer to 5 cmH₂O.

Changing the opening pressure of the *proGAV 2.0* requires the following steps:

1. Localisation

WARNING

The *proGAV 2.0 Compass* should be placed as centrally as possible over the valve, as the determined opening pressure may otherwise be incorrect.

NOTE

- ▶ The *proGAV 2.0 Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *proGAV 2.0 Adjustment Instrument* should not be placed in the immediate vicinity of the *proGAV 2.0 Compass* when determining the opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin may make adjustment difficult for a few days after surgery. If the valve setting cannot be checked conclusively using the *proGAV 2.0 Compass*, we recommend checking it with an imaging method.

Opening the *proGAV 2.0 Compass* reveals a template that can be used to localise the valve in the patient's head using the index finger (Fig. 11).


Fig. 11: Localising the valve

After this, the template of the *proGAV 2.0 Compass* is adjusted to the flow direction of the CSF and placed over the valve. The "proximal" and "distal" direction markings show the flow direction.

2. Test procedure

The pressure level is displayed automatically when the *proGAV 2.0 Compass* is folded down.


Fig. 12: Determination of the pressure level with the *proGAV 2.0 Compass*

3. Adjustment process

WARNING

When adjusting the differential pressure unit of the *proGAV 2.0*, care must be taken to change the opening pressure by a maximum of 8 cmH₂O per adjustment process; otherwise, errors can result. Example: The opening pressure is to be changed from 3 to 18 cmH₂O. The correct method is an adjustment in two stages: an initial adjustment from 3 to 11 cmH₂O and subsequently from 11 to 18 cmH₂O.

The *proGAV 2.0 Adjustment Instrument* is centred over the valve. The recess in the middle of the tool makes it easy to identify the valve by touch using the index finger to place the tool correctly (Fig. 13). When doing so, the required pressure level must be displayed on the scale in the direction of the valve inlet or the *Ventricular Catheter*.

Fig. 13: Adjustment using the *proGAV 2.0 Adjustment Instrument*

Lightly pressing the index finger on the adjustable unit releases the rotor brake and changes the pressure level of the *proGAV 2.0* (Fig. 14).


Fig. 14: Adjustment using the *proGAV 2.0 Adjustment Instrument*

The differential pressure unit of the *proGAV 2.0* is fitted with a feedback mechanism. Due to the valve housing design, targeted pressure on the valve produces an audible acoustic signal (a clicking sound) and/or palpable resistance as soon as the rotor brake has been released. The valve thus shows both acoustically and haptically when the pressure is sufficient for uncoupling. Once this pressure has been released, the rotor is once again adjustment-proof. Although the click caused by releasing the rotor brake is easily audible before implantation, it may be considerably reduced

after implantation and the filling of the valve depending on its position and the condition of the implant surroundings. Normally, however, it should be audible to the patient or by using a stethoscope.

Checking after adjustment

After adjusting the valve opening pressure, it is recommended to check the set pressure level. To do this, proceed as in Points 1 and 2. Should the measured value not correspond with the required pressure level, the adjustment process should be repeated. To do this, start again at Point 3.

proGAV Checkmate

The *proGAV Checkmate* (Fig. 15) is supplied sterile and can be re-sterilised. The *proGAV Checkmate* can be used to change and check the pressure level before and during valve implantation directly on the *proGAV 2.0*. To determine the pressure level, the *proGAV Checkmate* is placed centrally on the *proGAV 2.0*. The *proGAV Checkmate* automatically aligns itself over the valve. The pressure level can be read from the direction of the proximal catheter (leading towards the valve). For pressure level adjustment, the *proGAV Checkmate* is placed centrally on the *proGAV 2.0*. When doing so, the required pressure level must point towards the proximal catheter (leading towards the valve). Lightly pressing the *proGAV Checkmate* onto the valve releases the rotor brake in the *proGAV 2.0* and sets the pressure level.

Fig. 15: *proGAV Checkmate*, Pressure levels 0-20 cmH₂O

Checking and adjusting without unpacking

To protect the impregnation, the *proGAV 2.0 XABO* is wrapped in opaque packaging. However, it is still possible to check and adjust the setting of the *proGAV 2.0 XABO*.

The position of the *proGAV 2.0 XABO* is indicated by an arrow mark. The arrow indicates the flow direction.

For localisation, the *proGAV 2.0 Compass* is placed centrally on the marking of the valve that needs to be checked. Because of the packaging, the *proGAV 2.0 Compass* is slightly tilted on the packaging.

Checking and/or adjusting without unpacking retains the sterility of the product and the integrity of the packaging.


Fig. 16: Packaging of the *proGAV 2.0 XABO*: marking of the adjustable valve with an arrow


Fig. 17: Position of the *proGAV 2.0 Compass* on the packaging of the *proGAV 2.0 XABO*

3.14 SYSTEM COMPONENT

Combination with shunt components

The *proGAV 2.0 XABO* can be safely combined with our range of implantable shunt components. We recommend using the Christoph Miethke GmbH & Co. KG products in combination with the *proGAV 2.0 XABO*.

Reservoirs

The use of shunt systems with a reservoir provides options for the withdrawal of cerebrospinal fluid, administration of drugs and pressure control.

Thanks to an integrated check valve in the *CONTROL RESERVOIR* and the *SPRUNG RESERVOIR*, cerebrospinal fluid can be pumped towards the valve, thus making it possible to check the distal part of the shunt system as well as the *Ventricular Catheter*.

During the pump action, access to the *Ventricular Catheter* is closed. The use of this Reservoir

does not increase the opening pressure of the shunt system. Puncturing the Reservoir should be performed as perpendicular as possible to the reservoir surface with a maximum cannula diameter of 0.9 mm. A firm titanium floor prevents the bottom surface from being pierced. 30 punctures are possible without any restrictions.

WARNING

Frequent pumping of the Reservoir can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.

Burrhole Deflector

Because of the tight fit on the *Ventricular Catheter*, the *Burrhole Deflector* makes it possible to choose the length of catheter penetrating into the skull prior to implantation. The *Ventricular Catheter* is deflected at a right angle in the burr hole (see chapter 4.5.5).

Tube systems

The *proGAV 2.0 XABO* is supplied as a shunt system with integrated catheters (internal diameter 1.2 mm, external diameter 2.5 mm).

For a new connection of catheter and connector, the catheters have to be carefully fixed with a ligature to the valve's *Titanium Connectors*.

3.15 FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

- ▶ The medical device has been designed to work accurately and reliably over long periods of time. We guarantee that our medical devices are free from defects and fully functional at the time of dispatch. Excluded from this guarantee are cases in which the medical device has been exchanged or rather explanted for technical or medical reasons outside of our sphere of responsibility.
- ▶ *proGAV 2.0 XABO* together with the entire shunt system can safely withstand pressure of up to 100 cm H₂O occurring during and after surgery.

- Nuclear magnetic resonance examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. *proGAV 2.0 XABO* is MR Conditional. The supplied catheters are MR Safe. *Reservoirs, Bur-hole Deflectors and Titanium Connectors* are MR Conditional.

The conditions for MRI security of the products are found on our website:
<https://miethke.com/downloads/>

4 PROPERTIES OF THE *proGAV 2.0 XABO*

4.1 PRODUCT DESCRIPTION

4.1.1 *proGAV 2.0 XABO* VARIANTS

The *proGAV 2.0 XABO* is available in a range of different variants. They differ in the pre-set pressure of the gravitational unit.

| Adjustable differential pressure unit | Gravitational unit |
|---------------------------------------|-----------------------|
| 0–20 cmH ₂ O | without |
| 0–20 cmH ₂ O | 10 cmH ₂ O |
| 0–20 cmH ₂ O | 15 cmH ₂ O |
| 0–20 cmH ₂ O | 20 cmH ₂ O |
| 0–20 cmH ₂ O | 25 cmH ₂ O |
| 0–20 cmH ₂ O | 30 cmH ₂ O |
| 0–20 cmH ₂ O | 35 cmH ₂ O |

These variants are also available as shunt systems and can come with the following components: *XABO Ventricular Catheter, Prechambers, Reservoirs*.

4.1.2 SCOPE OF DELIVERY

| Box content | Number |
|------------------------------------------------------------|--------|
| Sterile packaging with <i>proGAV 2.0 XABO</i> shunt system | 1 |
| Instructions for use for <i>proGAV 2.0 XABO</i> | 1 |
| Patient card | 1 |
| Recommended pressure level | 1 |

4.1.3 STERILITY

WARNING

The product must not be used if the sterile packaging or the product is damaged or after the expiry date.

proGAV 2.0 XABO is sterilised under strictly controlled conditions using irradiation. The respective expiry date is printed on the packaging. The inner sterile packaging consists of a gas- and liquid-tight barrier to protect the antibiotics.

4.1.4 REPEATED USE AND RESTERILISATION

WARNING

The product must not be resterilised or reprocessed in any other way as the safe functioning and sterility of the product cannot be guaranteed.

Products that have already been implanted in a patient must not be reused either on the same or a different patient in order to minimise the risk of infection.

4.1.5 SINGLE-USE PRODUCT

This product is intended for single use. Reprocessing may lead to significant changes to the properties of the *proGAV 2.0 XABO*. No guarantee can be assumed for the functional safety of resterilised products.

4.1.6 PRODUCT CONFORMITY

The product meets current regulatory requirements.

The requirements stipulate the comprehensive documentation of the whereabouts of medical devices used in humans. The individual identification number of the medical device should therefore be recorded in the patient's medical records to ensure complete traceability.

4.2 IMPORTANT SAFETY INFORMATION

4.2.1 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

**WARNING**

- ▶ **The product must not be used if the sterile packaging or the product is damaged or after the expiry date.**
- ▶ **Due to the risk of injury resulting from incorrect use of the product, the instructions for use must be carefully read and understood before the product is used for the first time.**
- ▶ **Prior to use, it is essential to check the product for completeness and integrity.**

4.2.2 COMPLICATIONS, SIDE EFFECTS, PRECAUTIONS AND RESIDUAL RISKS

The following complications can occur in conjunction with the *proGAV 2.0 XABO*:

- ▶ Headaches, dizzy spells, mental confusion, vomiting in cases of possible leakage from the shunt system and shunt dysfunction
- ▶ Redness/irritation of the skin and tightness around the implantation site as an indication of a possible infection at the implant
- ▶ Occlusions caused by protein and/or blood in the cerebrospinal fluid
- ▶ Allergic reaction / intolerance to product materials
- ▶ Overdrainage/underdrainage
- ▶ Noise development

Violent external shocks (accident, fall, etc.) may put the integrity of the shunt system at risk.

As a precaution, a physician must be consulted immediately if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

The following residual risks exist when using the *proGAV 2.0 XABO*:

- ▶ Persistent headache
- ▶ Severe infection (e.g. sepsis, meningitis) / allergic shock
- ▶ Acute and chronic hygroma/subdural haematoma
- ▶ Cerebrospinal fluid accumulations
- ▶ Tissue damage/puncture
- ▶ Skin irritation

- ▶ Local shunt irritation
- ▶ Allergic reactions to catheter components, particularly to the active antibiotic substances rifampicin and clindamycin hydrochloride

4.2.3 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible state authority.

4.3 PATIENT EDUCATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product (Ch. 4.2, 4.5).

4.4 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place. *proGAV 2.0 XABO* is to be protected from direct sunlight. The products should not be removed from the packaging until they are required.

4.4.1 TRANSPORT

Transport conditions

| | |
|----------------------------|---------|
| Ambient temperature | ≤ 40 °C |
|----------------------------|---------|

4.4.2 STORAGE

Storage conditions

| | |
|--------------------------------------|---------|
| Temperature range for storage | ≤ 30 °C |
|--------------------------------------|---------|

4.5 USE OF THE PRODUCT

4.5.1 INTRODUCTION

The *proGAV 2.0 XABO* is a posture-dependent valve system with an adjustable differential pressure unit and a pre-set gravitational unit (SA 2.0) combined with antibiotic impregnated *XABO Catheters*.

proGAV 2.0 XABO is used for shunting of cerebrospinal fluid (CSF) in the treatment of hydro-

cephalus. Valves and Reservoirs are placed in suitable positions along the course of the shunt.

4.5.2 SAFETY NOTICES AND WARNINGS

WARNING

- ▶ *proGAV 2.0 XABO* must not be immersed in antibiotic solutions. The contact of the catheter with sterile water or sterile saline solution should be limited to a minimal period of time and should be done immediately prior to implantation. The solution may take on a faint orange colour.
- ▶ Frequent pumping of the *Reservoir* can result in excessive shunting and thus lead to pressure conditions outside the normal physiological range. The patient should be properly informed about this risk.
- ▶ The gravitational unit of the *proGAV 2.0 XABO* is posture-dependent. For that reason, care must be taken to implant the gravitational unit (*SA 2.0*) is implanted parallel to the body axis.
- ▶ The adjustable valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e. g. underneath heavily scarred tissue). If the implantation site is unfavourably chosen or the skin above the valve is too thick, it may no longer be possible to adjust the adjustable unit. The valve then works with unchangeable pressure levels.
- ▶ If a magnetic field is being applied and pressure is applied to the valve at the same time, thus triggering the brake mechanism, and adjustment of the valve cannot be ruled out.

CAUTION

- ▶ Silicone is extremely electrostatic. Care must be taken to avoid the catheters coming into contact with dry cloths, talcum powder or rough surfaces. Clinging particles could lead to tissue reactions.
- ▶ When using sharp instruments, care should be taken to avoid cuts and scratches in the silicone elastomer.

CAUTION

- ▶ It must be ensured that the ligature is not tightened excessively. Damage may result in a loss of integrity of the shunt and thus necessitate a revision.
- ▶ The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.
- ▶ For people using cardiac pacemakers: It is possible that the function of the pacemaker is influenced by the implantation of a *proGAV 2.0 XABO*.

NOTE

- ▶ In MRI imaging *proGAV 2.0 XABO* creates artefacts that are larger than the valve itself.

4.5.3 REQUIRED MATERIALS

The *proGAV 2.0 XABO* is designed so that it can be safely used with the shunt components described in chapter 3.14. Catheters with an inner diameter of 1.2 mm and an outer diameter of approx. 2.5 mm should be used for connection. In any case, catheters have to be carefully fixed with a ligature to the connectors of the shunt components. Any kinks in the catheter must be avoided.

4.5.4 PREPARING FOR IMPLANTATION

Checking the sterile packaging

Immediately before using the product, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Preoperative valve test

The *proGAV 2.0 XABO* should be vented before implantation and checked for permeability. The most careful way of filling the valve is by aspiration through a sterile single-use syringe attached to the distal end of the catheter. The distal end of the valve is connected and immersed in a sterile physiological salt solution. The valve is patent if saline solution can be extracted (Fig. 18).


Fig. 18: Patency test

**WARNING**

- ▶ **No antibiotic solution must be used for valve testing and/or venting, as this may result in a reaction with the active substances in the impregnation.**
- ▶ **Contaminants in the solution used for testing can impair the product's performance.**
- ▶ **Pressurisation with a single-use syringe should be avoided both at the proximal and the distal end. (Fig. 19).**


Fig. 19: Avoidance of pressurisation

4.5.5 PERFORMING THE IMPLANTATION**Positioning of the XABO Ventricular Catheter**

Several surgical techniques are available for the positioning of the XABO Ventricular Catheter. The required skin incision should be made in form of a lobule pedicled towards the shunting catheter or by a straight skin incision. If a *Burrhole Reservoir* - or a *SPRUNG RESERVOIR* - is used, the skin incision should not be located right above the reservoir. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole.

proGAV 2.0 XABO is available in a range of different configurations: If using a *Burrhole Reservoir* - or a *SPRUNG RESERVOIR* -, the XABO Ventricular Catheter is implanted first. Once the mandrin has been removed, the patency of the XABO Ventricular Catheter can be tested by checking if cerebrospinal fluid is dripping out. The catheter is shortened and the *Burrhole Reservoir* - or the *SPRUNG RESERVOIR* - con-

nected and the connection secured with a ligature.

When using a shunt system with a *CONTROL RESERVOIR*, a *Burrhole Deflector* is included. The *Burrhole Deflector* is used to adjust the length of the catheter to be implanted and to position it inside the ventricle. The *Ventricular Catheter* is deflected by 90° and the *CONTROL RESERVOIR* put into place. The position of the XABO Ventricular Catheter should be inspected after the procedure by imaging (such as CT or MRI).

Positioning of the valve system

A location behind the ear is suitable as an implantation position, whereby the implantation height has no influence on the function of the valve system.

The adjustable valve should be touching the bone or the periosteum since pressure must be exerted on the valve during any later adjustment.

A large arch-shaped or a small straight skin cut with a pocket for the valve system should be made. The catheter is pushed forward from the burr hole to the selected valve implantation location, shortened if necessary, and secured to the proGAV 2.0 XABO by ligation. The valve system should not be located directly under the skin incision. The valve unit has an arrow in the flow direction (arrow towards distal or downwards). The surface of the valve with the arrow markings points to the outside.

**WARNING**

The gravitational unit of the proGAV 2.0 XABO is posture-dependent. For that reason, care must be taken to implant the gravitational unit (SA 2.0) is implanted parallel to the body axis.

Positioning of the XABO Peritoneal Catheter

The access site for the XABO Peritoneal Catheter is left to the surgeon's discretion. For example, it can be used in a paraumbilical application or applied at the level of the epigastrium. Likewise, various surgical techniques are available for placing the XABO Peritoneal Catheter. The recommendation is to pull the XABO Peritoneal Catheter from the valve to the intended position using a subcutaneous *Tunneller*, if necessary with the aid of an auxiliary

incision. The *XABO Peritoneal Catheter*, usually securely attached to the valve, has an open distal end and no wall slits. Following the exposure of the peritoneum or with the aid of a trocar, the *XABO Peritoneal Catheter* (shortened if necessary) is pushed forward into the open space of the abdominal cavity.

4.5.6 POSTOPERATIVE VALVE TEST

Postoperative valve test

The *proGAV 2.0 XABO* has been constructed as a reliably functioning unit without pump or test function. The valve test can be performed by flushing, pressure measurement or pumping via a *Reservoir* or a *Prechamber*.

4.6 EXPLANTATION AND DISPOSAL

4.6.1 EXPLANTATION

The explantation of the *proGAV 2.0 XABO* should be performed according to the state of the art and in compliance with medical practice.

4.6.2 DISPOSAL

proGAV 2.0 XABO and shunt components

Products and product parts not used in the implantation or surgically removed must be disposed of correctly as potentially infectious material in accordance with medical practice as well as respective regional laws and regulations.


Explanted medical devices must not be reused.

4.7 TECHNICAL INFORMATION

4.7.1 TECHNICAL DATA

| | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|
| Manufacturer | Christoph Miethke GmbH & Co. KG |
| Product designation | <i>proGAV 2.0 XABO</i> |
| Medical Purpose | Shunting of cerebrospinal fluid (CSF) |
| Sterilisability | Cannot be resterilised |
| Storage | Store in a clean and dry place at ≤ 30 °C |
| For single use only | |
| Sketch with outer dimensions: | |
| <p>The sketch shows a horizontal device with four dimension lines. From left to right: a connector of 4.5 mm, a main body of 17 mm, a valve of 12 mm, and a valve connector of 4.2 mm. The main body is labeled 'proGAV 2.0' and the valve is labeled 'SA 2.0'.</p> | |

4.8 SYMBOLS USED FOR LABELLING

| Symbol | Explanation |
|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  | EU conformity marking, xxxx indicates the identifier of the responsible notified body |
|  | Medical device |
|  | Manufacturer |
|  | Date of manufacture |
|  | Use-by date |
|  | Batch code |
|  | Catalogue number |
|  | Serial number |
|  | Sterilised using irradiation |
|  | Do not re-sterilise |
|  | Do not reuse |
|  | Do not use if package is damaged and consult instructions for use |
|  | Keep dry |
|  | Upper limit of temperature |
|  | Keep away from sunlight |
|  | Consult instructions for use / electronic instructions for use |
|  | Caution |
|  | Contains a medical substance |

| Symbol | Explanation |
|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------|
|  | Non-pyrogenic |
|  | Free of natural rubber latex, latex-free |
|  | Indicates that in the USA, the product may only be issued to physicians. |
|  | MR Conditional |
|  | Patient identification |
|  | Date |
|  | Healthcare centre or doctor |
|  | Patient information website |
|  | Model number / European Medical Device Nomenclature Code |

5 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany