Titanium Connector

Ⓛ Gebrauchsanweisung | Ⓜ Instructions for use | Ⓧ Mode d'emploi
Ⓝ Instrucciones de manejo | ⓘ Istruzioni per l'uso

INDIKATION

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INDICATION

One possible treatment in cases of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain. This shunt system is comprised of a series of catheters and a pressure relief valve. Connectors are used to link these components. *Titanium connectors* are biocompatible, very stable and generate hardly any artefacts during the performance of imaging procedures.

TUBE SYSTEMS

The titanium connectors allow the application of catheters with an inner diameter of 1.2 mm to 1.5 mm. The outer diameter of the catheter should be approximately twice that of the inner diameter. For spinal drainage special titanium connectors, step down are available. In every case, the catheter must be carefully fastened to the *titanium connector* by means of a ligature. Kinks in the catheters have to be avoided.

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

Owing to the fact that the titanium connectors are comprised of titanium, a nonmagnetic material, the product will not interfere with nuclear magnetic resonance and computer tomography procedures. The titanium connectors are MRI Conditional (ASTM F2503-08).

SAFETY MEASURES

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflammation and tension

around the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation or vomiting often occur in cases of shunt dysfunction. In the event such symptoms occur or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

STERILIZATION

All products are carefully and thoroughly sterilized. Owing to the fact that the products are packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances. The functional safety and reliability of resterilized products cannot be guaranteed.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz
Michaela Funk-Neubarth

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
14469 Potsdam, Germany
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:

AESCLAP AG
Am Aesculap Platz
78532 Tuttlingen, Germany
Phone: +49 (0) 7461 95-0
Fax: +49 (0) 7461 95-26 00
e-mail: information@aesculap.de

Service address in the US

AESFULAP Inc.
Attn. AESFULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042
AESFULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada

AESFULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesfulapusa.com



NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

GENERAL INFORMATION


Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>Titanium Connector</i>
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation with its external dimensions:	
<p>normal:</p>  <p>1,9 mm $d_i = 0,9 \text{ mm}$</p>	
<p>step down:</p>  <p>1,4 mm $d_i = 0,6 \text{ mm}$</p>	

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany