Titanium Connectors

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 www.miethke.com

 This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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


MAGRAM Konnektor, dual right angle

TABLE OF CONTENTS

INDICATION	6
TUBE SYSTEMS	6
RE-IMPLANTATION	6
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	6
SAFETY MEASURES	6
STERILIZATION	6
MEDICAL DEVICES CONSULTANTS	6
REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)	6
GENERAL INFORMATION	7
VARIANTS	8

INDICATION

One possible treatment in cases of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain. This shunt system is comprised of a series of catheters and a pressure relief valve. Connectors are used to link these components. Titanium connectors are biocompatible, very stable and generate hardly any artefacts during the performance of imaging procedures.

TUBE SYSTEMS

The titanium connectors allow the application of catheters with an inner diameter of 1.2 mm to 1.5 mm. The outer diameter of the catheter should be approximately twice that of the inner diameter. For spinal drainage special titanium connectors, step down are available. In every case, the catheter must be carefully fastened to the titanium connector by means of a ligature. Kinks in the catheters have to be avoided.

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

Owing to the fact that the titanium connectors are comprised of titanium, a non-magnetic material, the product will not interfere with nuclear magnetic resonance and computer tomography procedures. The titanium connectors are MR Conditional (ASTM F2503-2020).

SAFETY MEASURES

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflammation and tension around the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation or vomiting often occur in cases of shunt dysfunction. In the event such symptoms occur or if there is any leakage in the

shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

STERILIZATION

All products are carefully and thoroughly sterilized. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

The functional safety and reliability of resterilized products cannot be guaranteed.

MEDICAL DEVICES CONSULTANTS

In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical products consultants as contacts for all product-related questions.



You can reach our medical devices consultants at:

Phone +49 331 62083-0
info@miethke.com

REQUIREMENTS OF THE MDD (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the whereabouts of medical devices used in humans, especially for implants. The individual identification number of the implanted valve should therefore be recorded in the patient's medical records and patient passport to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website (<https://www.miethke.com/en/products/downloads>).

GENERAL INFORMATION






Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>Titanium Connectors</i>
Intended use	Treatment of hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation with its external dimensions:	
normal:	
	
step down:	
	

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany