Burrhole Reservoir

SPRUNG RESERVOIR

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 www.miethke.com

US This Instructions for Use is NOT intended for United State users.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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

pädiatrisches Bohrlochport Set


TABLE OF CONTENTS

INDICATION	8
TECHNICAL DESCRIPTION	8
IMPLANTATION	8
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	9
PRECAUTIONARY MEASURES AND CONTRAINDICATIONS	10
STERILISATION	10
ADVERSE REACTIONS AND INTERACTIONS	10
FUNCTIONAL SAFETY	10
REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)	10
MEDICAL DEVICES CONSULTANTS	10
GENERAL INFORMATION	11
VARIATIONS	12

INDICATION

Accessories for the treatment of hydrocephalus

TECHNICAL DESCRIPTION

One possible treatment in cases of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain to an appropriate region of the body. Our recommendation is to use a reservoir in the shunt system is to drain to the abdominal cavity in accordance to the indication of use. This shunt system consists of a series of catheters and a differential pressure valve. In addition, the option of positioning the Burrhole Reservoir in the burrhole of the cranium also exists. It then becomes possible to measure intraventricular pressure, inject drugs and check the status of the valve. Its solid titanium base is highly puncture-resistant. A puncture of the Burrhole Reservoir or the *SPRUNG RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0.9 mm. 30 times of punctures are able without any restrictions. The volume per pump procedure is approx. 0.08 ml for the Burrhole Reservoir. With the *SPRUNG RESERVOIR* it is possible to flush CSF towards the valve due to a non-return valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The volume per pump procedure is approx. 0.2 ml for the *SPRUNG RESERVOIR*. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

WARNING

Frequent pumping of the *SPRUNG RESERVOIR* can lead to an excessive drainage and thus to unphysiological pressure situations. The patient should be informed about the risk.

Nothing may be passed through the Burrhole Reservoir. A puncture should be performed as perpendicular to the reservoir surface as possible with a cannula whose maximum size is G21. The Burrhole Reservoir is available in various product configurations (individually or as a set) and as a Burrhole Port. The Burrhole Port doesn't have a lateral connector for a distal catheter. A pediatric Burrhole Reservoir and pediatric Burrhole Port adapted to the size of children's skulls are also available.

IMPLANTATION

The Ventricular Catheter is implanted with the aid of a introducing stylet. After being connected to the Burrhole Reservoir the catheter is fixed with a ligature. The reservoirs are intended for use with catheters with an inner diameter of approx. 1.2 mm and an outer diameter of approx. 2.5 mm. When implanting a Burrhole Reservoir without a distal catheter, care must be taken to ensure that the drainage catheter is also fixed in place with the aid of a ligature.

Burrhole Reservoir

It is recommended that a Burrhole with a diameter of 10 mm be performed. Owing to the fact that the ventricular catheter practically disappears in the Burrhole Reservoir, the requisite ligature can be made in a very secure fashion (see fig. 1).


Fig. 1: Burrhole Reservoir

a) additional security

Pediatric Burrhole Reservoir

Unlike the standard Burrhole Reservoir, the pediatric Burrhole Reservoir does not have an additional reduction at the bottom of the reservoir. The outer diameter is much smaller (see variants).

SPRUNG RESERVOIR

As additional new feature of the *SPRUNG RESERVOIR* CSF can be flushed towards the valve due to a non-return valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR* (see fig. 2).


Fig. 2: *SPRUNG RESERVOIR*
a) nonreturn valve

McLANAHAN RESERVOIR

The *McLANAHAN RESERVOIR* functions similar to the *SPRUNG RESERVOIR*: By pressing down the membrane CSF-flow is produced in distal direction (see fig 3).


Fig. 3: *McLANAHAN RESERVOIR*: function similar to the *SPRUNG RESERVOIR*
a) springback-leverage

In addition to the above function, it is possible to inject specific drugs or fluids into either the ventricles or in a distal direction into the shunt system. To inject into the ventricles: A cannula is used to puncture the reservoir. The cannula is then used to depress the floor of the reservoir. By depressing the floor of the reservoir completely, the distal outlet is closed and the proximal outlet is opened, thus allowing drainage into the ventricles (fig 4).


Fig. 4: *McLANAHAN RESERVOIR* closed distally and open proximally
a) distal closed

When a cannula is used without depressing the reservoir floor, the proximal outlet is blocked and the distal outlet is opened. During the transition between the above mentioned conditions, both outlets are open.

NOTE

The patency of the reservoir should be tested prior to implantation. In case of uncertainty, distilled water should be aspirated carefully from the distal connector and by immersing the proximal part in sterile liquid.

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, as a validated decontamination process will compromise the functionality of the valve.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

Like the distal catheter (when applied), the burrhole reservoirs consist of the non-magnetic substances titanium, sapphire and silicone. Consequently, the reservoir and catheters will not interfere with either nuclear magnetic resonance or computer tomographic procedures. The burrhole reservoirs are MRI Conditional. The provided catheters are MRI Safe. Deflectors and connectors are MRI Conditional.

PRECAUTIONARY MEASURES AND CONTRAINDICATIONS

Patients must be carefully monitored after implantation. Reddening of skin or tightness in the area of the drained tissue may be indications of infections at the shunt system. Symptoms such as headache, dizziness, confusion or vomiting often occur in conjunction with shunt dysfunction. These symptoms and a leakage within the shunt system require the immediate replacement of the affected shunt component or the entire shunt system. The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e.g. meningitis, ventriculitis, peritonitis, bacteraemia, septicaemia) in the region affected by the implantation.

STERILISATION

The products are sterilised with steam under strictly controlled conditions. The double wrapping in sterile bags ensures sterility for a five-year period. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used in any circumstances. No guarantee can be given for the functional safety and reliability of resterilised products.

ADVERSE REACTIONS AND INTERACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases noise development. Violent shocks from the outside (accident, fall) may put the integrity of the shunt system at risk.

WARNING

The *Burrhole Reservoir / SPRUNG RESERVOIR* must not be used in conjunction with hydrostatic valves as this may result in increased ventricular pressure outside of the physiological range. In case of doubt, please contact the medical products consultants at Christoph Miethke GmbH & Co. KG.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the fate of medical devices used in humans, especially for implants. The individual implant identification number should therefore be recorded in the patient's medical records and patient passport to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website (<https://www.miethke.com/en/products/downloads/>).



MEDICAL DEVICES CONSULTANTS

In compliance with the European Directive concerning medical devices (Directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical devices consultants as contacts for all product-related questions.

You can contact our medical devices consultants via:

Tel. +49 331 62083-0
info@miethke.com


GENERAL INFORMATION

Hersteller	Christoph Miethke GmbH & Co. KG
Product name	see Variations
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the shunt with its external dimensions:	
standard:	
	
pediatric:	
	

VARIATIONS

All reservoirs are available with distal connector or distal catheter.


Pediatric Burrhole Reservoir

Pediatric Burrhole Port


Burrhole Reservoir


Burrhole Port

SPRUNG RESERVOIR

McLANAHAN RESERVOIR

ETOU RESERVOIR

connector can be adjusted to various positions

The reservoirs are also available as a set, including a Ventricular Catheter.

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

Produttore	Christoph Miethke GmbH & Co. KG
Schizzo della valvola con misure esterne:	
standard:	
 <p>Technical drawing of a standard valve. The drawing shows a dome-shaped valve with a central stem. Dimension lines indicate a diameter of $\varnothing 20\text{ mm}$ and a height of 5.65 mm.</p>	
pediatrico:	
 <p>Technical drawing of a pediatric valve. The drawing shows a smaller dome-shaped valve with a central stem. Dimension lines indicate a diameter of $\varnothing 14\text{ mm}$ and a height of 4.8 mm.</p>	

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany