Bohrlochumlenker/ Deflector

Ⓓ Gebrauchsanweisung | ⒸB Instructions for use | Ⓕ Mode d'emploi
Ⓔ Instrucciones de manejo | ⓘ Istruzioni per l'uso

Ⓔ USA This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000.

A paper copy will be provided to you upon request at no additional cost.

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CAUTION

Federal law restricts this device to sale by or on order of a physician!

INDICATION

The *deflector* is part of a shuntsystem for the treatment of hydrocephalus. It serves as a right-angular guide for the ventricular catheter coming out of the cranial borehole of a hydrocephalus patient.

DESCRIPTION

One possible treatment in cases of hydrocephalus is the implantation of a shuntsystem that achieves drainage of cerebrospinal fluid away from the ventricles of the brain into a suitable area of the body (in most cases the intestinal cavity).

This shuntsystem consists of a valve, a proximal and distal catheter and a reservoir. If the *borehole reservoir* is placed into the borehole, there is no need for a *deflector*. If the reservoir is a *prechamber*, the reservoir is in line with the valve and the *deflector* is part of the shuntsystem. The *deflector* is positioned in the borehole in the cranium and it is possible to select the length of the catheter to be inserted into the cranium prior to implantation. The *deflector* provides a turn (at a 90 degree angle) for the ventricular catheter in the borehole, which it seals off completely. Additionally to the standard deflector a small pediatric deflector is also available. Due to the eyelets on two sides of the standard deflector (and eyelets on three sides for the pediatric deflector), it is now possible to fix the deflector with ligatures to the tissue.

Fig. 1 Deflector with ventricular catheter

IMPLANTATION

The ventricular catheter is implanted with the aid of a mandrin. The length of catheter to be implanted in the cranium is selected in advance through the positioning of the *deflector* on the catheter. The *deflector* is designed and intended for use with a catheter with an inner diameter of ca. 1.2 mm and an outer diameter of ca. 2.5 mm. The lower crank of the standard variant has 9,5 mm in diameter (see fig. 2). Therefore a borehole with a diameter of 10 mm or more should be prepared.


Fig. 2 Deflector

REIMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, as a validated decontamination process will compromise the functionality of the valve.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

The *deflectors* and the catheters consist of the non-magnetic substances titanium and silicone. Consequently, the *deflector* and catheters will not interfere with either nuclear magnetic resonance or computer tomographic procedures. The *deflectors* are MRI Conditional (ASTM F2503-08). The provided catheters are MRI Safe.

SAFETY MEASURES

Following implantation, the patient's condition should be carefully and thoroughly monitored. Inflamed skin and tension in the area of the drainage tissues could possibly be a sign of

infection in the shunt system. Symptoms such as headache, dizziness, disorientation, or vomiting often occur in cases of shunt dysfunction. In the event such symptoms do occur, or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay. Particular care must be taken not to confuse the *deflector* with the *borehole reservoir*. The *borehole reservoir* is manufactured by Christoph Miethke GmbH & Co.KG. It allows both measurement of pressure and injection of medication to be carried out by means of a puncture.

STERILIZATION

All products are carefully and thoroughly steam sterilized. Owing to the fact that the product is packaged in two layers of sterile packaging, five years of sterility is guaranteed. The expiration date is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

RESTERILIZATION

The safe and reliable performance of resterilized products cannot be guaranteed.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted product (Charge/LOT/SN) are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
14469 Potsdam · Germany
Phone: +49(0) 7000 MIETHKE or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap-Platz
78532 Tuttlingen · Germany
Phone: +49 (0) 7461 95 0
Fax: +49 (0) 7461 95 26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>Deflector</i>
Intended use	Treatment of hydrocephalus
Intended for one-time use only (disposable)	
Store in a clean and dry place	
Schematic representation of the <i>deflector</i> with its external dimensions:	
<i>Deflector</i>	
<i>Pediatric deflector</i>	

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Deutschland