## *Burrhole Deflector*

ⒹE Gebrauchsanweisung | ⒺN Instructions for use | ⒻR Mode d'emploi | ⒺS Instrucciones de manejo

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 [www.miethke.com](http://www.miethke.com)

ⒺS This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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**TABLE OF CONTENTS**

INDICATION	8
DESCRIPTION	8
IMPLANTATION	8
REIMPLANTATION	8
INTERACTION WITH PRODUCTS FROM OTHER MANUFACTURERS	8
FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	8
SAFETY MEASURES AND CONTRAINDICATIONS	9
ADVERSE REACTIONS	9
STERILISATION	9
RESTERILISATION	9
REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)	9
MEDICAL DEVICES CONSULTANTS	9
GENERAL INFORMATION	10

## INDICATION

The *Burrhole Deflector* is part of a shuntsystem for the treatment of hydrocephalus. It serves as a right-angular guide for the *Ventricular Catheter* coming out of the cranial burrhole of a hydrocephalus patient.

## DESCRIPTION

One possible treatment in cases of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid away from the ventricles of the brain into a suitable area of the body (in most cases the intestinal cavity).

This shunt system consists of a valve, a proximal and distal catheter and a reservoir. If the *Burrhole Reservoir* is placed into the burrhole, there is no need for a *Burrhole Deflector*. If the reservoir is a prechamber, the reservoir is in line with the valve and the *Burrhole Deflector* is part of the shunt system. The *Burrhole Deflector* is positioned in the burrhole in the cranium and it is possible to select the length of the catheter to be inserted into the cranium prior to implantation. The *Burrhole Deflector* provides a turn (at a 90 degree angle) for the *Ventricular Catheter* in the burrhole, which it seals off completely. Additionally to the standard *Burrhole Deflector* a small *Pediatric Burrhole Deflector* is also available. Due to the eyelets on two sides of the standard *Burrhole Deflector* (and eyelets on three sides for the *Pediatric Burrhole Deflector*), it is now possible to fix the *Burrhole Deflector* with ligatures to the tissue.

Fig. 1: *Burrhole Deflector* with *Ventricular Catheter*

## IMPLANTATION

The *Ventricular Catheter* is implanted with the aid of a mandrin. The length of catheter to be implanted in the cranium is selected in advance through the positioning of the *Burrhole Deflector* on the catheter.

The *Burrhole Deflector* (including paediatric version) is designed and intended for use with a catheter with an inner diameter of approx. 1.2 mm and an outer diameter of approx. 2.5 mm. The lower crank of the standard variant has 9.5 mm in diameter (see fig. 2). Therefore a burrhole with a diameter of 10 mm or more should be prepared.


Fig. 2: *Burrhole Deflector*

For the paediatric version, the burrhole opening should be selected correspondingly smaller (fig. 3).


Fig. 3: *Pediatric Burrhole Deflector*

## REIMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, as a validated decontamination process will compromise the functionality of the valve.

## INTERACTION WITH PRODUCTS FROM OTHER MANUFACTURERS

When using the *Burrhole Deflector*, the use of catheters of the Christoph Miethke GmbH & Co. KG is recommended.

## FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. How-

ever, we cannot guarantee that these medical devices will require replacement for medical or technical reasons. These medical devices have to be stored in a clean and dry environment at all times.

The *Burrhole Deflectors* and the catheters consist of the non-magnetic substances titanium and silicone. Consequently, the *Burrhole Deflector* and catheters will not interfere with either nuclear magnetic resonance or computer tomographic procedures. The *Burrhole Deflectors* are MR Conditional. The provided catheters are MR Safe.

The conditions for MR compatibility of the products can be found on our website: <https://www.miethke.com/downloads/>

### **SAFETY MEASURES AND CONTRAINDICATIONS**

Following implantation, the patient's condition should be carefully and thoroughly monitored. Inflamed skin and tension in the area of the drainage tissues could possibly be a sign of infection in the shunt system. Symptoms such as headache, dizziness, disorientation, or vomiting often occur in cases of shunt dysfunction. In the event such symptoms do occur, or if there is any leakage in the shunt system, either individual shunt system components or the entire shunt system must be replaced without delay.

Particular care must be taken not to confuse the *Burrhole Deflector* with the *Burrhole Reservoir*. The *Burrhole Reservoir* is manufactured by Christoph Miethke GmbH & Co. KG. It allows both measurement of pressure and injection of medication to be carried out by means of a puncture.

The implantation of medical products is contraindicated if an infection (e. g. meningitis, ventriculitis, peritonitis, bacteremia, septicemia) or the suspicion of an infection is diagnosed in a region of the patient's body where the shunt system will be implanted.

### **ADVERSE REACTIONS**

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in very rare cases noise development. Violent

shocks for the outside (accident, fall etc.) may put the integrity of the shunt system at risk.

### **STERILISATION**

The products are sterilised with steam under strictly controlled conditions. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used in any circumstances. No guarantee can be given for the functional safety and reliability of resterilised products.

### **RESTERILISATION**

The safe and reliable performance of resterilized products cannot be guaranteed.

### **REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)**

The Medical Device Directive requires comprehensive documentation of the fate of medical devices used in humans, especially for implants. The individual implant identification number should therefore be recorded in the patient's medical records and patient card to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website: <https://www.miethke.com/downloads/>

### **MEDICAL DEVICES CONSULTANTS**

In compliance with the European Directive concerning medical devices (Directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical devices consultants as contacts for all product-related questions.

You can contact our medical devices consultants via:

**Tel. +49 331 62083-0**  
**info@miethke.com**

## GENERAL INFORMATION

<b>Manufacturer</b>	<b>Christoph Miethke GmbH &amp; Co. KG</b>
Product name	<i>Burrhole Deflector</i>
Intended use	Treatment of hydrocephalus

Intended for one-time use only (disposable)

Store in a clean and dry place

Schematic representation of the *Burrhole Deflector* with its external dimensions:


Fig. 4: Burrhole Deflector


Fig. 5: Pediatric Burrhole Deflector

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany