Prechamber **CONTROL RESERVOIR®**

ⓓ Gebrauchsanweisung | Ⓜ Instructions for use | ⓕ Mode d'emploi
ⓔ Instrucciones de manejo | ⓘ Istruzioni per l'uso

Ⓤ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
Maßstab 1:1

INDICATION

The *prechamber* and the *CONTROL RESERVOIR* are used in cases of hydrocephalus to remove cerebrospinal fluid and monitor pressure.

DESCRIPTION

One possible treatment of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid from the ventricles of the brain to an appropriate region of the body (in most instances to the abdominal cavity). This shunt system consists of a series of catheters and a differential pressure valve. In addition, the option a *prechamber* also exists. It then becomes possible to measure intraventricular pressure, inject medication and check the status of the valve. With the *CONTROL RESERVOIR* it is possible to flush CSF towards the valve due to a nonreturn valve in the proximal inlet connector of the reservoir. By this mechanism a flow in the direction of the proximal catheter is avoided during the pumping procedure. The volume per pump procedure is approx. 0.25 ml. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

Caution: Frequent pumping of the *CONTROL RESERVOIR* can lead to an excessive drainage and thus to unphysiological pressure situations. The patient should be informed about the risk.


Fig. 1: CONTROL RESERVOIR

Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the control reservoir should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0,9 mm. 30 times of punctures are possible without any restrictions.

The reservoirs are available in various product configurations (individually or as a set) and as a CSF Port. The CSF Port has only one lateral connector for the ventricular catheter.


Fig. 2: CSF Port

A *pediatric prechamber* is available also, it is adapted to the size of children's skulls.

IMPLANTATION

When implanting a reservoir, care must be taken to ensure catheters are secured by means of a ligature, if this has not already been done by the catheter manufacturer.

The reservoirs are intended for use with catheters with an inner diameter of approx. 1.2 mm and an outer diameter of approx. 2.5 mm. When implanting a reservoir set, the neurosurgeon should determine how the procedure is to be performed in accordance with experience. The ventricle catheter is implanted with the aid of a introducing stylet, after which it can be brought through the deflector into the 90°-angle. During this process (and prior to implantation of the ventricular catheter), the deflector can be precisely adjusted, and then used as a stop. The reservoir is used in conjunction, and is connected to a complete shunt system. A final check of the complete shunt system is recommended to ensure that it is not obstructed in any way.

REIMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

Like other shunt components of the Christoph Miethke GmbH & Co. KG (when applied), the *prechamber* and the *CONTROL RESERVOIR* consist of the nonmagnetic substances titanium, sapphire and silicone. Consequently, the reservoirs and catheters will not interfere with either nuclear magnetic resonance (MRI) or computer tomographic (CT) procedures. The *prechamber* and the *CONTROL RESERVOIR* are MR Conditional (ASTM F2503-13). The provided catheters are MRI Safe. All components are visible via X-ray.

SAFETY MEASURES

Following implantation of a shunt system, the patient's condition should be carefully and thoroughly monitored. Inflammation and tension around the drainage tissues could possibly be a sign of infection in the shunt system.

Symptoms such as headache, dizziness, disorientation or vomiting often occur in cases of shunt dysfunction.

In the event such symptoms occur or if there is any leakage in the shunt system, either the individual shunt system components or the entire shunt system must be replaced without delay.

STERILIZATION

All products are carefully and thoroughly sterilised. Owing to the fact that the products are packaged in two layers of sterile packaging, five years of sterility are guaranteed. The expiration date for each item is indicated on the package. If the packaging is damaged in any way, the product should not be used under any circumstances.

RESTERILIZATION

The functional safety and reliability of re-sterilised products cannot be guaranteed.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
14469 Potsdam · Germany
Phone: +49(0) 7000 MIETHKE or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com



Please address any enquiries to:
AESCULAP AG
Am Aesculap-Platz
78532 Tuttlingen · Germany
Phone: +49 (0) 7461 95 0
Fax: +49 (0) 7461 95 26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	see variations
Intended use	Treatment of hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the <i>prechamber</i> with its external dimensions:	
standard:	
pediatric:	

VARIATIONS

Prechamber

Prechamber with distal catheter

Pediatric prechamber

CONTROL RESERVOIR

CONTROL RESERVOIR with distal catheter

Prechamber for L-P derivation

CSF Port


CSF Port with ventricular catheter

Prechamber Set and CONTROL RESERVOIR Set

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany