*Prechamber*

*CONTROL RESERVOIR*

DE Gebrauchsanweisung | EN Instructions for use | FR Mode d'emploi | ES Instrucciones de manejo

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 [www.miethke.com](http://www.miethke.com)

US This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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**Tel. +49 331 62983-0**  
**info@miethke.com**

**TABLE OF CONTENTS**

INDICATION	8
TECHNICAL DESCRIPTION	8
IMPLANTATION	9
REIMPLANTATIONS	9
PRECAUTIONARY MEASURES AND CONTRAINDICATIONS	9
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	9
FUNCTIONAL SAFETY	9
ADVERSE REACTION	9
STERILISATION	9
REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)	9
MEDICAL DEVICE CONSULTANTS	10

## INDICATION

Accessories for the treatment of hydrocephalus

## TECHNICAL DESCRIPTION

One possible treatment of hydrocephalus is the implantation of a shunt system that achieves drainage of cerebrospinal fluid (CSF) from the ventricles of the brain to an appropriate region of the body (in most instances to the abdominal cavity). This shunt system consists of a series of catheters and a differential pressure valve. In addition, a reservoir (e.g. *Prechamber*) can be positioned on top of the skull. It allows for measurement of the intraventricular pressure, the injection of medication and the checking of the valve's status.


Fig. 1: Prechamber

With the *CONTROL RESERVOIR* it is possible to flush CSF towards the valve due to a non-return valve in the proximal inlet connector of the reservoir (Fig. 2). By this mechanism a flow in the direction of the proximal catheter is avoided during the pumping procedure. The volume per pump procedure is approx. 0.25 ml. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.


Fig. 2: CONTROL RESERVOIR with non-return valve  
a) non-return valve

A *Pediatric Prechamber* or a *Pediatric CONTROL RESERVOIR* is also available (Fig. 3). It is adapted to the size of children's skulls.


Fig. 3: Pediatric Prechamber,  
Pediatric CONTROL RESERVOIR

Its solid titanium base is highly puncture-resistant. A puncture of the *Prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0.9 mm. 30 punctures are possible without any restrictions. The reservoirs are available in various product configurations (individually or as a set) and as a *CSF-PORT* (Fig. 4). The *CSF-PORT* has only one lateral connector for the ventricular catheter.


Fig. 4: CSF-PORT  
b) proximal connector

## WARNING

**Frequent pumping can lead to an excessive drainage and thus to unphysiological pressure situations. The patient should be informed about this risk.**

## IMPLANTATION

When implanting a reservoir, care must be taken to ensure catheters are secured by means of a ligature, if this has not already been done by the catheter manufacturer. The reservoirs are intended for use with catheters with an inner diameter of approx. 1.2 mm and an outer diameter of approx. 2.5 mm. When implanting a reservoir set, the neurosurgeon should determine how the procedure is to be performed in accordance with experience. The ventricular catheter is implanted with the aid of an introducing stylet, after which it can be brought through the deflector into a 90° angle. During this process (and prior to implantation of the ventricular catheter), the deflector can be precisely adjusted and then used as a stop. The reservoir is used in conjunction, and is connected to a complete shunt system. A final check of the complete shunt system is recommended to ensure that it is not obstructed in any way.

## REIMPLANTATIONS

Under no circumstances should products that had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

## PRECAUTIONARY MEASURES AND CONTRAINDICATIONS

Patients must be carefully monitored after implantation. Reddening of skin or tightness in the area of the drained tissue may be indications of infections at the shunt system. Symptoms such as headache, dizziness, confusion or vomiting often occur in conjunction with shunt dysfunction. These symptoms and a leakage within the shunt system require the immediate replacement of the affected shunt component or the entire shunt system. The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e.g. meningitis, ventriculitis, peritonitis, bacteriaemia, septicaemia) in the region affected by the implantation.

## COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

Nuclear magnetic resonance (MRI) examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the patient. The *Prechamber* and the *CONTROL RESERVOIR* are MR Conditional. The supplied catheters are MR Safe. Reservoirs, deflectors and connectors are MR Conditional.

The conditions for MR compatibility of the products can be found on our website:

<https://www.miethke.com/downloads/>

## FUNCTIONAL SAFETY

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons. These medical devices are able to resist positive and negative pressures up to 100 cmH<sub>2</sub>O during and after implantation. These medical devices have to be stored in a clean and dry environment at all times.

## ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in very rare cases, noise development. Due to violent shocks from the outside (accident, fall etc.) the integrity of the shunt may be endangered.

## STERILISATION

The products are sterilised with steam under strictly controlled conditions. The expiry date is printed on the wrapping of each individual product. If the packaging is damaged, the product must not be used in any circumstances. No guarantee can be given for the functional safety and reliability of resterilised products.

## REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the fate of medical devices used in humans, especially for

implants. The individual implant identification number should therefore be recorded in the patient's medical records and patient card to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website:  
<https://www.miethke.com/downloads/>

### **MEDICAL DEVICE CONSULTANTS**

In compliance with the requirements of the European law MDD 93/42/ECC, Christoph

Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany