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
# GAV<sup>®</sup> 2.0 / GAV<sup>®</sup> 2.0 LP

## RECOMMENDED PRESSURE SETTINGS

Patient	Selection of pressure levels
Newborn 	5 / 25
Children from 6 months 	10 / 25
Children from 3 years 	10 / 30
Adults < 1,60 m  > 1,80 m	5 / 30 5 / 25 5 / 35
Adults from 65 years < 1,60 m  > 1,80 m	5 / 25 5 / 20 5 / 30

The GAV 2.0 comprises a differential pressure unit and a gravitational unit.

The opening pressure of the GAV 2.0 in the horizontal position is determined by the differential pressure unit.

The opening pressure, which works in the vertical position is the **sum of both**: differential pressure unit and gravitational unit.

Values apply for mobile patients. In case of less mobile patients or a high BMI the gravitational unit should be set lower than it is recommended here.

The Miethke App demonstrates the functionality of GAV 2.0 and the effect of patient posture on valve opening pressure.

Now available in the App and Google Play Store.

Our Strong Partner in Neurosurgery:

**B | BRAUN**  
SHARING EXPERTISE

www.bbraun.com

All indicated pressure levels are in cmH<sub>2</sub>O. This is a non-binding recommendation for the attending physician. According to his diagnosis, the physician decides each case independently, without instructions and individually. The stated values consider the current scientific knowledge up to June 2018.

www.miethke.com