paediSHUNTASSISTANT SHUNTASSISTANT

DE Gebrauchsanweisung | **EN** Instructions for use | **FR** Mode d'emploi
ES Instrucciones de manejo | **IT** Istruzioni per l'uso

 www.miethke.com

US This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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CONTENT

| | |
|---|----|
| INDICATION | 13 |
| TECHNICAL DESCRIPTION | 13 |
| PHYSICS BACKGROUND | 14 |
| FUNCTION OF THE VALVE | 14 |
| SELECTING THE APPROPRIATE VALVE | 15 |
| READING THE PRESSURE SETTING FROM AN X-RAY IMAGE | 16 |
| TUBE SYSTEMS | 16 |
| SURGICAL PROCEDURE | 16 |
| PRESSURE-FLOW CHARACTERISTICS | 17 |
| PREOPERATIVE VALVE TEST | 18 |
| INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS | 18 |
| RE-IMPLANTATION | 18 |
| SAFETY MEASURES | 18 |
| COMPATIBILITY WITH DIAGNOSTIC PROCEDURES | 18 |
| POSTOPERATIVE VALVE TEST | 18 |
| FUNCTIONAL SAFETY | 18 |
| ADVERSE REACTION | 19 |
| STERILISATION | 19 |
| RESTERILISATION | 19 |
| REQUIREMENTS OF THE MDD 93/42/EEC | 19 |
| NOTE ON THE INSTRUCTIONS FOR USE | 19 |
| MEDICAL DEVICES CONSULTANT | 19 |
| GENERAL INFORMATION | 19 |
| VARIANTS | 20 |

INDICATION

The *paediSHUNTASSISTANT/SHUNTASSISTANT* is intended to shunt cerebrospinal fluid (CSF) from the lateral ventricles of the brain into the peritoneum.

TECHNICAL DESCRIPTION

The *paediSHUNTASSISTANT/SHUNTASSISTANT* is intended for use as a supplementary valve, in order to allow for monitoring and control of the siphon effect in conjunction with either a conventional or adjustable differential pressure valve. Particular attention was paid to develop a device that can forestall negative effects such as changes in subcutaneous pressure (e.g. owing to tissue capsules or an unfavorable patient position).

The *paediSHUNTASSISTANT/SHUNTASSISTANT* utilises the same mechanism as the *SHUNTASSISTANT*. However, it is considerably smaller and its functions are intended to meet the special needs of pediatric hydrocephalus.

The *paediSHUNTASSISTANT/SHUNTASSISTANT* consists of a robust titanium casing (1), a weight-ball (2), as well as a small sealing ball (3) that rests on the ball seat (4). The cerebrospinal fluid flows in through the inlet port (5) and flows out by way of the outlet port (6).

- 1 titanium casing
- 2 tantalum ball
- 3 sapphire ball
- 4 ball seat
- 5 inlet port
- 6 outlet port


Fig. 1: Schematic cross section of the paedISHUNTASSISTANT/SHUNTASSISTANT

PHYSICS BACKGROUND

With the body in a horizontal position, the intraventricular pressure in healthy individuals is positive. In order to regulate this pressure through a differential-pressure valve, only the pressure in the abdominal cavity needs to be taken into account.

In standing position, the ventricular pressure should be slightly negative. To adjust this pressure through valve drainage, the *paediSHUNTASSISTANT/SHUNTASSISTANT* is implanted in addition to a differential-pressure valve. Only in this way the hydrostatic pressure minus the pressure in the abdominal cavity can be compensated for.

FUNCTION OF THE VALVE

With the *paediSHUNTASSISTANT/SHUNTASSISTANT* in horizontal position, the valve sealing ball (3) and the weight-ball (2) are not in the valve seat. Hence, in this configuration the *paediSHUNTASSISTANT/SHUNTASSISTANT* is open, and the force brought to bear by the weight of the ball does not counteract the fluid pressure. The ventricular pressure is now controlled exclusively by the differential-pressure valve installed.


Fig. 3: horizontal position of the *paediSHUNTASSISTANT/SHUNTASSISTANT*

$$\text{horizontal: } IVP = P_L + P_B$$

$$\text{vertical: } IVP = P_H - P_S - P_B$$

| | |
|-------|--|
| IVP | <i>Intraventricular pressure</i> |
| P_L | <i>Opening pressure in horizontal position (adjustable DP-unit only)</i> |
| P_S | <i>Opening pressure in vertical position (adjustable DP-unit + gravitational unit)</i> |
| P_B | <i>Pressure in the abdominal cavity</i> |
| P_H | <i>Hydrostatic pressure</i> |


Fig. 2: Calculating the intraventricular pressure

In the vertical position, the *paediSHUNTASSISTANT/SHUNTASSISTANT* and the differential pressure valve work together. If the pressure at the ball seat exceeds the force exerted by the weight of the tantalum ball, the valve sealing ball and weight ball are raised, thereby allowing the lateral flow towards distal of fluid through the outlet port. This prevents any further increase in pressure. In order to achieve optimal pressure regulation for each individual patient, the *SHUNTASSISTANT* is available with five different pressure ratings, in addition to the *paediSHUNTASSISTANT*, which was developed especially for pediatric hydrocephalus.

It is recommended that the *paediSHUNTASSISTANT/SHUNTASSISTANT* is selected for which the sum total of the opening pressure of the *paediSHUNTASSISTANT/SHUNTASSISTANT* and of the differential pressure valve – measured in terms of water column height – equals approximately 5 cm less than the distance from the peritoneum (diaphragm) to the third ventricle (Foramen Monroi) (see recommendation of pressure settings at <https://www.miethke.com/downloads/>).

Suggested procedure (example):

- 1) Measure the distance between the peritoneum (diaphragm) and the third ventricle (Foramen Monroi) (e.g. 42 cm).
- 2) Select the differential pressure valve to be implanted as a supplement to the *paediSHUNTASSISTANT/SHUNTASSISTANT*, and take a reading of the valve pressure rating (e.g. 13 cmH₂O). If that value is available in mmHg only, the pressure has to be converted into cmH₂O.
1 mmHg = 1.36 cmH₂O (e.g. 9.6 mmHg = 13 cmH₂O)
- 3) Subtract the pressure of the differential pressure valve from the distance measured between the peritoneum and the third ventricle (e.g. 42 cmH₂O - 13 cmH₂O = 29 cmH₂O)
- 4) Subtract 5 cmH₂O from the computed figure (e.g. 29 cmH₂O - 5 cmH₂O = 24 cmH₂O).
- 5) Choose the *paediSHUNTASSISTANT/SHUNTASSISTANT* specified for the pressure rating closest to this value. For the example given here, the 25 cmH₂O pressure rating would be selected. The ventricular pressure in standing position would then be:

$25 \text{ cmH}_2\text{O} + 13 \text{ cmH}_2\text{O} - 42 \text{ cmH}_2\text{O} = -4 \text{ cmH}_2\text{O}$ thus inside the physiological range. Conflicting statements are found in the literature as to what constitutes an ideal ventricular pressure in hydrocephalus patients in a standing position. In the example above, the assumption was made that a slightly negative pressure comes closest to fulfilling the patient's needs. A definite standard value cannot be given here, however. Applicable values will vary according to patient and anamnesis.


Fig. 4: vertical position of the *paediSHUNTASSISTANT/SHUNTASSISTANT* a) closed b) open

SELECTING THE APPROPRIATE VALVE

The choice of an appropriate *paediSHUNTASSISTANT/SHUNTASSISTANT* is determined by the ventricular pressure that needs to be achieved for the standing position. The sum total of the opening pressures of the differential pressure valve and the *paediSHUNTASSISTANT/SHUNTASSISTANT* used determines the ventricular pressure for the patient in standing position.


READING THE PRESSURE SETTING FROM AN X-RAY IMAGE

Each *paediSHUNTASSISTANT/SHUNTASSISTANT* is calibrated under conditions of strict quality control. The *paediSHUNTASSISTANT* is available with a pressure rating of 10 cmH₂O; the *SHUNTASSISTANT* can be supplied with the following pressure ratings:


| Pressure rating (cmH ₂ O) | Ring code |
|--------------------------------------|---|
| 10 paediSHUNTASSISTANT |  |
| 15 SHUNTASSISTANT |  |
| 20 SHUNTASSISTANT |  |
| 25 SHUNTASSISTANT |  |
| 30 SHUNTASSISTANT |  |
| 35 SHUNTASSISTANT |  |

Post-operatively, the pressure rating selected can be read from x-ray images.

The functioning of the *paediSHUNTASSISTANT/ SHUNTASSISTANT* too, can be checked through postoperative X-rays. Radiography performed with the patient in a horizontal position should reveal the tantalum ball to be situated at a slightly greater distance from the ball seat than in radiography carried out with the patient standing up.

a)

b)

Fig. 5:

a) X-ray image of the *paediSHUNTASSISTANT* pressure setting 10 cmH₂O

b) X-ray image of the *SHUNTASSISTANT* pressure setting 20 cmH₂O

TUBE SYSTEMS

The *paediSHUNTASSISTANT/SHUNTASSISTANT* is designed and built to ensure optimal ventricular pressure when applied in conjunction with a supplementary conventional or adjustable differential pressure valve. For connections, catheters with an internal diameter of 1.2 mm and an external diameter of 2.5 mm should be used preferably. In any case, the catheter must be carefully fixed to the valve connector by means of a ligature. Kinks in the catheters should be avoided.

The provided catheters have virtually no effect on the pressure-flow characteristics.

SURGICAL PROCEDURE

Depending on the anatomic situation, it is recommended to implant the *SHUNTASSISTANT* in the thoracic region. Because of growth issues, the *paediSHUNTASSISTANT* should be implanted retroauricularly and subcutaneously. In any case, it has to be taken care that the *paediSHUNTASSISTANT/SHUNTASSISTANT* is fixated in parallel to the body axis. Where a *paediSHUNTASSISTANT/SHUNTASSISTANT* is used without integrated catheters, the catheters (proximal and distal) have to be fastened with a ligature. The *paediSHUNTASSISTANT/SHUNTASSISTANT* should be positioned in such a way that the skin incision does not lie directly over the valve. The peritoneal catheter can be used to pull the implant to the desired implantation height. In doing this, the catheter should not be stretched excessively. Furthermore, care must be taken that no parts (e.g. tissue fragments) get into the *paediSHUNTASSISTANT/SHUNTASSISTANT* during implantation. An arrow on the housing indicates the direction of flow. The position of the ventricular catheter should be inspected again by post-operative CT or MR imaging. The implantation of the conventional or (adjustable) differential-pressure valve is carried out according to the instructions given by the manufacturer of the valve selected.







CAUTION

The catheters should only be blocked with a sheathed clamp and not directly behind the valve as they might be damaged otherwise.

PRESSURE-FLOW CHARACTERISTICS

The diagrams below show the pressure-flow characteristics for the pressure ratings in which the paedISHUNTASSISTANT/SHUNTASSISTANT is available.

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher

PREOPERATIVE VALVE TEST


Fig. 6: Patency test

The paedISHUNTASSISTANT/SHUNTASSISTANT can be filled most gently by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see fig. 6).

WARNING

Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end. Contaminations in the solution used for the test can impair the product's performance.

INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS

The paedISHUNTASSISTANT/SHUNTASSISTANT should not be used under any circumstances in conjunction with hydrostatic valves, as this can bring about an unphysiologically high ventricular pressure. Hydrostatic valves allow for changes in hydrostatic pressure in the drainage system caused by changes in position. If in doubt, please contact the medical product consultants at Christoph Miethke GmbH & CO. KG.

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently re-implanted in another, as a validated decontamination process will com-

promise the functionality of the valve.

SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt.

The paedISHUNTASSISTANT/SHUNTASSISTANT is MR Conditional. All components are visible via X-ray. The provided catheters are MRI Safe. Reservoirs, deflectors and connectors are MRI Conditional.

The conditions for MR compatibility of the products can be found on our website: <https://www.miethke.com/downloads/>.

POSTOPERATIVE VALVE TEST

The paedISHUNTASSISTANT/SHUNTASSISTANT has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the unit if a shunt system with a *prechamber* or a *burrhole reservoir* is used. Valve tests can be carried out by flushing or pressure measurements.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons.

The valve and the valve system are able to

resist positive and negative pressure up to 100 cmH₂O during and after implantation.

ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

STERILISATION

The products are sterilized with steam under closely monitored conditions. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILISATION

The functional safety and reliability of resterilized products cannot be guaranteed, therefore resterilisation is not recommended.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

Translations of these instructions for use into additional languages can be found on our website: <https://www.miethke.com/downloads>

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

MEDICAL DEVICES CONSULTANT



In compliance with European directive on medical devices (directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical products consultants as contacts for all product-related questions.

You can reach our medical devices consultants at:

Phone +49 331 62083-0

info@miethke.com

GENERAL INFORMATION

| Manufacturer | Christoph Miethke GmbH & Co. KG |
|---|---------------------------------|
| Product name | SHUNTASSISTANT |
| Intended use | Treatment of hydrocephalus |
| Intended for one-time use (disposable) | |
| Store in a clean, dry place | |
| Schematic representation of the shunt with its external dimensions: | |
|  | |
| SHUNTASSISTANT | |
|  | |
| paediSHUNTASSISTANT | |

VARIANTS

The *paediSHUNTASSISTANT/SHUNTASSISTANT* is available as a single valve or with an integrated catheter.

paediSHUNTASSISTANT

SHUNTASSISTANT

SHUNTASSISTANT with distal catheter

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany